



**OSSTEM<sup>®</sup>**  
IMPLANT

# US SYSTEM CATALOG



**OSSTEM<sup>®</sup>**  
IMPLANT



# US SYSTEM CATALOG

**Osstem Implant**  
**2018-19 Comprehensive Catalog**

**Overall Planning/Editing** PR Department Design Team  
**Supervision** Implant Lab, Marketing PM  
**Production/Distribution** Marketing & Planning Team  
**Date of Publication** 2018.08  
**Publisher** Osstem Implant  
8th FL, World Meridian II, 123, Gasan digital 2-ro,  
Geumcheon-gu, Seoul, Korea  
**Phone** +82.2.2016.7000  
**Fax** +82.2.2016.7001  
www.osstem.com

<b>003</b>	INTRODUCTION
<b>012</b>	CONTENTS
<b>016</b>	US SYSTEM
<b>072</b>	USER MANUAL

### Providing cutting edge technology and superior quality

Making products that dentists want to use,  
trust, and are satisfied with :  
This is our mission at **OSSEM IMPLANT**

**We are forever grateful  
to all the dentists  
who have given  
unwavering support to  
OSSTEM IMPLANT**


Thank you for using Osstem Implant. Osstem, Korea's first implant manufacturer, has secured world-class implant competitiveness through continuous R&D investment and quality innovation. It has grown to become Asia-Pacific No.1 and World No.5 Implant Company. In addition to dental implants and treatment tools, we are leading the development of products that are essential for dentists, including dental equipment, dental materials, and dental IT, and contribute to the development of the dental industry. The comprehensive catalog of the 2018-19 product series published here shows Osstem's technology-rich products. We have focused on catalog structure so that it is convenient to browse and order products. In particular, in the case of fixtures, abutments, and surgical tools, we introduced the diameter, length, and functions in detail.

GBR products are also easy to order by type, size and capacity. In addition, the product release date and time are displayed so that customers can understand when the existing product is released and what the newly released product is. We also introduced the CAD/CAM product in terms of preparing the digital dentistry, a major trend in the dentistry. In terms of design, we also implemented high-quality images of representative products by specification. By applying representative colors for each product system, it is easy to sort by category. We hope this will help you effectively find and purchase the products you need from the dental clinic of 2018-19. Osstem Implant will continue to develop products that the dentist can trust. We will work to create greater customer value. Thank you.

CEO of OSSTEM IMPLANT  
Choi Kyu-ok (DDS.Ph.D)



# Worldwide & History



## 1997

- 01 Established 'Osstem Co., Ltd.'
- 12 Released 'Doobunae' (health insurance claim application software program)

## 2001

- 01 Obtained CE-0434 certification
- 03 Established AIC training center

## 2006

- 03 Changed the company name to Osstem Implant Co., Ltd
- 04 Obtained GOST-R certification (russia)
- 12 Established 12 overseas branches (first round)

## 2008

- 01 Established osstem bone science research center
- 12 Selected as a managing organization for the national strategic technology development project

## 2010

- 03 Launched TSIII SA line
- 06 Launched TSIII HA line

## 2012

- 06 Launched TSIII CA line
- 07 Established osstem dental equipment research institute

## 2014

- 05 Selected as 'World Class 300'
- 05 Released 'HyFlex', an impression material
- 08 Released 'BeauTis' whitening material

## 2016

- 01 Established Vussen Co., Ltd.
- 03 Acquired Cardiotec Co., Ltd.
- 08 Acquired Hubit Co., Ltd.
- 11 Launched OneGuide system

## 2000

- 06 Released 'Hanaro' (dentistry management software)
- 10 Acquired sumin comprehensive dental materials

## 2002

- 01 Established Osstem Implant R&D center
- 08 Obtained FDA certification, launched USII line
- 10 Launched SSII line

## 2007

- 02 Listed on KOSDAQ and began trading publicly
- 06 Selected as No.1 products for the next generation and obtained TGA certification (australia)

## 2009

- 10 Obtained approval for medical device manufacturing and sale from the ministry of health, labor and welfare, japan

## 2011

- 06 Osstem Implant R&D center was selected as ATC (advanced technology center)
- 07 Selected as 'World Champ' business
- 12 Launched 'K2 unit chair', which was selected as a 'World Class Product'

## 2013

- 01 Launched osstem xenograft material 'A-Oss'
- 09 Launched 'K3 unit chair'
- 10 Selected as a 'Hidden Champion' company

## 2015

- 03 Established Osstem BioPharma Co., Ltd.
- 12 Awarded 'USD 50 Million Export Tower'

## 2017

- 12 2017 presidential commendation for job creation


## 2018

- 01 TS exceeded 10 million production



# OSSTEM<sup>®</sup> Implant Design feature

OSSTEM IMPLANT has revolutionized implant dentistry in South Korea. With a focus on aggressive R&D, a commitment to education and a dedication to manufacturing the best products, Osstem Implant's ultimate goal is to become the global leader in implant dentistry.



Each implant system has its own unique color code

## Submerged type implant with an internal hex and 11 tapered connection

- Internal connection type - Mini / Regular
- Excellent initial stability in soft bone due to smaller threads in the upper section
- Corkscrew thread with cutting edges
  - Strong self-threading effect for easy fixture path
  - Higher initial stability and consistent insertion torque
- Different body types to properly match the patient's bone quality and clinical condition
  - TSII (straight body) : easy to adjust depth
  - TSIII (1.5° tapered body) : excellent initial stability necessary for immediate loading, even in soft bone
  - TSIV (6° tapered body) : specifically designed for the maxillary sinus and soft bone, excellent initial stability
- Available surface types - SA / CA / HA / BA / SOI

## Non-submerged type implant with an internal octa and 8 tapered connection

- Internal connection type - Regular / Wide
- Corkscrew thread with cutting edges
  - Strong self-threading effect for easy fixture path
  - Higher initial stability and consistent insertion torque
- Different body types to properly match the patient's bone quality and clinical condition
  - SSII (straight body) : easy to adjust the insertion depth
  - SSIII (1.5° tapered body) : excellent initial stability necessary for immediate loading, even in soft bone
- Available surface types - SA / CA / HA / BA

## Submerged type implant with an external hex connection structure

- Internal connection type - Mini / Regular / Wide / Wide PS
- Corkscrew thread with cutting edges
  - Strong self-threading effect for easy fixture path
  - Higher initial stability and consistent insertion torque
- Different body types to properly match the patient's bone quality and clinical condition
  - USII (straight body) : easy to adjust the insertion depth
  - USIII (1.5° tapered body) : excellent initial stability necessary for immediate loading, even in soft bone
  - USIV (6° tapered body) : specifically designed for the maxillary sinus and soft bone, excellent initial stability
- Available surface types - SA / CA



# OSSTEM<sup>®</sup> Implant Surface feature

The key factor in providing implant treatment safely and efficiently is surface technology.

OSSTEM IMPLANT is proud of its cutting-edge surface technology.

SA CA HA BA SOI

## Acid Treated Optimized Surface

- Ra 2.5~3.0 $\mu$ m surface roughness (note : the upper 0.5mm part of the implant has Ra 0.5~0.6 $\mu$ m)
- Consistent surface micro pits between 1 to 3 $\mu$ m
- Surface area is increased by 46 percent compared to RBM treated implants

## In-vitro & In-vivo Bone Response

- 20% improvement in osteoblast separation and ossification compared to RBM
- Initial bone reaction performance in animal model (mini-pig)
  - 48% improvement in initial stability (RT, 4 weeks) compared to RBM
  - 20% improvement in ossification (BIC, 4 weeks) compared to RBM

## Super-hydrophilic SA surface suspended in a calcium solution

- Same SA surface morphology
- Optimizing surface reaction by suspension in a calcium (CaCl<sub>2</sub>) solution
- Increased new bone formation area due to the excellent blood wettability
- Bone response improved in early osseointegration stage compared to standard SA surface

## In-vitro & In-vivo Bone Response

- Protein and cellular adhesion tripled compared to SA surfaces
- Initial cellular differentiation by 19 percent compared to SA surfaces (7 days)
- Initial stability increased by 34 percent compared to SA surfaces (RT at 4 weeks)
- Ossification rate Increased by 26 percent compared to SA surfaces (BIC at 4 weeks)

## Premium high-crystalline HA-coated surface

- 30 to 60 $\mu$ m thick high-crystalline HA coating
- HA coated onto a RBM surface (Ra 3.0 to 3.5 $\mu$ m)
- High HA crystalline over 98 percent
- Solved the problem with low-crystalline HA resorption

## In-vitro & In-vivo Bone Response

- Excellent biocompatibility in HA that is similar to bone
- Initial ossification by osteoblasts doubled compared to SA surfaces (5 days)
- 40% improvement in initial stability (RT, 4 weeks) in animal models compared to SA
- Suitable for poor bone quality, tooth extraction sites or immediate implant insertion

## Premium low crystalline nano-HA coated SA surface

- SA surface (Ra 2.5 to 3.0 $\mu$ m) coated with HA
- 10nm ultra-thin HA coating
- Dual function between titanium and HA
  - HA is naturally resorbed during ossification

## In-vitro & In-vivo Bone Response

- Advantages of both SA and HA surfaces
  - SA's ability to maintain an optimal surface
  - HA's ability to form high quality initial bone, even in a poor bone quality
- 40% improvement in ossification (BIC) compared to SA
- It is applicable to all types of bone quality

## Next-generation surface coated with special material (K material)

- Activation of blood clot formation
- Avoid carbon adsorption in air
- Coating of K material on SA surface (Ra 2,0~3,0 $\mu$ m)
- Superior blood wettability with super hydrophilic surface.

## In-vitro & In-vivo Bone Response

- Protein and cellular adhesion 130 times increase compared to SA surface
- Initial stability increased by 57 percent compared to SA surfaces (RT at 4 weeks)
- Surface with the shortest duration of surgery



US SYSTEM Contents



016  
USII SA  
Fixture



018  
USII CA  
Fixture



020  
USIII SA  
Fixture




022  
USIII CA  
Fixture




024  
USIV SA  
Fixture




056  
Esthetic-  
low  
Abutment




058  
Esthetic-  
low Healing  
Cap




058  
Esthetic-  
low Gold  
Cylinder




058  
Esthetic-  
low  
Plastic  
Cylinder




059  
Esthetic-  
low  
Temporary  
Cylinder  
Standard Type




026  
Simple  
Mount




026  
Cover  
Screw




026  
Headless  
Cover  
Screw




027  
Healing  
Abutment




029  
Cement  
Abutment




059  
Esthetic-  
low  
Temporary  
Cylinder  
Narrow Type




060  
Esthetic-  
low Pick-up  
Impression  
Coping




060  
Esthetic-  
low  
Transfer  
Impression  
Coping



060  
Esthetic-  
low Lab  
Analog



061  
Esthetic-  
low  
Polishing  
Protector



032  
Cement ID  
Abutment



035  
Angled  
Abutment



037  
UCLA  
Gold  
Abutment



038  
UCLA  
NP-Cast  
Abutment



039  
UCLA  
Plastic  
Abutment



062  
Multi  
Angled  
Abutment



065  
O-ring  
Abutment



066  
O-ring  
Retainer  
Cap Set




066  
O-ring  
Retainer  
Set



066  
O-ring  
Set



040  
UCLA  
Temporary  
Abutment



041  
Fixture  
Pick-Up  
Impression  
Coping




042  
Fixture  
Transfer  
Impression  
Coping



043  
Fixture  
Lab  
Analog



043  
UCLA  
Polishing  
Protector



066  
O-ring  
Lab  
Analog



067  
Locator®  
Abutment



068  
Locator®  
Male  
Processing  
Kit



068  
Locator®  
Replacement  
Male



068  
Locator®  
Extended  
Replacement  
Male



044  
OneFit  
Abutment



044  
Scan  
Boby



045  
Pre-Milled  
Abutment



046  
ZioCera  
Abutment



047  
ZioCera  
Angled  
Abutment



069  
Locator®  
Black  
Processing  
Male



069  
Locator®  
Block Out  
Spacers



069  
Locator®  
Impression  
Coping




069  
Locator®  
Lab Analog




070  
Locator®  
Core Tool




048  
Safe  
Abutment




052  
Esthetic  
Abutment




053  
Esthetic  
Healing  
Cap




053  
Esthetic  
Gold  
Cylinder




053  
Esthetic  
Plastic  
Cylinder




070  
Locator®  
Torque  
Driver




071  
OneSeal




054  
Esthetic  
Temporary  
Cylinder




054  
Esthetic  
Pick-up  
Impression  
Coping




055  
Esthetic  
Transfer  
Impression  
Coping



055  
Esthetic  
Lab  
Analog



055  
Esthetic  
Polishing  
Protector





# US SYSTEM

OSSTEM<sup>®</sup>  
IMPLANT

## FIXTURE

- 016** USII SA Fixture
- 018** USII CA Fixture
- 020** USIII SA Fixture
- 022** USIII CA Fixture
- 024** USIV SA Fixture
- 026** Simple Mount
- 026** Cover Screw
- 026** Headless Cover Screw
- 027** Healing Abutment

## COMPONENTS


- 028** PROSTHETIC FLOW DIAGRAM 1
- 029** Cement Abutment
- 035** Angled Abutment
- 037** UCLA Gold Abutment
- 038** UCLA NP-Cast Abutment
- 039** UCLA Plastic Abutment
- 040** UCLA Temporary Abutment
- 044** OneFit Abutment
- 045** Pre-Milled Abutment
- 046** ZioCera Abutment
- 047** ZioCera Angled Abutment
- 048** Safe Abutment
- 050** PROSTHETIC FLOW DIAGRAM 2
- 052** Esthetic Abutment
- 056** Esthetic-low Abutment
- 062** Multi Angled Abutment
- 064** PROSTHETIC FLOW DIAGRAM 3
- 065** O-ring Abutment
- 067** Locator<sup>®</sup> Abutment
- 071** OneSeal




USII SA Fixture

- Submerged type implant with external hex connection
- Optimized screw thread design with the ideal SA surface
- Straight body design allows easy insertion depth adjustments
- Corkscrew threading with excellent self-threading effect
- Recommended insertion torque : ≤40 Ncm
- ※ Fixtures with a diameter of 4.5mm or more are recommended for the posterior area


NoMount fixture order code  
: fixture product code (ex : US2R4010S)  
Pre-Mounted fixture (fixture + mount + cover screw) order code  
: A + fixture product code (ex : AUS2R4010S)




D Ø5.0  
P Ø5.1  
Hex 3.4  
W




D Ø5.0  
P Ø5.0  
Hex 2.7  
W<sup>PS</sup>




D Ø3.5  
P Ø3.5  
Hex 2.4  
M



D Ø4.0  
P Ø4.1  
Hex 2.7  
R



D Ø4.5  
P Ø4.1  
Hex 2.7  
R



Nominal and actual diameters may slightly differ  
Caution For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.

USII CA Fixture




- Submerged type implant with an external hex connection
- Super-hydrophilic SA surface suspended in a calcium solution
- Straight body design allows easy insertion depth adjustments
- Corkscrew threading with excellent self-threading effect


**Ultra-wide**

- Ideal for an extracted tooth site in the posterior area, for immediate placement, or for replacing a failed implant
- Recommended insertion torque : ≤40 Ncm
- ※ Fixtures with a diameter of 4.5mm or more are recommended for the posterior area

**Pre-Mounted fixture** (fixture + mount + cover screw) **order code**  
: **A** + fixture product code (ex : **AUS2R4010C**)




**D Ø5.0**  
**P Ø5.1**  
**Hex 3.4**  
**W**




Ultra-wide

**D Ø6.0**  
**P Ø5.1**  
**Hex 3.4**  
**W**



**D Ø7.0**  
**P Ø5.1**  
**Hex 3.4**  
**W**



**Nominal and actual diameters may slightly differ**  
**Caution** For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.


US SYSTEM

US SYSTEM


018

019


**D Ø3.5**  
**P Ø3.5**  
**Hex 2.4**  
**M**



**D Ø4.0**  
**P Ø4.1**  
**Hex 2.7**  
**R**



**D Ø4.5**  
**P Ø4.1**  
**Hex 2.7**  
**R**





USIII SA Fixture

- A submerged type Implant with external hex connection
- Optimal screw thread design for optimal SA surface implementation
- Taper body design with excellent initial fixation
- Powerful self-threading effect with corkscrew thread
- Acquired initial fixation force for immediate loading in soft bone

Ultra-wide

- Useful for Posterior extraction and immediate implant placement case, and for replacement of failed implant case
- Optimized apex design ensures stable initial fixation force even at the extraction and bottom 3mm
- Recommended implant placement torque : below 40Ncm


※ In posterior single case, fixture with diameter of 4.5mm or more is recommended

NoMount fixture order code


: fixture product code (ex : US3R4010S)

Pre-Mounted fixture (fixture + mount + cover screw) order code

: A + fixture product code (ex : AUS3R4010S)




D Ø5.0  
P Ø5.1  
Hex 3.4



US3W5006S US3W5007S US3W5008S US3W5010S US3W5011S US3W5013S


D Ø5.0  
P Ø5.0  
Hex 2.7



US3P5006S US3P5007S US3P5008S US3P5010S US3P5011S US3P5013S


Ultra-wide

D Ø6.0  
P Ø5.1  
Hex 3.4



US3W6006S US3W6007S US3W6008S US3W6010S US3W6011S US3W6013S

D Ø7.0  
P Ø5.1  
Hex 3.4




US3W7006S US3W7007S US3W7008S US3W7010S US3W7011S US3W7013S

Nominal and actual diameters may slightly differ


Caution For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.

D Ø3.5  
P Ø3.5  
Hex 2.4




US3M3508S US3M3510S US3M3511S US3M3513S

D Ø4.0  
P Ø4.1  
Hex 2.7



US3R4007S US3R4008S US3R4010S US3R4011S US3R4013S

D Ø4.5  
P Ø4.1  
Hex 2.7




US3R4507S US3R4508S US3R4510S US3R4511S US3R4513S


USIII CA Fixture

- A submerged type Implant with external hex connection
  - Super-hydrophilic SA surface in calcium solution.
  - Straight body design that can be easily adjusted in depth
  - Powerful self-threading effect with corkscrew thread
  - Acquired initial fixation force for immediate loading in soft bone
- Ultra-wide**
- Useful for Posterior extraction and immediate implant placement case, and for replacement of failed implant case
  - Recommended implant placement torque : below 40Ncm
- ※ In posterior single case, fixture with diameter of 4.5mm or more is recommended

**Pre-Mounted fixture** (fixture + mount + cover screw) **order code**  
: **A** + fixture product code (ex : **AUS3R4010S**)




**D Ø5.0**  
**P Ø5.1**  
**Hex 3.4**  
**W**




Ultra-wide

**D Ø6.0**  
**P Ø5.1**  
**Hex 3.4**  
**W**




**D Ø7.0**  
**P Ø5.1**  
**Hex 3.4**  
**W**




Nominal and actual diameters may slightly differ  
**Caution** For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.


**D Ø3.5**  
**P Ø3.5**  
**Hex 2.4**  
**M**



**D Ø4.0**  
**P Ø4.1**  
**Hex 2.7**  
**R**



**D Ø4.5**  
**P Ø4.1**  
**Hex 2.7**  
**R**






USIV SA Fixture




- A submerged type implant with external hex connection
- Optimal screw thread design for optimal SA surface implementation
- Fixture for maxillary sinus and soft bone
- Powerful self-threading effect with corkscrew thread
- With a sharp apex design, it can be placed after D4 bone 2.0/3.0mm drilling
- Recommended implant placement torque : below 40Ncm
- ※ In posterior single case, fixture with diameter of 4.5mm or more is recommended
- ※ USIV fixture is recommended to be lowered to 15rpm or less because the pitch of the thread is large and the implant placement speed is fast

NoMount fixture (fixture + mount + cover screw) order code  
: fixture product code (ex : US4R4010S)

Pre-Mounted fixture (fixture + mount + cover screw) order code  
: A + fixture product code (ex : AUS4R4010S)




D Ø5.0  
P Ø5.0  
Hex 2.7  
W<sup>PS</sup>




US4P5007S US4P5008S US4P5010S US4P5011S US4P5013S

D Ø4.0  
P Ø4.1  
Hex 2.7




US4R4007S US4R4008S US4R4010S US4R4011S US4R4013S

D Ø4.5  
P Ø4.1  
Hex 2.7



US4R4507S US4R4508S US4R4510S US4R4511S US4R4513S

D Ø5.0  
P Ø5.1  
Hex 3.4



US4W5007S US4W5008S US4W5010S US4W5011S US4W5013S


Nominal and actual diameters may slightly differ

Mount & Screw

Simple Mount

- Select according to the fixture platform
- Use of 1.2 hex driver with hand force
- ※ Disposable, Do not reuse
- P = Platform


- M Mini
- R Regular
- W Wide



Cover Screw

- Select according to the fixture platform
- Use of 1.2 hex driver with hand force
- P = Platform


- M Mini
- R Regular
- W Wide



Headless Cover Screw

- Used when the soft tissue of the suture part is insufficient
- 0.9 hex(only mini), Use of 1.2 hex driver with hand force
- P = Platform


- M Mini
- R Regular
- W Wide



Healing Abutment

- Use a 1.2 hex driver
- Recommended tightening torque : 5~8Ncm

- M Mini
- R Regular
- W Wide



D \ H	2.0	3.0	4.0	5.5	7.0
Ø 4.0	-	AIHM403	-	AIHM405	-
Ø 5.0	-	AIHM503	-	AIHM505	-
D \ H	2.0	3.0	4.0	5.5	7.0
Ø 4.1 One Piece	-	AIOHR403	-	AIOHR405	AIOHR407
Ø 4.1 Two Piece	-	AIHR403	-	AIHR405	AIHR407
Ø 5.0	AIHR502	AIHR503	AIHR504	AIHR505	AIHR507
Ø 6.0	AIHR602	AIHR603	AIHR604	AIHR605	AIHR607
D \ H	2.0	3.0	4.0	5.5	7.0
Ø 5.1 One Piece	-	AIOHW503	-	AIOHW505	-
Ø 5.1 Two Piece	-	AIHW503	-	AIHW505	-
Ø 6.0	AIHW602	AIHW603	AIHW604	AIHW605	-
Ø 7.0	AIHW702	AIHW703	AIHW704	AIHW705	-
Ø 6.0 PS	-	TIHW603	-	TIHW605	-





Cement Abutment

D Ø5.0

R

Ti screw  
: ASR200

H \ G/H	1.0	2.0	3.0	4.0	1.0	2.0	3.0	4.0
Type	Hex				Non-Hex			
4.0	CAR514	CAR524	CAR534	CAR544	CAR514N	CAR524N	CAR534N	CAR544N
5.5	CAR515	CAR525	CAR535	CAR545	CAR515N	CAR525N	CAR535N	CAR545N
7.0	CAR517	CAR527	CAR537	CAR547	CAR517N	CAR527N	CAR537N	CAR547N

D Ø7.0

W

Ti screw  
: ASW200

H \ G/H	1.0	2.0	3.0	4.0	1.0	2.0	3.0	4.0
Type	Hex				Non-Hex			
5.5	CAW715	CAW725	CAW735	CAW745	CAW715N	CAW725N	CAW735N	CAW745N

D Ø6.0

R

Ti screw  
: ASR200

H \ G/H	1.0	2.0	3.0	4.0	1.0	2.0	3.0	4.0
Type	Hex				Non-Hex			
5.5	CAR615	CAR625	CAR635	CAR645	CAR615N	CAR625N	CAR635N	CAR645N

D Ø6.0

W<sup>PS</sup>

Ti screw  
: ASR200

H \ G/H	1.0	2.0	3.0	4.0	1.0	2.0	3.0	4.0
Type	Hex				Non-Hex			
7.0	-	TCAW627	-	TCAW647	-	TCAW627N	-	TCAW647N

D Ø5.1

W

Ti screw  
: ASW200

H \ G/H	1.0	2.0	3.0	4.0	1.0	2.0	3.0	4.0
Type	Hex				Non-Hex			
7.0	-	-	CAW537		-	-	CAW537N	-

D Ø6.0

W

Ti screw  
: ASW200

H \ G/H	1.0	2.0	3.0	4.0	1.0	2.0	3.0	4.0
Type	Hex				Non-Hex			
4.0	CAW614	CAW624	CAW634	CAW644	CAW614N	CAW624N	CAW634N	CAW644N
5.5	CAW615	CAW625	CAW635	CAW645	CAW615N	CAW625N	CAW635N	CAW645N
7.0	CAW617	CAW627	CAW637	CAW647	CAW617N	CAW627N	CAW637N	CAW647N




Cement ID Abutment



• Non-insured cement abutment

Abutment + Ti screw order code  
: product code + TH (ex : BCAR525TH)



D Ø5.0



Ti screw  
: ASR200

H \ G/H Type	1.0	2.0	3.0	4.0					1.0	2.0	3.0	4.0
	Hex				Non-Hex							
4.0												
5.5												
7.0												

D Ø6.0



Ti screw  
: ASR200

H \ G/H Type	1.0	2.0	3.0	4.0					1.0	2.0	3.0	4.0
	Hex				Non-Hex							
5.5												

D Ø4.0



Ti screw  
: USABSMT

H \ G/H Type	1.0	2.0	3.0	4.0					1.0	2.0	3.0	4.0
	Hex				Non-Hex							
7.0	-		-		-		-					

D Ø4.1



Ti screw  
: ASR200

H \ G/H Type	1.0	2.0	3.0	4.0					1.0	2.0	3.0	4.0
	Hex				Non-Hex							
7.0	-	-		-	-	-		-				

D Ø5.1



Ti screw  
: ASW200

H \ G/H Type	1.0	2.0	3.0	4.0					1.0	2.0	3.0	4.0
	Hex				Non-Hex							
7.0	-	-		-	-	-		-				

D Ø6.0



Ti screw  
: ASW200






H \ G/H Type	1.0	2.0	3.0	4.0					1.0	2.0	3.0	4.0
	Hex				Non-Hex							
4.0												
5.5												
7.0												

Cement ID Abutment

D Ø7.0

W

Ti screw : ASW200

H \ G/H	1.0	2.0	3.0	4.0	1.0	2.0	3.0	4.0
Type	Hex				Non-Hex			
5.5								
	BCAW715	BCAW725	BCAW735	BCAW745	BCAW715N	BCAW725N	BCAW735N	BCAW745N

D Ø6.0

W PS

Ti screw : ASR200


H \ G/H	1.0	2.0	3.0	4.0	1.0	2.0	3.0	4.0
Type	Hex				Non-Hex			
7.0	-		-		-		-	
	-	BTCAW627	-	BTCAW647	-	BTCAW627N	-	BTCAW647N

Angled Abutment



- Cement/combination-retained prosthesis
- Angle compensation of 15/25°
- Double hex(dodecagon) connection allows 12 positions
- Fixture level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw





Abutment + Ti screw order code  
: product code + TH (ex : AAR5152CTH)



D Ø4.0

M





Ti screw : USABSMT

G/H	2.0	4.0	2.0	4.0
Angle	15°		25°	
				
	AAM4152C	AAM4154C	AAM4252C	AAM4254C

D Ø5.0

R





Ti screw : ASR200

G/H	2.0	4.0	2.0	4.0
Angle	15°		25°	
				
	AAR5152C	AAR5154C	AAR5252C	AAR5254C

D Ø6.0

W

Ti screw : ASW200

G/H	2.0	4.0	2.0	4.0
Angle	15°		25°	
				
	AAW6152C	AAW6154C	AAW6252C	AAW6254C



Angled Abutment


<b>D Ø6.0</b> <b>W</b> <sup>PS</sup> Ti screw : ASR200	G/H	2.0	4.0	2.0	4.0
	Angle	15°		25°	
					
		TAAW6152C	TAAW6154C	TAAW6252C	TAAW6254C

UCLA Gold Abutment



- Cement/combination/screw abutment for prosthesis manufacturing
- Used to make customized prosthesis by casting with gold alloy
- Melting temperature of Abutment : 1,400~1,450°C
- Fixture level impression
- Use of 1.2 hex driver
- Recommended torque of tightening screw : 30Ncm
- Packing unit : abutment + Ti screw



Abutment + Ti screw order code  
: product code + TH (ex : GCR200TH)



<b>D Ø4.0</b> <b>M</b> Ti screw : USABSMT	Type	Hex	Non-Hex
			
		GCM200	GCM100

<b>D Ø4.5</b> <b>R</b> Ti screw : ASR200	Type	Hex	Non-Hex
			
		GCR200	GCR100

<b>D Ø5.5</b> <b>W</b> Ti screw : ASW200	Type	Hex	Non-Hex
			
		GCW200	GCW100

<b>D Ø5.5</b> <b>W</b> <sup>PS</sup> Ti screw : ASR200	Type	Hex	Non-Hex
			
		TGCW200	TGCW100


# UCLA NP-Cast Abutment

# UCLA Plastic Abutment




- Cement/combination/screw-retained prosthesis
- Customized prosthesis cast with non precious alloys
- Abutment melting point : 1400~1550°C
- Fixture level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw

Abutment + Ti screw order code  
: product code + TH (ex : NCR200TH)



- Cement/combination/screw-retained prosthesis
- Customized prosthesis, cast with gold alloys up to the joint
- Fixture level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw

Abutment + Ti screw order code  
: product code + TH (ex : PSR200TH)



US SYSTEM

038

D Ø4.0

M

Ti screw  
: USABSMT

Type	Hex	Non-Hex
	NCM200	NCM100

D Ø4.5

R

Ti Screw  
: ASR200

Type	Hex	Non-Hex
	NCR200	NCR100

D Ø5.5

W

Ti screw  
: ASR200

Type	Hex	Non-Hex
	NCW200	NCW100

D Ø5.5

W<sup>PS</sup>

Ti screw  
: ASR200

Type	Hex	Non-Hex
	TNCW200	TNCW100

US SYSTEM

039

D Ø4.0

M

Ti screw  
: USABSMT

Type	Hex	Non-Hex
	PSM200	PSM100

D Ø4.5

R

Ti screw  
: ASR200

Type	Hex	Non-Hex
	PSR200	PSR100

D Ø5.5

W

Ti screw  
: ASW200

Type	Hex	Non-Hex
	PSW200	PSW100

D Ø5.5

W<sup>PS</sup>

Ti screw  
: ASR200

Type	Hex	Non-Hex
	TPSW200	TPSW100




UCLA Temporary Abutment

UCLA Abutment Components

- Cement/screw-retained prosthesis
- A trim able provisional prosthesis (made of Ti Gr-3)
- Fixture level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw


Abutment + Ti screw order code  
: product code + TH (ex : TAR200TH)










Fixture Pick-Up Impression Coping

- Components for fixture level impression taking
- For open tray impressions
- Unique design that is fixed position in the impression material
- Use a 1.2 hex driver (torque manually)
- Packing unit : impression coping body + guide pin(\*)

- M** Mini (Yellow)
- R** Regular (Green)
- W** Wide (Blue)



D \ L	7				12			Guide Pin		
	Type	Hex	Non-Hex	Hex	Non-Hex	10	15	17		
										
Ø 4.0	-	-	ICFM400	ICFM400N	-	CSM150	-			
Ø 5.0	ICSR500	ICSR500N	ICFR500	ICFR500N	CSR100	CSR150*	CSR170			
Ø 6.0	-	-	ICFR600	ICFR600N	-	-	-			
Ø 6.0	ICSW600	ICSW600N	ICFW600	ICFW600N	CSW100	CSW150*	-			
Ø 6.0 <sub>PS</sub>	-	-	TICFW600	TICFW600N	-	TCSW150	-			

D Ø 4.0	Type	Hex	Non-Hex
<b>M</b>			
Ti screw : USABSMT		TAM200	TAM100

D Ø 4.5	Type	Hex	Non-Hex
<b>R</b>			
Ti screw : ASR200		TAR200	TAR100

D Ø 5.5	Type	Hex	Non-Hex
<b>W</b>			
Ti screw : ASW200		TAW200	TAW100


D Ø 5.5 <sup>PS</sup>	Type	Hex	Non-Hex
<b>W<sup>PS</sup></b>			
Ti screw : ASR200		TTAW200	TTAW100

UCLA Abutment Components

Fixture Transfer Impression Coping

- Components for fixture level impression taking
- For closed tray impressions
- Triangular arc ensures precise placement
- Use a 1.2 hex driver (torque manually)
- Packing unit
  - Hex : impression coping body + guide pin
  - Non-hex : impression coping

- M** Mini (Yellow)
- R** Regular (Green)
- W** Wide (Blue)



D \ L	10.5		13.5	
	Hex	Non-Hex	Hex	Non-Hex
Ø 4.0	ICPM402S	ICPM401S	ICPM402L	ICPM401L
Ø 5.0	ICPR502S	ICPR501S	ICPR502L	ICPR501L
Ø 6.0	ICPW602S	ICPW601S	ICPW602L	ICPW601L
Ø 6.0 PS	-	-	TICPW602	TICPW601

UCLA Abutment Components

Fixture Lab Analog

- A lab analog for fixture level impressions
- Packing unit : lab analog

- M** Mini
- R** Regular
- W** Wide

C	
<b>M</b>	FAM300
<b>R</b>	FAR300
<b>W</b>	FAW300
<b>W</b> <sup>PS</sup>	TFAW300
<b>W</b> <sup>R-type</sup>	RFAW300

UCLA Polishing Protector

- Protects GoldCast/plastic cylinder joints during polishing process
- Use a 1.2 hex driver (torque manually)

- M** Mini
- R** Regular
- W** Wide

C	
<b>M</b>	UPCM100
<b>R</b>	UPCR100
<b>W</b>	UPCW100
<b>W</b> <sup>PS</sup>	TUPCW100
<b>W</b> <sup>R-type</sup>	RUPCW100

# OneFit Abutment

# Pre-Milled Abutment

- Cement/combination abutment for prosthesis manufacturing
- Custom abutment manufactured by CAD/CAM
- Fixture level impression
- Abutment level impression when scan healing abutment is used
- Production period (based on working day)
  - Titanium : 5 days
  - Titanium + gold color : 7 days
- Use of 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw



- Manufacturing custom abutment with dental milling equipment
- Superior tightening accuracy compared to non-genuine
- Packing unit : pre-milled abutment + Ti

Pre-milled abutment + screw order code  
: product code + TH (ex : USPM10AGRTH)



## Scan Body

- Scan body for titanium custom abutment production
- Convenient specification classification through screw color
- Use of 1.2 hex driver with hand force
- Packing unit : scan body + Ti screw

Scan body + screw order code  
: product code + TH (ex : USSBMTH)

- M** Mini
- R** Regular
- W** Wide



Equipment	Implant	D	Specifications		Code
DooWon ARUM	Osstem US	Ø 10	Mini	Hex	USPM10AGMTH
			Mini	Non-hex	USPM10AGMNTH
			Regular	Hex	USPM10AGRTH
			Regular	Non-hex	USPM10AGRNTH
			Wide	Hex	USPM10AGWTH
			Wide	Non-hex	USPM10AGWNTH




# ZioCera Abutment



- Cement/combination/screw-retained prosthesis
- Zirconium abutment for an all ceramic prostheses ideal for anterior - esthetic zone
- Fixture level impression
- Use a 1.2 hex driver
- Abutment screw included
- Recommended tightening torque : 30Ncm(regular)
- Packing unit : abutment + Ti screw

Abutment + Ti screw order code  
: product code + TH (ex : ZAR537TH)




# ZioCera Angled Abutment



- Cement/combination/screw-retained prosthesis
- Zirconium abutment for all ceramic prostheses ideal for the anterior esthetic zone
- Angle compensation up to 23° without the need for additional adjustment
- Fixture level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm(regular)
- Packing unit : abutment + Ti screw

Abutment + Ti screw order code  
: product code + TH (ex : ZAAR5173TH)



US SYSTEM

046

D Ø5.0	H \ G/H	3.0	5.0
M			
Ti screw : ASM200	7.0	ZAM537	ZAM557

D Ø5.0	H \ G/H	3.0	5.0
R			
Ti screw : ASR200	7.0	ZAR537	ZAR557

D Ø6.0	H \ G/H	3.0	5.0
R			
Ti screw : ASR200	7.0	ZAR637	ZAR657

US SYSTEM

047


D Ø5.5	H \ G/H	3.0
R		
Ti screw : ASR200	9.0	ZAAR5173

D Ø6.5	H \ G/H	4.0
R		
Ti screw : ASR200	9.0	ZAAR6174

Safe Abutment



- Cement-retained prosthesis
- Designed to minimize screw loosening
- No modification or trimming is required
- Fixture/abutment level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw + carrier cap + protect cap



OSSTEM<sup>®</sup>  
IMPLANT

D Ø4.8



H \ G/H	1.0	2.5	4.0
4.0	SFAR514SC	SFAR524SC	SFAR544SC
5.5	SFAR515SC	SFAR525SC	SFAR545SC
7.0	SFAR517SC	SFAR527SC	SFAR547SC

D Ø6.0




H \ G/H	1.0	2.5	4.0
4.0	SFAW614SC	SFAW624SC	SFAW644SC
5.5	SFAW615SC	SFAW625SC	SFAW645SC

Esthetic / Esthetic-low / Multi Angled

Abutment Level Impression


US SYSTEM

050



US SYSTEM

051






# Esthetic Abutment



- Screw-retained prosthesis in multiple scenarios
- Connection is located near the gum line
- Abutment-level impression
- Angle compensation of up to 30°
- Tightening with exclusive outer driver (code : TIHD20L / TIHD20S)
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw


Abutment + Ti screw order code  
: product code + **TH** (ex : EAR200**TH**)



D Ø4.8



G/H 1.0 2.0 3.0 4.0



# Esthetic Abutment Components

## Esthetic Healing Cap

- Protective cap
- Use a 1,2 hex driver (torque manually)

**R** Regular

D \ H 6.0

Ø4.8



EHC100

## Esthetic Gold Cylinder

- Screw-retained prosthesis
- Customizable prosthesis, cast with gold alloy
- Cylinder melting point : 1400~1450°C
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + EbonyGold screw order Code  
: product code + **TH** (ex : EGC200**TH**)

**R** Regular

D \ Type Hex Non-Hex

Ø4.8

Ti screw  
: TS200



EGC200



EGC100



## Esthetic Plastic Cylinder

- Screw-retained prosthesis
- Customizable prosthesis, cast with non precious alloys
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti Screw order code  
: product code + **TH** (ex : ETT200**TH**)

**R** Regular

D \ Type Hex Non-Hex

Ø4.8

Ti screw  
: TS200



EPS200



EPS100




Esthetic Abutment Components

Esthetic Temporary Cylinder

- Used in fabricating a temporary prosthesis (Ti Gr-3)
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti screw order code  
: product code + TH (ex : ETT200TH)

**R** Regular

D \ Type	Hex	Non-Hex
		
	ETT200	ETT100

Ø 4.8

Ti screw : TS200

Esthetic Transfer Impression Coping

- Transfer Impression coping for esthetic abutment
- Use a 1.2 hex driver (torque manually)

**R** Regular


D \ H	8.0
	
	ETR100







Ø 4.8

Esthetic Pick-up Impression Coping

- Pick up impression coping for esthetic abutment
- Use a 1.2 hex driver (fastened manually)
- Packing unit : impression coping body + guide pin(\*)

**R** Regular



D \ L	Hex	Non-Hex	5	10	12	15
						
	ESR200	ESR100	GP100	GP150*	GP170	GP200

Guide Pin

Ø 4.8

Esthetic Lab Analog

- Lab analog for esthetic abutment
- Use a 1.2 hex driver (torque manually)

**R** Regular

D	
	
	ERR300

Ø 4.8

Esthetic Polishing Protector

- Protects GoldCast/plastic cylinder joints during polishing process
- Use a 1.2 hex driver (torque manually)

**R** Regular

D	
	
	EPCR100


Ø 4.8

# Esthetic-low Abutment



- Screw-retained prosthesis for multiple case scenarios
- The prosthetic connection matches the crest of the soft tissue
- Abutment level impression
- Angle compensation of up to 48°
- Torque using a dedicated outer driver
  - Regular : 2.0 internal hex driver (code : TIHD20S/TIHD20L)
  - Wide : 2.7 internal hex driver (code : TIHD27)
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw

Abutment + Ti screw order code  
: product code + **TH** (ex : MER200**TH**)



**D Ø5.5**  
**W**

G/H	1.0	2.0	3.0	4.0
	MEW100	MEW200	MEW300	MEW400

**D Ø5.5**  
**W<sup>PS</sup>**

G/H	1.0	2.0	3.0	4.0
	TMEW100	TMEW200	TMEW300	TMEW400

G/H	1.0	2.0	3.0	4.0
	MEM100	MEM200	MEM300	MEM400

G/H	1.0	2.0	3.0	4.0
	MER100	MER200	MER300	MER400

US SYSTEM

US SYSTEM

056

057

**D Ø4.8**  
**M**

**D Ø4.8**  
**R**





Esthetic-low Abutment Components

Esthetic-low Healing Cap

- A protective cap
- Use a 1.2 hex driver (torque manually)

- M Mini
- R Regular
- W Wide

D \ H	6.0
	
	MHCR100 MHCW100



Ø 4.8 / Ø 4.8  
Ø 5.5 / Ø 5.5 PS

Esthetic-low Gold Cylinder

- Screw-retained prosthesis
- Cast with gold alloys
- Cylinder melting point : 1400~1450°C
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti screw order code  
: product code + TH (ex : MGR200TH)

- M Mini
- R Regular
- W Wide

D \ Type	Hex	Non-Hex
		
	MGR200 MGW200	MGR100 MGW100

Ø 4.8 / Ø 4.8  
Ø 5.5 / Ø 5.5 PS



Ti screw  
: MTS200 (Ø 4.8 / Ø 4.8)  
: WTS200 (Ø 5.5 / Ø 5.5PS)

Esthetic-low Plastic Cylinder

- Screw-retained prosthesis
- Cast with non-precious metal alloys
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti screw order code  
: Product code + TH (ex : MEPR200TH)

- M Mini
- R Regular
- W Wide

D \ Type	Hex	Non-Hex
		
	MEPR200 MEPW200	MEPR100 MEPW100

Ø 4.8 / Ø 4.8  
Ø 5.5 / Ø 5.5 PS

Ti screw  
: MTS200 (Ø 4.8 / Ø 4.8)  
: WTS200 (Ø 5.5 / Ø 5.5PS)



Esthetic-low Temporary Cylinder

Standard Type

- Used in fabricating a temporary prosthesis (Ti Gr-3)
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti screw order code  
: product code + TH (ex : MTR200TH)

- M Mini
- R Regular
- W Wide

D \ Type	Hex	Non-Hex
		
	MTR200 MTW200	MTR100 MTW100

Ø 4.8 / Ø 4.8  
Ø 5.5 / Ø 5.5 PS



Ti screw  
: MTS200 (Ø 4.8 / Ø 4.8)  
: WTS200 (Ø 5.5 / Ø 5.5PS)

Narrow Type

- Used in fabricating a temporary prosthesis (Ti Gr-3)
- Ideal for overdenture cases due to its smaller diameter
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti screw order code  
: product code + TH (ex : NMTR200TH)

- M Mini
- R Regular
- W Wide

D \ Type	Hex	Non-Hex
		
	NMTR200 NMTW200	NMTR100 NMTW100

Ø 4.8 / Ø 4.8  
Ø 5.5 / Ø 5.5 PS


Ti screw  
: MTS200 (Ø 4.8 / Ø 4.8)  
: WTS200 (Ø 5.5 / Ø 5.5PS)

# Esthetic-low Abutment Components

## Esthetic-low Pick-up Impression Coping

- Components for fixture level impression taking
- Pick up impression coping for esthetic-low abutment
- Use a 1.2 hex driver (torque manually)
- Packing unit : impression coping body + guide pin(\*)

- M Mini
- R Regular
- W Wide



D \ L	Hex		Non-Hex		Guide Pin			
					5	10	12	15
		MSR200	MSR100		GP100	GP150*	GP170	GP200
		MSW200	MSW100		GPW100	GPW150*	-	-

Ø 4.8 / Ø 4.8  
Ø 5.5 / Ø 5.5 PS

## Esthetic-low Transfer Impression Coping

- Transfer impression coping for esthetic-low abutment
- Use a 1.2 hex driver (fastened manually)

- M Mini
- R Regular
- W Wide

D \ H	8.0
	MTTR100
	MTTW100

Ø 4.8 / Ø 4.8  
Ø 5.5 / Ø 5.5 PS

## Esthetic-low Lab Analog

- Lab analog for esthetic-low abutment
- Use a 1.2 hex driver (fastened manually)

- M Mini
- R Regular
- W Wide

D	
	MERR300
	MERW300

Ø 4.8 / Ø 4.8  
Ø 5.5 / Ø 5.5 PS

## Esthetic-low Polishing Protector

- Protects GoldCast/plastic cylinder joints during polishing process
- Use a 1.2 hex driver (torque manually)


- M Mini
- R Regular
- W Wide

D	
	MPCR100
	MPCW100

Ø 4.8 / Ø 4.8  
Ø 5.5 / Ø 5.5 PS



# Multi Angled Abutment







- Used in the manufacture of screw retained prosthesis in multiple cases
- It compensates the angle of the fixture up to 108°
- The same platform as esthetic low abutment
- Prosthetics production with US esthetic low cylinder (regular/non-hex)
- Use dedicated abutment screw
- Use of 1.2 hex driver
- Recommended tightening torque : 20Ncm(mini), 30Ncm(regular)
- Packing unit : abutment + Ti screw




Abutment + Ti screw order code  
: product code + TH (ex : US17MAR4830TH)

OSSTEM<sup>®</sup>  
IMPLANT

D Ø4.8 M Ti screw : USMABSM	H/G	2.0	3.0	4.0	3.0	4.0	5.0
	Angle		17°			30°	
				-	-	-	-
		US17MAM4820	US17MAM4830	-	-	-	-


D Ø4.8 R Ti screw : USMABSR	H/G	2.0	3.0	4.0	3.0	4.0	5.0
	Angle		17°			30°	
							
		US17MAR4820	US17MAR4830	US17MAR4840	US30MAR4830	US30MAR4840	US30MAR4850





# O-ring Abutment

- Retains overdentures with an o-ring system
- Angle compensation of up to 20°
- Torque using a dedicated outer driver (code : AORD)
- Recommended tightening torque : 30Ncm



**D Ø5.0**  
**R**

G/H	2.0	3.0	4.0	5.0	6.0
	OAA200	OAA300	OAA400	OAA500	OAA600

**D Ø5.6**  
**W**

G/H	2.0	3.0	4.0	5.0	6.0
	OAAW200	OAAW300	OAAW400	OAAW500	OAAW600

# O-ring Abutment Components

## O-ring Retainer Cap Set

- O-ring housing
- Place appropriate o-ring in the metal housing before connecting to the abutment
- Packing unit : retainer cap + o-ring



## O-ring Retainer Set

- Used when vertical dimension is shorter than the retainer cap
- Packing unit : retainer cap + o-ring



## O-ring Set

- O-ring set
- Packing unit : o-ring x 5ea




## O-ring Lab Analog

- Lab analog for o-ring abutment








# Locator® Abutment

- Genuine zest anchors abutment
- Angle compensation of up to 40°
- 1.5mm lower profile, attachment with various and stable retention forces
- Torque using a dedicated outer driver (code : TWLDLK/TWLDLSK)
- Recommended tightening torque : 30Ncm








D Ø3.5

M

G/H	1.0	2.0	3.0	4.0	5.0
					
	HULCA3510M	HULCA3520M	HULCA3530M	HULCA3540M	HULCA3550M






D Ø4.1

R

G/H	1.0	2.0	3.0	4.0	5.0
					
	HULCA4010R	HULCA4020R	HULCA4030R	HULCA4040R	HULCA4050R

D Ø5.0

W

G/H	1.0	2.0	3.0	4.0	5.0
					
	HULCA5010W	HULCA5020W	HULCA5030W	HULCA5040W	HULCA5050W

# Locator® Abutment Components



## Locator® Male Processing Kit

- **Components**
  - Block out spacer / denture cap connected to black processing male
  - Replacement male blue/pink/clear
- A full range of retentive males are included with each denture cap to allow personalized retention for each specific patient
- Locator core tool places and removes nylon retentive males
- Packing unit : 2set



## Locator® Replacement Male

- Retention force : approx. 6N
- Angle compensation of up to 20°
- Packing unit : 4ea



LRM06S

- Retention force : approx. 12N
- Angle compensation of up to 20°
- Packing unit : 4ea



LRM12S

- Retention force : approx. 22N
- Angle compensation of up to 20°
- Packing unit : 4ea



LRM22S

## Locator® Extended Replacement Male

- Retention force : approx. 6N
- Angle compensation from 20~40°
- Packing unit : 4ea



LEM06S

- Retention force : approx. 12N
- Angle compensation from 20~40°
- Packing unit : 4ea



LEM12S

## Locator® Black Processing Male

- A nylon male used in prosthesis fabrication process
- Packing unit : 4ea



LBPS

## Locator® Block Out Spacers

- Place block-out spacers on the heads of the locator abutments. Position denture cap with integrated black processing onto the locator abutments. If necessary, add additional block-out spacers until no gap is visible between female, block-out spacer and gum.
- Packing unit : 20ea



LBSS

## Locator® Impression Coping

- A pick up impression coping
- Closed tray is used
- Packing unit : 4ea



LICS

## Locator® Lab Analog

- A lab analog for locator abutment
- Packing unit : 4ea



LAL40S



Locator<sup>®</sup> Abutment Components

Locator<sup>®</sup> Core Tool

- Places and removes nylon retentive males in the denture cap
- Separated into three different tools, includes a hand driver for locator abutment



Locator<sup>®</sup> Torque Driver


- Torque driver



OneSeal

OneSeal

- Disposable medical devices for internal filling of superstructure
- Cut to desired length (medical silicone)
- Packing unit : long 5ea



# Instructions for Use (AUG. 2017, Ver. 5.5)

## Description of Osstem implant system

Osstem Implant is a brand for implant materials for dental practices, and the fixture is made mainly of titanium. The abutment, prosthetic components and tools for the Osstem Implant system are compatible with the Osstem Implant fixture only. Using this product in combination with products from other manufacturers may cause various problems including loosening and fracture due to incomplete locking and compatibility issues. Refer to the manual or the catalogue or our website (www.osstem.com) for details. See the product label for the product code, specifications, manufacturing date, and expiration date.

## Sterility

The fixture, cover screw, and healing abutment are cleansed and sterilized with gamma radiation. This product is a disposable sterilized medical device intended for one-time use. In order to prevent contamination or infection of the product or operated site, the product must be used using a sterilized instrument in a sterilized environment. Damaged products, products with open packaging, or expired products must be discarded due to potential risks of contamination, infection, or osseointegration failure. Re-sterilization or re-use of the product may result in infection, osseointegration failure, or implant damage due to reduced accuracy.

## Storage condition

Keep the product in a dry place at room temperature(1~30℃). Keep away from direct sunlight.

## General precautions

The surgical technology of dental implant involves an expert, complex procedure. Formal training is required to perform implant surgery. Careful considerations must be made before the operation in case of bone disorders (osteoporosis, osteomalacia) or metabolic disorders of the bone.

## Precautions

Determine the local anatomy and suitability of the available bone for implant placement. Prepare the implant considering the expected situations and cautions. Excessive occlusal load may cause loosening or fracture of an implant. In order to avoid this condition, the implant must be placed in accurate location and direction considering the relationship between the implant and opposing dentition. Visual inspection as well as panoramic and periapical radiographs are essential to determine anatomical landmarks, occlusal conditions, periodontal status, and the adequacy of the bone. Adequate radiographs, direct palpation, and visual inspection of the implant site are necessary prior to implant surgery.

## Procedural precautions

Osstem Implant System is for single and two stage surgical procedures. As much as possible, try to minimize damage to the cell tissue and surgical trauma, pay special attention to maintaining the temperature at the implant site and removal of the source of contamination and infection. All drills and taps must be sufficiently and continuously irrigated for cooling during use. Implant placement should be accomplished at very low speed (25-30 rpm) or manually. Excessive torque (greater than 55Ncm) in the fixture placement can have adverse effects such as partial fracture or necrosis of the bone. Placing an implant tilted by 30° or higher is not recommended due to possible fracture of implant. Immediate loading to the fixture right after the surgery should be avoided. The bone quality and initial stability after fixture placement are important elements in determining the appropriate loading time. Mini-diameter implant or implant with diameter of 4.0 or less and which integrates with angled abutment may be fractured due to limitations of structural rigidity. They are not recommended for use in a posterior area. The Ultra-Wide fixtures are intended to be used only to replace molar teeth and

that angled abutments are not to be used with the Ultra-Wide fixtures. Evaluate the quantity of bone and radiographs to assess any potential anatomical contraindications to use of the Ultra-Wide fixture. For the placement of the Short Implant (diameter is 5mm or more and length is shorter than 7mm) which is used on the molar region only, clinicians should closely examine the patients for any of the following conditions: 1) periimplant bone loss, 2) changes to implant's response to percussion, 3) radiographic changes in bone to implant contact along the implant's length. If a short implant shows mobility or greater than 50% bone loss, the implant should be considered for possible removal. And clinicians should consider a two-stage surgical approach, splinting a short implant to an additional implant, and placement of the widest possible fixture. Allow longer healing periods for osseointegration before fabrication of the prosthesis and avoid immediate loading. Products with diameter of 3.25mm or less must be used exclusively for mandibular anterior teeth in order to prevent fracture due to excessive occlusal load. It is recommended that you should avoid applying HA coated fixture to hard bone, and the insertion torque of the implant should be less than 35Ncm, because cracks or damages might occur in the coated layer during implant placement. The surfaces of CA and SOI have the same physical shape as the SA surface made through blasting and etching treatments. After the SA surface treatment, to prevent the products' exposure to the atmosphere, CA is stored in solution, whereas SOI is stored in water-film coating form; it is designed to maintain the chemically activated state of the SA surface. Thus, CA or SOI products should be implanted in the target region at least within 15 minutes of taking them out of the container.

## Warning

The selection of inappropriate patients and surgical methods can cause implant failure or loss of bone supporting the implant. Osstem implants must not be used for purposes other than the recommended use and must not be remodeled. Implant mobility, bone loss, and chronic infection can result in failure of the implant surgery.

## Indications for use

The Osstem Implant System is an artificial dental root that has been designed for use in dental implant treatment in order to recover lost teeth. The system is implanted via a surgical method in maxillary or mandibular bone to replace natural dental root. The Osstem Implant System is indicated for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-units restorations including; cemented retained, screw retained, or overdenture restorations, and final or temporary abutment support for fixed bridgework. It is intended for delayed loading. Products with diameter of 3.25mm or less must be used exclusively for mandibular anterior teeth in order to prevent fracture due to excessive occlusal load.

## Side effects

A few problems may occur after the operation (loss of implant stability, damage of prosthesis, etc.). Deficient quality and quantity of the remaining bone, infection, allergic reaction, inferior oral hygiene or uncooperativeness of patient, implant mobility, partial deterioration of tissue, and improper position or arrangement of implants may cause the above mentioned problems.

## Contraindications

- Contraindications include the following, but are not limited to:
- Patients with hemophilia or difficulties related to bone or wound treatment
  - Patients with uncontrollable diabetes, heavy smoker or alcoholic
  - Patients whose immunity system is inactive due to chemical therapy or radiation therapy
  - Patients with oral infection or inflammation (improper oral hygiene, bruxism)
  - Patients with untreatable occlusion/joint disorder, insufficient dental arch space
  - Any patient who is not suitable for an surgery



Manufacturer : Osstem Implant Co., Ltd.  
203, Geonje-daero, Yeonje-gu, Busan, Korea  
TEL 82-51-850-2500 FAX 82-51-861-4693



DEUTSCHE OSSTEM GmbH.  
Mergenthalerallee 25  
65760 Eschborn, Germany  
+49-(0)6196-777-550

Storage condition  
Dry place at room temperature

## Rx only

For USA only : Federal law restricts this device to sale by or on the order of a dentist



2460



Sterilized using irradiation



Use by



Manufacture



Do not reuse



Date of manufacture



Keep away from sunlight



Catalogue number



Non-Sterile



Keep dry



Batch code



Do not resterilize



Caution, Consult accompanying documents



**OSSTEM<sup>®</sup>**  
**IMPLANT**

