

OSSTEM[®]
IMPLANT

SS SYSTEM CATALOG



OSSTEM[®]
IMPLANT



SS SYSTEM CATALOG



Osstem Implant 2018-19 Comprehensive Catalog

Overall Planning/Editing PR Department Design Team

Supervision Implant Lab, Marketing PM

Production/Distribution Marketing & Planning Team

Date of Publication 2018.08

Publisher Osstem Implant

8th FL, World Meridian II, 123, Gasan digital 2-ro,

Geumcheon-gu, Seoul, Korea

Phone +82.2.2016.7000

Fax +82.2.2016.7001

www.osstem.com

2018 Printed in Seoul, Korea

003	INTRODUCTION
012	CONTENTS
016	SS SYSTEM
076	USER MANUAL

**Providing cutting
edge technology and
superior quality**

Making products that dentists want to use,
trust, and are satisfied with :
This is our mission at **OSSEM IMPLANT**

**We are forever grateful
to all the dentists
who have given
unwavering support to
OSSTEM IMPLANT**


Thank you for using Osstem Implant. Osstem, Korea's first implant manufacturer, has secured world-class implant competitiveness through continuous R&D investment and quality innovation. It has grown to become Asia-Pacific No.1 and World No.5 Implant Company. In addition to dental implants and treatment tools, we are leading the development of products that are essential for dentists, including dental equipment, dental materials, and dental IT, and contribute to the development of the dental industry. The comprehensive catalog of the 2018-19 product series published here shows Osstem's technology-rich products. We have focused on catalog structure so that it is convenient to browse and order products. In particular, in the case of fixtures, abutments, and surgical tools, we introduced the diameter, length, and functions in detail.

GBR products are also easy to order by type, size and capacity. In addition, the product release date and time are displayed so that customers can understand when the existing product is released and what the newly released product is. We also introduced the CAD/CAM product in terms of preparing the digital dentistry, a major trend in the dentistry. In terms of design, we also implemented high-quality images of representative products by specification. By applying representative colors for each product system, it is easy to sort by category. We hope this will help you effectively find and purchase the products you need from the dental clinic of 2018-19. Osstem Implant will continue to develop products that the dentist can trust. We will work to create greater customer value. Thank you.

CEO of OSSTEM IMPLANT
Choi Kyu-ok (DDS.Ph.D)



Worldwide & History



1997

- 01 Established 'Osstem Co., Ltd.'
- 12 Released 'Doobunae' (health insurance claim application software program)

2001

- 01 Obtained CE-0434 certification
- 03 Established AIC training center

2006

- 03 Changed the company name to Osstem Implant Co., Ltd
- 04 Obtained GOST-R certification (russia)
- 12 Established 12 overseas branches (first round)

2008

- 01 Established osstem bone science research center
- 12 Selected as a managing organization for the national strategic technology development project

2010

- 03 Launched TSIII SA line
- 06 Launched TSIII HA line

2012

- 06 Launched TSIII CA line
- 07 Established osstem dental equipment research institute

2014

- 05 Selected as 'World Class 300'
- 05 Released 'HyFlex', an impression material
- 08 Released 'BeauTis' whitening material

2016

- 01 Established Vussen Co., Ltd.
- 03 Acquired Cardiotec Co., Ltd.
- 08 Acquired Hubit Co., Ltd.
- 11 Launched OneGuide system

2000

- 06 Released 'Hanaro' (dentistry management software)
- 10 Acquired sumin comprehensive dental materials

2002

- 01 Established Osstem Implant R&D center
- 08 Obtained FDA certification, launched USII line
- 10 Launched SSII line

2007

- 02 Listed on KOSDAQ and began trading publicly
- 06 Selected as No.1 products for the next generation and obtained TGA certification (australia)

2009

- 10 Obtained approval for medical device manufacturing and sale from the ministry of health, labor and welfare, japan

2011

- 06 Osstem Implant R&D center was selected as ATC (advanced technology center)
- 07 Selected as 'World Champ' business
- 12 Launched 'K2 unit chair', which was selected as a 'World Class Product'

2013

- 01 Launched osstem xenograft material 'A-Oss'
- 09 Launched 'K3 unit chair'
- 10 Selected as a 'Hidden Champion' company

2015

- 03 Established Osstem BioPharma Co., Ltd.
- 12 Awarded 'USD 50 Million Export Tower'

2017

- 12 2017 presidential commendation for job creation


2018

- 01 TS exceeded 10 million production

OSSTEM[®] Implant Design feature



OSSTEM IMPLANT has revolutionized implant dentistry in South Korea. With a focus on aggressive R&D, a commitment to education and a dedication to manufacturing the best products, Osstem Implant's ultimate goal is to become the global leader in implant dentistry.



TS

SS

US



Each implant system has its own unique color code

Submerged type implant with an internal hex and 11 tapered connection

- Internal connection type - Mini / Regular
- Excellent initial stability in soft bone due to smaller threads in the upper section
- Corkscrew thread with cutting edges
 - Strong self-threading effect for easy fixture path
 - Higher initial stability and consistent insertion torque
- Different body types to properly match the patient's bone quality and clinical condition
 - TSII (straight body) : easy to adjust depth
 - TSIII (1.5° tapered body) : excellent initial stability necessary for immediate loading, even in soft bone
 - TSIV (6° tapered body) : specifically designed for the maxillary sinus and soft bone, excellent initial stability
- Available surface types - SA / CA / HA / BA / SOI

Non-submerged type implant with an internal octa and 8 tapered connection

- Internal connection type - Regular / Wide
- Corkscrew thread with cutting edges
 - Strong self-threading effect for easy fixture path
 - Higher initial stability and consistent insertion torque
- Different body types to properly match the patient's bone quality and clinical condition
 - SSII (straight body) : easy to adjust the insertion depth
 - SSIII (1.5° tapered body) : excellent initial stability necessary for immediate loading, even in soft bone
- Available surface types - SA / CA / HA / BA

Submerged type implant with an external hex connection structure

- Internal connection type - Mini / Regular / Wide / Wide PS
- Corkscrew thread with cutting edges
 - Strong self-threading effect for easy fixture path
 - Higher initial stability and consistent insertion torque
- Different body types to properly match the patient's bone quality and clinical condition
 - USII (straight body) : easy to adjust the insertion depth
 - USIII (1.5° tapered body) : excellent initial stability necessary for immediate loading, even in soft bone
 - USIV (6° tapered body) : specifically designed for the maxillary sinus and soft bone, excellent initial stability
- Available surface types - SA / CA

OSSTEM[®] Implant Surface feature

The key factor in providing implant treatment safely and efficiently is surface technology.

OSSTEM IMPLANT is proud of its cutting-edge surface technology.

SA CA HA BA SOI

Acid Treated Optimized Surface

- Ra 2.5~3.0 μ m surface roughness (note : the upper 0.5mm part of the implant has Ra 0.5~0.6 μ m)
- Consistent surface micro pits between 1 to 3 μ m
- Surface area is increased by 46 percent compared to RBM treated implants

In-vitro & In-vivo Bone Response

- 20% improvement in osteoblast separation and ossification compared to RBM
- Initial bone reaction performance in animal model (mini-pig)
- 48% improvement in initial stability (RT, 4 weeks) compared to RBM
- 20% improvement in ossification (BIC, 4 weeks) compared to RBM

Super-hydrophilic SA surface suspended in a calcium solution

- Same SA surface morphology
- Optimizing surface reaction by suspension in a calcium (CaCl₂) solution
- Increased new bone formation area due to the excellent blood wettability
- Bone response improved in early osseointegration stage compared to standard SA surface

In-vitro & In-vivo Bone Response

- Protein and cellular adhesion tripled compared to SA surfaces
- Initial cellular differentiation by 19 percent compared to SA surfaces (7 days)
- Initial stability increased by 34 percent compared to SA surfaces (RT at 4 weeks)
- Ossification rate increased by 26 percent compared to SA surfaces (BIC at 4 weeks)

Premium high-crystalline HA-coated surface

- 30 to 60 μ m thick high-crystalline HA coating
- HA coated onto a RBM surface (Ra 3.0 to 3.5 μ m)
- High HA crystalline over 98 percent
- Solved the problem with low-crystalline HA resorption

In-vitro & In-vivo Bone Response

- Excellent biocompatibility in HA that is similar to bone
- Initial ossification by osteoblasts doubled compared to SA surfaces (5 days)
- 40% improvement in initial stability (RT, 4 weeks) in animal models compared to SA
- Suitable for poor bone quality, tooth extraction sites or immediate implant insertion

Premium low crystalline nano-HA coated SA surface

- SA surface (Ra 2.5 to 3.0 μ m) coated with HA
- 10nm ultra-thin HA coating
- Dual function between titanium and HA
- HA is naturally resorbed during ossification

In-vitro & In-vivo Bone Response

- Advantages of both SA and HA surfaces
- SA's ability to maintain an optimal surface
- HA's ability to form high quality initial bone, even in a poor bone quality
- 40% improvement in ossification (BIC) compared to SA
- It is applicable to all types of bone quality

Next-generation surface coated with special material (K material)

- Activation of blood clot formation
- Avoid carbon adsorption in air
- Coating of K material on SA surface (Ra 2.0~3.0 μ m)
- Superior blood wettability with super hydrophilic surface.

In-vitro & In-vivo Bone Response






- Protein and cellular adhesion 130 times increase compared to SA surface
- Initial stability increased by 57 percent compared to SA surfaces (RT at 4 weeks)
- Surface with the shortest duration of surgery






SS SYSTEM Contents













016 SSII SA Fixture  018 SSII CA Fixture  020 SSII BA Fixture  022 SSIII SA Fixture  024 SSIII CA Fixture 






053 ComOcta Gold Abutment  054 ComOcta NP-Cast Abutment  055 ComOcta Temporary Abutment  056 OneFit Abutment  056 Scan Baby 






026 SSIII BA Fixture  028 SSIII HA Fixture  030 Simple Mount  030 Cover Screw  031 Closing Screw 



057 Pre-Milled Abutment  058 ComOcta Angled Abutment  059 ComOcta Retraction Cap  059 Fixture Pick-up Impression Coping  060 Fixture Transfer Impression Coping 






032 Healing Abutment  035 Solid Abutment  036 Solid Protect Cap  036 Solid Retraction Cap  037 Solid Impression Coping 





060 Fixture Lab Analog  063 Octa Abutment  064 Octa Protect Cap  064 Octa Gold Cylinder  064 Octa Combination Cylinder 

037 Solid Lab Analog  037 Solid Burn-out Cylinder  038 Solid Impression Cap  038 Solid Shoulder Analog  038 Solid Shoulder Analog Pin 

065 Octa Temporary Cylinder  065 Octa Plastic Cylinder  066 Octa Pick-up Impression Coping  066 Octa Transfer Impression Coping  066 Octa Lab Analog 

039 Excellent Solid Abutment  040 Excellent Solid Protect Cap  040 Excellent Solid Retraction Cap  041 Excellent Solid Impression Coping  041 Excellent Solid Lab Analog 

069 O-ring Abutment  070 O-ring Retainer Cap Set  070 O-ring Retainer Set  070 O-ring Set  070 O-ring Lab Analog 

041 Excellent Solid Burn-out Cylinder  042 Excellent Solid Impression Cap  042 Excellent Solid Shoulder Analog  042 Excellent Solid Shoulder Analog Pin  045 ComOcta Abutment 

071 Locator® Abutment  072 Locator® Male Processing Kit  072 Locator® Replacement Male  072 Locator® Extended Replacement Male  073 Locator® Black Processing Male 

046 ComOcta ID Abutment  047 ComOcta Plus Abutment  049 ComOcta Plus ID Abutment  051 ComOcta Milling Abutment  052 ComOcta Milling ID Abutment 

073 Locator® Impression Coping  073 Locator® Lab Analog  074 Locator® Core Tool  074 Locator® Torque Driver  075 OneSeal 



SS SYSTEM

OSSTEM[®]
IMPLANT

FIXTURE

- 016** SSII SA Fixture
- 018** SSII CA Fixture
- 020** SSII BA Fixture
- 022** SSIII SA Fixture
- 024** SSIII CA Fixture
- 026** SSIII BA Fixture
- 028** SSIII HA Fixture
- 030** Simple Mount
- 030** Cover Screw
- 031** Closing Screw
- 032** Healing Abutment

COMPONENTS

- 034** PROSTHETIC FLOW DIAGRAM 1
- 035** Solid Abutment
- 039** Excellent Solid Abutment
- 044** PROSTHETIC FLOW DIAGRAM 2
- 045** ComOcta Abutment
- 047** ComOcta Plus Abutment
- 051** ComOcta Milling Abutment
- 053** ComOcta Gold Abutment
- 054** ComOcta NP-Cast Abutment
- 055** ComOcta Temporary Abutment
- 056** OneFit Abutment
- 057** Pre-Milled Abutment
- 058** ComOcta Angled Abutment
- 062** PROSTHETIC FLOW DIAGRAM 3
- 063** Octa Abutment
- 068** PROSTHETIC FLOW DIAGRAM 4
- 069** O-ring Abutment
- 071** Locator[®] Abutment
- 075** OneSeal

SSII SA Fixture


- Non-submerged type implant with an internal octa and 8° tapered connection
- Optimized screw thread design with the ideal SA surface
- Straight body design allows easy insertion depth adjustments
- Excellent initial stability in soft bone due to small threads in the upper section
- Corkscrew threading with excellent self-threading effect
- Recommended insertion torque : ≤40 Ncm
- ※ Fixtures with a diameter of 4.5mm or more are recommended for the posterior area

NoMount fixture order code

: fixture product code (ex : SS2R4011S18)

Pre-Mounted fixture (fixture + simple mount + cover screw) order code

: **A** + fixture product code (ex : ASS2R4011S18)



D Ø5.0
P Ø6.0
W

G/H \ L	6	7	8.5	10	11.5	13
Short						
2.0	SS2W5006S20	SS2W5007S20	SS2W5008S20	SS2W5010S20	SS2W5011S20	SS2W5013S20

016

D Ø4.0
P Ø4.8
R

G/H \ L	7	8.5	10	11.5	13
1.8	SS2R4007S18	SS2R4008S18	SS2R4010S18	SS2R4011S18	SS2R4013S18
2.8		SS2R4008S28	SS2R4010S28	SS2R4011S28	SS2R4013S28

D Ø4.5
P Ø4.8
R

G/H \ L	7	8.5	10	11.5	13
1.8	SS2R4507S18	SS2R4508S18	SS2R4510S18	SS2R4511S18	SS2R4513S18
2.8		SS2R4508S28	SS2R4510S28	SS2R4511S28	SS2R4513S28

D Ø4.5
P Ø6.0
W

G/H \ L	7	8.5	10	11.5	13
2.0	SS2W4507S20	SS2W4508S20	SS2W4510S20	SS2W4511S20	SS2W4513S20

017

SS SYSTEM

Nominal and actual diameters may slightly differ

Caution For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.

SSII CA Fixture

- Non-submerged type implant with an internal octa and 8° tapered connection
- Super-hydrophilic SA surface suspended in a calcium solution
- Straight body design allows easy insertion depth adjustments
- Excellent initial stability in soft bone due to small threads in the upper section
- Corkscrew threading with excellent self-threading effect


Ultra-wide

- Ideal for an extracted tooth site in the posterior area, for immediate placement, or for replacing a failed implant
- Apex is designed specially for excellent initial stability in an extracted tooth site
- Recommended insertion torque : ≤40 Ncm


※ Fixtures with a diameter of 4.5mm or more are recommended for the posterior area

NoMount fixture order code

: fixture product code (ex : SS2R4010C18)




D Ø5.0
P Ø6.0




Ultra-wide

D Ø6.0
P Ø6.0




D Ø7.0
P Ø6.0




Nominal and actual diameters may slightly differ

Caution For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.


D Ø4.0
P Ø4.8



D Ø4.5
P Ø4.8



D Ø4.5
P Ø6.0



SSII BA Fixture



- Non-submerged type implant with an internal octa and 8° tapered connection
- SA surface coated with premium low crystalline nano-HA
- No risk of fracturing or flaking off with Bio-absorbable coating
- Straight body design allows easy insertion depth adjustments
- Excellent self-threading effect with corkscrew threading

Ultra-wide


- Ideal for an extracted tooth site in the posterior area, for immediate placement, or for replacing a failed implant
- Apex is designed specially for excellent initial stability in an extracted tooth site
- Recommended insertion torque : ≤ 40 Ncm
- ※ Fixtures with a diameter of 4.5mm or more are recommended for the posterior single case

NoMount fixture order code


: fixture product code (ex : SS2R4011B18)

Pre-Mounted fixture (fixture + simple mount + cover screw) order code

: A + fixture product code (ex : ASS2R4011B18)




D $\varnothing 5.0$
P $\varnothing 6.0$
W




Ultra-wide

D $\varnothing 6.0$
P $\varnothing 6.0$
W




D $\varnothing 7.0$
P $\varnothing 6.0$
W




Nominal and actual diameters may slightly differ

Caution For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.


D $\varnothing 4.0$
P $\varnothing 4.8$
R



D $\varnothing 4.5$
P $\varnothing 4.8$
R



D $\varnothing 4.5$
P $\varnothing 6.0$
W



SSIII SA Fixture

- Non-submerged type implant with an internal octa and 8° tapered connection
- Optimized screw thread design with the ideal SA surface
- Tapered body design with high initial stability
- Excellent initial stability in soft bone due to the small thread on the upper part
- Corkscrew threading with excellent self-threading effect
- Excellent initial stability necessary for immediate loading, even in soft bone

Ultra-wide


- Ideal for an extracted tooth site in the posterior area, for immediate placement, or for replacing a failed implant
- Apex is specifically design for excellent initial stability in an extracted tooth site
- Recommended insertion torque : ≤ 40 Ncm
- ※ Fixtures with a diameter of 4.5mm or more are recommended for the posterior area

NoMount fixture order code


: fixture product code (ex : SS3R4011S18)

Pre-Mounted fixture (fixture + simple mount + cover screw) order code


: A + fixture product code (ex : ASS3R4011S18)



D $\varnothing 4.5$
P $\varnothing 6.0$




D $\varnothing 5.0$
P $\varnothing 6.0$




Ultra-wide

D $\varnothing 6.0$
P $\varnothing 6.0$




D $\varnothing 7.0$
P $\varnothing 6.0$




Nominal and actual diameters may slightly differ

Caution For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.


D $\varnothing 3.5$
P $\varnothing 4.8$



D $\varnothing 4.0$
P $\varnothing 4.8$



D $\varnothing 4.5$
P $\varnothing 4.8$



SSIII CA Fixture

- Non-submerged type implant with an internal octa and 8° tapered connection
- Tapered body design with high initial stability
- Corkscrew threading with excellent self-threading effect
- Excellent initial stability necessary for immediate loading, even in soft bone


Ultra-wide

- Ideal for an extracted tooth site in the posterior area, for immediate placement, or for replacing a failed implant
- Apex is designed specially for excellent initial stability in an extracted tooth site
- Recommended insertion torque : ≤40 Ncm

※ Fixtures with a diameter of 4.5mm or more are recommended for the posterior area

NoMount fixture order code

: fixture product code (ex : SS3R4011C18)



D Ø4.5
P Ø6.0



2.0 SS3W4507C20 SS3W4508C20 SS3W4510C20 SS3W4511C20 SS3W4513C20


D Ø5.0
P Ø6.0



2.0 SS3W5006C20 SS3W5007C20 SS3W5008C20 SS3W5010C20 SS3W5011C20 SS3W5013C20

Ultra-wide

D Ø6.0
P Ø6.0



2.0 SS3WB6006C20 SS3WB6007C20 SS3WB6008C20 SS3WB6010C20 SS3WB6011C20 SS3WB6013C20

D Ø7.0
P Ø6.0



2.0 SS3WB7006C20 SS3WB7007C20 SS3WB7008C20 SS3WB7010C20 SS3WB7011C20 SS3WB7013C20

Nominal and actual diameters may slightly differ

Caution For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.

024


SS SYSTEM

D Ø3.5
P Ø4.8




1.8 SS3R3508C18 SS3R3510C18 SS3R3511C18 SS3R3513C18
2.8 SS3R3508C28 SS3R3510C28 SS3R3511C28 SS3R3513C28

D Ø4.0
P Ø4.8



1.8 SS3R4007C18 SS3R4008C18 SS3R4010C18 SS3R4011C18 SS3R4013C18
2.8 SS3R4008C28 SS3R4010C28 SS3R4011C28 SS3R4013C28

D Ø4.5
P Ø4.8



1.8 SS3R4507C18 SS3R4508C18 SS3R4510C18 SS3R4511C18 SS3R4513C18
2.8 SS3R4508C28 SS3R4510C28 SS3R4511C28 SS3R4513C28

025

SS SYSTEM


SSIII BA Fixture



- Non-submerged type implant with an internal octa and 8° tapered connection
- SA surface coated with premium low crystalline nano-HA
- No risk of fracturing or flaking off with Bio-absorbable coating
- Tapered body design with high initial stability
- Excellent self-threading effect with corkscrew threading
- Excellent initial stability necessary for immediate loading, even in soft bone

Ultra-wide

- Ideal for an extracted tooth site in the posterior area, for immediate placement, or for replacing a failed implant
 - Apex is designed specially for excellent initial stability in an extracted tooth site
 - Recommended insertion torque : ≤ 40 Ncm
- ※ Fixtures with a diameter of 4.5mm or more are recommended for the posterior single case



NoMount fixture order code

: fixture product code (ex : SS3R4011B18)

Pre-Mounted fixture (fixture + simple mount + cover screw) order code

: A + fixture product code (ex : ASS3R4011B18)

D \varnothing 3.5
P \varnothing 4.8



G/H \ L	8.5	10	11.5	13
1.8	SS3R3508B18	SS3R3510B18	SS3R3511B18	SS3R3513B18
2.8	SS3R3508B28	SS3R3510B28	SS3R3511B28	SS3R3513B28

D \varnothing 4.0
P \varnothing 4.8



G/H \ L	7	8.5	10	11.5	13
1.8	SS3R4007B18	SS3R4008B18	SS3R4010B18	SS3R4011B18	SS3R4013B18
2.8	SS3R4008B28	SS3R4010B28	SS3R4011B28	SS3R4013B28	

D \varnothing 4.5
P \varnothing 4.8



G/H \ L	7	8.5	10	11.5	13
1.8	SS3R4507B18	SS3R4508B18	SS3R4510B18	SS3R4511B18	SS3R4513B18
2.8	SS3R4508B28	SS3R4510B28	SS3R4511B28	SS3R4513B28	

D \varnothing 4.5
P \varnothing 6.0



G/H \ L	7	8.5	10	11.5	13
2.0	SS3W4507B20	SS3W4508B20	SS3W4510B20	SS3W4511B20	SS3W4513B20

D \varnothing 5.0
P \varnothing 6.0



G/H \ L	6	7	8.5	10	11.5	13
2.0	Short SS3W5006B20	SS3W5007B20	SS3W5008B20	SS3W5010B20	SS3W5011B20	SS3W5013B20

Ultra-wide

D \varnothing 6.0
P \varnothing 6.0



G/H \ L	6	7	8.5	10	11.5	13
2.0	Short SS3W6006B20	SS3W6007B20	SS3W6008B20	SS3W6010B20	SS3W6011B20	SS3W6013B20

D \varnothing 7.0
P \varnothing 6.0



G/H \ L	6	7	8.5	10	11.5	13
2.0	Short SS3W7006B20	SS3W7007B20	SS3W7008B20	SS3W7010B20	SS3W7011B20	SS3W7013B20


Nominal and actual diameters may slightly differ

Caution For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.

SSIII HA Fixture

- Non-submerged type implant with an internal octa and 8° tapered connection
- Premium high-crystalline HA-coated surface
- Tapered body design with high initial stability
- Corkscrew threading with excellent self-threading effect
- Excellent initial stability necessary for immediate loading, even in soft bone
- Recommended insertion torque : ≤ 40 Ncm
- ※ Fixtures with a diameter of 4.5mm or more are recommended for the posterior area
- ※ Recommended to use carefully, as the high-crystalline HA coating has a risk of fracturing or flaking off

Pre-Mounted fixture (fixture + simple mount + cover screw) **order code**
: **A** + fixture product code (ex : **ASS3R4011H18**)



D $\varnothing 5.0$
P $\varnothing 6.0$
W

G/H \ L	6	7	8.5	10	11.5	13
	Short					
2.0	SS3W5006H20	SS3W5007H20	SS3W5008H20	SS3W5010H20	SS3W5011H20	SS2W5013H20

028

D $\varnothing 4.0$
P $\varnothing 4.8$

R

G/H \ L	7	8.5	10	11.5	13
1.8	SS3R4007H18	SS3R4008H18	SS3R4010H18	SS3R4011H18	SS3R4013H18
2.8		SS3R4008H28	SS3R4010H28	SS3R4011H28	SS3R4013H28

D $\varnothing 4.5$
P $\varnothing 4.8$

R

G/H \ L	7	8.5	10	11.5	13
1.8	SS3R4507H18	SS3R4508H18	SS3R4510H18	SS3R4511H18	SS3R4513H18
2.8		SS3R4508H28	SS3R4510H28	SS3R4511H28	SS3R4513H28

D $\varnothing 4.5$
P $\varnothing 6.0$

W

G/H \ L	7	8.5	10	11.5	13
2.0	SS3W4507H20	SS3W4508H20	SS3W4510H20	SS3W4511H20	SS3W4513H20

029

SS SYSTEM

Nominal and actual diameters may slightly differ

Caution For a short implant, a sufficient healing period is strongly recommended. A short implant should be splinted with another implant when considering prosthetic options.


Mount & Screw

Simple Mount

- Select appropriate mount according to the fixture platform
- Use a 1.2 hex driver (fasten manually)
- ※ Disposable; do not re-use
- P = Platform

R Regular

W Wide



Closing Screw

- Used when the soft tissue of the suture part is insufficient
- Use of 1.2 hex driver with hand force
- P = Platform

R Regular

W Wide




Cover Screw

- Select appropriate mount according to the fixture platform
- Use a 1.2 hex driver (fasten manually)
- P = Platform

R Regular

W Wide



Healing Abutment

- Select appropriate mount according to the fixture platform
- Use a 1.2 hex driver (fasten manually)
- P = Platform

R Regular
W Wide



OSSTEM[®]
 IMPLANT

032

P \ H	2.0	3.0	4.0	5.0
Ø 4.8	SSH482	SSH483	SSH484	SSH485

P \ H	2.0	3.0	4.0	5.0
Ø 6.0	-	SSH603	SSH604	SSH605

033

SS SYSTEM

SS SYSTEM

Solid / Excellent Solid


Abutment Level Impression




Solid Abutment

- Cement-retained prosthesis
- Abutment level impression
- $\varnothing 4.8$: connect using a solid abutment driver (code : SDSL/SDSS)
- $\varnothing 6.0$: connect using a 1.2 hex driver or solid abutment driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + protect cap

Abutment + protect cap order code
: product code + **P** (ex : SSS485**P**)



034



D $\varnothing 4.8$



035

SS SYSTEM

D $\varnothing 6.0$



Solid Abutment Components

Solid Protect Cap

- Protects the solid abutment and minimizes patient irritation
- Can be used as the base for a provisional crown

R Regular

W Wide

D \ H	4.0	5.5	7.0
Ø 4.8	SSC484	SSC485	SSC487
Ø 6.0	SSC604	SSC605	SSC607

Solid Impression Coping

- Components for solid abutment impression
- Possibility of precise prosthesis using lab analog
- Color coded by abutment height

R Regular

W Wide

D \ H	4.0	5.5	7.0
Ø 4.8	SSIC484	SSIC485	SSIC487
Ø 6.0	SSIC604	SSIC605	SSIC607

Solid Retraction Cap

- Possible to get clear margin by pushing out gingiva around margin in solid abutment direct impressions
- Used for accurate margin reproduction when taking a direct impression
- Can be used as the base for a provisional crown

R Regular

W Wide

D \ H	4.0	5.5	7.0
Ø 4.8	SSSRC484	SSSRC485	SSSRC487
Ø 6.0	SSSRC604	SSSRC605	SSSRC607

Solid Lab Analog

- Components that replace resin caps before wax up using solid abutments
- Used in the same color as solid impression coping

R Regular

W Wide

D \ H	4.0	5.5	7.0
Ø 4.8	SSSA484	SSSA485	SSSA487
Ø 6.0	SSSA604	SSSA605	SSSA607

Solid Burn-out Cylinder

- Solid abutment components that reproduce them on the model after impression taking
- Sophisticated prosthesis can be produced inside
- After casting, remove the lower part of the margin holding part

R Regular

W Wide

D \ H	Single	Bridge
Ø 4.8	SSSP480S	SSSP480B
Ø 6.0	SSSP600S	SSSP600B

Solid Abutment Components



Solid Impression Cap

- An impression cap used when the solid abutment is trimmed
- Used with a solid shoulder analog and analog pin


R Regular

W Wide



Excellent Solid Abutment

- Cement-retained prosthesis
- Ideal for molar cases due to its larger volume (compared to the solid abutment), trim as needed
- Abutment level impression
- Ø 4.8 : connect using a 1.2 hex driver or an excellent solid abutment driver (code : ESDSS/ESDSL)
- Ø 6.0 : connect using a 1.2 hex driver or an excellent solid abutment driver (code : ESD60S)
- Recommended tightening torque: 30Ncm
- Packing unit : abutment + protect cap




Abutment + protect cap order code
: product code + **P** (ex : SSE485**P**)

Solid Shoulder Analog

- Impression product used when removing solid abutment
- Reproducing the platform part of the fixture in the work model
- Used with solid impression cap and shoulder analog pin


R Regular

W Wide



D Ø 4.8

R




Solid Shoulder Analog Pin

- An impression coping component used when the solid abutment is trimmed
- Reinforces the narrow part of the abutment
- Used with a solid shoulder analog and impression cap


R Regular

W Wide



D Ø 6.0

W









Excellent Solid Abutment Components

Excellent Solid Protect Cap

- Protects the solid abutment and minimizes patient irritation
- Can be used as the base for a provisional crown

R Regular

W Wide







D \ H	4.0	5.5	7.0
Ø 4.8	 SSEC484	 SSEC485	 SSEC487
Ø 6.0	 SSEC604	 SSEC605	 SSEC607

Excellent Solid Impression Coping

- Impression components for Excellent solid abutment
- Possibility of precise prosthesis using lab analog
- Color coded by abutment height

R Regular

W Wide







D \ H	4.0	5.5	7.0
Ø 4.8	 SSEIC484	 SSEIC485	 SSEIC487
Ø 6.0	 SSEIC604	 SSEIC605	 SSEIC607

Excellent Solid Retraction Cap

- Used for accurate margin reproduction when taking a direct impression
- Possible to get clear margin in excellent solid abutment direct impressions
- Can be used as the base for a provisional crown

R Regular

W Wide

D \ H	4.0	5.5	7.0
Ø 4.8	 SSERC484	 SSERC485	 SSERC487
Ø 6.0	 SSERC604	 SSERC605	 SSERC607

Excellent Solid Lab Analog

- Components that replace resin caps before wax up using excellent solid abutments
- Connect to the appropriate color coded with solid impression coping

R Regular

W Wide

D \ H	4.0	5.5	7.0
Ø 4.8	 SSEA484	 SSEA485	 SSEA487
Ø 6.0	 SSEA604	 SSEA605	 SSEA607

Excellent Solid Burn-out Cylinder

- Excellent solid abutment components that reproduce this on the model after impression taking
- Sophisticated prosthesis can be produced inside
- After casting, remove the lower part of the margin holding part

R Regular

W Wide

D \ H	Single	Bridge
Ø 4.8	 SSEP480S	 SSEP480B
Ø 6.0	 SSEP600S	 SSEP600B

Excellent Solid Abutment Components

Excellent Solid Impression Cap

- An impression cap used when the solid abutment is trimmed
- Used with an excellent solid shoulder analog and an analog pin

R Regular

W Wide

D \

Ø 4.8
Ø 6.0



SSEIP480
SSEIP600

Excellent Solid Shoulder Analog

- Impression product used when removing excellent solid abutment
- Reproducing the platform part of the fixture in the work model
- Used with excellent solid impression cap and shoulder analog pin

R Regular

W Wide

D \

Ø 4.8
Ø 6.0



SSSLA480
SSSLA600

Excellent Solid Shoulder Analog Pin

- An impression coping component used when the solid abutment is trimmed
- Reinforces the narrow part of the abutment
- Used with an excellent solid shoulder analog and an impression cap

R Regular

W Wide

D \

Ø 4.8
Ø 6.0




SSSAP480

OSSTEM[®]
IMPLANT


ComOcta / OneFit

Fixture Level Impression






ComOcta Abutment

- Cement/combination-retained prosthesis
- Fixture level impression
- Able to take abutment level impression using a retraction cap
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw



Abutment + Ti screw order code
: product code + TH (ex : SSCA485TH)



D Ø4.8	H 4.0			H 5.5			H 7.0		
	Type			Type			Type		
R	Octa			Octa			Non-Octa		
	SSCA484	SSCA485	SSCA487	SSCA484N	SSCA485N	SSCA487N			

D Ø6.0	H 4.0			H 5.5			H 7.0		
	Type			Type			Type		
W	Octa			Octa			Non-Octa		
	SSCA604	SSCA605	SSCA607	SSCA604N	SSCA605N	SSCA607N			

044

SS SYSTEM

045

SS SYSTEM


ComOcta ID Abutment

ComOcta Plus Abutment




- Non-insured ComOcta abutment
- Easy to identify abutment platform by applying color coding

Abutment + Ti screw order code
: product code + **TH** (ex : BSSCA485**TH**)



- Cement / combination abutment for prosthesis manufacturing
- Used when the gingiva is thick or the fixture is deeply implanted
- Abutment - fixture interlock is 45° platform contact
- Fixture level impression
- Use of 1.2 hex driver
- Recommended torque of tightening screw : 30Ncm
- Packing unit : abutment + Ti screw

Abutment + Ti screw order code
: product code + **TH** (ex : SSCAP4826**CTH**)



P Ø4.8



Ti screw : ASR200

H	4.0	5.5	7.0	4.0	5.5	7.0
Type	Octa			Non-Octa		
	BSSCA484	BSSCA485	BSSCA487	BSSCA484N	BSSCA485N	BSSCA487N

P Ø4.8



Ti screw : ASR200

G/H	1.0	2.0	3.0	4.0
D	Ø5.5	Ø5.5	Ø5.5	Ø5.5
Octa	SSCAP4816C	SSCAP4826C	SSCAP4836C	SSCAP4846C
Non-Octa	SSCAP4816CN	SSCAP4826CN	SSCAP4836CN	SSCAP4846CN

047

P Ø6.0



Ti screw : ASR200

H	4.0	5.5	7.0	4.0	5.5	7.0
Type	Octa			Non-Octa		
	BSSCA604	BSSCA605	BSSCA607	BSSCA604N	BSSCA605N	BSSCA607N

G/H	2.0	3.0	4.0
D	Ø6.0	Ø6.5	Ø7.0
Octa	SSCAP4826E	SSCAP4836E	SSCAP4846E
Non-Octa	SSCAP4826EN	SSCAP4836EN	SSCAP4846EN

ComOcta Plus Abutment

P Ø6.0



Ti screw
: ASR200


	G/H	1.0	2.0	3.0	4.0
D		Ø6.5	Ø6.5	Ø6.5	Ø6.5
Octa					
Non-Octa		SSCAP6016CN	SSCAP6026CN	SSCAP6036CN	SSCAP6046CN

	G/H	2.0	3.0	4.0
D		Ø6.8	Ø7.2	Ø7.6
Octa				
Non-Octa		SSCAP6026EN	SSCAP6036EN	SSCAP6046EN

ComOcta Plus ID Abutment

• Non-insured ComOcta abutment

Abutment + Ti screw order code
: product code + TH (ex : BSSCAP4826C**TH**)



P Ø4.8



Ti screw
: ASR200

	G/H	1.0	2.0	3.0	4.0
D		Ø5.5	Ø5.5	Ø5.5	Ø5.5
Octa					
Non-Octa		BSSCAP4816CN	BSSCAP4826CN	BSSCAP4836CN	BSSCAP4846CN

	G/H	2.0	3.0	4.0
D		Ø6.0	Ø6.5	Ø7.0
Octa				
Non-Octa		BSSCAP4826EN	BSSCAP4836EN	BSSCAP4846EN

ComOcta Plus ID Abutment

P Ø6.0



Ti screw : ASR200


	G/H	1.0	2.0	3.0	4.0
D		Ø6.5	Ø6.5	Ø6.5	Ø6.5
Octa					
Non-Octa		BSSCAP6016CN	BSSCAP6026CN	BSSCAP6036CN	BSSCAP6046CN

	G/H	2.0	3.0	4.0
D		Ø6.8	Ø7.2	Ø7.6
Octa				
Non-Octa		BSSCAP6026EN	BSSCAP6036EN	BSSCAP6046EN

ComOcta Milling Abutment

- Cement/combination-retained prosthesis
- Used when the abutment's margin shape needs corrections
- Abutment fastens to the platform at a 45° angle
- Fixture level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw

Abutment + Ti screw order code
: product code + TH (ex : SSCMA4830TH)



D Ø4.8



Ti screw : ASR200

G/H

3.0



SSCMA4830

D Ø6.0



Ti screw : ASR200

G/H

3.0



SSCMA6030

050

SS SYSTEM

051


SS SYSTEM

ComOcta Milling ID Abutment

ComOcta Gold Abutment




• Non-insured ComOcta milling abutment
Abutment + Ti screw order code
 : product code + **TH** (ex : BSSCMA4830**TH**)



- Cement/combination/screw-retained prosthesis
- Customized prosthesis cast with gold alloy
- Abutment fastens to the platform at a 45° angle
- Abutment melting point : 1400-1450°C
- Fixture level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw

Abutment + Ti screw order code
 : product code + **TH** (ex : COG480**STH**)



P Ø4.8 \ G/H **3.0**

R

Ti screw : ASR200



BSSCMA4830

P Ø6.0 \ G/H **3.0**

W

Ti screw : ASR200



BSSCMA6030

P Ø4.8 \ Type **Octa** **Non-Octa**

R

Ti screw : ASR200






COG480S **COG480B**

P Ø6.0 \ Type **Octa** **Non-Octa**

W

Ti screw : ASR200





COG600S **COG600B**

ComOcta NP-Cast Abutment

- Cement/combination/screw-retained prosthesis
- Customized prosthesis cast with non-precious alloys
- Abutment fastens to the platform at a 45° angle
- Abutment melting point : 1400~1450°C
- Fixture level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw


Abutment + Ti screw order code
: product code + TH (ex : CON480STH)



ComOcta Temporary Abutment

- Cement/screw-retained prosthesis
- A trim able provisional prosthesis (made of Ti Gr-3)
- Fixture level impression
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : abutment + Ti screw

Abutment + Ti screw order code
: product code + TH (ex : SSTA0480TH)



054

D Ø4.8 \ Type Octa Non-Octa

R

Ti screw : ASR200

CON480S CON480B

D Ø6.0 \ Type Octa Non-Octa

W

Ti screw : ASR200

CON600S CON600B

055

D Ø4.8 \ G/H Type 0 2.0 \ 0 2.0

R

Ti screw : ASR200

Octa Non-Octa

SSTA0480 SSTA0482 SSTAN480 SSTAN482

D Ø6.0 \ G/H Type 0 2.0 \ 0 2.0

W

Ti screw : ASR200

Octa Non-Octa

SSTA0600 SSTA0602 SSTAN600 SSTAN602

SS SYSTEM

SS SYSTEM

OneFit Abutment


- Cement/combination-retained prosthesis
- CAD/CAM designed and milled customized abutments
- Fixture level impression
- Possibility of abutment level impression when using Scan healing abutment
- Lead time (by working days)
 - Titanium : 5 days
 - Titanium + gold color : 7 days
- Use a 1.2 hex driver
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + Ti screw



Pre-Milled Abutment

- Manufacturing custom abutment with dental milling equipment
- Excellent tightening precision compared to unauthentic
- Packing unit : abutment + Ti screw

Pre-milled abutment + screw order code
: product code + **TH** (ex : SSPM10AGR**TH**)



Scan Body

- Scan body for manufacturing a titanium SmartFit abutment
- Use a 1.2 hex driver (fastened manually)
- Packing unit : scan body + Ti screw

Scan body + screw order code
: product code + **TH** (ex : SSSB**TH**)

- R** Regular
- W** Wide



Green color screw
: USSBRS



SSSBR



Green color screw
: USSBRS




SSSBW

Equipment	Implant	D	Specifications		Code
DooWon ARUM	Osstem SS	Ø10	Regular	Octa	SSPM10AGRTH
			Regular	Non-octa	SSPM10AGRNTH
			Wide	Octa	SSPM10AGWTH
			Wide	Non-octa	SSPM10AGWNTH

ComOcta Angled Abutment



- Cement/combination-retained prosthesis
- Angle compensation between 15°/20°
- Use dedicated abutment screw
- Fixture level impression
- Fastened using a 1.2 hex driver
- Recommended tightening torque: 30Ncm
- Packing unit : abutment + Ti screw(only angled)



Abutment + Ti screw order code
: product code + **TH** (ex : SSA4815**TH**)

D Ø4.8



Ti screw : ASS200

Angle Type	15°	20°	15°	20°
Type	Octa		Non-Octa	
	SSA4815	SSA4820	SSA4815N	SSA4820N

D Ø6.0



Ti screw : ASS200

Angle Type	15°	20°	15°	20°
Type	Octa		Non-Octa	
	SSA6015	SSA6020	SSA6015N	SSA6020N

ComOcta Abutment Components

ComOcta Retraction Cap

- Used for accurate margin reproduction when taking a direct impression
- Can be used as the base for a provisional crown

R Regular

W Wide


D \ H	4.0	5.5	7.0
Ø 4.8	 SSCRC484	 SSCRC485	 SSCRC487
Ø 6.0	 SSCRC604	 SSCRC605	 SSCRC607

Fixture Pick-up Impression Coping

- Components for fixture level impression taking
- For open tray impressions
- Unique design that is fixed position in the impression material
- Use a 1.2 hex driver (torque manually)
- Packing unit : impression coping body + guide pin(*)

R Regular (Green)

W Wide (Blue)




P \ L	5		10		Guide Pin		
Type	Octa	Non-Octa	Octa	Non-Octa	10	15	17
Ø 4.8							
Ø 6.0							

ComOcta Abutment Components

Fixture Transfer Impression Coping

- Components for fixture level impression taking
- For closed tray impressions
- Triangular arc enabling precise placement
- Use a 1.2 hex driver (torque manually)
- Packing unit
 - Octa : impression coping body + guide pin
 - Non-octa : impression coping



- R** Regular (Green)
- W** Wide (Blue)

P \ L	Type			
	9.5		12.5	
	Octa		Non-Octa	
	SSCTIS480	SSCTIL480	SSCTIS480N	SSCTIL480N
	SSCTIS600	SSCTIL600	SSCTIS600N	SSCTIL600N

Ø 4.8
Ø 6.0

OSSTEM[®]
IMPLANT

061

SS SYSTEM

Fixture Lab Analog

- A lab analog for fixture level impression
- Select an appropriate fixture platform; Ø 4.8/6.0

- R** Regular (Silver)
- W** Wide (Blue)

D \



Ø 4.8
Ø 6.0


SSFA480
SSFA600

060

SS SYSTEM


Octa

Abutment Level Impression



Octa Abutment

- Screw-retained prosthesis for multiple prosthetic options
- Angle compensation of up to 60°
- Use a dedicated outer driver (code : ODSL/ODSS)
- Recommended tightening torque : 30Ncm



D Ø4.8



SSOA480

D Ø6.0



SSOA600

063

SS SYSTEM

062

SS SYSTEM

Octa Abutment Components

Octa Protect Cap

- Protective cap
- Use a 1.2 hex driver (fastened manually)
- Packing unit : protect cap + Ti screw

Protect cap + Ti screw order code
: product code + **TH** (ex : SSHC480**TH**)

- R** Regular
- W** Wide




Octa Temporary Cylinder

- Provisional prosthesis (Ti Gr-3)
- Fastened using a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti screw order code
: product code + **TH** (ex : SSTCO480**TH**)

- R** Regular
- W** Wide




Octa Gold Cylinder

- Screw-retained prosthesis
- Customized prosthesis cast with gold alloy
- Cylinder melting point : 1400~1450°C
- Fastened using a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti screw order code
: product code + **TH** (ex : SSGCO480**TH**)

- R** Regular
- W** Wide




Octa Plastic Cylinder

- Screw-retained prosthesis
- Customized prosthesis cast with non-precious alloys
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti screw order code
: product code + **TH** (ex : SSPSO480**TH**)

- R** Regular
- W** Wide




Octa Combination Cylinder

- Combination-retained prosthesis
- Compatible with bone octa/non-octa specs
- Use a 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti cylinder screw

Cylinder + Ti screw order code
: product code + **TH** (ex : SSOCC480**TH**)

- R** Regular
- W** Wide




Octa Abutment Components

Octa Pick-up Impression Coping

- A pick up impression coping for octa abutment
- Use a 1.2 hex driver (fastened manually)
- Packing unit : impression coping body + guide pin(*)

- R** Regular (Silver)
- W** Wide (Blue)



D \ L	Coping		Guide Pin	
	Octa	Non-Octa	10	15
Ø 4.8	SSICO480	SSICN480N	SSGS100	SSGS150*
Ø 6.0	SSICO600	SSICN600N		



Octa Transfer Impression Coping

- A pick-up impression coping
- Use a 1.2 hex driver (fastened manually)
- Packing unit : Impression coping body + guide pin

- R** Regular (Silver)
- W** Wide (Blue)

D	Transfer Coping
Ø 4.8	SSOTI480
Ø 6.0	SSOTI600

Octa Lab Analog

- A lab analog
- Use a 1.2 hex driver (fastened manually)


- R** Regular (Silver)
- W** Wide (Blue)

D	Lab Analog
Ø 4.8	SSLA480
Ø 6.0	SSLA600


O-ring / Locator®

Overdenture

O-ring Abutment



- Retains overdenture with o-ring system
- Angle compensation of up to 20°
- Torque using an outer driver (code : AORD)
- Recommended tightening torque : 30Ncm



D Ø4.8

R



068

SS SYSTEM

069

SS SYSTEM

O-ring Abutment Components



O-ring Retainer Cap Set

- O-ring housing
- Place an appropriate o-ring in the metal housing before connecting to the abutment
- Packing unit : retainer cap + o-ring



O-ring Retainer Set

- Used when vertical dimension is shorter than the retainer cap
- Packing unit : retainer cap + o-ring



O-ring Set

- O-ring set
- Packing unit : o-ring x 5ea




O-ring Lab Analog

- A lab analog for o-ring abutment



Locator[®] Abutment

- Genuine zest anchors abutment
- Angle compensation of up to 40°
- 1.5mm lower profile, attachment with various and stable retention forces
- Torque using a dedicated outer driver (code : TWLDLK/TWLDLSK)
- Recommended tightening torque : 30Ncm



D Ø4.8

R

G/H	0.7	2.0	3.0	4.0
	HSLCA4810R	HSLCA4820R	HSLCA4830R	HSLCA4840R

Locator® Abutment Components

Locator® Male Processing Kit

- Components
 - Block out spacer / denture cap connected black processing male
 - Replacement male blue/pink/clear
- A full range of retentive males are included with each denture cap to allow personalized retention for each specific patient
- Locator core tool places and removes nylon retentive males
- Packing unit : 2set



Locator® Black Processing Male

- A nylon male used in prosthesis fabrication process
- Packing unit : 4ea



Locator® Replacement Male

- Retention force : approx. 6N
 - Angle compensation up to 20°
 - Packing unit : 4ea
- Retention force : approx. 12N
 - Angle compensation up to 20°
 - Packing unit : 4ea
- Retention force : approx. 22N
 - Angle compensation up to 20°
 - Packing unit : 4ea



Locator® Block Out Spacers

- Place block-out spacers on the heads of the locator abutments. Position denture cap with integrated black processing onto the locator abutments. If necessary, add additional block-out spacers until no gap is visible between female, block-out spacer and gum.
- Packing unit : 20ea




Locator® Impression Coping

- A pick up impression coping
- Closed tray
- Packing unit : 4ea



Locator® Extended Replacement Male

- Retention force : approx. 6N
 - Angle compensation up to 20~40°
 - Packing unit : 4ea
- Retention force : approx. 12N
 - Angle compensation up to 20~40°
 - Packing unit : 4ea



Locator® Lab Analog

- A lab analog for locator abutment
- Packing unit : 4ea



Locator® Abutment Components

Locator® Core Tool

- Places and removes nylon retentive males in the denture cap
- Separated into three different tools, includes a hand driver for locator abutment



Locator® Torque Driver


- A torque driver for locator abutment



OneSeal

OneSeal

- Disposable medical devices for internal filling of superstructure
- Cut to desired length (medical silicone)
- Packing unit : long 5ea



Instructions for Use (AUG. 2017, Ver. 5.5)

Description of Osstem implant system

Osstem Implant is a brand for implant materials for dental practices, and the fixture is made mainly of titanium. The abutment, prosthetic components and tools for the Osstem Implant system are compatible with the Osstem Implant fixture only. Using this product in combination with products from other manufacturers may cause various problems including loosening and fracture due to incomplete locking and compatibility issues. Refer to the manual or the catalogue or our website (www.osstem.com) for details. See the product label for the product code, specifications, manufacturing date, and expiration date.

Sterility

The fixture, cover screw, and healing abutment are cleansed and sterilized with gamma radiation. This product is a disposable sterilized medical device intended for one-time use. In order to prevent contamination or infection of the product or operated site, the product must be used using a sterilized instrument in a sterilized environment. Damaged products, products with open packaging, or expired products must be discarded due to potential risks of contamination, infection, or osseointegration failure. Re-sterilization or re-use of the product may result in infection, osseointegration failure, or implant damage due to reduced accuracy.

Storage condition

Keep the product in a dry place at room temperature(1~30°C). Keep away from direct sunlight.

General precautions

The surgical technology of dental implant involves an expert, complex procedure. Formal training is required to perform implant surgery. Careful considerations must be made before the operation in case of bone disorders (osteoporosis, osteomalacia) or metabolic disorders of the bone.

Precautions

Determine the local anatomy and suitability of the available bone for implant placement. Prepare the implant considering the expected situations and cautions. Excessive occlusal load may cause loosening or fracture of an implant. In order to avoid this condition, the implant must be placed in accurate location and direction considering the relationship between the implant and opposing dentition. Visual inspection as well as panoramic and periapical radiographs are essential to determine anatomical landmarks, occlusal conditions, periodontal status, and the adequacy of the bone. Adequate radiographs, direct palpation, and visual inspection of the implant site are necessary prior to implant surgery.

Procedural precautions

Osstem Implant System is for single and two stage surgical procedures. As much as possible, try to minimize damage to the cell tissue and surgical trauma, pay special attention to maintaining the temperature at the implant site and removal of the source of contamination and infection. All drills and taps must be sufficiently and continuously irrigated for cooling during use. Implant placement should be accomplished at very low speed (25-30 rpm) or manually. Excessive torque (greater than 55Ncm) in the fixture placement can have adverse effects such as partial fracture or necrosis of the bone. Placing an implant tilted by 30° or higher is not recommended due to possible fracture of implant. Immediate loading to the fixture right after the surgery should be avoided. The bone quality and initial stability after fixture placement are important elements in determining the appropriate loading time. Mini-diameter implant or implant with diameter of 4.0 or less and which integrates with angled abutment may be fractured due to limitations of structural rigidity. They are not recommended for use in a posterior area. The Ultra-Wide fixtures are intended to be used only to replace molar teeth and

that angled abutments are not to be used with the Ultra-Wide fixtures. Evaluate the quantity of bone and radiographs to assess any potential anatomical contraindications to use of the Ultra-Wide fixture. For the placement of the Short Implant (diameter is 5mm or more and length is shorter than 7mm) which is used on the molar region only, clinicians should closely examine the patients for any of the following conditions: 1) periimplant bone loss, 2) changes to implant's response to percussion, 3) radiographic changes in bone to implant contact along the implant's length. If a short implant shows mobility or greater than 50% bone loss, the implant should be considered for possible removal. And clinicians should consider a two-stage surgical approach, splinting a short implant to an additional implant, and placement of the widest possible fixture. Allow longer healing periods for osseointegration before fabrication of the prosthesis and avoid immediate loading. Products with diameter of 3.25mm or less must be used exclusively for mandibular anterior teeth in order to prevent fracture due to excessive occlusal load. It is recommended that you should avoid applying HA coated fixture to hard bone, and the insertion torque of the implant should be less than 35Ncm, because cracks or damages might occur in the coated layer during implant placement. The surfaces of CA and SOI have the same physical shape as the SA surface made through blasting and etching treatments. After the SA surface treatment, to prevent the products' exposure to the atmosphere, CA is stored in solution, whereas SOI is stored in water-film coating form; it is designed to maintain the chemically activated state of the SA surface. Thus, CA or SOI products should be implanted in the target region at least within 15 minutes of taking them out of the container.

Warning

The selection of inappropriate patients and surgical methods can cause implant failure or loss of bone supporting the implant. Osstem implants must not be used for purposes other than the recommended use and must not be remodeled. Implant mobility, bone loss, and chronic infection can result in failure of the implant surgery.

Indications for use

The Osstem Implant System is an artificial dental root that has been designed for use in dental implant treatment in order to recover lost teeth. The system is implanted via a surgical method in maxillary or mandibular bone to replace natural dental root. The Osstem Implant System is indicated for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-units restorations including; cemented retained, screw retained, or overdenture restorations, and final or temporary abutment support for fixed bridgework. It is intended for delayed loading. Products with diameter of 3.25mm or less must be used exclusively for mandibular anterior teeth in order to prevent fracture due to excessive occlusal load.

Side effects

A few problems may occur after the operation (loss of implant stability, damage of prosthesis, etc.). Deficient quality and quantity of the remaining bone, infection, allergic reaction, inferior oral hygiene or uncooperativeness of patient, implant mobility, partial deterioration of tissue, and improper position or arrangement of implants may cause the above mentioned problems.

Contraindications

- Contraindications include the following, but are not limited to:
- Patients with hemophilia or difficulties related to bone or wound treatment
 - Patients with uncontrollable diabetes, heavy smoker or alcoholic
 - Patients whose immunity system is inactive due to chemical therapy or radiation therapy
 - Patients with oral infection or inflammation (improper oral hygiene, bruxism)
 - Patients with untreatable occlusion/joint disorder, insufficient dental arch space
 - Any patient who is not suitable for an surgery

OSSTEM[®]
IMPLANT

Manufacturer : Osstem Implant Co., Ltd.
203, Geoje-daero, Yeonje-gu, Busan, Korea
TEL 82-51-850-2500 FAX 82-51-861-4693



2460



Sterilized using irradiation



Use by



Manufacture



DEUTSCHE OSSTEM GmbH.
Mergenthalerallee 25
65760 Eschborn, Germany
+49-(0)6196-777-550



Do not reuse



Date of manufacture



Keep away from sunlight



Catalogue number



Non-Sterile



Keep dry

Rx only

For USA only : Federal law restricts this device to sale by or on the order of a dentist



Batch code



Do not resterilize



Caution, Consult accompanying documents

