


Hands-on Guide

Osstem Implant KIT

OSSTEM[®]
IMPLANT

OSSTEM[®]
IMPLANT






Hands-on Guide

- | | |
|---|---|
| 1 | Short Implant Solution
485 KIT |
| 2 | Simple Drilling Solution
122 Taper KIT |
| 3 | Implant Aligning Solution
Positioning Guide KIT |
| 4 | Sinus Surgery Solution
CAS KIT |
| 5 | Sinus Surgery Solution
LAS KIT |
| 6 | Narrow Ridge Solution
ESSET KIT |

Short Implant Solution

The KIT is composed of specially designed drills to avoid nerve damage and place short-body implants. (4mm ~ 8.5mm)



Procedure

- 01
KIT & Case
Introduction
- 02
Normal
Drilling
: 1st Model
- 03
485 Drilling
: 2nd Model
- ★ 04
Comparing
2 Models
(Normal vs 485)
- 05
Implant
Placement

- PREPARATION**
- Set the drill speed to 1,200 rpm for drilling.
 - Use 2 models (1st model : normal drill sequence / 2nd model : 485 drill sequence).
 - Prepare normal option drills (122 Taper Drill).

★ **SPECIAL NOTE**

The practice model was designed to show how the 485 KIT is safer for alveolar nerve than the standard KIT by comparing how transparent the nerve is at the drilled hole.

* **Two models are required for a single hands-on**
(One for 122 taper drilling / the other for 485 drilling).

Do the hands-on

485 Drill

VS

Normal Drill

No exposure of nerve

Exposure of nerve

Materials

KIT

485 KIT

Drill

122TPD3504	122TPD3505
122TPD3506	122TPD3507
122TPD5004	122TPD5005
122TPD5006	122TPD5007

Model

485 Model X2

*A Custom order is required.


Implant

5.0 X 4mm	5.0 X 5mm
5.0 X 6mm	5.0 X 7mm

Short Implant Solution

485 KIT

- 485 Model *A Custom order is required.
- RPM: 1,200
- Normal Bone



Drilling Sequence for The Normal Drill Ø 5.0mm

Use 1st Model

Bone Quality

Soft	▶	▶	▶	▶	▶
Normal	▶	▶	▶	▶	▶
Hard	▶	▶	▶	▶	▶

Implant placement

Drilling Sequence for The 485 Drill Ø 5.0mm

Use 2nd Model

Bone Quality

Normal	▶	▶	▶	▶
Hard	▶	▶	▶	▶


Implant placement

Caution 1

This Hands-on case was designed to demonstrate the superior safety of 485 KIT comparing the transparencies of alveolar nerves at the drilled holes. Please note that a 2mm distance must be secured between drilling depth (implant size) and alveolar bone in actual clinical case.

Caution 2

The stopper was design to be stopped at +1mm of its size. Please drill to the laser mark.



Simple Drilling Solution

By using 122 Taper KIT, it is able to place an implant with drilling only once for soft bone and twice for normal bone & hard bone.



Procedure

- 01 KIT & Case Introduction
- 02 Drilling : Bone Block
- 03 Drilling : Model

PREPARATION

1. Set the drill speed to 1,200 rpm for drilling.
2. Drill the bone block using three different drill sequences : Soft / Normal / Hard bone.
3. Drill the practice model using the normal bone drill sequence.

★ **SPECIAL NOTE**

Operators can first practice the basic drilling steps on the bone block to experience how simple they are. Then, they can practice the drilling sequence on the practice model.

*Please explain the 122 drilling sequence at the bone block practice stage.

Materials

KIT

122 Taper KIT

Model

SUMBB01 X1 GL07 X1 BL07 X1


Implant

3.5X10 mm X1 4.5X10 mm X2

Simple Drilling Solution

122 Taper KIT

- Model: GL07 / BL07
- RPM: 1,200
- Normal Bone



Drilling Sequence for Mini Implant Ø 3.5mm #41



Bone Quality	Sidecut Drill	122 Taper Drill (F3.5)	122 Taper Drill (F4.0)	Ø 3.5 Implant
Soft	▶			
Normal	▶	▶		Implant placement
Hard	▶	▶	▶	

Drilling Sequence for Regular Implant Ø 4.5mm #35,36



Bone Quality	Sidecut Drill	122 Taper Drill (F3.5)	122 Taper Drill (F4.0)	122 Taper Drill (F4.5)	122 Taper Drill (F5.0)	Ø 3.5 Implant
Soft	▶		▶			
Normal	▶	▶	▶	▶		Implant placement
Hard	▶	▶			▶	

Implant Aligning Solution

The KIT helps to find both right spot and angle for implant placement without using a surgical template.
Also, it can be used to place implants on the center of teeth with the consideration of final prosthesis positions.

○ Good Case

✗ Bad Case



Implant placement on the center of teeth

Incorrect Implant Position

Implant placed too close to the adjacent teeth

Implant placed too close to each other

Procedure

01
KIT & Case Introduction

02
Case 1 : Single

03
Case 2 : Bridge

04
Case 3 : Free End

★
05
(Implant Placement) : 122 Taper KIT
* Optional

PREPARATION

1. Set the drill speed to 1,200 rpm for drilling.
2. Prepare the 122 Taper KIT for implant placement.
3. Prepare Implants for the implant placement.

★
SPECIAL NOTE

The Positioning Guide KIT provides guidance of initial drilling. Please follow the instruction of 122 Taper drill for the implant placement.

* For the practice of the implant placement, 122 Taper KIT is required.

Positioning Guide KIT



122 Taper KIT



Materials

KIT



Positioning Guide KIT

122 Taper KIT
* Optional

Model



OR



GL07 X1

BL07 X1

Implant



3.5X10 mm X1

4.5X10 mm X4

* Optional

* Optional

Implant Aligning Solution

Positioning Guide KIT

- Model: GL07 / BL07
- RPM: 1,200
- Normal Bone

Case 3

R 4.5X10 mm

R 4.5X10 mm


Case 2

R 4.5X10 mm

R 4.5X10 mm

Case 1


M 3.5X10 mm



Case 1 Single (Mini)

Case 2 Bridge (Regular)

Case 3 Free End Case



Single Guide X1

Ø 5.0

Single Guide X2

Ø 9.0

Ø 10.0

Single X1

Bridge X1

Ø 9.0

9-10.5mm

Drilling Sequence for Single Guide

Case 1

Case 2

1

2

3

Select Guide

Connect guide and perform drilling

Check path and crown width

Drilling Sequence for Single & Bridge Guide

Case 3

1

2

3

4

Select single guide and perform drilling

Connect bridge guided

Perform drilling

Check path and crown width, then place implant

Caution 1

- When using the single guide, drill the guide towards the adjacent teeth.
- Check Buccal & lingual part carefully and drill the guide at the center of them.


Caution 2

When connecting the single guide to the drill, make sure that the open side is facing up.



Sinus Surgery Solution

Crestal Approach Sinus Surgery KIT is designed to perform sinus lifting without causing membrane perforation. The KIT is composed of patented conical inverse drills and secures easier and safer sinus membrane elevation.



Procedure

- 01
KIT & Case
Introduction
- 02
Drilling
- 03
Membrane
Lifting
- ★ 04
Check
Perforation
- 05
Bone
Graft
- 06
Implant
Placement

* Optional

- PREPARATION**
- Set the drill speed to 800 rpm for drilling.
 - Fill 3cc syringe with 1.5 cc of water.
 - Prepare bone graft materials.

★
**SPECIAL
NOTE**

The practice model was designed to check the condition of sinus membrane from the open window. An operator can check the sinus perforation during the hands-on practice in real time.

* After lifting the sinus membrane, check the sinus perforation.


Materials

KIT	Model	Syringe	Bone	Implant
CAS KIT	BU05 X1	3cc Syringe X1	A-oss 1cc X1 * Alternative	4.5X8.5mm X1

Sinus Surgery Solution

CAS KIT

- Model: BU05
- RPM: 800
- Normal Bone



Adjust the Depth by Fastening the Stopper

Stopper Drill

Exclusive Design for Sinus Surgery

Initial Drilling

Drilling

Drilling (CAS Drill)

*** Lifting**

*** Bone Graft**

Depth Check

Connect a stopper to the depth gauge to check the condition of a drilled hole whether it is open and ready to be lifted.

Membrane Lifting


Fill the 3cc Syringe with 1.5cc of saline
Repeat injecting and drawing saline in the following volumes. : 0.5cc – 1 cc – 1,5 cc

Bone Graft

Herry Y and Lee DY, 2005	
Lift Heigh	Volume of Bone Matrix
3mm	0.36cc
4mm	0.5cc
5mm	0.7cc
6mm	0.9cc

Sinus Surgery Solution

LAS KIT is composed with dome drills, used to form a lateral window for sinus surgery. It is safer and less time consuming to form a lateral window, compared to cut the plate by using a round bur. Different size of stoppers can be connected to the drill to control the drilling depth.



Procedure


- 01 KIT & Case Introduction
- 02 Drill Choice : Core or Dome Drill
- 03 Drilling
- 04 Check Perforation
- 05 (Implant Placement) : 122 Taper KIT
* Optional

PREPARATION

1. Set the drill speed to 1,200 rpm for drilling.
2. Give the two options : the core drill and the dome drill
3. Follow the drill sequence for the selected drill.

★ **SPECIAL NOTE**

When using a **Dome Drill**,
Please explain that a large lateral window can be secured by sliding the dome drill.



When using a **Core Drill**,
Please explain that the bone lid generated by the core drill can be used to cover the window.



Materials

KIT	Model	Bone	Implant
 LAS KIT	 BU05 X1	 A-oss 3cc X1 * Alternative	 4.5X8.5mm X2 * Optional
 122 Taper KIT * Optional			

Sinus Surgery Solution

LAS KIT

- Model: BU05
- RPM: 1,200



Prevents the Membrane Damage using Special Drills



Dome Drill

OR



Core Drill

Adjust the Depth by connecting the Stopper to the Drill



Stopper

Smooth cutting and minimal unseating

Safely creation of the bone lid

* Use side wall drill after sequence (1,500 rpm)


Drilling Sequence for The Dome Drill




Approaching the Sinus Membrane



Widening the Window




Drilling Sequence for The Core Drill



Approaching the Sinus Membrane




Separating the Bone Lid



Narrow Ridge Solution

By using ESSET KIT, ridge split and expansion can be much less traumatic for patients.
No need to use a chisel or a mallet.



Procedure

- 01
KIT & Case
Introduction
- 02
Ridge
Removal
- ★ 03
Sawing
- ★ 04
Expansion
- ★ 05
Implant
Placement

PREPARATION

1. Set the drill speed to 1,200 and 35rpm for drilling and expansion.
2. Prepare a saw protector. (*Optional)

★
**SPECIAL
NOTE**

To prevent buccal plate fracture, Follow the guided process in the order of middle, distal and mesial.

Sawing Expansion Implant Placement

Materials

KIT

ESSET KIT

Model

RS-01 X1


Implant

4.0X10mm X1 4.5X10mm X2

Narrow Ridge Solution

ESSET KIT

- Model: RS-01
- RPM: 1,200 / 35



Approaching Sinus	Cut Alveolar Bone	Expand Alveolar Bone	Place Implant
1 Mesial to Distal RPM: 1,200 Crestal Remover * Move the crestal remover in a counter drill wise to prevent irregular vibrations	2 RPM: 1,200 Twist Drill & Saw	3 RPM: 35 Expansion Drill	4 4010 4510 Implant
 Crestal Remover Secure 4mm in bone width	 Twist Drill Mark the positions for expansion drills	 Expansion Drill ①,② Expand alveolar bone (Middle-Distal-Mesial)	 4.5X10mm Remove expansion drills and place implants
 Crestal Remover Make a flat area for position marking	 Saw Saw alveolar bone (Middle-Distal-Mesial)	 Expansion Drill ③,④ Maintain each of the last expansion drills after the final expansion	 4.0X10mm Place implants (Middle-Distal-Mesial)