



## Minimal Incision

# Implant Finder


- Probe with explorer after detecting precise implant inserted area
- Shorten the healing period by minimal incision

**OSSTEM<sup>®</sup>**  
IMPLANT

### Detecting Precise Implant Inserted Area

- A number of lights vary according to number of fixtures implanted and its depth
- Shorten the healing period by minimal incision as it is able to detect precise inserted area

Flap open



General incision

Implant finder incision

### Optimize Surgical Efficiency by Unique Body Design

- Easiness of treatment and excellence in grip feeling by streamlined design minimizing interproximal contact



**Reset Button**  
Reset Sensor



**Implant Detecting LED**  
More number of led lights turned on as finder gets closer to implant inserted area.




**Power LED**  
A single light turned on when power on

### Excellent Sensor Performance


- Sense widely up to 5mm radius(height & depth)
- Protect patient's inner oral by silicon cover
- 10 sets of sensors included

### Surgical Sequence


① Probe mesial, distal direction




② Probe lingual, buccal direction



③ Probe with explore



④ Incision with tissue punch (blade)







## One Touch, Automatic Pressure Adjustment


# SlowJec

Minimize patient pain with automatic pressure adjustment function,  
excellent usability with one-touch operation


**OSSTEM<sup>®</sup>**  
IMPLANT

### Automatic Pressure Adjustment by Sensor

- Pause when over pressure is detected by measuring the pressure loaded on the needle
- When the pressure decreases, automatically inject the remaining amount



Pause when over pressure is detected



Automatic injection when the pressure is decreased

### Normal Syringe vs SlowJec

	Normal syringe	SlowJec
Injection time	With in 10seconds	<b>1min 4seconds</b>
Pressure	High	<b>Low</b>
Pain	Strong	<b>Very weak</b>

※ Criteria : Inject 100% 1.8ml, in one touch mode

**LCD Screen**  
Check Injection quantity,  
pressure quantity,  
injection speed

**Voice Support**  
Beeps to notify

**Aspiration Function**  
Check the penetration of the vessel before injection

**One Touch Button**


**Compatible with General Types of Needles**  
Wide range of drugs can be injected

**Replaceable Battery**  
Removable

**Various Injection Volume Control and Speed Adjustment**  
5-level injection volume controls

### Easy One Touch Operation

- Press once to complete the injection
- No need of gripping the syringe tight for the whole time



※ Check Injection quantity, pressure quantity, injection speed

### Convenient Pen Type & Cordless

- Easy to use like a pen · Compact size
- Light weight (150g including a battery)





Machine Screw Driver Handle Attached to The Low-Speed Handpiece

# Torque Handle

**OSSTEM<sup>®</sup>**  
IMPLANT

## Convenient for The Following Applications

- Connect cover screw and healing abutment
- Fasten abutment screw
- Insert corrective screw
- Removal surrounding bone with Bone profiler
- Less cumbersome than a high-cost torque driver

## Compatible with All Low Speed Handpieces

KaVo, NSK, W&H, etc. connectable to all low-speed handpieces  
(Economical and inexpensive, compared to other products which require an exclusive headpiece)



## Tools

- Machine screw driver
- Corrective driver
- Bone profiler



## Good Grip Feeling and Stable Torque

- Long facing and linear design, light(91g) with good grip feeling
- Stable torque by big diameter handle at the end

## Ease access to Narrow Mouth

Once Torque handle is connected to handpiece, it is being screwed and unscrewed out of the oral cavity, allowing safe procedure in molars



## Enables Implantation of Orthodontic Screw, with No Surgical Engine

For dental clinics specialized in orthodontic, there is no need to purchase an expensive surgical engine, as orthodontic screw insertion can be performed with the torque handle.



## Reduced Chair Time

>50% reduction in chair time compared to using a general driver. It is increasing satisfaction for both patient and doctor

