





Gingival Retraction Cords



- #101 AtriaPak #000 (Black)
- #102 AtriaPak #00 (Yellow)
- #103 AtriaPak #0 (Purple)
- #104 AtriaPak #1 (Blue)
- #105 AtriaPak #2 (Green)
- #106 AtriaPak #3 (Red)



AtriaPak Knitted Retraction Cords Non-Medicated

AtriaPak cords 254 cm (100 inches). AtriaPak knitted cords provide an excellent medium for gel absorption a placement. AtriaPak cords are used in the retraction of gingival tissue around tooth. The unique knitted design exerts a gentle, continues outward force following placement as the knitted loops seek to open.



Cord Stand (4 holes)

- #101AS AtriaPak Aluminium Sulfate #000 (Black)
- #102AS AtriaPak Aluminium Sulfate #00 (Yellow)
- #103AS AtriaPak Aluminium Sulfate #0 (Purple)
- #104AS AtriaPak Aluminium Sulfate #1 (Blue)
- #105AS AtriaPak Aluminium Sulfate #2 (Green)
- #106AS AtriaPak Aluminium Sulfate #3 (Red)

AtriaPak Knitted Retraction with Aluminium Sulfate


AtriaPak cords 254 cm (100 inches). AtriaPak cords with aluminium sulfate are mad of 100% cotton and impregnated with a solution of aluminium sulfatem formulated to contain 0.35 ± 0.10 mg per cm, then dried by a special process, assuring both consistent aluminium sulfate content and hemostatic retraction properties. AtriaPak cords are used in the retraction of gingival tissue around tooth. The unique knitted design exerts a gentle, continues outward force following placement as the knitted loops seek to open.



Mouth Mirrors High quality mouth mirrors for best reflection and image






The Thinner, the better !



Simple Stem (SS type)



Cone Socket (CS type)



Rear surface mirrors



12 pcs. / box


- * The glass is very thin only 1.10 mm
- * Rear surface mirrors are silver coated on the back space.
- * The full silicone seal to prevent water from entering autoclave sterilization.

Front surface mirrors




12 pcs. / box

- * Front surface mirrors will not give a double image and so are ideal for endodontic procedures.



Cone Socket - 3CS, 4CS, 5CS
Simple Stem - 3SS, 4SS, 5SS



Cone Socket - F3CS, F4CS, F5CS
Simple Stem - F3SS, F4SS, F5SS

Mirror Holders



MHO CS (Hexagonal, All stainless steel material)



MHO SS (Hexagonal, All stainless steel material)



MHR CS (Round, All stainless steel material)



MHR SS (Round, All stainless steel material)



MHCS (Plastic, Blue, Autoclaveable)





MHSS (Plastic, Green, Autoclaveable)



Modeling Waxes (Paraffin, Base plate wax)

For modeling denture bases and acrylic saddle supporting the teeth
For vestibular parts of wax models of removable dentures, pre-forms of jacket crowns, bite blocks and manufacturing of spacers for custom impression trays




WP 17 (Pink) - 1 kg, 170 x 95 x 1.7 mm
WR 17 (Red) - 1 kg, 170 x 95 x 1.7 mm
WP 500 (Pink) - 500 g, 170 x 95 x 1.5 mm
WR 500 (Red) - 500 g, 170 x 95 x 1.5 mm
WP 225 (Pink) - 225 g, 144 x 74 x 1.6 mm
WR 225 (Red) - 225 g, 144 x 74 x 1.6 mm

Utility Waxes

Round strip is used for preserving the periphery roll and extending and building up impression trays and post-damming.

Red normal - 320 g (5 x 300 mm)
Red soft - 320 g (5 x 300 mm)
White normal - 320 g (5 x 300 mm)



Orthodontic Waxes

Patient relief wax sticks help soothe and prevent irritation to the gums caused by Braces and other dental appliances.


Clear Soft - 320 g (4.7 x 5.2 x 300 mm)
OWC - 5 pcs./pack (4.7 x 5.2 x 50 mm)



Wax profiles (Sprue wax, Line wax)

Wax profile is used for the production of wax sprue, exhaust passage, palate, such as a stent wax and then embedding, dewaxing and casting metal.

Round - R7, R10, R15, R20, R25, R35, R40, R45, R50
Half - H15, H20, H25, H30, H35, H40
Half Pear - L35, L40, L45
Assortments (Net 160 g / pack)
- R10, R15, R25, L40, H15, H20, H30, H40



Sheet Waxes

Sheet waxes can be used as a tracing wax for transferring illustrated patterns and perfect for making models for the lost wax casting process

Color : Blue, Pink (170 x 95 mm)
Thickness - 0.20, 0.25, 0.30, 0.35, 0.40, 0.45, 0.50, 0.55, 0.60, 0.65, 0.70, 0.75



Stippled Waxes

Color : Blue (170 x 95 mm)
Thickness - 0.20, 0.25, 0.30, 0.35, 0.40, 0.45, 0.50, 0.55, 0.60, 0.65, 0.70, 0.75




A-Sleeve (A551600)

Disposal Drapes, prevent cross-contamination suitable for exposed cable, cords and hoses reduces clean-up time.
Protects and gives handpiece or suction that disinfection does not become when operate.

Disposable drape for covering dental and medical handpieces and hoses
55 mm x 1600 mm (24 pcs. / box)



Plastic Cassettes



AT-10

For 10 instruments
Size - 206 x 136 x 22 mm

Color - Red, Yellow, Black, Blue, Green
It made of an exceptionally durable special plastic.
It can sterilize max 180 °C / 360 °F at autoclave.

AT-5

For 5 instruments
Size - 206 x 77 x 22 mm



Dental Anesthetic Syringes / Stainless Steel



AS-100SS

Not aspirating



AS-200SS

Aspirating



AS-300SS

Aspirating



AS-400SS

Self aspiration



Bone Crusher with Titanium



BC-100

All Titanium materials, diameter : 40 mm



Surgical Mallet

SM-210

190 mm, 210 g
With interchangeable plastic autoclaveable disks

Amalgam Band Retainer



Tofflemire
TR-600

Secure lock
Safe band retention

Amalgam Bands



Tofflemire

AT01-0015

0.0015" (100 ea / pack)

AT01-002

0.002" (100 ea / pack)

Smooth Stainless Steel material, Fine finished burr free edges
Band will be rigid against the existing tooth structure and prevent gingival overhang

Cord Packers



Packing 1 (Non serrated)



CSI-1 (serrated)



CP 113 (serrated)



CP 170 (45° serrated large tips)



CP 171 (45° serrated small tips)



CP 172 (90° serrated large tips)



CP 174 (90° serrated small tips)



Tissue Punches for contra-angle



- TP-35** (Inside 3.5 mm)
- TP-40** (Inside 4.0 mm)
- TP-45** (Inside 4.5 mm)
- TP-50** (Inside 5.0 mm)
- TP-60** (Inside 6.0 mm)



Tray (for 6 pcs.)


Crown Remover



It is using spring without screw to connect tips into a body.

CR-400

Metal Suction tips



For vacuum valve For saliva ejector

Dental Cotton Container



2 Chambers
It is autoclaveable

Instruments Tray



218 X 118 mm
Stainless Steel for autoclave



Irrigation Tubing Set

Able to reuse as it is safe to autoclave



For Implant Motor

- * NI - for NSK
- * KI - for KAVO
- * WI - for W&H
- * GI - for GOOF

For Piezo Surgery

- * PI - for PiezoSurgery
- * SI - for SurgyBone
- * EI - for Esacrom
- * AI - for PiezoArt

Tips for Crown Tractor



Green Tips : Hard
40 pcs. / pack

Yellow Tips : Soft
40 pcs. / pack

Plastic-Polymer tips.
It is autoclaveable.

Color Code Rings

Standard size

Autoclave to 134°C or 275°F / 18 minutes

Dry heat to 160°C or 350°F / 2 hours



Set (140 pcs.)
7 colors x 20 pcs.



Refill (50 pcs.)

Color : Red, Orange, Yellow, White, Green, Blue, Black