



www.befrest.com

PRODUCT CATALOG

HIOSSEN

003	INTRODUCTION
006	CONTENTS
014	ET SYSTEM
070	HIOSEN KIT
122	GBR
138	REFERENCE




HIOSSEN Design features

HiOssen Implant,

We stand out with our passion for strategic R&D and innovative products, creating globally trend-setting implants.

www.befrest.com



Submerged type implant with an Internal hex 11° taper connection structure

- Connection type and color - Mini / Regular
- Highest initial stability in soft bone by using upper-section small thread
- Corkscrew thread & cutting edge
 - Easy path adjustment through a superior self-threading effect
 - Acquires insertion torque with an increase in soft bone initial stability and without deviation according to the drill diameters
- The various body shape options are available according to the bone and patient's clinical condition
 - ETII (Straight Body) : Easily adjustable insertion depth
 - ETIII (1.5 Taper Body) : Strong initial stability
 - TSIV (6° Taper Body): Superior initial stability (recommended for soft bone and maxillary sinus)
- Applied Surface - SA / CA / BA / HA

HIOSSEN Surface features

HiOssen Implant,

provides world-class surface treatment technologies for fast and safe procedures



www.befrest.com

SA BA

Provides optimum surface through acid treatment

- Provides Ra 2.5~3.0 μm surface roughness
However, upper section 0.5mm area is Ra 0.5~0.6 μm
- Achieved uniform micro-pit 1.3 μm in size
- 46% greater surface area compared to RBM

SA Surface treatment performance (in-vitro and in-vivo)

- 20% improvement in osteoblast separation and ossification compared to RBM
- 48% improvement in initial stability (RT, 4 weeks) compared to RBM
- 20% improvement in ossification (BIC, 4 weeks) compared to RBM

Surface coated with low crystalline Nano-HA in SA

- Ultra-thin film with HA coating and 10nm or lower thickness
- HA coating on SA surface (Ra 2.5~3.0 μm)
- Dual function of titanium and HA
- HA is naturally removed during ossification process

BA Surface treatment performance (in-vitro and in-vivo)

- Fused surface having advantages of both SA and HA
- Maintains advantage of SA optimum surface formation
- Superior early ossification of the HA in soft bone condition
- 30% improvement in ossification (BIC) compared to SA



ET SYSTEM Contents

014
ETIII SA
Fixture



016
ETII SA
Fixture



018
ETIII BA
Fixture



020
Cover
Screw



021
Healing
Abutment



022
Custom
Healing
Abutment



026
SmartFit
Abutment



027
Transfer
Abutment



029
GoldCast
Abutment



030
NP-Cast
Abutment



031
Temporary
Abutment



032
Quick Temp
Abutment



033
Angled
Abutment



035
Angled
Abutment
Selector



036
ZioCera
Abutment



038
ZioCera
Angled
Abutment



039
FreeForm
ST
Abutment



041
Fixture
Transfer
Impression
Coping



042
Fixture
Pick-up
Impression
Coping



043
Fixture
Lab
Analog



045
Rigid
Abutment



048
Rigid
Protect
Cap



048
Rigid
Impression
Coping



048
Rigid
Burn-out
Cylinder



049
Rigid
Lab
Analog



051
Locator®
Abutment



052
Locator®
Male
Processing
Kit



052
Locator®
Replacement
Male



052
Locator®
Extended
Replacement
Male



053
Locator®
Black
Processing
Male



053
Locator®
Block Out
Spacers



053
Locator®
Impression
Coping



053
Locator®
Lab Analog



054
Locator®
Core Tool



054
Locator®
Torque
Driver





055
Stud
Abutment



056
O-ring
Retainer
Cap Set



056
O-ring
Retainer
Set



056
O-ring
Set



056
O-ring
Set



059
Multi
Abutment



060
Multi
Abutment
Outer Driver



060
Multi
Angled
Abutment
Screw



060
Multi
Angled
Abutment
Carrier



061
Esthetic-
low Healing
Cap



061
Esthetic-
low Gold
Cylinder



061
Esthetic-
low Plastic
Cylinder



062
Esthetic-
low
Temporary
Cylinder



062
Esthetic-
low Pick-up
Impression
Coping



063
Esthetic-
low Transfer
Impression
Coping



063
Esthetic-
low Lab
Analog



063
Esthetic-
low
Polishing
Protector



064
Multi
Angled
Abutment



065
EbonyGlad
Abutment
Screw



065
Titanium
Abutment
Screw



066
Laboratory
Screw





KIT Contents

www.befrest.com

008

CONTENTS

070

Taper
KIT



071

Standard
KIT



072

Taper
Mini
KIT



073

Taper
Ultra
KIT



074

Ultra KIT



075

123 KIT



076

Lance Drill -
Guide Drill



076

123 Guide
Drill



076

Sidecut
Drill



076

Drill
Extension



077

Twist Drill
- Stopper
Drill



077

Twist Drill -
Non Stopper
Drill



078

123 Twist
Drill



078

123 Drill
Stopper



079

Direct Drill



079

Taper Drill



080

Taper Ultra
Drill



080

Long
Shank
Pilot Drill



080

Cortical
Drill 3
for Taper
Fixture



081

Taper
Cortical
Drill for
Taper
Fixture



081

Cortical
Drill for
Ultra-Wide



081

Parallel
Pin



082

Depth
Gauge



082

Depth
Gauge



082

Torque
Wrench -
Bar Type



082

NoMount
Driver for
ET



083

NoMount
Torque
Driver for
ET



083

Fixture
Driver for
ET



083

Torque
Extension



083

Simple
Mount
Driver



084

Simple
Mount
Extension



084

Removal
Tool for
Fixture
Mount



085

Tissue
Punch



085

Bone
Profiler



086

Trephine
Drill





086
Bone
Mill



087
Prosthetic
KIT



088
Hand
Driver



088
Machine
Screw
Driver



089
Torque
Driver



089
O-ring
Abutment
Driver



089
Rigid Outer
Driver



090
Solid
Abutment
Driver



090
Excellent
Solid
Abutment
Driver



091
Octa
Abutment
Driver



091
Driver
Handle



091
Finishing
Reamer
Set



092
Reamer
Bite



093
CAS-KIT



094
CAS-
Drill



094
Guide
Drill



094
Ø 2.2
Twist Drill



094
Hydraulic
Membrane
Lifter Set



095
Stopper



095
Bone
Carrier



095
Bone
Carrier
Head



095
Bone
Condenser



096
Hydraulic
Membrane
Lifter Tube



096
Depth
Gauge



097
LAS-KIT



098
Dome
Drill



098
Core
Drill



098
Side Wall
Drill



099
Bone
Separator



099
Stopper



100
ESSET
KIT



101
Crest
Remover



101
Twist Drill



101
Saw



102
SET Drill





KIT Contents (continued)

www.befrest.com

010

CONTENTS

102
Mount
Extension



102
Torque
Wrench



102
Depth
Gauge



103
Osstem
Guide
KIT



104
Osstem
Guide
Twist Drill



105
Osstem
Guide
Mount



105
Osstem
Guide
Cylinder
(Pin)



106
MS KIT



107
ESR KIT



108
Reverse
Driver



108
Screw
Removal Drill
(SR Drill)



108
Guide



109
Torque
Handle



109
Screw
Removal
Tip
(SR Tip)



109
Screw
Holder



109
Re-tap



110
ESR
Handle



110
Abutment
Removal
Tip



110
Transfer
Abutment
Separate
Tool



110
Slot Driver



111
EFR KIT



112
Remover
Screw



112
Screw
Driver



112
Remover
Body



113
Torque
Extension



113
Torque
Wrench



113
Fixture
Wrench



114
Osteo
KIT



115
Sinus KIT



116
Bone
Spreader
KIT



117
Ridge
Split KIT -
Straight



118
Ridge
Split KIT -
Offset



GBR Contents



123
SmartBuilder
SB1



124
SmartBuilder
SB2



125
SmartBuilder
SB3



126
Healing
Cap



126
Cover
Cap



126
SB Anchor
for ET



127
Tenting
Screw



127
Cover
Cap Driver
Ejector Set



127
Defect
Gauge



128
AutoBone
Collector



128
Stopper



128
Bone
Ejector



129
Allogeneic
bone
SureOss
(FDBA)



129
Allogeneic
bone
SureOss-D
(DFDBA)



130
Allogeneic
bone
OsteOss
(FDBA)



131
Allogeneic
bone
SureFuse II



131
Allogeneic
bone
ExFuse II



132
Xenograft
A-Oss



133
Synthetic
bone
Q-Oss



133
Synthetic
bone
Q-Oss+



133
Synthetic
bone
Bongros



134
Absorbent
membrane
Cytoplast
RTM
Collagen



135
Absorbent
membrane
SureDerm



136
Impermeable
membrane
Cytoplast
TXT-200



136
Impermeable
membrane
Cytoplast
TI-250





HROSSEN



ET SYSTEM

FIXTURE

- 014** ETIII SA Fixture
- 016** ETII SA Fixture
- 018** ETIII BA Fixture
- 020** Cover Screw
- 021** Healing Abutment
- 022** Custom Healing Abutment

COMPONENTS

- 024** PROSTHETIC FLOW DIAGRAM 1
- 025** PROSTHETIC FLOW DIAGRAM 2
- 026** SmartFit Abutment
- 027** Transfer Abutment
- 029** GoldCast Abutment
- 030** NP-Cast Abutment
- 031** Temporary Abutment
- 032** Quick Temporary Abutment
- 033** Angled Abutment
- 036** Ziocera (Angled) Abutment
- 039** FreeForm ST Abutment
- 044** PROSTHETIC FLOW DIAGRAM 3
- 045** Rigid Abutment
- 050** PROSTHETIC FLOW DIAGRAM 4
- 051** Locator® Abutment
- 055** Stud Abutment
- 058** PROSTHETIC FLOW DIAGRAM 5
- 059** Multi (Angled) Abutment
- 065** Abutment Screw

ETIII SA Fixture



- Submerged type implant with 11° Morse taper internal hex connection
- Taper body with optimized screw thread design for superior initial stability
- World class SA surface treatment initializes osseointegration efficiently

Ultra-wide

- Ideal for posterior immediate placement after extraction or failed implant
- Optimized apex design for obtaining initial stability
- Recommended insertion torque : 40Ncm or less


※ For a single posterior implant case, using 4.5mm or bigger diameter is recommended.

NoMount fixture order code


: **C** + fixture product code (ex: **CET3R4010S**)

Pre-Mounted fixture order code (fixture + mount + cover screw)


: **A** + fixture product code (ex : **AET3R4010S**)




D Ø3.5
Hex 2.1



D Ø4.0
Hex 2.5



D Ø4.5
Hex 2.5





D Ø5.0
Hex 2.5



L	6	7	8.5	10	11.5	13	15
							
	Short						
	ET3R5006S Pre-Mount Only	ET3R5007S	ET3R5008S	ET3R5010S	ET3R5011S	ET3R5013S	ET3R5015S

Ultra-Wide

D Ø6.0
Hex 2.5



L	6	7	8.5	10	11.5	13
						
	Short					
	ET3R6006S Pre-Mount Only	ET3R6007S	ET3R6008S	ET3R6010S	ET3R6011S	ET3R6013S

D Ø7.0
Hex 2.5



L	6	7	8.5	10	11.5	13
						
	Short					
	ET3R7006S Pre-Mount Only	ET3R7007S	ET3R7008S	ET3R7010S	ET3R7011S	ET3R7013S



ETII SA Fixture

- Submerged type implant with 11° Morse taper internal hex connection
- World class SA surface treatment initializes osseointegration efficiently
- Straight body design for easy adjustment of insertion depth
- Unique thread design has enhanced initial stability in soft bone and powerful self threading capabilities
- Recommended insertion torque : 40Ncm or less


※ For a single posterior implant case, using 4.5mm or bigger diameter is recommended.

NoMount fixture order code

: C + fixture product code (ex:CET2R4010S)

Pre-Mounted fixture order code (fixture + mount + cover screw)

: A + fixture product code (ex : AET2R4010S)



D Ø3.5

Hex 2.1



L	8.5	10	11.5	13	15
	ET2M3508S	ET2M3510S	ET2M3511S	ET2M3513S	ET2M3515S

D Ø4.0

Hex 2.5



L	7	8.5	10	11.5	13	15
	ET2R4007S	ET2R4008S	ET2R4010S	ET2R4011S	ET2R4013S	ET2R4015S

D Ø4.5

Hex 2.5



L	7	8.5	10	11.5	13	15
	ET2R4507S	ET2R4508S	ET2R4510S	ET2R4511S	ET2R4513S	ET2R4515S



D Ø5.0
Hex 2.5



L	6	7	8.5	10	11.5	13	15
							
	Short						
	ET2R5006S Pre-Mount Only	ET2R5007S	ET2R5008S	ET2R5010S	ET2R5011S	ET2R5013S	ET2R5015S



ETIII BA Fixture

- Submerged type implant with 11° Morse taper internal hex connection
- Taper body with optimized screw thread design for superior initial stability
- SA surface coated with low crystalline Nano-HA
- Recommended insertion torque : 40Ncm or less


※ For a single posterior implant case, using 4.5mm or bigger diameter is recommended.

NoMount fixture order code


: C + fixture product code (ex: CET3R4010B)

Pre-Mounted fixture order code (fixture + mount + cover screw)


: A + fixture product code (ex : AET3R4010B)




D Ø3.5
Hex 2.1



D Ø4.0
Hex 2.5




D Ø4.5
Hex 2.5





D Ø5.0
Hex 2.5





Screw

Cover Screw

- 1.2 Hex driver
- Recommended tightening torque : 5~8Ncm
- C = Connection

M Mini

R Regular

C \ H 0.4

M

Fixture level



GSCS35

C \ H 0.4

R

Fixture level



GSCS40S-B

Healing Abutment




- 1.2 Hex driver
- Tightening torque : 5~8Ncm

M Mini
R Regular

Selection table

Healing Abutment	H	3.0	4.0	5.0	7.0
Abutment	G/H	1.0	2.0 or 3.0	3.0 or 4.0	5.0
Impression coping	Type	Short	Short	Long	Long



D \ H 3.0 4.0 5.0 7.0

Ø4.0
Ø4.5



ETHLA4003M
ETHLA4503M



ETHLA4004M
ETHLA4504M



ETHLA4005M
ETHLA4505M



ETHLA4007M
ETHLA4507M

D \ H 3.0 4.0 5.0 7.0

Ø4.0
Ø4.5
Ø5.0
Ø6.0
Ø7.0



ETHLA4003R
ETHLA4503R
ETHLA5003R
ETHLA6003R
ETHLA7003R



ETHLA4004R
ETHLA4504R
ETHLA5004R
ETHLA6004R
ETHLA7004R



ETHLA4005R
ETHLA4505R
ETHLA5005R
ETHLA6005R
ETHLA7005R



ETHLA4007R
ETHLA4507R
ETHLA5007R
ETHLA6007R
ETHLA7007R




Custom Healing Abutment

Available Soon

- Material : medical PEEK
- Customizable : Easy to adjust to achieve ideal emergence profile
- 2.1 Hex driver
- Recommended tightening torque : 5~8Ncm
- Packing unit : abutment + Ti screw
- C = Connection

M Mini

R Regular



C \ H

5.0

7.0

9.0

M

Ti screw
: GSCHABSMT

TSCHAP**M5**

TSCHAP**M7**

TSCHAP**M9**

R

Ti screw
: GSCHABSST

TSCHAP**R5**

TSCHAP**R7**

TSCHAP**R9**



HROSSEN



PROSTHETIC FLOW DIAGRAM 1


Transfer / Angled / FreeForm ST / GoldCast / NP-Cast

Fixture Level Impression

www.befrest.com

ET SYSTEM

024






PROSTHETIC FLOW DIAGRAM 2

SmartFit / ZioCera / ZioCera Angled Temporary / Quick Temporary

Fixture Level Impression





SmartFit Abutment

- CAD/CAM Custom Abutment
- 1.2 hex driver
- Recommended tightening torque : 20Ncm(mini), 30Ncm(regular)
- Recommended clinical case
 - Case where implant insertion area or angle is not ideal (max correction 30°)
 - Multiple cases requiring consistent path and stable guide
 - Anterior case where aesthetic design is required
 - Irregular or exceedingly deep gingiva case

How to Order


- Fill out an order sheet (*Contact your local rep for the order sheet)
- Send necessary items for each case to SmartFit Center
- Working time : 5~7days



Transfer Abutment



- Cement-retained prosthetic
- Gold colored gingiva region for improved aesthetics
- 1.2 hex driver
- Recommended tightening torque : 20Ncm(mini), 30Ncm(regular)
- Packing unit : abutment + EbonyGold screw



ET SYSTEM

027

D Ø4.5



EbonyGold screw
: HGATSM

Hex
Non-Hex **5.5**

H \ G/H	1.0	2.0	3.0	4.0	5.0
					
	HGTFA4515MHW	HGTFA4525MHW	HGTFA4535MHW	HGTFA4545MHW	HGTFA4555MHW
	HGTFA4515MNW	HGTFA4525MNW	HGTFA4535MNW	HGTFA4545MNW	HGTFA4555MNW

D Ø5.0



EbonyGold screw
: HGATSM

Hex
Non-Hex **5.5**

H \ G/H	1.0	2.0	3.0	4.0	5.0
					
	HGTFA5015MHW	HGTFA5025MHW	HGTFA5035MHW	HGTFA5045MHW	-
	HGTFA5015MNW	HGTFA5025MNW	HGTFA5035MNW	HGTFA5045MNW	-



Transfer Abutment

D Ø4.5



EbonyGold screw
: HGATSS

	H \ G/H	1.0 2.0 3.0 4.0 5.0				
						
Hex	5.5	HGTFA4515SHW	HGTFA4525SHW	HGTFA4535SHW	HGTFA4545SHW	HGTFA4555SHW
Non-Hex		HGTFA4515SNW	HGTFA4525SNW	HGTFA4535SNW	HGTFA4545SNW	HGTFA4555SNW

D Ø5.0



EbonyGold screw
: HGATSS

	H \ G/H	1.0 2.0 3.0 4.0 5.0				
						
Hex	4.0	HGTFA5014SHW	HGTFA5024SHW	HGTFA5034SHW	HGTFA5044SHW	HGTFA5054SHW
	5.5	HGTFA5015SHW	HGTFA5025SHW	HGTFA5035SHW	HGTFA5045SHW	HGTFA5055SHW
	7.0	HGTFA5017SHW	HGTFA5027SHW	HGTFA5037SHW	HGTFA5047SHW	HGTFA5057SHW
Non-Hex	4.0	HGTFA5014SNW	HGTFA5024SNW	HGTFA5034SNW	HGTFA5044SNW	HGTFA5054SNW
	5.5	HGTFA5015SNW	HGTFA5025SNW	HGTFA5035SNW	HGTFA5045SNW	HGTFA5055SNW
	7.0	HGTFA5017SNW	HGTFA5027SNW	HGTFA5037SNW	HGTFA5047SNW	HGTFA5057SNW

D Ø6.0




EbonyGold screw
: HGATSS

	H \ G/H	1.0 2.0 3.0 4.0 5.0				
						
Hex	4.0	HGTFA6014SHW	HGTFA6024SHW	HGTFA6034SHW	HGTFA6044SHW	HGTFA6054SHW
	5.5	HGTFA6015SHW	HGTFA6025SHW	HGTFA6035SHW	HGTFA6045SHW	HGTFA6055SHW
	7.0	HGTFA6017SHW	HGTFA6027SHW	HGTFA6037SHW	HGTFA6047SHW	HGTFA6057SHW
Non-Hex	4.0	HGTFA6014SNW	HGTFA6024SNW	HGTFA6034SNW	HGTFA6044SNW	HGTFA6054SNW
	5.5	HGTFA6015SNW	HGTFA6025SNW	HGTFA6035SNW	HGTFA6045SNW	HGTFA6055SNW
	7.0	HGTFA6017SNW	HGTFA6027SNW	HGTFA6037SNW	HGTFA6047SNW	HGTFA6057SNW

GoldCast Abutment



- Screw-retained prosthesis
- Used for cases with path, aesthetic, and/or spatial constraints
- Prosthesis must be cast with dental grade gold alloy
- Melting Point : 1400 °C~1450 °C
(cast of non-precious metal alloy is prohibited)
- 1.2 hex driver
- Recommended tightening torque : 20Ncm(mini), 30Ncm(regular)
- Packing unit : abutment + EbonyGold screw



D Ø4.0



EbonyGold screw
: HGATSM

G/H Type	1.0		3.0	
	Hex	Non-Hex	Hex	Non-Hex
	HGGCA401MHW	HGGCA401MNW	HGGCA403MHW	HGGCA403MNW

D Ø4.5




EbonyGold screw
: HGATSS

G/H Type	1.0		3.0	
	Hex	Non-Hex	Hex	Non-Hex
	HGGCA451SHW	HGGCA451SNW	HGGCA453SHW	HGGCA453SNW



NP-Cast Abutment

- Screw-retained prosthesis
- Used for cases with path, aesthetic, and/or spatial constraints
- Prostheses must be cast with dental grade non-precious metal alloy
- Melting Point : 1400°C~1650°C
- 1.2 hex driver
- Recommended tightening torque : 20Ncm(mini), 30Ncm(regular)
- Packing unit : abutment + EbonyGold screw
- **Contains nickel:** Not recommended for patients with nickel allergy



D Ø4.0



EbonyGold screw
: HGATSM

G/H Type	1.0		3.0	
	Hex	Non-Hex	Hex	Non-Hex
	HGNCA401MHW	HGNCA401MNW	HGNCA403MHW	HGNCA403MNW

D Ø4.5




EbonyGold screw
: HGATSS

G/H Type	1.0		3.0	
	Hex	Non-Hex	Hex	Non-Hex
	HGNCA451SHW	HGNCA451SNW	HGNCA453SHW	HGNCA453SNW

Temporary Abutment



- For temporary prostheses (Material: Ti Gr-3)
- Easily customizable
- 1.2 hex driver
- Recommended tightening torque : 20Ncm(mini/regular)
- Packing unit : abutment + Ti screw



D Ø4.0



Ti screw
: HGATSMT

G/H Type	1.0		3.0	
	Hex	Non-Hex	Hex	Non-Hex
	HGTPA401MHT	HGTPA401MNT	HGTPA403MHT	HGTPA403MNT

D Ø4.5




Ti screw
: HGATSST

G/H Type	1.0		3.0	
	Hex	Non-Hex	Hex	Non-Hex
	HGTPA451SHT	HGTPA451SNT	HGTPA453SHT	HGTPA453SNT



Quick Temporary Abutment

- Material : Medical grade PEEK
- Used to fabricate a temporary prostheses for immediate loading
- Easy to modify
- Titanium base provides a stable connection to the fixture
- Use up to 180 days in the oral cavity
- Recommended tightening torque : 20Ncm(mini/regular)
- Packing unit : abutment + Ti screw



D Ø4.5



Ti screw
: HGATSMT

G/H
Type

5.0

Hex

Non-Hex



HGQTA455MHT

HGQTA455MNT

D Ø5.5



Ti screw
: HGATSST

G/H
Type

5.0

Hex

Non-Hex




HGQTA555SHT

HGQTA555SNT

Angled Abutment



- For path adjustment up to 17° axial angle
- Gold colored for aesthetics
- 12 orientation choices with two types of hexes, A and B
- Select appropriate hex type by using an Angled Abutment Selector
- 1.2 hex driver
- Recommended tightening torque : 20Ncm(mini), 30Ncm(regular)
- Packing unit : abutment + EbonyGold screw



ET SYSTEM

033

D Ø4.3



EbonyGold screw
: HGATSM

G/H Type	2.0			4.0		
	Hex A	Hex B	Non-Hex	Hex A	Hex B	Non-Hex
	HGAGA432MAW	HGAGA432MBW	HGAGA432MNW	HGAGA434MAW	HGAGA434MBW	HGAGA434MNW

D Ø4.5



EbonyGold screw
: HGATSS

G/H Type	2.0			4.0		
	Hex A	Hex B	Non-Hex	Hex A	Hex B	Non-Hex
	HGAGA452SAW	HGAGA452SBW	HGAGA452SNW	HGAGA454SAW	HGAGA454BW	HGAGA454SNW



Angled Abutment

D Ø5.5



EbonyGold screw
: HGATSS

G/H Type	2.0			4.0		
	Hex A	Hex B	Non-Hex	Hex A	Hex B	Non-Hex
						
	HGAGA552SAW	HGAGA552SBW	HGAGA552SNW	HGAGA554SAW	HGAGA554SBW	HGAGA554SNW



Angled Abutment Components

Angled Abutment Selector

- Used for selecting Angled Abutment A or B type, diameter, and G/H when in oral cavity or in a working model



Mini




Regular

D	G/H Type	2.0		4.0	
		Hex A	Hex B	Hex A	Hex B
Ø 4.3					
		HGAAS432MA	HGAAS432MB	HGAAS434MA	HGAAS434MB
Ø 4.5					
		HGAAS452SA	HGAAS452SB	HGAAS454SA	HGAAS454SB
Ø 5.5					
		HGAAS552SA	HGAAS552SB	HGAAS554SA	HGAAS554SB

ZioCera Abutment



- Zirconia material with superior biocompatibility
- Mimics color of natural tooth
- 1.2 hex driver
- Recommended tightening torque : 20Ncm(mini), 30Ncm(regular)
- Packing unit : abutment + EbonyGold screw



D Ø4.5



EbonyGold screw
: HGABSM

G/H Type	3.5		5.0	
	Hex	Non-Hex	Hex	Non-Hex
				
	HGZCA453MHW	HGZCA453MNW	HGZCA455MHW	HGZCA455MNW

D Ø4.5



EbonyGold screw
: HGABSS

G/H Type	3.5		5.0	
	Hex	Non-Hex	Hex	Non-Hex
				
	HGZCA453SHW	HGZCA453SNW	HGZCA455SHW	HGZCA455SNW



D Ø5.5



EbonyGold screw
: HGABSS

G/H Type	3.5		5.0	
	Hex	Non-Hex	Hex	Non-Hex
				
	HGZCA553SHW	HGZCA553SNW	HGZCA555SHW	HGZCA555SNW

D Ø6.5




EbonyGold screw
: HGABSS

G/H Type	3.5		5.0	
	Hex	Non-Hex	Hex	Non-Hex
				
	HGZCA653SHW	HGZCA653SNW	HGZCA655SHW	HGZCA655SNW



ZioCera Angled Abutment

- For path adjustment up to 17° axial angle
- Mimics color of natural tooth
- 1.2 hex driver
- Recommended tightening torque : 30Ncm(regular)
- Packing unit : abutment + EbonyGold screw



D Ø5.5



EbonyGold screw
: HGABSS

G/H 3.0
Type Hex Non-Hex



HGZAA553SHW



HGZAA553SNW

D Ø6.5



EbonyGold screw
: HGABSS

G/H 4.0
Type Hex Non-Hex



HGZAA654SHW




HGZAA654SNW

FreeForm ST Abutment



- Used to alter the path of an abutment or to customize the margin area
- Gold colored gingiva region for aesthetics
- Recommended tightening torque : 20Ncm(mini), 30Ncm(regular)
- Packing unit : abutment + EbonyGold screw



D Ø4.0



EbonyGold screw
: HGATSM

G/H Type	1.5		3.0	
	Hex	Non-Hex	Hex	Non-Hex
	HGFS401MHW	HGFS401MNW	HGFS403MHW	HGFS403MNW

D Ø4.0



EbonyGold screw
: HGATSS

G/H Type	1.5		3.0	
	Hex	Non-Hex	Hex	Non-Hex
	HGFS401SHW	HGFS401SNW	HGFS403SHW	HGFS403SNW



FreeForm ST Abutment

D Ø5.5



EbonyGold screw
: HGATSS

G/H Type	1.5		3.0	
	Hex	Non-Hex	Hex	Non-Hex
				
	HGFA551SHW	HGFA551SNW	HGFA553SHW	HGFA553SNW

D Ø7.0



EbonyGold screw
: HGATSS

G/H Type	1.5		3.0	
	Hex	Non-Hex	Hex	Non-Hex
				
	HGFA701SHW	HGFA701SNW	HGFA703SHW	HGFA703SNW

Impression Components



Fixture Transfer Impression Coping

- For closed tray impressions
- Accurate impression with gemstone-shaped structure ()
- 1.2 hex driver
- Packing unit : hex - impression coping + guide pin
non-hex - impression coping(pre-assembled)

M

Mini

R

Regular

D \ L	11		14	
	Hex	Non-Hex	Hex	Non-Hex
				
Ø 4.0	HGFTI4011MH	HGFTI4011MN	HGFTI4014MH	HGFTI4014MN
Ø 4.0	HGFTI4011SH	HGFTI4011SN	HGFTI4014SH	HGFTI4014SN
Ø 5.0	HGFTI5011SH	HGFTI5011SN	HGFTI5014SH	HGFTI5014SN
Ø 6.0	HGFTI6011SH	HGFTI6011SN	HGFTI6014SH	HGFTI6014SN



Impression Components

Fixture Pick-up Impression Coping

- For open tray impressions
- Extend flag with a hole adds superior impression stability
- 1.2 hex driver
- Packing unit : Impression coping body + guide pin

M Mini

R Regular

D \ L		11		Guide Pin		
Type		Hex	Non-Hex	0	5.0	
Ø 4.0 Ø 4.0 Ø 5.0 Ø 6.0						
		HGFPI4011MH	HGFPI4011MN	-	HGFPG00M	HGFPG05M*
		HGFPI4011SH	HGFPI4011SN	]	HGFPG00S	HGFPG05S*
		HGFPI5011SH	HGFPI5011SN			
		HGFPI6011SH	HGFPI6011SN			

D \ L		15		Guide Pin		
Type		Hex	Non-Hex	0	5.0	
Ø 4.0 Ø 4.0 Ø 5.0 Ø 6.0						
		HGFPI4015MH	HGFPI4015MN	-	HGFPG00ML	HGFPG05ML*
		HGFPI4015SH	HGFPI4015SN	]	HGFPG00SL	HGFPG05SL*
		HGFPI5015SH	HGFPI5015SN			
		HGFPI6015SH	HGFPI6015SN			



Fixture Lab Analog

- Replicates fixture's position in a working model



Mini



Regular



HGFLAM



HGFLAS



PROSTHETIC FLOW DIAGRAM 3


Rigid / Transfer

Abutment Level Impression

www.befrest.com

ET SYSTEM


044



Rigid Abutment



- Cement-retained prosthetic
- Gold colored gingiva region for improved aesthetics
- Outer driver : $\varnothing 4.0$
- Outer driver or 1.2 hex driver : $\varnothing 4.5$, $\varnothing 5.0$, $\varnothing 6.0$
- 1.2 hex driver : $\varnothing 7.0$
- Recommended tightening torque : 30Ncm
- Packing unit : abutment + protect cap



ET SYSTEM

045

D $\varnothing 4.0$

M

H \ G/H	1.0	2.0	3.0	4.0
				
4.0	HGRGA4014MP	HGRGA4024MP	HGRGA4034MP	HGRGA4044MP
5.5	HGRGA4015MP	HGRGA4025MP	HGRGA4035MP	HGRGA4045MP
7.0	HGRGA4017MP	HGRGA4027MP	HGRGA4037MP	HGRGA4047MP

D $\varnothing 4.5$

M

H \ G/H	1.0	2.0	3.0	4.0
				
4.0	HGRGA4514MP	HGRGA4524MP	HGRGA4534MP	HGRGA4544MP
5.5	HGRGA4515MP	HGRGA4525MP	HGRGA4535MP	HGRGA4545MP
7.0	HGRGA4517MP	HGRGA4527MP	HGRGA4537MP	HGRGA4547MP



Rigid Abutment

D Ø4.0



H \ G/H	1.0	2.0	3.0	4.0
				
4.0	HGRGA4014SP	HGRGA4024SP	HGRGA4034SP	HGRGA4044SP
5.5	HGRGA4015SP	HGRGA4025SP	HGRGA4035SP	HGRGA4045SP
7.0	HGRGA4017SP	HGRGA4027SP	HGRGA4037SP	HGRGA4047SP

D Ø4.5



H \ G/H	1.0	2.0	3.0	4.0
				
4.0	HGRGA4514SP	HGRGA4524SP	HGRGA4534SP	HGRGA4534SP
5.5	HGRGA4515SP	HGRGA4525SP	HGRGA4535SP	HGRGA4535SP
7.0	HGRGA4517SP	HGRGA4527SP	HGRGA4537SP	HGRGA4537SP

D Ø5.0



H \ G/H	1.0	2.0	3.0	4.0	5.0
					
4.0	HGRGA5014SP	HGRGA5024SP	HGRGA5034SP	HGRGA5044SP	HGRGA5054SP
5.5	HGRGA5015SP	HGRGA5025SP	HGRGA5035SP	HGRGA5045SP	HGRGA5055SP
7.0	HGRGA5017SP	HGRGA5027SP	HGRGA5037SP	HGRGA5047SP	HGRGA5057SP



D Ø6.0



H \ G/H	1.0	2.0	3.0	4.0	5.0
					
4.0	HGRGA6014SP	HGRGA6024SP	HGRGA6034SP	HGRGA6044SP	HGRGA6054SP
5.5	HGRGA6015SP	HGRGA6025SP	HGRGA6035SP	HGRGA6045SP	HGRGA6055SP
7.0	HGRGA6017SP	HGRGA6027SP	HGRGA6037SP	HGRGA6047SP	HGRGA6057SP

D Ø7.0



H \ G/H	1.0	2.0	3.0	4.0	5.0
					
5.5	HGRGA7015SP	HGRGA7025SP	HGRGA7035SP	HGRGA7045SP	HGRGA7055SP



Rigid Abutment Components

Rigid Protect Cap

- Used to protect rigid abutments in the oral cavity and minimize foreign body sensation by the patient
- Can be used as a substructure for temporary prosthesis

M Mini
R Regular

D \ H	4.0	5.5	7.0
			
Ø 4.0 / Ø 4.0	HGRPC4004	HGRPC4005	HGRPC4007
Ø 4.5 / Ø 4.5	HGRPC4504	HGRPC4505	HGRPC4507
Ø 5.0	HGRPC5004S	HGRPC5005S	HGRPC5007S
Ø 6.0	HGRPC6004S	HGRPC6005S	HGRPC6007S
Ø 7.0	-	HGRPC7005S	-

Rigid Impression Coping

- Used in taking impressions of rigid abutments

M Mini
R Regular

D \ H	4.0	5.5	7.0
			
Ø 4.0 / Ø 4.0	HGRIC4004S	HGRIC4005S	HGRIC4007S
Ø 4.5 / Ø 4.5	HGRIC4504S	HGRIC4505S	HGRIC4507S
Ø 5.0	HGRIC5004SS	HGRIC5005SS	HGRIC5007SS
Ø 6.0	HGRIC6004SS	HGRIC6005SS	HGRIC6007SS
Ø 7.0	-	HGRIC7005SS	-

Rigid Burn-out Cylinder

- Used in attaching a rigid lab analog in framework for prosthetics
- After casting a prosthetic, margin area is adjusted using specialized reamer

M Mini
R Regular

D \ Type	Single	Bridge
		
Ø 4.0 / Ø 4.0	HGRBC40S	HGRBC40B
Ø 4.5 / Ø 4.5	HGRBC45S	HGRBC45B
Ø 5.0	HGRBC50SS	HGRBC50SB
Ø 6.0	HGRBC60SS	HGRBC60SB
Ø 7.0	HGRBC70SS	HGRBC70SB



Rigid Lab Analog

- Replicates shape of rigid abutment in a working model

M Mini
R Regular

D \ H	4.0	5.5	7.0
			
Ø 4.0 / Ø 4.0	HGRLA4004	HGRLA4005	HGRLA4007
Ø 4.5 / Ø 4.5	HGRLA4504	HGRLA4505	HGRLA4507
Ø 5.0	HGRLA5004S	HGRLA5005S	HGRLA5007S
Ø 6.0	HGRLA6004S	HGRLA6005S	HGRLA6007S
Ø 7.0	-	HGRLA7005S	-



PROSTHETIC FLOW DIAGRAM 4


Stud / Locator®

Overdenture

www.befrest.com

ET SYSTEM


050



Locator[®] Abutment



- Low vertical dimension, stability, and various attachments with retention
- Accommodates up to 40° divergence between two implants
- Tighten by using a locator torque driver
- Recommended tightening torque : 30Ncm



D Ø3.7

M

G/H	1.0	2.0	3.0	4.0	5.0
					
	HGLCA3510M	HGLCA3520M	HGLCA3530M	HGLCA3540M	HGLCA3550M

D Ø3.7

R

G/H	1.0	2.0	3.0	4.0	5.0
					
	HGLCA4010S	HGLCA4020S	HGLCA4030S	HGLCA4040S	HGLCA4050S



Locator[®] Abutment Components

Locator[®] Male Processing Kit

- Component
 - Block-out spacer / denture cap connected black processing male
 - Replacement male: blue/pink/clear
- Used after selecting retention males that are appropriate for the case
- Use Core tool to replace males
- Packing unit : 2 sets



Locator[®] Replacement Male

- Retention: Approximately 6N
- Allows up to 20 degrees of divergence between implants
- Packing unit : blue replacement male (4 ea)



- Retention: Approximately 12N
- Allows up to 20 degrees of divergence between implants
- Packing unit : pink replacement male (4 ea)




- Retention: Approximately 22N
- Allows up to 20 degrees of divergence between implants
- Packing unit : clear replacement male (4 ea)



Locator[®] Extended Replacement Male

- Retention: Approximately 6N
- Allows up to 40 degrees of divergence between implants
- Packing unit : red extended replacement male (4 ea)



- Retention: Approximately 12N
- Allows up to 40 degrees of divergence between implants
- Packing unit : green extended replacement male (4 ea)





Locator® Black Processing Male

- Used in lab. process
- Packing unit : 4 ea



LBPS

Locator® Block Out Spacers

- Seals gap between denture cap and abutment
- Packing unit : 20 ea



LBSS

Locator® Impression Coping

Used to take impressions after attaching locator abutment
Packing unit : 4 ea



LICS

Locator® Lab Analog

- Replicates Locator Abutment in a working model
- Packing unit : 4 ea



LAL40S



Locator[®] Abutment Components

Locator[®] Core Tool

- Used in attaching and changing replacement males



LOCT

Locator[®] Torque Driver

- Used to tighten locator abutment

Type	Short	Long
------	-------	------



TWLDSK




TWLDLK

Stud Abutment



- Stud type overdenture prostheses
- Compensates the path up to 20°
- O-ring abutment driver (AORD)
- Recommended tightening torque : 30Ncm



D Ø3.5

M

G/H	1.0	2.0	3.0	4.0	5.0	6.0
						
	HGSAO351M	HGSAO352M	HGSAO353M	HGSAO354M	HGSAO355M	HGSAO356M

D Ø3.5

R

G/H	1.0	2.0	3.0	4.0	5.0	6.0
						
	HGSAO351S	HGSAO352S	HGSAO353S	HGSAO354S	HGSAO355S	HGSAO356S



Stud Abutment Components

O-ring Retainer Cap Set

- Retaining system for stud type overdenture prostheses
- Packing unit : retainer cap + o-ring



OARCS

O-ring Retainer Set

- For low occlusal clearance
- Packing unit : retainer + o-ring



OARS

O-ring Set

Packing unit : 5 ea



OAON01S

O-ring Lab Analog

Replicates O-ring abutment in a working model



OAL




HROSSEN



PROSTHETIC FLOW DIAGRAM 5

Multi / Multi Angled


Abutment Level Impression



Multi Abutment



- Uses specialized outer driver (MAOD)
- Used in screw type prostheses for multiple cases
- Prostheses is fabriacted with US aesthetic-low cylinders
- Recommended tightening torque : 30Ncm



D Ø4.8



G/H	1.0	2.0	3.0	4.0	5.0
					
	ETMTA501M	ETMTA502M	ETMTA503M	ETMTA504M	ETMTA505M

D Ø4.8



G/H	1.0	2.0	3.0	4.0	5.0
					
	ETMTA501R	ETMTA502R	ETMTA503R	ETMTA504R	ETMTA505R



Multi Abutment Components

Multi Abutment Outer Driver

- Torque driver specialized for multi-abutments



MAOD

Multi Angled Abutment Screw

- 1.2 hex driver
- For Multi Angled abutments only

M Mini

R Regular

M



HGMATSM

R



HGMATSS

Multi Angled Abutment Carrier

- Provides aid in carrying and connecting multi angled abutment to the fixture.



HGMAC



Esthetic-low Healing Cap

- Used when protecting aesthetic-low abutment in the oral cavity and minimizing foreign body sensation by the patient
- 1.2 hex driver
- Recommended tightening torque : 20Ncm

R Regular

D \ H

6.0

Ø 4.8



MHCR100

Esthetic-low Gold Cylinder

- Cast with dental-grade gold alloy
- Melting Point : 1400°C ~ 1450°C
(cast using non-precious metal alloy is prohibited)
- 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti screw


R Regular

D \ Type

Non-Hex

Ø 4.8

Ti screw
: MTS200 (Ø 4.8)



MGR100

Esthetic-low Plastic Cylinder

- Cast with dental-grade alloy (gold, non-precious metals) after customization
- Connection area slightly different compared to gold cylinder
- 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti screw


R Regular

D \ Type

Non-Hex

Ø 4.8

Ti screw
: MTS200 (Ø 4.8)



MEPR100



Multi Abutment Components

Esthetic-low Temporary Cylinder

- Temporary prostheses (Material: Ti Gr-3)
- 1.2 hex driver
- Recommended tightening torque : 20Ncm
- Packing unit : cylinder + Ti screw


R Regular

D \ Type

Non-Hex

Ø 4.8

Ti screw
: MTS200 (Ø 4.8)



MTR200TH

Esthetic-low Pick-up Impression Coping

- Open tray impression
- Holed flag has superior impression stability
- Packing unit : impression coping body + guide pin

R Regular

D \ L

Non-Hex

Guide Pin
15

Ø 4.8



MSR100



GP150



Esthetic-low Transfer Impression Coping

- Closed tray impression

R Regular

D \ H

8.0

Ø 4.8



MTTR100

Esthetic-low Lab Analog

- Replicates abutment's position in a working model

R Regular

D

Ø 4.8



MER300

Esthetic-low Polishing Protector

- Used to prevent damage to the connection area of cylinder when polishing.

R Regular

D

Ø 4.8




MPCR100



Multi Angled Abutment

- UScrewed retained prostheses
- Compensation path in multiple implant case
- Up to 60° path compensation (two implant standard)
- Prosthetic fabrication on US aesthetic-low cylinder(non-hex)
- Recommended tightening torque : 20Ncm(mini), 30Ncm(regular)
- Packing unit : multi angled abutment + EbonyGold screw



D Ø4.8



EbonyGold screw
: GSMABSM

Angle \ G/H	2.5	3.0	4.0
17°	 HGMAA217MHW	 HGMAA317MHW	 HGMAA417MHW
Angle \ G/H	3.5	4.0	5.0
30°	 HGMAA330MHW	 HGMAA430MHW	 HGMAA530MHW

D Ø4.8



EbonyGold screw
: GSMABSS

Angle \ G/H	2.5	3.0	4.0
17°	 HGMAA217SHW	 HGMAA317SHW	 HGMAA417SHW
Angle \ G/H	3.5	4.0	5.0
30°	 HGMAA330SHW	 HGMAA430SHW	 HGMAA530SHW

Abutment Screw



EbonyGold Abutment Screw

- TypeA - Transfer Abutment
Angled Abutment
FreeForm ST Abutment
GoldCast Abutment
- Type B - ZioCera Abutment

M Mini

R Regular

Type A

Type B

M



HGATSM



HGABSM

R



HGATSS



HGABSS

Titanium Abutment Screw

- Applicable abutments
- Temporary abutment

M Mini

R Regular

M



HGATSMT



HGABSM

R



HGATSST



HGABSST



Lab Screw

Laboratory Screw

- Lab Screw : Use for laboratory work instead of abutment screw.
- Waxing Screw : Use for making a screw hole in a transfer jig or wax-up part.
- Applicable abutments
 - Transfer abutment
 - Angled abutment
 - FreeForm ST abutment
 - GoldCast abutment

M Mini

R Regular

Lab Screw

Waxing Screw

M



HGATSMML



HGATSMW

R



HGATSSL



HGATSSW

HROSSEN



HROSSEN





Hiossen KITS

070	Taper KIT	103	OsstemGuide KIT
071	Standard KIT	106	MS KIT
072	Taper Mini KIT	107	ESR KIT
073	Taper Ultra KIT	111	EFR KIT
074	Ultra KIT	114	Osteo KIT
075	123 KIT	115	Sinus KIT
087	Prosthetic KIT	116	Bone Spreader KIT
093	CAS-KIT	117	Ridge Split KIT - Straight
097	LAS-KIT	118	Ridge Split KIT - Offset
100	ESSET KIT		



Taper KIT (HTAPEK)

Available for use

ETIII


SSIII

USIII

www.befrest.com

070

HIOSEN KITS



Standard KIT (HSTDDK)




Available for use

ETII / III TSII / III

SSII / III

USII / III



071

HIOSEN KITS



Taper Mini KIT (HTAPEMK)


Available for use

ETII / III TSII / III

SSII / III

USII / III

Ultra-wide



Taper Ultra KIT (HULTPK)



Base components

Open Wrench
ASOW




Ratchet Wrench
CITQW-1185A



Available for use

Ultra-wide

www.befrest.com



Ø6.0 Taper Ultra Drill



TPD3C6006, TPD3C6007
TPD3C6008, TPD3C6010
TPD3C6011, TPD3C6013

Ø7.0 Taper Ultra Drill



TPD3C7006, TPD3C7007
TPD3C7008, TPD3C7010
TPD3C7011, TPD3C7013

073

HIOSEN KITS

Ultra KIT (HULTRK)



Base components

Open Wrench
ASOW




Ratchet Wrench
CITQW-1185A



Available for use

Ultra-wide



123 KIT (O123K)




Available for use

ETII / III TSII / III

SSII / III

USII / III



075

HIOSEN KITS




Surgical Instruments

Lance Drill - Guide Drill

Mark the initial osteotomy position

- Bone density can be judged through drilling

L	Short	Long
	AGDSC	AGDLC



123 Guide Drill

- Drill that forms a hole in the bone to make initial drilling easy
- Easy depth control with stoppers


D	Ø2.0
	OGD2027L



Sidecut Drill

- Correct the angle or position of the osteotomy cutting edge

L	D1 / D2	Ø1.5 / 2.0	Ø2.0 / 2.5	Ø3.0 / 3.5
13		OSLMDS		-
16.5		-	-	OSLMD30L
20		OSLMDL	OSLMD20L	-



Drill Extension

- Tool used to extend the length of drills and other miscellaneous hand piece tools
- It may bend or break if excessive force is applied
- Drill length is extended by 16.9mm with drill extension


ODE





Twist Drill - Stopper Drill


- Long stopper(6mm)
- Built-in stopper
- Procedure possible without drill extension when performing posterior procedures
- Color coded stopper indicates the drill length



L \ TL	D	Ø2.2	Ø3.0	Ø3.3	Ø3.6	Ø3.8	Ø4.1	Ø4.3	Ø4.6
6.0	30.5	2D2206LC	3D3006LC	-	-	3D3806LC	-	-	-
7.0	31.5	2D2207LC	3D3007LC01	-	-	3D3807LC01	-	-	-
8.5	33	2D2208LC	3D3008LC01	-	-	3D3808LC01	-	-	-
10.0	34.5	2D2210LC	3D3010LC01	-	-	3D3810LC01	-	-	-
11.5	34.5	2D2211LC	3D3011LC01	3D3311LC01	3D3611LC01	3D3811LC01	3D4111LC01	3D4311LC01	3D4611LC01
13.0	36	2D2213LC	3D3013LC01	-	-	3D3813LC01	-	-	-
Y-Dim.		0.6	0.9	1.0	1.0	1.0	1.0	1.0	1.0

Twist Drill - Non Stopper Drill

- A drill that can be used when stopper drill's accessibility in the oral cavity is limited
- Marking drill with short and long specifications
- Refer to the image for marking measurements



TL \ D	Ø1.5	Ø2.0	Ø2.2	Ø2.7	Ø3.0	Ø3.3
33	2D1518FNLC	2D2018FNLC	2D2218FNLC	3D2718FNLC	3D3018FNLC	3D3318FNLC
41	-	-	2D2215FNLC	3D2715FNLC01	3D3015FNLC01	3D3315FNLC01


TL \ D	Ø3.6	Ø3.8	Ø4.1	Ø4.3	Ø4.6
33	3D3618FNLC	3D3818FNLC	3D4118FNLC	3D4318FNLC	3D4618FNLC
41	3D3615FNLC01	3D3815FNLC01	3D4115FNLC01	3D4315FNLC01	3D4615FNLC01



Surgical Instruments

123 Twist Drill

- Color coded with corresponding fixture sizes
- Easy depth control with drill stopper
- Designed to be used with drill stopper
- F = Fixture



TL	D1 / D2 / D3			
	F3.5(Ø2.2 / 3.0)	F4.0(Ø3.0 / 3.6)	F4.5(Ø3.0 / 3.6 / 4.1)	F5.0(Ø3.0 / 4.1 / 4.6)
34	2D2230FNS	2D3036FNS	2D3041FNS	2D3046FNS
40.4	2D2230FNL	2D3036FNL	2D3041FNL	2D3046FNL
Color	Yellow	Green	Blue	Red

078

HIOSEN KITS

123 Drill Stopper

Available Soon


- The length of the Stopper indicates the drilling depth
- Color coded for each length

L	6.2	7	8	9.5	11	12.5	14	16
	ODST05	ODST06	ODST07	ODST08	ODST10	ODST11	ODST13	ODST15



Direct Drill


- Direct drill: Dual-gear drill that has both pilot drill and twist drill capabilities
- Final drilling can be done right away without pilot drilling
- To prepare the apex area of the osteotomy after extraction socket



D1 / D2	Ø4.6 / 5.2	Ø4.6 / 5.5	Ø5.5 / 6.2	Ø5.5 / 6.5
	3D5213FNLC	3D5513FNLC	3D6213FNLC	3D6513FNLC

Taper Drill

- Taper drill specialized for taper (III type) fixture with built-in stopper
- Stopper drill with 1mm margin for countersink
- Color coded grip part indicates fixture diameter
- F = Fixture



L	TL	F3.5	F4.0	F4.5	F5.0
6.0	30.5	TPTD3C3506T	TPTD3C4006T	TPTD3C4506T	TPTD3C5006T
7.0	31.5	TPTD3C3507T	TPTD3C4007T	TPTD3C4507T	TPTD3C5007T
8.5	33	TPTD3C3508T	TPTD3C4008T	TPTD3C4508T	TPTD3C5008T
10.0	34.5	TPTD3C3510T	TPTD3C4010T	TPTD3C4510T	TPTD3C5010T
11.5	34.5	TPTD3C3511T	TPTD3C4011T	TPTD3C4511T	TPTD3C5011T
13.0	36	TPTD3C3513T	TPTD3C4013T	TPTD3C4513T	TPTD3C5013T
15.0	38	TPTD3C3515T	TPTD3C4015T	TPTD3C4515T	TPTD3C5015T
Y-Dim.		0.8	0.9	1.0	1.0
Color		Yellow	Green	Blue	Red




Surgical Instruments

Taper Ultra Drill

- Taper drill specialized for taper ultra-wide fixture with built-in stopper
- Stopper drill with 1mm margin for countersink
- Color coded grip part indicates fixture diameter
- F = Fixture


L	F6.0	F7.0
6	TPTD3C6006T	TPTD3C7006T
7	TPTD3C6007T	TPTD3C7007T
8.5	TPTD3C6008T	TPTD3C7008T
10	TPTD3C6010T	TPTD3C7010T
11.5	TPTD3C6011T	TPTD3C7011T
13	TPTD3C6013T	TPTD3C7013T
Color	Green	Blue



Long Shank Pilot Drill

- Guide drill for the next drill sequence
- Can be used to alter the path

D1 / D2	Ø2.0 / 2.7	Ø2.0 / 3.0	Ø3.0 / 3.8	Ø3.0 / 4.1
	APD270C	APD300C	APD380C	APD410C



Cortical Drill 3 for Taper Fixture (ETIII, SSIII, USIII)

- Drill used to remove cortical bone in hard bone (for III type)
- Bottom band is for normal bone, and the top is for hard bone
- Recommend drilling to the bottom of the band
- F = Fixture

F3.5	F4.0	F4.5	F5.0
CD4C35	CD4C40	CD4C45	CD4C50





Taper Cortical Drill for Taper Fixture (ETIII, SSIII, USIII)

- Drill used to remove cortical bone in hard bone (used immediately after taper drill)
- Bottom of the band is for fixture insertion less than 8.5mm
Top of the band is for fixture insertion more than 10mm
- Recommend drilling to the bottom of the band
- F = Fixture


F3.5	F4.0	F4.5	F5.0
TCD4C35	TCD4C40	TCD4C45	TCD4C50



Cortical Drill for Ultra-Wide

- Drill used to remove cortical bone in hard bone (for ultra-wide)
- For use for different fixture diameters
- Recommend drilling to the bottom of the band
- F = Fixture

F6.0	F7.0
CD4C60	CD4C70




081

HIOSEN KITS

Parallel Pin

- Check the position and orientation of osteotomy after the initial drilling sequence

D	Ø4.0	Ø5.0	Ø6.0	Full Set
APP400	APP500	APP600	APPS	






Surgical Instruments


Depth Gauge

- A: to measure drill depth (7~15mm)
- B: to measure soft tissue height



Depth Gauge (NEW)

- Combination of depth gauge and an open wrench



Torque Wrench - Bar Type


- Used to adjust the implant location or tighten abutment, screw, etc.
- Apply torque by pulling bar to the line indicating the torque value



NoMount Driver for ET

- Driver that can directly attach to fixture when inserting, using hand piece for procedure
- C = Connection

L \ C	Mini	Regular
Short	TSNMDMS	TSNMDRS
Long	TSNMDML	TSNMDRL






NoMount Torque Driver for ET

- Directly attaches to the fixture for insertion with wrench
- Check for correct and complete fit before insertion.
- Incorrect usage may cause it to break.
- C = Connection


L \ C	Mini	Regular
Short	GSNMT32S	GSNMT35S
Long	GSNMT32L	GSNMT35L



Fixture Driver for ET

- Directly attached to fixture, used to adjust final depth of insertion using wrench
- C = Connection

L \ C	Mini	Regular
Short	GSMFDS	GSRFDS
Long	GSMFDL	GSRFDL



Torque Extension


- Extends the length of instruments by 10mm, connects to wrench



Simple Mount Driver

- Used to insert fixture with mount using hand piece

L \ C	
Short	ASMDS
Long	ASMDL




Surgical Instruments



Simple Mount Extension

- Used for manual torque application during fixture placement

L	
Short	ASMES
Long	ASMEL



Removal Tool for Fixture Mount

- Used for removing the mount if the mount is stuck. Remove the mount screw before use.
- Used with Driver Handle or Torque Wrench
- Insert into mount and turning clockwise
- App = Application




App	Mini (ET,US)	Regular (ET,SS,US) / Wide (SS)	Wide (US)
	ERFM	HRFR	ERFW



Tissue Punch

- Tool used for flapless surgery procedure
 - Laser marking at 2mm interval, easily measures gingiva height
 - Packing unit : tissue punch + guide pin
- ※ Recommend using use tissue punch with a 0.7~1.5mm smaller diameter than the healing abutment



D	Ø3.3	Ø3.8	Ø4.3	Ø4.8	Ø5.3
	OSTP33	OSTP38	OSTP43	OSTP48	OSTP53
ET	Ø 4.0/4.5	Ø 4.5/5.0	Ø 5.0	Ø 6.0	Ø 6.0
SS	-	Ø 4.8	-	Ø 6.0	Ø 6.0
US	Ø 4.0	Ø 5.0	Ø 5.0	Ø 6.0	Ø 6.0

The application of Healing abutment

085

HIOSEN KITS

Bone Profiler

- Used to remove bone around fixture after first or second stage surgery
- Connect Guide Screw to the fixture and remove excess bone
- Guide screw protects fixture shoulder
- Packing unit : bone profiler + guide screw
- C = Connection




C \ Healing Abutment D	Ø4.5	Ø5.5	Ø6.5/7.5
Mini / Regular	GSBP45	GSBP55	GSBP75
	Mini + Regular guide screw	Mini + Regular guide screw	Regular guide screw

Surgical Instruments



Trephine Drill


- Used to collect bone or remove a failed fixture
- Used to remove septal bone
- Can be used as an initial drill when inserting an ultra fixture



L \ D (Inner / Outer)	3.7 / 4.5	4.2 / 5.0	4.7 / 5.5	5.2 / 6.0	5.7 / 6.5	6.2 / 7.0
Short	TD37S	TD42S	TD47S	TD52S	TD57S	TD62S

Bone Mill

- Grinds collected autogenous bone



Prosthetic KIT (HPRSTKA)




Available for use

ETII / III TSII / III

SSII / III

USII / III

Ultra-wide



087


HIOSEN KITS

Prosthetic Instruments




Hand Driver

- Manual driver
- Internal hex type L is 11mm



0.9 hex




1.2 hex

L \ Type	0.9 Hex	1.2 Hex
Ex.Short (8)	AHD09MSH	AHD12MSH
Short (13)	AHD09SH	AHD12SH
Middle (15)	-	AHD12MH
Long (18)	AHD09LH	AHD12LH
Ex.Long (25)	-	AHD12EH

Machine Screw Driver

- 1.2 hex driver for engine hand piece
- Internal hex type L is 8mm



L \ Type	0.9 Hex	1.2 Hex
Short (5.6)	AMSD09S	AMSD12S
Long (11.6)	AMSD09L	AMSD12L
Ex.Long (17.6)	-	AMSD12E

Application

Driver applicable product

(common for hand, machine screw, and torque driver)

Cover screw (US mini)

Healing abutment, UCLA, Cemented abutment screw, Mount screw


Esthetic abutment screw regular, Esthetic-low abutment screw, standard

Wide esthetic-low abutment screw



Torque Driver


- Driver for torque wrench
- May bend or break if excessive torque is applied
- Damage is possible even at low torque if not fully engaged.
- When applying torque, apply vertically pressure (do not tilt)
- If the tip is bent or stripped due replace immediately



L \ Type	0.5 Slot	0.9 Hex	1.2 Hex
Ex.Short (8)	-	-	TRHD12MS
Short (13)	TRSD05S	TRHD09S	TRHD12S
Middle (15)	-	-	TRHD12M
Long (20)	TRSD05L	TRHD09L	TRHD12L
Ex.Long (25)	TRSD05E	-	TRHD12E

O-ring Abutment Driver

- Driver specific for O-ring abutments



Rigid Outer Driver

- Driver specific for rigid abutments
- Recommended tightening torque : 30Ncm

L \ Abutment D	Ø4.0	Ø4.5	Ø5.0	Ø6.0
Short (16.5)	ORDMS	ORD45S	ORDRS	ORDWS
Long (21.5)	ORDML	ORD45L	ORDRL	ORDWL






Prosthetic Instruments

Solid Abutment Driver

- Driver specific for solid abutments
- AConnect to the Solid Abutment by matching up the groove with the triangle indicator
- Recommended tightening torque : 30Ncm

Regular

L \ Type	Square	Round
Short (6)	 SDSS	 SDRS
Long (12)	 SDSL	 SDRL



Wide

L \ Type	Square
Short (10)	 SD60S

Excellent Solid Abutment Driver

- Driver specific for excellent solid abutments
- Connect to Excellent Solid Abutment by matching up the groove with the triangle indicator
- Recommended tightening torque : 30Ncm

Regular

L \ Type	Square	Round
Short (6)	 ESDSS	 ESDRS
Long (12)	 ESDSL	 ESDRLC



Wide


L \ Type	Square
Short (10)	 ESD60S



Octa Abutment Driver


- Driver specific for octa abutment
- Recommended tightening torque : 30Ncm

L \ Type	Square	Round
Short	ODSS	ODRSC
Long	ODSL	ODRL



Driver Handle

- Connects to torque driver




091

HIOSEN KITS


Finishing Reamer Set

- Tool used to remove the excess cast after plastic coping is cast



Reamer user guide

1. Select a Reamer Tip of the same size as the abutment and connect to burn-out cylinder after casting.
2. Hold casting body and turn reamer bite with consistent force
3. Ream until cutting stops/ends





Prosthetic Instruments

Reamer Bite

- Cutting edge removes excess cast after plastic coping is cast



FRBC

CAS-KIT (HCRSNK)



Base components

Bone Carrier Head
SNBCH30



Bone Carrier
SNBCS35



Bone Condenser
SNBC1114




Available for use

ETII / III TSII / III

SSII / III

USII / III



093

HIOSEN KITS




CAS-KIT Surgical Instruments

CAS-Drill

- Crestal Approach Sinus Drill forms a conical bone lid and lifts membrane safely
- Superior bone removing capability from low speed to high speed, collects autogenous bone at low speeds
- Safely advance to the floor of sinus using stoppers (1mm increment)
- Final drill diameter based on bone density, regardless of straight or tapered fixture type


D	Ø2.8	Ø3.1	Ø3.3	Ø3.6	Ø3.8	Ø4.1
	SNDR2813TL	SNDR3113TL	SNDR3313TL	SNDR3613TL	SNDR3813TL	SNDR4113TL



Guide Drill

- Drill used to mark fixture insertion location
- Used to modify side wall of extraction site with its side blade
- Marker on apex at 2mm


D	Ø2.0 / 2.7
	SNGD2027TL



Ø2.2 Twist Drill

- Recommend drilling 2mm below the floor of sinus
- Use with stopper for safety
- Apex tip 0.6mm


D	Ø2.2
	SNTD2213TL



Hydraulic Membrane Lifter Set

- Lift sinus membrane using water pressure
- Used for all CAS-drill Ø2.8~Ø4.1 drilled holes
- Contraindication:
 - Not recommended for patients with inflammation of maxillary sinus (Sinusitis)
 - Not recommended for patients with severe or complex sinus floor morphologies

D	Ø2.6 / 5.0
	SNMLS





Stopper

- Number on each stopper indicates drilling depth length.
- Color coded per length
- Can be used up to 50 times before replacement

L	2	3	4	5	6	7	8	9	10	11	12
											
Color	Purple	Blue	Yellow	Purple	Blue	Yellow	Purple	Blue	Yellow	Purple	Blue

Bone Carrier

- Used to deliver bone material into the sinus cavity
- Adjust/Fixes head by tightening the back of body part
- Head(SNBCH30 or SNBCH35) can be replaced




SNBCS35

Bone Carrier Head

- Used to deliver bone material into the sinus cavity
- SNBCH30: Use after drilling with CAS-drill $\phi 3.1/\phi 3.3$
- SNBCH35: Use after drilling with CAS-drill $\phi 3.6/\phi 3.8/\phi 4.1$
- Fill bone material upto the bottom of the band, push material through gradually with bone condenser into the sinus cavity completely, and repeat the procedure

D	$\phi 3.1$	$\phi 3.6$
	SNBCH30	SNBCH35



Bone Condenser

- Tool to push bone material into the sinus cavity
- SNBCH30: Uses $\phi 1.1$ / SNBCH35 : Uses $\phi 1.4$

D	$\phi 1.1 / 1.4$
	SNBC1114





CAS-KIT Surgical Instruments


Hydraulic Membrane Lifter Tube

- Attach to hydraulic membrane lifter



Depth Gauge (CAS)


- Check the opening of the sinus floor and measure the depth of remaining bone



LAS-KIT (HLRSNK)



- Lateral Approach - Sinus KIT (LAS-KIT) is optimized for a lateral approach sinus lift
- Dome drill and core drill, safely form lateral window and come in diameters of $\varnothing 5.5$ and $\varnothing 7.0$
- LAS drill uses a stopper system (0.5mm increment), which adjusts the depth, and safely forms a window without perforating the membrane



097

HIOSEN KITS





LAS-KIT Surgical Instruments

Dome Drill

- Forms window and collects autogenous bone
- Optimized cutting force due to the combination of macro cutting edge and micro cutting edge
- Stopper attachment controls depth
- Cutting speed : 1,200~1,500rpm
- ※ Excessive drilling will damage the membrane


L \ D	Ø5.5	Ø7.0	Wide Ø7.0
25	LSDR554TD	LSDR74TD	LSDR74WTD



Core Drill


- Creates a bone lid when forming the lateral window
- Same cutting edge design as the CAS-KIT enhances safety
- Cutting speed : 1,200~1,500rpm
- ※ Excessive drilling will damage the membrane

L \ D	Ø5.5	Ø7.0
25	LSDR554TC	LSDR74TC



Side Wall Drill

- Expands window and trims the rough edges around the window
- Cutting speed : 1,500rpm



Height of side cutting part (mm)	1.0	2.0	3.0	4.0	5.0
CAS-KIT stopper(mm)	8.0	9.0	10	11	12
Side wall drill + CAS-KIT stopper					




※ Used for all CAS-KIT stoppers and can adjust depth



Bone Separator

- Removes bone lid inside the core drill



Stopper

- Number on each stopper indicates length of the drilling depth
- Color coded length
- Can be used up to 50 times before replacement



						
	LSNST0.5	LSNST1.0	LSNST1.5	LSNST2.0	LSNST2.5	LSNST3.0
Color	Yellow	Purple	Blue	Yellow	Purple	Blue

ESSET KIT (HESEK)



Base components

Torque Wrench
TQWCB / TW30B



Depth Gauge
ODG



Available for use

ETII / III TSII / III

SSII / III

USII / III

Ultra-wide

Saw

RA231DC070
RA231DC100
58231DC204130



Ø5.0 Crest Remover

CERM50A
CERM50S



CERM70A


Ø7.0 Crest Remover



2D1808LC01
2D1810LC01
2D1811LC01



Twist Drill




ESSET KIT Surgical Instruments



Crest Remover

- Horizontal modification of the crest bone
- Recommended speed for angled type : 1,200~1,500rpm
- Recommended speed for straight type : 15,000~30,000rpm


L \ D	Ø5.0	Ø7.0
29	CERM50A	CERM70A
45	CERM50S	-



Twist Drill

- Marks the fixture location
- Controls depth with built-in stopper
- Recommended speed : 1,200~1,500rpm


L \ TL \ D	Ø1.8
8.5 33	2D1808LC01
10 34.5	2D1810LC01
11 36	2D1811LC01



Saw

- For ridge modification and splitting
- Cut vertically and incise the whole part from mesial → distal direction
- Recommended speed : 1,200~1,500rpm
- T = Thickness

T \	Ø7.0	Ø10.0	Ø13.0
0.3	RA231DC070	RA231DC100	58231DC204130






ESSET KIT Surgical Instruments

SET Drill


- Gradually expands bone width
- Use SET drill sequentially according to the diameter of fixture
i.e. F4.0 : I → II → III / F4.5 : I → II → III → IV
- Recommended speed : 25~35rpm

L \ Type	I	II	III	IV
D1 / D2	Ø1.6 / 2.8	Ø2.2 / 3.6	Ø2.7 / 4.1	Ø3.1 / 4.5
8.5	SET162808	SET223608	SET274108	SET314508
10	SET162810	SET223610	SET274110	SET314510
11.5	SET162811	SET223611	SET274111	SET314511



Mount Extension

- Used to apply manual torque during the process of inserting / removing SET drill into alveolar bone



Torque Wrench

- Used to apply torque to SET drill



Depth Gauge (ESSET)

- Combination of Open wrench and depth gauge



OsstemGuide KIT (OGDK)



Base components

Kit Steel bowl

AKB*

Ratchet Wrench

CITQW-1185A*




Open Wrench

ASOW*



- Optimized for Tapered implants with internal connection.
- ※ 7mm and 15mm implants require additional drills (not included in the kit)






OsstemGuide KIT Components

OsstemGuide Twist Drill

- All drills feature stoppers.
 - Drills are designed to minimize soft tissue damage during flapless surgery.
 - Drills are 10 mm longer than regular surgical drill
- ※ 7mm and 15mm drills are available through Special Order Only; Not included in the Kit.



L(TL) \ D	Ø2.0	Ø3.0	Ø3.3	Ø3.8	Ø4.3
*7.0 (32.4)	QGTD2007	QGTD3007	QGTD3307	QGTD3807	QGTD4307
8.5 (33.9)	QGTD2008	QGTD3008	QGTD3308	QGTD3808	QGTD4308
10 (35.4)	QGTD2010	QGTD3010	QGTD3310	QGTD3810	QGTD4310
11.5 (36.9)	QGTD2011	QGTD3011	QGTD3311	QGTD3811	QGTD4311
13 (38.4)	QGTD2013	QGTD3013	QGTD3313	QGTD3813	QGTD4313
*15 (40.4)	QGTD2015	QGTD3015	QGTD3315	QGTD3815	QGTD4315



OsstemGuide Mount

- Special mount for OsstemGuide.
- Color-coded to match the surgical guide sleeve

ET System



OsstemGuide Cylinder (Pin)

- For prosthetic fabrication.
- Attach lab analog through surgical guide for casting
- Color-coded to match the surgical guide sleeve

ET System





MS KIT (OMSK)

Available for use


MS

Base components

Depth Gauge (MS)
MSDG



Ratchet Wrench
CITQW-1185A



ESR KIT Easy Screw Removal KIT (OESRK) NEW



Base components

Abutment Removal Tip


OARTML

OARTRL



Slot Driver

OTSD07






ESR KIT Surgical Instruments

Reverse Driver

- Tool used to remove a fractured screw
- Must be used with the Guide that fits the fixture
- Insert until the red band is in the Guide and turn counterclock wise to remove the screw
- *Use the Screw Holder to remove fractured screws
- Use manually / Rotate counterclock wise / Number of use: 10 times


L \ F	Mini	Regular / Wide
Short	-	ORVDRS
Long	ORVDMML	ORVDRL



Screw Removal Drill (SR Drill)

- Used to create a hole in the fractured screw
- Shavings must be suction and irrigate with water throughly
- Short and long specifications suitable for different intermaxillary spaces
- Drill until the red band cannot be seen
- Recommended speed: 1,200~1,500rpm in reverse
- Number of uses: 5
- ※ Must be used with the guide / Do not apply excessive vertical force / Do not immerse in hydrogen peroxide
- F = Fixture


L \ F	Mini(ET/GS)	Regular / Wide(ET/GS/SS/US)
Short	OSRD08S	OSRD10S
Long	OSRD08L	OSRD10L



Guide

- Guide used for center and prevent wobbling of the reverse driver, screw removal drill (SR drill), and re-tap
- Short and long used according to intermaxillary spaces
- Used in combination with the ESR handle
- F = Fixture

F \ Type	ET/TS (Hex)		SS (Octa)		US (Hex)	
	Short	Long	Short	Long	Short	Long
Mini	OGTMS	OGTML	OGUMS	OGUML	OGUMS	OGUML
Regular	OGTRS	OGTRL	OGSRS	OGSRL	OGURS	OGURL
Wide	-	-	OGSRS	OGSRL	OGUWS	OGUWL





Torque Handle

- Handle for manual torquing




Screw Removal Tip (SR Tip)

- Removes fractured screw. Insert into hole drill by Screw Removal Drill.
- Rotate counterclock wise



Short Long



Screw Holder

- If part of fractured screw is protruding, remove with screw holder
- Color coded for easy selection
- F = Fixture



Mini Regular Wide




Re-tap

- Tool that rethreads the fixture's internal connection
- Used with the torque wrench or ratchet wrench, forms thread by hand
- F = Fixture



Mini Regular Wide



ESR KIT Surgical Instruments



ESR Handle

- Stabilizes connecting guide

OARH



Abutment Removal Tip

- Removes fractured or jammed abutments and mounts
- Insert into fractured abutment hole, turn counterclockwise. Rock back and forth to loosen and remove with forceps
- Can be used to remove Mini abutment screw with a stripped hex.
- After attaching to hex turn counterclockwise
- F = Fixture

F Mini Regular
OARTML OARTRL



Mini



Regular

Transfer Abutment Separate Tool

- Remove jammed abutment of non-hex type transfer abutment
- The tip is for mini abutments; the next step for regular abutments
- Remove abutment screw, insert Separate Tool Body into the abutment, tighten clockwise with Driver and remove the abutment. However, if there is difficulty separating the abutment, attach a ratchet wrench for extra torque

Driver Body Set
TASD TASB TAST



Driver



Body

Slot Driver

- Used after grinding a slot with a $\varnothing 0.8$ bur when the hex of healing abutment, cover screw, and abutment screw is stripped

OTSD07



EFR KIT Easy Fixture Removal KIT (OSFRK) NEW



Base components

Fixture Wrench

FRDFE



Torque Wrench

TW400B



Available for use

ETII / III TSII / III

SSII / III

USII / III

Ultra-wide



Remover Screw




Remover Screw



Screw Driver

Remover Screw

FRSM35
FRSR40
FRSW50



Remover Body



FRBM35S
FRBM35L
FRBR40S
FRBR40L
FRBW50S
FRBW50L

Torque Extension

OTE



111

HIOSEN KITS



EFR KIT Surgical Instruments

Remover Screw

- When securely fastened to the fixture, it serves as supporting structure for the Remover Body
- Use the proper type that matches the diameter of the fixture to be removed (ET/SS/US, normal/fracture)
- Use "Fracture" when removing fractured fixture
- Recommend tightening torque : 100Ncm
- F = Fixture
- P = Platform



Type \ F / P	Mini Ø3.5 / -	Regular Ø4.0~4.5 / P4.8	Wide Ø5.0 / P6.0
ET/SS Normal	FRSM35	FRSR40	FRSW50
Fracture	FRSM35F	FRSR40F	FRSW50F
US	FRSM35US	FRSR40US	FRSW50US

Screw Driver

- Driver for Remover Screw
- Recommend tightening torque : 100Ncm
- F = Fixture



F	Mini	Regular	Wide
	FRSDM23	FRSDR25	FRSDW30

Remover Body

- Connects to Remover Screw to apply torque to the fixture.
- Use the proper type that matches the diameter of the fixture to be removed
- F = Fixture



F	Mini	Regular	Wide
Short	FRBM35S	FRBR40S	FRBW50S
Long	FRBM35L	FRBR40L	FRBW50L



Torque Extension

- Extends the length (10mm) of screw driver and remover body



Torque Wrench

- Tightens Screw Driver and remove fixture using Remover Body
- Maximum torque : 400Ncm (increments of 100Ncm)
- After use, store after washing and sterilizing



Fixture Wrench


- Wrench used to separate the fixture from the Remover Body





Osteo KIT (OSTK)


- Concave : Instruments for maxillary sinus floor augmentation to vertically increase the amount of alveolar bone that can be used within maxillary molar area
- Expanding : Instrument for increasing initial stability of an implant by preserving bone and densifying trabecular bone of poor bone quality
- Stopper system to adjust depth



Sinus KIT (ASLK)



- Instruments used for maxillary sinus membrane elevation and grafting procedure
- Sinus Curette for lateral approach sinus lift
- Components
 - Freer elevator : OFE
 - Bone graft carrier : OBGC
 - Membrane separator (circle type) : OMSC
 - Sinus curette-short : OSCS
 - Sinus curette-long : OSCL





115

HIOSEN KITS




Bone Spreader KIT (HBNSPK)

- Instruments for expanding a narrow alveolar ridge
- Offset type for more convenient procedures
- Components
 - OBSO22F, OBSO28F, OBSO35F, OBSO35R



- Use for alveolar bone expansion
- Offset type for easy operation
- Markers corresponds to the implant length.

(Unit : mm)






Code		Tip length	7	8.5	10	11.5
Spec.						
ORSO22F	Thickness		1.15	1.3	1.45	1.6
	Width		2.1	2.2	2.2	2.2
ORSO28F	Thickness		1.15	1.3	1.45	1.6
	Width		2.65	2.8	2.8	2.8
ORSO35F	Thickness		1.3	1.45	1.6	1.8
	Width		3.3	3.5	3.5	3.5
ORSO35R (round type)	Thickness		1.85	2.1	2.3	2.55
	Width		3.3	3.5	3.5	3.5



Ridge Split KIT- Straight (ORSSK)

- Chisel: Instruments for expanding a narrow alveolar ridge
- Blade holder: When it is difficult to incise bone using a bur in poor bone, malletting can be done by attaching #15 blade
- Components
 - Ridge split chisel : ORSS15, ORSS20, ORSS25, ORSS30
 - Blade holder : ORSBH






(Unit : mm)

Code	Spec.	Tip length			
		7	8.5	10	11.5
ORSS15	Thickness	1.1	1.27	1.5	1.5
	Width	4	4	4	4
ORSS20	Thickness	1.45	1.7	2.0	2.0
	Width	4	4	4	4
ORSS25	Thickness	1.8	2.15	2.5	2.5
	Width	4	4	4	4
ORSS30	Thickness	2.15	2.5	3.0	3.0
	Width	4	4	4	4



Ridge Split KIT- Offset (HRGSPOK)

- Chisel: Instruments for expanding a narrow aleolar ridge
- Blade holder: When it is difficult to incise bone using a bur in poor bone, malletting can be done by attaching #15 blade
- Components
 - Ridge split chisel : ORSO15, ORSO20, ORSO25,ORSO30
 - Blade holder : ORSBH



		Tip length			
Code	Spec.	7	8.5	10	11.5
ORSO15	Thickness	1.1	1.27	1.5	1.5
	Width	4	4	4	4
ORSO20	Thickness	1.45	1.7	2.0	2.0
	Width	4	4	4	4
ORSO25	Thickness	1.8	2.15	2.5	2.5
	Width	4	4	4	4
ORSO30	Thickness	2.15	2.5	3.0	3.0
	Width	4	4	4	4

(Unit : mm)



HIOSSSEN



HROSSEN



GBR

122 SmartBuilder

128 AutoBone Collector

129 Allogeneic bone

132 Xenograft

133 Synthetic bone

134 Absorbent membrane

136 Impermeable membrane



SmartBuilder - Titanium


SmartBuilder™

- **SmartBuilder™** is a titanium mesh membrane that provides space maintenance necessary for bone regeneration in alveolar bone defects. Optimized pore sizes provide mechanical strength while still being flexible, prevents shifting or migration of grafting material and allows ideal diffusion of blood and growth factors required for the promotion of healing and bone regeneration.
- **Smart Design** Smooth edges and preformed shaped minimize potential exposure through soft tissue.
- **Smart Coverage** Enhanced space-maintaining properties, providing stable coverage and securing bone material in place.
- **Smart Protection** Rigid shape withstands deformative forces - the membrane does not collapse into the bone graft
- **Smart Application** Optimal preformed shape eliminates the need for additional modifications. No bending or cutting is necessary, and no need for tacks! Conveniently affix the SmartBuilder to a dental implant or to an Anchor Screw
- **Smart Healing** Optimized pore patterns allow for diffusion of fluids and growth factors, while creating a barrier for cell occlusions between soft and hard tissue.

SB1

Papilla Builder


Preformed titanium membrane secures bone grafts to vertically augment alveolar bone up to 1.5mm, restoring papilla height for esthetically sensitive cases.



SB2

Lateral Builder


A semi-dome shaped titanium mesh membrane ideal for simple socket extractions, fenestrations and dehiscence cases that require horizontal and/or vertical bone augmentation of the alveolar ridge.



SB3

Jaw Builder

A clamshell shaped titanium mesh membrane secures bone grafts for horizontal and vertical ridge augmentation of up to 5 - 10mm.





SB1 Papilla Builder

P = Proximal
BL = Buccal Length
BD = Buccal Distance
LD = Lingual Distance

P			P	BL	BD	LD	
7mm			Scallop shape				
			7	3	5	5	SB1P 73A
			7	5	5	5	SB1P 75A
10mm			Scallop shape				
			10	3	5.5	5.5	SB1P 103A
			10	5	5.5	5.5	SB1P 105A
12mm			Scallop shape				
			12	3	6.5	6.5	SB1P 123A
			12	5	6.5	6.5	SB1P 125A




SmartBuilder - Titanium

P = Proximal
BW = Buccal Width
BL = Buccal Length
BD = Buccal Distance

SB2 Lateral Builder

Augmentation



1 Wall



P	BW	BL	BD	
4	8	7	5.5	SM1W487SB
4	10	7	5.5	SM1W4107SB
4	10	9	5.5	SM1W4109SB

2 Wall



Buccal-Proximal



7	9	7	5.5	SM2W797SB
7	9	9	5.5	SM2W799SB
10	12	7	5.5	SM2W10127SB
10	12	9	5.5	SM2W10129SB
12	12	7	5.5	SM2W12127SB
12	12	9	5.5	SM2W12129SB



2 Wall

Buccal-Lingual



5	9	7	5.5	SM2W597SE
5	9	9	5.5	SM2W599SE
6	12	7	5.5	SM2W6127SE
6	12	9	5.5	SM2W6129SE
7	12	7	5.5	SM2W7127SE
7	12	9	5.5	SM2W7129SE

3 Wall



7	9	7	5.5	SM3W797SB
7	9	9	5.5	SM3W799SB
10	12	7	5.5	SM3W10127SB
10	12	9	5.5	SM3W10129SB
12	12	7	5.5	SM3W12127SB
12	12	9	5.5	SM3W12129SB





BW = Buccal Width
BL = Buccal Length
LL = Lingual Length
BD = Buccal Distance
LD = Lingual Distance

SB3 Jaw Builder Available Soon

Augmentation



BW BL LL BD LD

Horizontal





10	7	3.5	5.5	3.7	SB3H107B
10	9	4.5	5.5	3.7	SB3H109B
10	11	6	5.5	3.7	SB3H1011B

Horizontal





20	7	3.5	5.5	3.7	SB3H207B
20	9	4.5	5.5	3.7	SB3H209B
20	11	6	5.5	3.7	SB3H2011B

Vertical



10	7	7	5.5	5.5	SB3V107B
10	9	9	5.5	5.5	SB3V109B
10	11	11	5.5	5.5	SB3V1011B

Vertical




20	7	7	5.5	5.5	SB3V207B
20	9	9	5.5	5.5	SB3V209B
20	11	11	5.5	5.5	SB3V2011B



SmartBuilder Components

Healing Cap


- For non-submerged closure
- Compatible with SB1, SB2 & SB3
- NOTE: 5.0 is not compatible with SB1 or SB3
- Connect using a 1.2 hex hand driver
- Recommended tightening torque: 5~8Ncm



D \ H	3.0	4.0
Ø4.0	SMHA443R	SMHA444R
Ø5.0	SMHA553R	SMHA554R

Cover Cap


- For submerged closure
- Compatible with SB1, SB2 & SB3
- Connect using a 0.9 hex hand driver
- Recommended tightening torque: 5~8Ncm



D \ H	1.5
Ø4.0	SMCC415

SB Anchor for ETIII

- For use with ETIII fixture (shoulder contact)
- Compatible with SB1, SB2 & SB3
- Connect using a 1.2 hex hand driver
- Recommended tightening torque: 12~15Ncm




D \ H	0.5	1.0	1.5	2.0	2.5	3.0
Ø3.5	SMHI305TSM	SMHI310TSM	SMHI315TSM	SMHI320TSM	SMHI325TSM	SMHI330TSM
Ø4.0	SMHI405TSR	SMHI410TSR	SMHI415TSR	SMHI420TSR	SMHI425TSR	SMHI430TSR



Anchor Screw Available Soon

- For use without placing an implant
- Compatible with SB1, SB2 & SB3
- Insertion depth: 3~5mm
- Connect using a 1.2 hex hand driver

L	8.5	10	11.5	13
	SBS2008E	SBS2010E	SBS2011E	SBS2013E



Cover Cap Driver Ejector Set

- Use to screw and unscrew the cover cap
 - Friction type driver
 - After unscrewing the Cover Cap, use the Cap ejector to release the Cover Cap from the driver
 - Packing unit : cap driver + cap ejector
- ※ For the posterior region, it is recommended to use the Short version

Short	Long
SMCDESS	SMCDES



Cover cap driver Cap ejector

Defect Gauge

- Laser markings at every 1mm
 - Thick markings at 4-5mm, 9-10mm, 14-15mm
- ※ Use the Defect Gauge to measure the bone defect's buccal width, length and distance in order to select the proper SmartBuilder shape and size.

SMDG






AutoBone Collector

AutoBone Collector

- $\varnothing 5.0$ and $\varnothing 6.0$ specifications include drill + stopper set
- Recommended speed : 300 ~ 600rpm
- Can be used up to 50 times before replacement
- ※ Before drilling, slide the Stopper up the drill until it engages with the first notch.
Commence collecting autogenous bone.
 - Stop drilling and pull back when the Stopper reaches the hilt of the drill.

L \ D	$\varnothing 5.0$	$\varnothing 6.0$
Short	ABC504S	ABC604S
Long	ABC504L	ABC604L



Stopper

- Prevents over drilling and collects autogenous bone

L \ D	$\varnothing 5.0$	$\varnothing 6.0$
Short	ABC2ST504S	ABC2ST604S
Long	ABC2ST504L	ABC2ST604L




Bone Ejector

- Used for dislodging collected bone inside the Stopper
- Can be used with both 5.0 and 6.0 Stoppers

D	$\varnothing 5.0 / 6.0$
	ABBE52L



Stopper Attachment Guide -2 Stage Connection System



※ Refer to videos for more detailed instructions

Allogeneic Bone - Powder/Chip




SureOss (FDBA)

- FDBA (Freeze Dried Bone)
- 100% allograft cortical bone
- Superior biocompatibility
- Particle size - powder 200~850 μ m
- chip 850~1000 μ m
- Manufacturer: Hansbiomed Co., Ltd.
Korea

CC Powder

0.25	SOP25UH
0.5	SOP50UH
1.0	SOP100UH



CC Chip

0.25	SOC25UH
0.5	SOC50UH
1.0	SOC100UH



SureOss-D (DFDBA)

- DFDBA (Demineralized Freeze Dried Bone)
- 100% demineralized allograft bone
- Superior biocompatibility
- Particle size - powder 200~850 μ m
- chip 850~1000 μ m
- Manufacturer: Hansbiomed Co., Ltd.
Korea


CC Powder

0.25	DSOP25UH
0.5	DSOP50UH
1.0	DSOP100UH



CC Chip

0.25	DSOC25UH
0.5	DSOC50UH
1.0	DSOC100UH





Allogeneic Bone - Powder/Chip

OsteOss (FDBA)

- FDBA
(Freeze Dried Bone)
- 100% allograft bone
(cortical + cancellous bone)
- Particle size - powder 200~850 μ m
- chip 850~1000 μ m
- Manufacturer: Hansbiomed Co., Ltd.
Korea

CC Powder

0.25	CCP 25UH
0.5	CCP 50UH
1.0	CCP 100UH



CC Chip

0.25	CCC 25UH
0.5	CCC 50UH
1.0	CCC 100UH




Allogeneic Bone - DBM



SureFuse II

- DBM (allograft)
- Bone conduction capabilities
- Can be manipulated to form any shape & size
- Superior biocompatibility
- Manufacturer: Hansbiomed Co., Ltd.
Korea

CC	Gel
0.3	SG03UH
0.5	SG03UH
1.0	SG1UH
3.0	SG3UH
5.0	SG5UH



CC	Putty
0.3	SP03UH
0.5	SP03UH
1.0	SP1UH
3.0	SP3UH
5.0	SP5UH



ExFuse II

- DBM (allograft)
- Demineralized allograft + cancellous bone
- Bone conduction capabilities
- Can be manipulated to form any shape & size
- Superior biocompatibility
- Manufacturer: Hansbiomed Co., Ltd.
Korea

CC	Gel
0.3	EG03UH
0.5	EG05UH
1.0	EG1UH
3.0	EG3UH
5.0	EG5UH



CC	Putty
0.3	EP03UH
0.5	EP05UH
1.0	EP1UH
3.0	EP3UH
5.0	EP5UH





Xenograft

A-Oss Available Soon


- DBB
(Deproteinized Bovine Bone)
- Osteoconduction capability
- Excellent volume maintenance
- Superior biocompatibility
- Manufacturer:
Osstem Implant Co., Ltd. Korea

g \ Particle Size 0.25~1.0mm

0.1 (0.2cc)	BBS01
0.25 (0.5cc)	BBS02
0.5 (1.0cc)	BBS05
1.0 (2.0cc)	BBS10
2.0 (4.0cc)	BBS20

g \ Particle Size 1.0~2.0mm

0.1 (0.3cc)	BBL01
0.25 (0.75cc)	BBL02
0.5 (1.5cc)	BBL05
1.0 (3.0cc)	BBL10
2.0 (6.0cc)	BBL20



Synthetic Bone



Q-Oss Available Soon


- Synthetic bone composed of HA 75% + β -TCP 25%
- Porous structure with interconnected micropores
- Excellent osteoconductivity
- Long-term volume maintenance
- Manufacturer: Osstem Implant Co., Ltd. Korea

g 0.5~1.0mm

0.1 (0.2cc)	BQS01
0.25 (0.4cc)	BQS02
0.5 (0.8cc)	BQS05
1.0 (1.5cc)	BQS10
2.0 (3.0cc)	BQS20

g 1.0~2.0mm

0.1 (0.2cc)	BQL01
0.25 (0.5cc)	BQL02
0.5 (1.0cc)	BQL05
1.0 (2.0cc)	BQL10
2.0 (4.0cc)	BQL20



Q-Oss+ Available Soon


- Synthetic bone composed of HA 20% + β -TCP 80%
- Porous structure with interconnected micropores
- Excellent osteoconductivity
- Gradual resorbence
- Manufacturer: Osstem Implant Co., Ltd. Korea

g 0.5~1.0mm

0.1 (0.2cc)	BQ+S01
0.25 (0.4cc)	BQ+S02
0.5 (0.8cc)	BQ+S05
1.0 (1.5cc)	BQ+S10
2.0 (3.0cc)	BQ+S20

g 1.0~2.0mm

0.1 (0.2cc)	BQ+L01
0.25 (0.5cc)	BQ+L02
0.5 (1.0cc)	BQ+L05
1.0 (2.0cc)	BQ+L10
2.0 (4.0cc)	BQ+L20



Bongros

- Synthetic bone
- Carbon-apatite ingredient
- Bone conduction capability
- Manufacturer: CGBio. Korea

g

0.25	HAGM 0.25g
0.5	HAGM 0.5g
1.0	HAGL 1.0g





Absorbent Membrane - Collagen

Cytoplast RTM Collagen

- Bovine collagen
- RTM = Resorbable Tissue Matrix
- Longer, predictable resorption time
(26 ~ 38 weeks)
- Packing unit : 2ea
- Manufacturer: Osteogenics. USA
- T = Thickness

mm \ T 0.3mm

15 × 20

RTM1520

20 × 30

RTM2030

30 × 40

RTM3040



Absorbent Membrane - Allogeneic Dermal




SureDerm

- Human dermis
(Allogeneic acellular dermal)
- Reconstructing effect on
membrane and gingiva
- Freeze drying and rehydration required
- Excellent biocompatibility
- Manufacturer:
Hansbiomed Co., Ltd. Korea
- T = Thickness


mm \ T 0.25~0.59mm

10 × 20	G102002UH
10 × 40	G102004UH
20 × 20	G102104UH
20 × 40	G102008UH




mm \ T 0.60~0.99mm

10 × 20	G202002UH
10 × 40	G202004UH
20 × 20	G202104UH
20 × 40	G202008UH



mm \ T 1.00~1.39mm

10 × 20	I302002UH
10 × 40	I302004UH
20 × 20	I302104UH
20 × 40	I302008UH



GBR

135



Impermeable Membrane

Cytoplast TXT-200


- Non resorbable membrane (non resorbable membrane)
- 100% medical grade PTFE
- Minimizes the risk of bacterial infection
 - Membrane pore size 0.3mm
 - Bacteria size 2~5mm
- Superior biocompatibility
- Manufacturer: Osteogenics. USA

mm \ TXT-200 Singles

12 × 20 TXT1220-1

mm \ TXT-200

25 × 30 TXT2530-1



Cytoplast TI-250

- Non resorbable membrane + titanium reinforced
- Increased space retaining capability with titanium reinforcement
- Minimizes the risk of bacterial infection (100% medical grade PTEE material)
 - Membrane pore size 0.3mm
 - Bacteria size 2~5mm
- Superior biocompatibility
- Manufacturer: Osteogenics. USA

mm \ Ti-250 Anterior Narrow

12 × 20 TI250AN-1

mm \ Ti-250 Anterior Singles


14 × 24 TI250AS-1

mm \ Ti-250 Posterior Singles

20 × 25 TI250PS-1

mm \ TXT-200 Posterior Large

25 × 30 TI250PL-1





HROSSEN

User Manual 2013.02 ver.4.0 “Disposable, re-use prohibited, medical appliance”



Hiossen/Osstem Implant product information

Osstem Implant products are manufactured using medical grade titanium.

Osstem Implant abutments, denture materials and surgical tools are only compatible with Osstem fixtures. For more details about each product, please refer to the user manuals, catalogs, or please visit our corporate website (www.osstem.com). Please check all product labels for product codes, specifications, manufactured dates and expiration dates.

Sterility

Fixtures, cover screws and healing abutments are cleansed and gamma-sterilized. These products are disposable sterile medical appliances and must be used in a sterile field. If the package is damaged or has expired, it must not be used. If the product package has been opened but not used, there is a risk of contamination and it is not recommended that the product be resterilized.

Storage conditions

Store all products in a dry place at room temperature (30°C). Avoid direct sunlight.

General precautions

Dental implant surgery requires formal training and education.

Cautions before dental surgery

A thorough review of the patient health history, oral and radiographic examinations must be completed to determine bone quality and proper treatment planning.

Cautions during dental implant surgery

Osstem Implant System can be placed in a single or two stage dental implant procedure. In order to minimize damage to the patient's tissue, special attention to temperature, surgical lesions and eliminating all sources of contamination and infections are required. Any deviation from the standard surgical protocol increases the risk of failure. When inserting the dental implant, supply sufficient irrigation (with either water or saline) and avoid excessive torque force (greater than 55Ncm). Placing dental implants using a force greater than 300Ncm will result in fracturing the dental implant. Direct pressure to the fixture should be avoided right after surgery. Immediate or delayed loading of the fixture must be determined after a thorough examination of the patient's bone condition and dental implant initial stability.

Dental implants with a diameter less than 4.0mm are not recommended for the posterior region. Wide body dental implants (diameter greater than 5.0mm) are recommended for the posterior region but should not be used with angled abutments. If considering a wide body implant, proper radiographic evaluation must be done to determine the bone quality and potential anatomical restrictions. Wide diameter implants should be performed as a two stage procedure. Sufficient healing period must pass before splinting with another implant or when loading. Immediate loading is not recommended. Short body implants (diameter greater than 5.0mm and shorter than 7mm) are only to be used in the posterior region.

The clinician must thoroughly evaluate the patient's condition and recognize the following: 1) bone loss due to peri-implantitis, 2) changes to the dental implant condition, 3) proper osseointegration determined radiographically. If there is movement or if there is bone loss more than 50%, immediately remove the implant.

Dental implants with a HA coating are prone to cracking or fracturing under high torque, therefore placing in hard bone should be avoided and use less than 35Ncm of force. CA and SOSI coated dental implants are suspended in a solution to prevent the chemically treated surface from reacting to air. After removing the CA or SOSI dental implant, place the implant within 15 minutes to minimize surface degradation.

Warning

Improper patient selection and treatment planning may result in dental implant failure or loss of bone. Osstem Implants must not be used for any other purpose other than what has been prescribed and must not be alter in any shape or form. Implant movement, bone loss or chronic infections can result in implant failure.

Indications

Osstem Implant Systems are designed to replace a patient's tooth or teeth.

They can be placed in both the maxillary and submaxillary alveolar bone and after full osseointegration can be restored.

Osstem Implant Systems offer both temporary and final prosthetic components which can be cement or screw retained, overdenture or fixed bridge.

Side effects

There are possible side effects after implant surgery (loss of implant stability, damage to dentures, etc...). These issues can be the result of lack of bone volume or poor bone quality, an infection, patient's poor hygiene, non compliance with post op procedures, improper implant placement, improper restoration.

Contraindications

Patients with the following conditions are not eligible for dental implants:

- Blood clotting or wound healing issues
- Diabetes
- Smoking or drinking excessively
- Compromised immune system due to disease, chemo or radiation therapy
- Chronic oral infection or inflammation due to improper oral hygiene or teeth grinding
- Incurable malocclusion/arthritis and insufficient arch space

Manufacturer : Osstem Implant Co., Ltd.
203, Geojje-daero, Yeonje-gu, Busan, Korea
TEL 82-51-850-2500 FAX 82-51-861-4693



DEUTSCHE OSTEM GmbH.
Mergenthalerallee 25
65760 Eschborn, Germany
+49-(0)6196-777-550

Storage condition

Dry place at room temperature

Rx only

For USA only : Federal law restricts this device to sale by or on the order of a dentist



0434



Sterilized using irradiation



Use by



Manufacture



Do not reuse



Date of manufacture



Keep away from sunlight



Catalogue number



Non-Sterile



Keep dry



Batch code




Do not resterilize



Caution, Consult accompanying documents



How to Maintain Surgical KITS



- ① During the surgery, be sure to keep the used tools in saline or distilled water.
- ② After the surgery, wash all tools used during surgery in alcohol.
Caution : Do NOT use hydrogen peroxide.
Exposure to hydrogen peroxide may cause discoloration of the laser marking and/or TiN coating.
- ③ Wash the tool with distilled water or under running water until all blood stains and/or foreign objects are removed.
- ④ Remove moisture completely with dry cloth or a warm fan.
- ⑤ Place the dried tools inside the Kit case.
(Refer to the color-coding for easy placement.)
- ⑥ After drying the Kit in the Autoclave for 15 minutes at 132°C, store the Kit at room temperature.

Precautions

Separate, wash and store all tools used immediately after the surgery.

It is advised to disinfect the Hiossen Surgical KIT again prior to the surgery (at 132°C for 15 minutes).

Although the Hiossen Surgical KIT is covered under the product warranty for one year after opening the Kit, all drills and drivers may be used up to 50 times only.

HROSSEN



www.befrest.com



HROSSEN

