



# EQ-Condenser

Versatility; Enhance performances



## I Superior flexibility with NiTi.

Flexible NiTi tip will let you compact apical third part more thoroughly



## I Flat End & Smooth Edge

Prevent from being stuck within dentin and allow you to condense exactly



Stainless Steel tip

## I Laser marking (5mm on the tip)

This laser marking will let you save time for measuring

## I Systematic color;

Matching with EQ-V packtip color

## I Size

40 NiTi/80 SS  
50 NiTi/100 SS  
60 NiTi/120 SS



## I Resistant to chemical & Autoclave

Stress free from chemical & heat damage