



EQ-S

Item List

Item Code	Description	Remark
026-210	EQ-S Gray Kit	Tips and Protective Sleeves included
026-220	EQ-S Pink Kit	Tips and Protective Sleeves included
126-230	EQ-S Tip 15/02	30ea/1set (Color: White)
126-240	EQ-S Tip 25/02	30ea/1set (Color: Red)
126-250	EQ-S Tip 35/02	30ea/1set (Color: Green)
126-070	Protective sleeve	100ea/1set

Specifications

		Remark
Dimensions	W24 x H25 x L183	Including tip
Weight	43g	
Battery	1.5 V [Alkaline(AA)]	Not included in the Kit
Operation Time	9 hours (continuous operation)	



Meta Biomed Co.,Ltd - Head Office

270 Osongsaengmyeong 1-ro,
Osong-eup, Heungdeok-gu,
Cheongju-si, Chungbuk, Korea

Phone : +82 43 218 1981
dental@meta-biomed.com
www.meta-biomed.com

Meta Biomed Inc.

3015 Advance Lane Colmar,
PA 18915, USA

Phone : +1 267 282 5893
america@metabiomed-inc.com
www.metabiomedamericas.com

Meta Biomed Europe GmbH

Wiesenstr. 35,
45473 Mülheim an der Ruhr,
Germany

Phone: +49 208 309 91910
europe@meta-europe.com
www.meta-europe.com

DO-056-1



Let it Flow!



EQ-S

The most effective, efficient and safe

www.facebook.com/metabiomed.Endodontics



EQ-S

The most effective, efficient and safe

- ✓ Effective
- ✓ Efficient
- ✓ Safe



Effective 3D Irrigation


Multi-direction strongly and quickly sends irrigant throughout the canal

Active Flow™ Technology

Rotational & Vertical Movement


Indentation along the tip, combined with counter clockwise motion, ensures 3D Active Flow™. The unique shape and design of the EQ-S tip create stronger force, which enhances flow in small canals.

Active Flow™ Polymer Tip



- Less pressure to the apical foramen
- Easier access to curved canals
- Perfect interaction with irrigant
- Color coded

Powerful Cleaning



After the usage of EQ-S with MD-Cleanser

Testimonial



My clinical feeling of EQ-S, as a multidirectional sonic device and the novel design of the activating tips, when working inside of the root canal is quite good. It is also useful in my clinical practice for removing medicaments, spreading bioceramic sealers inside of the root canal previous to the placement of the master cone and cements previous to the placement of fiber posts

Dr Jenner Argueta, DDS, M.Sc.



EQ-S is an innovative device, versatile, of easy clinical use, offering something that so far had not been seen in a sonic device. My clinical experience is the multidirectional movement of the tip causing more fluid dynamics, creating continuous waves and having better cleaning, disinfection and removal of the smear layer. High power can be used for irrigation or removing Calcium hydroxide and low can be used for spreading bioceramic sealer before placing the master cone.

Dr Johnny Onori, DDS, M.Sc.