

Ezfil

Nano Hybrid Composite



Feel the way It fills

Ezfil can make your daily restorations easier.





It maintains strength and aesthetics

Experience our new improved composite Ezfil

Nanoscale fillers of various sizes enhance both strength and naturalness of restoration. Our own optimized filler combination makes your treatment easier than ever. It works for both anterior and posterior teeth (Class I, II, III, IV, V)

Non-Sticky, Non-Slumping





<Competitor>



<Ezfil>

Ezfil has 74.36% of optimal filler content which makes it non-sticky and non-slumping and very easy to manipulate before placement.




<SEM images of Fillers: B.G & Aerosil>

Low Shrinkage

Low shrinkage is an important property to avoid Post-Op sensitivity.

Ezfil has 2% of shrinkage rate for superior marginal integrity



Various Shades (Matches with VITA® Shades)

A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, OP, TL

Components

- Ezfil Kit
- 4g Syringe each of : A1, A2, A3, A3.5, B1, B2, OP
- Meta P&Bond 1bottle × 5g (with Microbrush 10EA)
- Meta Etchant 1syringe × 3g (with Disposable Tips 20EA)
- Ezfil Refill : 4g Syringe 1EA

Clinical Cases




Before



After (Class I)

Courtesy of Dr. Jenner O. Argueta (Guatemala)



Before



After (Class II)

Courtesy of Dr. Johnny Onori (Spain)

Meta Biomed Co., Ltd — Head Office
 270 Osongsangmyeong1-ro, Osong-eup
 Heungdeok-gu, Cheongju-si,
 Chungbuk, Korea
 Phone: +82 43 216 0433
dental@meta-biomed.com
www.meta-biomed.com

Meta Biomed Inc.
 3015 Advance Lane
 Colmar, PA 18915,
 USA
 Phone: +1 267 282 5893
america@metabiomed-inc.com
www.metabiomedamericas.com

Meta Biomed Europe GmbH
 Wiesenstr. 35,
 45473 Mülheim an der Ruhr,
 Germany
 Phone: + 49 208 309 91910
europe@meta-europe.com
www.meta-europe.com

DO-061-0

