



Faster Dissociation

Faster medicinal effect & easily removed

Antibacterial with high pH: 12.15 (1day) / 12.35 (14day)

Increasing success probability of root canal treatment due to the high antibacterial effect

Premixed Syringe type

Quick and Easy treatment with minimal waste & clean up

PPG Free

Reducing Post-op sensitivity, ensuring more stable treatment

Improved flowability: 19-20mm (ISO 6876:2012)

Convenient and accurate placement to the apex enhanced

Advantages

Applications

- Weeping Canal
- Formation of a hard tissue barrier (Apexification)
- Temporary root canal filling
- Internal & external root resorption

Composition

- Calcium Hydroxide, Barium Sulfate, Purified Water

Shelf life

- Two years from the production date

Components

- 2g syringe 1EA, Disposable Tips 20EA
- 2g syringe 2EA, Disposable Tips 20EA



**Plus
Metapex**
Calcium Hydroxide Temporary Filling Material



Premixed Calcium Hydroxide Paste with Iodoform

Metapex Plus

Antibacterial with high pH: 12.00 (1day) / 12.4 (14day)

Increasing success probability of root canal treatment due to the high antibacterial effect

Improved flowability: 19-21mm (ISO 6876:2012)

Convenient and accurate placement to the apex enhanced

Premixed syringe type

Quick and Easy treatment with minimal waste & clean up

High radiopacity: 5.4 Almm

Accurate diagnosis available

Advantages

Applications

- Weeping Canal
- Formation of a hard tissue barrier (Apexification)
- Root canal filling for deciduous teeth
- Internal & external root resorption

Composition

- Calcium Hydroxide, Iodoform, Silicon oil

Shelf life

- Two years from the production date

Components

- 2.2g syringe 1EA, Disposable Tips 20EA
- 2.2g syringe 2EA, Disposable Tips 20EA

Meta Biomed Co., Ltd - Head Office

270 Osongsaengmyeong1-ro, Osong-eup,
Heungdeok-gu, Cheongju-si, Chungbuk, Korea

Phone: +82-43 216 0433
dental@mata-biomed.com
www.meta-biomed.com

Meta Biomed Inc.

3015 Advance Lane Colmar, PA 18915, USA

Phone: +1 267 282 5893
america@metabiomed-inc.com
www.metabiomedamericas.com

Meta Biomed Europe GmbH

Wiesenstr. 35, 45473 Mülheim an der Ruhr,
Germany

Phone: +49 208 309 91910
europe@mata-europe.com
www.meta-europe.com