



## Trantemp

Esthetic Temporization Material  
Bis-Acrylate-Composite



### Indications

- Crowns and Bridges
- Veneers, Inlays and Onlays
- Also indicated for long-term temporization

### Benefits

- Stronger temporary with long-term durability
- Quick intraoral set
- Highest fracture resistance
- Four natural looking shades


### Package

- Trantemp Automix  
Cartridge (65g X1ea)  
(Shade : A1, A3, A3, B2)  
Mixing tips (15ea)

## T STAT

Paste Retraction System

15% Aluminium Chloride Topical Gigival Retraction Paste Hemostasis and Retraction



### Indications

- Gingival retraction
- Stop all minor intraoral bleeding

### Benefits

- Excellent handling Non slumping
  - Quick effective retraction: 2 minutes
  - Safe and effective
- :15% aluminium chloride paste is highly effective at controlling bleeding without staining, Discoloring surrounding tissue

### Application

Before a CAD/CAM Scan



T STAT is applied  
around prepared tooth



After 2 minutes,  
T STAT easy rinses away

### Crown and Bridge



Before



Place



After 2 minutes,



A full arch retracted view

- Package** rinse away
- T STAT (0.9g X2syringes.)
  - Applicable tips

## T Hole-in

Light Curing Semi-gel Abutment Hole Filling

Implant Abutment Hole Filling



### Indications

- Implant-screw hole filling
- Temporary filling for inlay, onlay


### Benefits

- Excellent sealing ability
- Semi-gel type temporary restoration

### Package

- T Hole In (3g X3syringes.) Blue, Yellow
- Applicable tips

### Application



Implant abutment hole filling



# Endodontic materials



## T SEAL

Root Canal Sealer



T SEAL is a epoxy resin based root canal sealer which is a paste-paste type of dual syringe.

### Indications

- Permanent root canal sealing
- Thermal obturation

### Package

- T SEAL (12g X1 syringes.)
- Spatula
- Mixing pad

### Benefits

- Excellent biocompatibility
- Easy to mix paste-paste
- Good sealing ability
- Insoluble in tissue fluids

## T PEX

Calcium Hydroxide with Iodoform  
Temporary Root Canal Filling Material



### Indications

- Weeping canal
- Formation of a hard tissue barrier
- Root canal filling for deciduous teeth
- Internal & external root resorption

### Benefits

- Excellent antibacterial effect and radiopacity
- Premixed paste
- Syringe typ for easy delivery
- Disposable tips

### Package

- T PEX (0.9g X3 syringes.)
- Applicable tips

## T Paste

Calcium Hydroxide with Barium sulfate  
Water Soluble Temporary Root Canal Filling Material



### Indications

- Weeping canal
- Formation of a hard tissue barrier
- Temporary root canal filling
- Internal & external root resorption

### Benefits

- Easy cleaning and removing with good water solubility
- Excellent antibacterial effect and radiopacity
- Premixed paste
- Disposable tips

### Package

- T Paste (0.9g X3 syringes.)
- Applicable tips

27, Sanseong-ro 92beon-gil, Sangdang-gu, Cheongju-si, Chungcheongbuk-do, Korea  
Tel. 82-43.298.9775 Fax. 82-43.298.9777  
Email. info@nexobio.com, nexobio@gmail.com Homepage. [Http://www.nexobio.com](http://www.nexobio.com)