




BeauTis 15

Stronger whitening effect, less toothache symptoms

Professional whitening material with a high density of 15% Hydrogen Peroxide

OSSTEM
IMPLANT

Excellent whitening effect achieved through 'special activator' and 'high intensity curing light'




High intensity
Blue LED
curing light



Maximizing "free-radical"
creation which removes colored substances

Stronger effect, less toothache symptom

- BeauTis 15 is a safe whitening material with high density certified by the KFDA
- BeauTis, which contains anti-hyperesthesia, reduces toothache symptoms during the treatment



Whitening System that can be applied to all cases

15% Hydrogen Peroxide
Professional whitening material
BeauTis15

Whitening effect maximized by
special activator



3% Hydrogen Peroxide whitening material
for patients with sensitive hyperesthesia
BeauTis3

3% of Hydrogen Peroxide density creates the same
whitening effect as 15% Hydrogen Peroxide



Teeth whitening KIT
BTS Procedure KIT

Convenient treatment for teeth whitening



Curing Light
BeauTis Light

Safe power blue LED that maximizes
the whitening effect






Tooth Whitening System that is Easy to Use in Daily Life

BEAUTIS WHITE

- Maintain whitening effect after expert whitening treatment
- Colored teeth whitening after orthodontic treatment
- Safe whitening for sensitive patients


OSSTEM[®]
IMPLANT

Excellent Self-Whitener




Excellent Whitening Effect

A hypoallergenic tooth whitener for removing colored matter



Excellent Application

Silicone head that can meticulously paste the whitener



Easy Amount Control


Easy to control amount of whitener by pressing syringe method

How to Use

Set of Whitener + Light Transmitter (5minutes 3times a day, 1week use)



Paste



Light transmitting (5min)



Whitened



Beautis White Light Gel



Beautis White Light



Small
Large

Whitener Only (30minutes 2time a day, 2week use)



Paste and maintain (30minutes)




Whitened



Beautis White Gel

LED Light Scanner to Enhance Whitening Effect

- Light energy evenly dispersed from the anterior to the oral cavity
- Quickly maximize whitening effect



Wearable Light without Time and Space Constraints

- Necklace-type light transmitter that can be used anywhere at any time
- Simple and light fit

Dr. Cho's instrument Kit

Order Code : DCHO KIT

Optimal Implant Surgery Instrument Kit with Long-term Clinical Know-how



OSSTEM[®]
IMPLANT


1. Periosteal Elevator

- Comparing to Molt 9, with narrow width, fine and precise surgery possible
- The instrument lifting the mucosal periosteum after cutting the gingival tissue




2. Minesota Retractor

- Improving gripping area for long time use, which reduces hand fatigue (patent)
- The device that can secure the visibility by pulling the mouth, cheeks, etc.



3. Periosteal Elevator Selden

- The device that lifts and fixes the mucosal periosteum after cutting the gingival tissue




4. Needle Holder, Crile-Wood, TC

- Preventing the situation of being caught intermittently by bevel processing in the middle of the needle holder
- Tungsten carbide treatment in the Beak area enables retaining fine sutures, thereby increasing the convenience of the procedure.




5. Ochsenbein Chisel

- Useful when collecting a small amount of autogenous bone
- The blade is toward the front




6. Tissue Forcep, ADSON

- The instrument catching the gingival without hurting the soft tissue due to no protrusion on the inner surface.



7. Titanium Suction Tip

- High durability with titanium tip (inner diameter 3mm)



8. Dr.Cho's Instrument Pouch

- Using for the storage and sterilization of instruments
- Size 470X400mm

9. Extension Hose

- Hose for extension, connecting to unit chair suction
- Possible for autoclave sterilization (silicone material)
- Cleaning in transparent, Easy in dealing with

10. Extension Hose Adapter

- Adaptor to connect to unit chair suction



Osstem Basic instrument Kit

Product Code : OB KIT

The Kit with 25 Kinds of Instruments Composed of Commonly used Surgical Instruments



Periosteal Elevator



Selden



MOLT9



Prichard



P24G

Mirror



Mirror

Scissor



Tissue scissor



Scissor, LaGrange

Curette



Gracey



Surgical Curette, CM10



Surgical Curette, CM11

Hemostat



Mosquito

Needle Holder



Mayo-Hegar

Towel Clamp



Backhaus

Tissue Forcep



ADSON (ToothO)

Tweezer




Wide Tweezer

Scalpel Handle



Scalpel Handle (straight)



Scalpel Handle (Flat)

Caliper



Castroviejo

Retractor



Minesota Retractor

Chisel



Periodontal Chisel

Rongeurs



Bone Rongeurs, Friedman

Mallet



Mallet

Bone Well



Bone Well

Suction Tip



Titanium Suction Tip

Instrument Pouch



NO.	Code	Description
1	EP9-OB	Periosteal Elevator Molt9
2	EPPR3-OB	Periosteal Elevator Prichard
3	EP23-OB	Periosteal Elevator Selden
4	HTM130C-OB	Hemostat 12cm (Curved)
5	PT42-OB	Tissue Plier, ADSON(SAW TYPE)12cm
6	SCTC115-OB	Iris Scissor T/C
7	SCLC115-OB	La grange Scissors 11cm
8	RNGF140-OB	Bone Rongeur 14cm
9	ML25-OB	Mallet ML-01
10	CGR11-12-OB	Gracey Curette
11	RTCRM-OB	Minesota Retractor
12	NH160TC-OB	Needle Holder T/C 16cm
13	SHS-OB	Scalpel Handle (Straight)
14	SHF-OB	Scalpel Handle (Flat)
15	CHCO2-OB	Periodontal Chisel Ochsenbein & Fedi Curved
16	SN3TI-OB	Titanium Suction Tip
17	PCW150-OB	Wide Tweezer
18	DMCS4-OB	Mirror Set
19	LPC90-OB	Castro Gauge (Curved)
20	EP24G-OB	Periosteal Elevator P24G
21	URCM10-OB	Surgical Curette Curved CM10
22	URCM11-OB	Surgical Curette Curved CM11
23	CPTC135-OB	Towel Clamp
24	BWSUS1-OB	Bone Well (Small)
25	WPA-OB	Instrument Pouch