


ABSOLUTELY HARMONIZED

EsCom[®] 250

Light-Cured Restorative
Nano-Hybrid Composite Resin




Chameleon Effect

“Good Shade match with surrounding teeth.”



✓ Nano-Hybrid Fillers & Zirconium silicate filler


EsCom250 is formulated with advanced **nano-hybrid fillers** and addition of **zirconium silicate**(ZrSiO₄) enhances compressive strength and flexural strength, which provide best support for indirect restorations.



✓ Outstanding Wear Resistance

As nano fillers are dispersed uniformly, it has very low surface roughness, which does not interrupt with the opposite tooth.


Additionally, with smaller particles, it has exceptional wear resistance.




✓ Superior Shade Matching

One shade can invisibly blend with several tooth shades and shows superior shade match.

Perfect shade match with just 2 shades (EsCom250 A2 / TW)





*Source by Parkseongwon Dental Clinic



This innovative composition also contributes to :



- **Good Handling** – Smooth consistency for easy application and control
- **Sculptable** – Excellent adaptability for precise anatomical shaping
- **Easy to Polish** – Achieves a natural, long-lasting gloss with minimal effort
- **Good Shade Match (Chameleon Effect)** – Blends seamlessly with surrounding tooth structure for a natural look
- **Extended working time** – Low light sensitivity, allowing for a working time of up to 90 seconds

✓ Excellent Mechanical Properties



EsCom250 offers exceptional compressive and flexural strength, making it ideal for withstanding the daily stresses of the oral environment.


Its nano-hybrid filler technology ensures superior wear resistance and long-term performance, even in high-load bearing areas.



EsCom250 has very low polymerization shrinkage and has no worries of marginal leakage and secondary carries. It also minimized light sensitivity and gives you more time to contour.

Clinical case

Anterior restoration *Source by Jung. J. K. Dental Clinic



Anterior restoration(Diastema) *Source by Parkseongwon Dental Clinic




Posterior restoration *Source by Kyunghee Univ.




Posterior restoration *Source by Dr. Tatyana Kiyamova



Posterior restoration *Source by Dr. Tatyana Kiyamova



☑ Shade Guide



☑ Package



Available	Shade	Contents
EsCom250	TW, A1, A2, A3, A3.5, A4, B1, B2, C2, A02, A03, D2, G1	4g x 1 syringes
EsCom250 System Kit	A1, A2, A3, A3.5, B2	4g x 5 syringes, EsBond, FineEtch
EsCom250 Syringe Kit (4s)	A1, A2 x 2, A3	4g x 4 syringes, K-Bond Universal, FineEtch
EsCom250 Syringe Kit	A1, A2 x 2, A3 x 2, A3.5 B2, A03	4g x 8 syringes
EsCom250 Capsules	TW, A1, A2, A3, A3.5, A4, B1, B2, C2, A02, A03, D2, G1	0.25g x 20 capsules



VITA Classical shade

EsCom250's Shade Blendability

with A2 & A3 shade



SPIDENT CO., LTD.

203&312, Korea Industrial Complex,
722 Gojan-dong, Namdong-gu, Incheon, Korea, 21692
T +82-32-821-0071 F +82-32-821-0074
E synergy@spident.co.kr export@spident.co.kr