



dental component catalog

General Ordering Information

When ordering Forest products, the following information should be included:

- Quantity
- Model and/or part number
- Description
- Date required
- Purchase order number
- Method of shipping

Terms

All prices unless otherwise listed F.O.B. Hillsboro, Oregon USA.

Payments are net 30 days of invoice.

A service charge of 1½% per month (annual rate 18%) will be charged on past due accounts, \$1.00 minimum.

Minimum order \$50.00 retail.

All prices are subject to change without notice.

Return Merchandise

Dealers wishing to return merchandise for credit consideration must obtain prior factory authorization. Authorization is valid for an extended period of 90 days. After that time, if a product is not returned, authorization becomes invalid. A minimum handling fee of 10% will be charged on returned items approved for credit. Returned items must be shipped freight prepaid.

Special orders deviating from our standard product line as well as designer friendly equipment cannot be returned for credit consideration.

Customer Service

Forest Customer Service is available by phone for order entry or technical assistance Monday through Friday, 6:00 am to 4:30 pm (PST). Fax is available 24 hours a day. Voicemail is available 4:30 pm to 6:00 am (PST).

Phone: (800) 423 - 3555
Phone: (503) 640 - 3012
Fax: (503) 693 - 9715
Email: sales@forestmed.com
Website: www.forestdental.com

Warranty

Forest warrants its products to be free from defects in material and workmanship only. No other warranties are expressed or implied. This warranty shall be extended for five years after the purchase date unless otherwise stated. Written notice of breach must be given to Forest within this time period. Buyer's remedy for breach of this warranty is limited to repair parts or replacement of items by Forest. This warranty is void if items are carelessly used or improperly maintained or installed. This warranty does not cover damage that results from using cleaning, disinfecting or sterilization chemicals and processes. This warranty does not cover light bulbs. Handpiece illumination and heated syringe tubings are warranted for six months to be free from defects in material and workmanship only. Factory repairs are covered for a period of 90 days. No claim for labor or consequential damages will be allowed.

Forest Dental Products, Inc.
P.O. Box 989
Hillsboro, Oregon, 97123, USA

Shipping Address
6200 NE Campus Court
Hillsboro, Oregon 97124 USA



light service
items 1

handpiece
controls 2

arm
systems 3

syringes 4

foot
controls 5

vacuums 6

holders 7

valves 8

accessories 9

fittings 10

tubing 11

service items 12

light service items

visual identification:



Models:
S90



Models:
9040-9066



Models:
S2000/9000/9010

Note:

Part numbers with the description "9040-9066" function with the following light models:

- 9040
- 9041
- 9045
- 9046
- 9050
- 9051
- 9060
- 9061
- 9065
- 9066

face shield



0014-148



1110-816

part number	description	price
0014-148	faro, knob style, S2000/9000/9010	
1110-816	forest, 9040-9066	

face shield retainer



part number	description	price
1110-970	S2000	

on/off switch



1110-649



1110-656

part number	description	price
1110-649	2 prong, S90, S2000	
1110-656	4 prong, S2000/ 9040-9066	

reflector



1110-657



1110-808

part number	description	price
1110-657	S2000, S90	
1110-808	9040-9066	

bulb



0014-143




0014-300

part number	description	price
0014-143	S90/S2000/9000/9010	
0014-300	24V, 150 W, S90/S2000/9000/9010	

light service items

1

bulb socket



part number	description	price
0014-254	S90/S2000	
0014-301	9040-9066	

bulb shield



part number	description	price
0014-117	S90/S2000	
1110-815	9040-9066	

wire harness



part number	description	price
0014-099	yoke, S2000, low voltage, socket to switch	



0014-107	transformer to switch, arm, S90/S2000, high voltage, 3-wire	
-----------------	---	--



0014-314	head to transformer, yoke, low voltage, S90, 2 wire	
-----------------	---	--



0014-315	drop light, S2000, ceiling mount, 3-wire	
-----------------	--	--



0014-316	drop light, S90, ceiling mount, 2-wire	
-----------------	--	--



handpiece controls

2

air and water accessory panel



Standard features include:

- $\frac{3}{8}$ " female air quick disconnect
- Toggled shut off for air quick disconnect
- $\frac{1}{4}$ " female water quick disconnect
- Flow adjustment for water quick disconnect
- 4" x 3" flush mount panel
- White powder coated finish

part number	description	price
1300-030	air and water accessory panel	



LT 10 lab handpiece control kit

Includes control for one HP.

Standard features:

- Manual control for one HP, one dry
- Master on/off with indicator
- HP tubing- 2-hole Midwest, straight sterling
- Disc type foot control
- Air pressure regulator
- Umbilical, 12'



part number	description	price
2000-010	handpice control kit	

LT 100 lab handpiece control

Includes control for one HP.

Standard features include:

- Manual control for one HP, one dry
- Master on/off with indicator
- HP tubing- 2 hole Midwest, straight sterling
- Disc type foot control
- Air pressure regulator
- Umbilical, 12'



part number	description	price
2000-011	HP control	

arm systems

2" diameter ring mount flexible arm systems



Flexible counterbalanced arm system with pneumatic brake. Provides 27" vertical travel using the 24" flexible arm. Accommodates arm mounted HP control systems with tubing bundles concealed and protected within the arm. Adjustable for weights from 3 to 13½ lbs. Designed for use with 2" diameter post.

part number	description	price
1110-768	22" rigid w/ 18" flexible arm, horizontal reach 38"	
1110-769	22" rigid w/ 24" flexible arm, horizontal reach 44"	
1110-770	26" rigid w/ 18" flexible arm, horizontal reach 40"	
1110-771	26" rigid w/ 24" flexible arm, horizontal reach 66"	

top post mount flexible arm systems



Flexible counterbalanced arm system with pneumatic brake. Provides 27" vertical travel using the 24" flexible arm. Accommodates arm mounted HP control systems with tubing bundles concealed and protected within the arm. Adjustable for weights from 3 to 13½ lbs. Designed for use with 2" diameter control post with 1½" ID bearing or Arm Mount III (see page 6, p/n 1110-376).

part number	description	price
1110-760	12" rigid w /18" flexible arm, horizontal reach 26"	
1110-761	12" rigid w /24" flexible arm, horizontal reach 32"	
1110-762	22" rigid w /18" flexible arm, horizontal reach 38"	
1110-763	22" rigid w /24" flexible arm, horizontal reach 44"	
1110-764	26" rigid w /18" flexible arm, horizontal reach 40"	
1110-765	26" rigid w /24" flexible arm, horizontal reach 66"	

flexible arm



Flexible Counterbalanced Arm System with pneumatic brake. Provides 27" vertical travel. Accommodates arm mounted HP control systems with tubing bundles concealed and protected within the arm. Adjustable for weights from 3 to 13½ lbs. Designed for use with either top post or ring mounted rigid arm.

part number	description	price
1110-898	18" flexible arm	
1110-899	24" flexible arm	

arm systems

3

2' diameter ring mount rigid arm



Accommodates all lengths of flexible arms with tubing bundles concealed within the arm. Designed for use with 2" diameter post.

part number	description	price
1110-536	22" ring mount rigid arm	
1110-537	26" ring mount	

top post mount rigid arm



Accommodates all lengths of flexible arms with tubing bundles concealed within the arm. Designed for use with 2" diameter control post with 1½" inner diameter bearing or Arm Mount III (p/n 1110-376, see below).

part number	description	price
1110-525	12" top post mount	
1110-526	22" top post mount	
1110-527	26" top post mount	

bent rigid arm



Accommodates all lengths of Flexible Arms with tubing bundles concealed within the arm. Designed for use with 2" diameter control post with 1½" inner diameter bearing or Arm Mount III (p/n 1110-376, see below).

part number	description	price
1110-464	26" top post mount	

wall mount



Designed for mounting all top mount arms to any flat vertical surface and provides leveling adjustment. Mounting plate 6" x 6".

part number	description	price
1110-376	arm mount III	

wall mount



Designed for clamping a 2" diameter post to adapt to a flat vertical surface and provides leveling adjustment. Mounting plate 6" x 6".

part number	description	price
1110-386	post mount IV	

arm systems

telescoping arm system



Telescoping arm system provides 28⁵/₈" to 43¹/₂" overall horizontal reach. Designed for mounting assistants accessories and accepts 1/2" outside diameter pin.

Includes:

- Telescoping arm, 21 1/8" to 36" reach
- Rigid arm 7 1/2"
- 2" diameter post mount
- Safety release

(Order holder bar separately)

part number	description	price
1110-459	post mount	

folding arm rigid section



Accepts 1/2" outside diameter pin. Recommended: Three sections maximum.

part number	description	price
1110-446	16" section	
1110-447	12" section	
1110-450	8" section	

assistant's holder bar




For use with folding or telescoping arm systems. Accepts standard vacuum and universal holders.

Includes:

- Holder bar, 8"
- 1/2" Outside diameter pin

part number	description	price
0024-066	holder bar	

assistant's holder bar



Universal CV holder for use with folding or telescoping arm systems. Includes 1/2" outside diameter pin.


part number	description	price
0024-123	tri (3 universal)	
0024-124	quad (4 universal)	



arm systems

3


assistant's holder bar



Standard CV holder for use with folding or telescoping arm systems. Includes 1/2" outside diameter pin. Holds one standard CV valve.

part number	description	price
0024-133	tri (2 universal, 1 standard)	
0024-134	quad (3 universal, 1 standard)	

assistant's holder bar



Standard CV holder for use with folding or telescoping arm systems. Includes 1/2" outside diameter pin. Holds two standard CV valves.

part number	description	price
0024-135	tri (1 universal, 2 standard)	
0024-136	quad (2 universal, 2 standard)	

utility shelf assembly




Utility Shelf pivots on mounting bracket to provide easy access. 12" x 12" storage surface. Powder coated finish.

Includes:

- Utility shelf, 12x12"
- Mounting bracket

part number	description	price
1110-652	2" diameter post mount (includes 0023-406)	
1110-653	counter top mount (includes 1110-377)	
1110-654	wall mount (includes 1110-378)	
1110-655	2" diameter mobile stand mount (15° angle)	

tub & tray holder



For use with flexible or telescoping arm systems, the Tub & Tray Holder pivots on a rigid arm to provide easy access. Powder coated finish.

Includes:

- Utility shelf, 8³/₈" x 12 1/2"
- 12" rigid arm
- Tray pivots

part number	description	price
1110-744	utility shelf	

arm systems

arm mounted tray I



Designed for use with flexible or folding arm systems. May be used as a holder for plastic or stainless steel trays.

Includes:

- Stainless steel tray holder, 9 3/4" x 13 1/2"
- Wrap around holder bar, 1/2" diameter
- Maximum weight 15 lbs.

part number	description	price
1110-603	for folding arms	
1110-604	for forest flex	

3

arm mounted tray II



Designed for use with flexible or folding arm systems. May be used as a holder for plastic or stainless steel trays.

Includes:

- Stainless steel tray holder, 9 3/4" x 13 1/2"
- Mounting bracket with 1/2" outside diameter pin
- Wrap around holder bar, 1/2" diameter

part number	description	price
1110-606	without holder bar	
1110-608	with holder bar	

stainless steel tray



Stainless steel.

part number	description	price
0010-091	size b tray, 9 3/4" x 13 1/2"	
0010-191	tray, 10 5/8" x 15 1/8"	

plastic molded tray




Molded plastic with nonslip bottom.

part number	description	price
0010-200	size b tray, 9 5/8" x 13 3/8", white	
0010-200-B	size b tray, 9 5/8" x 13 3/8", black	



arm systems




folding/telescoping arm mounts

Designed for mounting folding and telescoping arm systems to 2" diameter post. Accepts 1/2" outside diameter pin.

part number	description	price
0023-406	2" post mount	

3

counter top mount



Designed for mounting folding and telescoping arm systems to any flat surface. Mounting plate 4" x 4". Accepts 1/2" outside diameter pin.

part number	description	price
1110-377	counter top mount	


wall mount



Designed for mounting folding and telescoping arm systems to any flat vertical surface. Mounting plate 4" x 4". Accepts 1/2" outside diameter pin.

part number	description	price
1110-378	folding/telescoping arm mount, wall mount	

mounting bracket III



Designed for mounting any counter surface on folding and telescoping arm systems. Mounting plate 4" x 4". Accepts 1/2" outside diameter pin.

part number	description	price
1110-380	mounting bracket III	

syringes

forestar-3 standard syringe assembly

Provides variable control for air, water and spray.

Includes:

- Universal autoclavable quick disconnect tip
- Tubing
- Service Tool



part number	description	price
-------------	-------------	-------

syringe only without tubing:

1107-086

with flow adjustment block:

1107-218	straight sterling tubing
1107-186	sheathed sterling tubing
1107-208	coiled sterling tubing
1107-217	straight gray tubing
1107-185	sheathed gray tubing
1107-207	coiled gray tubing
1107-213	straight black tubing
1107-286	sheathed black tubing
1107-203	coiled black tubing
1107-918	straight shadow tubing
1107-986	sheathed shadow tubing
1107-908	coiled shadow tubing

without flow adjustment block:

1107-238	straight sterling tubing
1107-196	sheathed sterling tubing
1107-228	coiled sterling tubing
1107-237	straight gray tubing
1107-195	sheathed gray tubing
1107-227	coiled gray tubing
1107-233	straight black tubing
1107-196-BLK	sheathed black tubing
1107-223	coiled black tubing
1107-938	straight shadow tubing
1107-196-SH	sheathed shadow tubing
1107-928	coiled shadow tubing

4



syringes

forestar-3 straight syringe assembly

Provides fingertip control for air, water and spray.



Includes:

- Universal autoclavable quick disconnect tip
- Tubing
- Service tool



part number	description	price
-------------	-------------	-------

syringe only without tubing:

1107-096

with flow adjustment block:

1107-246 sheathed sterling tubing

1107-245 sheathed gray tubing

1107-243 sheathed black tubing

1107-946 sheathed shadow tubing

without flow adjustment block:

1107-296 sheathed sterling tubing

1107-295 sheathed gray tubing

1107-244 sheathed black tubing

1107-995 sheathed shadow tubing



syringes

forestar-3 standard heated syringe: air actuated on/off



The Forestar-3 Standard Heated Syringe provides instant warm water at the syringe tip during normal operations. A heated element within the syringe tubing tempers the water as it advances through the system.

Includes:

- Standard syringe with universal autoclavable quick disconnect tip
- Power source with high/low and on/off controls
- Tubing
- Low voltage wire, 10'
- 110 VAC



part number	description	price
1107-418	sheathed sterling tubing	
1107-408	coiled sterling tubing	
1107-417	sheathed gray tubing	
1107-407	coiled gray tubing	
1107-413	sheathed black tubing	
1107-403	coiled black tubing	
1107-916	sheathed shadow tubing	
1107-906	coiled shadow tubing	

forestar-3 straight heated syringe



The Forestar-3 Straight Heated Syringe provides instant warm water at the syringe tip during normal operations. A heated element within the syringe tubing tempers the water as it advances through the system.

Includes:

- Straight syringe with universal autoclavable quick disconnect tip
- Power source with high/low and on/off controls
- Tubing
- Low voltage wire, 10'
- 110 VAC



part number	description	price
1107-468	sheathed sterling tubing	
1107-467	sheathed gray tubing	
1107-463	sheathed black tubing	
1107-969	sheathed shadow tubing	



syringes

heated syringe conversion kit: air actuated, normally open

Adapts to existing 3-way syringes to provide instant warm water during normal operations.



Note: Replacement tubing: see pages 72-73.

Includes:

- Power source with variable temperature and on/off controls (1/2 amp fuse, volts 10 1/2-12 with load applied)
- Size 2 3/4" x 5" x 2 1/2"
- Tubing
- Low voltage wire, 10'
- 110 VAC

part number	description	price
1107-448	sheathed sterling tubing	
1107-444	coiled sterling tubing	
1107-447	sheathed gray tubing	
1107-437	coiled gray tubing	
1107-443	sheathed black tubing	
1107-433	coiled black tubing	
1107-942	sheathed shadow tubing	
1107-947	coiled shadow tubing	
0014-004	power source only	

syringe tip kit, universal

For all Forestar-3 syringes.

Includes:

- One autoclavable syringe tip, nut, adapter and service tool
- Brown Viton o-ring



part number	description	price
1107-695	syringe tip kit	

locking syringe tip kit, forest

Locking syringe tip for Forestar-3 syringe.

Includes nut, cone, o-rings, adapter and allen wrench



part number	description	price
0015-096	o-ring kit	
1107-062	locking syringe tip kit	

syringes

autoclavable syringe tip kit (325°F), forest

Syringe tip for Forestar-3 syringe (prior to 1998).



part number	description	price
1107-020	autoclavable syringe tip (1)	
1107-063	syringe tip (package 5)	

autoclavable syringe tip kit (325°F), universal

Syringe tip for Forestar-3 & Asepsis-3 syringe (after to 1998).



part number	description	price
1107-691	universal (package 5)	
1107-693	universal (package 50)	

syringe replacement cartridge

Fits Forestar-3 syringe.



part number	description	price
1107-070	non-autoclavable with black o-rings	
1107-330	autoclavable with brown o-rings	

syringe replacement cartridge

Fits Forestar-3 syringe.



part number	description	price
1107-333	autoclavable with brown o-rings	
1107-770	non-autoclavable with black o-rings	

syringe button assembly

Syringe button assembly for Asepsis-3, DCI and Parts Warehouse syringe.



part number	description	price
1107-721	spool valve	

syringes




syringe service tool

Used for removal of cartridge, tip adapter and barbs. For Forestar-3 syringe.

part number	description	price
1107-057	syringe service tool	

4



flow adjustment block

Provides air/water metering to syringe.

Includes $\frac{1}{16}$ " & $\frac{1}{8}$ " barbs, washers and sleeves

part number	description	price
0016-092	flow adjustment block	



syringe o-ring kit, forest

Includes all o-rings and service tool. For Forestar-3 syringe.

part number	description	price
0015-089	autoclavable	
0015-099	non-autoclavable	



asepsis-3 syringe repair kit

Syringe repair kit for Asepsis-3, DCI and Parts Warehouse syringe.

part number	description	price
1107-694	spool valve	

foot controls

standard disc



Features:

- Die cast base
- Chrome plated cover
- Light touch, variable speed
- High Volume
- Vinyl tubing, 7'

Note: Service kits, see page 83 .

part number	description	price
1105-114	without relay, sterling	
1105-124	with relay, sterling	
1105-113	without relay, gray	
1105-123	with relay, gray	
1105-356	without relay, black	
1105-357	with relay, black	
1105-576	without relay, shadow	
1105-577	with relay, shadow	

5

with chip blower button

Includes all the features of the standard foot control above, plus toe operated momentary valve for chip blower.



part number	description	price
1105-134	without relay, sterling	
1105-144	with relay, sterling	
1105-133	without relay, gray	
1105-143	with relay, gray	
1105-359	without relay, black	
1105-358	with relay, black	
1105-579	without relay, shadow	
1105-578	with relay, shadow	



foot controls

wet/dry with on/off toggle



Includes all the features of the standard type (see page 19) plus toggle valve for control of water coolant.

part number	description	price
1105-174	without relay, sterling	
1105-184	with relay, sterling	
1105-173	without relay, gray	
1105-183	with relay, gray	
1105-361	without relay, black	
1105-360	with relay, black	
1105-581	without relay, shadow	
1105-580	with relay, shadow	

5

wet/dry toggle and chip blower button



Includes all the features of the standard type (see page 19) plus toggle valve for control of water coolant and toe operated momentary valve for chip blower.

part number	description	price
1105-224	without relay, sterling	
1105-234	with relay, sterling	
1105-223	without relay, gray	
1105-233	with relay, gray	
1105-363	without relay, black	
1105-362	with relay, black	
1105-583	without relay, sterling	
1105-582	with relay, sterling	

foot controls

support shelf



Provides storage for foot control while allowing cart mobility. For 2" diameter port only.

part number	description	price
0023-665	foot control support shelf	

chip blower button kit



Designed for adding chip blower feature to existing foot controls (not for use with standard disc foot control). Includes valve and mounting hardware.

Note: Order foot control tubing separately, see page 72.

part number	description	price
1105-079	2-way, chip blower	
1105-080	3-way, sealer	

wet/dry kit

Designed for adding wet/dry toggle feature to existing foot controls (not for use with standard disc foot control). Includes valve and mounting hardware.

Note: Order foot control tubing separately, see page 72.



part number	description	price
1105-082	wet/dry toggle kit	
1105-005-G	gray, toggle only (package 10)	
1105-005-BLK	black, toggle only (package 10)	
1105-006-G	gray, toggle (package 2)	
1105-006-BLK	black, toggle (package 2)	

vacuums



universal assistant kit (325°F)

Used with CV.

Includes:

- CV canister with disposable solids collector
- Autoclavable universal CV lever valve with tubing, 5'
- Autoclavable SE tip valve with tubing, 5'
- Forestar-3 standard syringe with tubing, 5'
- Three position holder
- Mounting hardware

part number				mount page #	description	price
Sterling	Gray	Black	Shadow			
1400-600	1400-500	1400-800	1400-900	36	fixed position	
1400-601	1400-501	1400-801	1400-901	33	horizontal swing	
1400-602	1400-502	1400-802	1400-902	33	vertical swing	
1400-603	1400-503	1400-803	1400-903	-	without holder	

standard assistant kit (325°F)



Used with CV.

Includes:

- CV canister with disposable solids collector
- Autoclavable standard CV lever valve with tubing, 5'
- Autoclavable SE tip valve with tubing, 5'
- Forestar-3 standard syringe with tubing, 5'
- Three position holder
- Mounting hardware

part number				mount page #	description	price
Sterling	Gray	Black	Shadow			
1400-604	1400-504	1400-804	1400-904	36	fixed position	
1400-605	1400-505	1400-805	1400-905	35	horizontal swing	
1400-606	1400-506	1400-806	1400-906	35	vertical swing	
1400-607	1400-507	1400-807	1400-907	-	without holder	

vacuums



autoclavable universal central vacuum kit, 325°F

Used with CV.

Includes:

- CV canister with disposable solids collector
- Autoclavable CV lever valves with tubing, 5'
- Autoclavable SE tip valve with tubing, 5'
- Three position holder
- Mounting hardware

part number				mount page #	description	price
Sterling	Gray	Black	Shadow			
1400-278	1400-078	1400-878	1400-978	36	fixed position	
1400-276	1400-177	1400-877	1400-977	33	horizontal swing	
1400-267	1400-167	1400-876	1400-976	33	vertical swing	
1400-274	1400-074	1400-866	1400-975	-	without holder	

6

autoclavable standard central vacuum kit (325°F)



Used with CV.

Includes:

- CV canister with disposable solids collector
- Autoclavable CV lever valves with tubing, 5'
- Autoclavable SE tip valve with tubing, 5'
- Three position holder
- Mounting hardware

part number				mount page #	description	price
Sterling	Gray	Black	Shadow			
1400-277	1400-077	1400-874	1400-974	36	fixed position	
1400-263	1400-163	1400-873	1400-973	35	horizontal swing	
1400-264	1400-164	1400-872	1400-972	35	vertical swing	
1400-282	1400-082	1400-871	1400-971	-	without holder	



vacuums

autoclavable high volume evacuator (325°F)



Autoclavable HVE features finger tip on/off and flow control with unique lever operated design. Swivel quick disconnect allows 360 degrees rotation of entire HVE and installs on any 1/2" inside diameter CV tubing. After HVE is installed on tubing, it can be removed at the swivel to allow autoclaving of HVE assembly (lever, body, and o-rings). Universal valve has one piece adapter added for use with universal HP holders.

HVE's accept 7/16" outside diameter tips.

Note: Replacement o-rings p/n 0015-313 (#13), see page 81.

part number	description	price
1100-342	standard	
1100-343	universal	

autoclavable saliva ejector (325°F)



Features finger tip on/off and flow control with unique lever operated design. Swivel quick disconnect works with universal holders allowing for 360 degrees rotation of entire SE and adapts to any 3/8" outside diameter SE tubing. After SE is installed, it can be removed at the base to allow autoclaving of SE assembly (lever, body and o-rings). HVE's accept 7/16" outside diameter tips.

Note: Replacement o-rings p/n 0015-308 (#8), see page 81.

part number	description	price
1100-341	SE	
1100-341-BLK	SE, black	

saliva ejector tip



For use with SE or assemblies.

Autoclavable SE tip 325°.

part number	description	price
1100-108	autoclavable, black	
1100-308	autoclavable, gray	

saliva ejector tip



For use with SE assemblies.

Note:

SE tip is not compatible with p/n's 1100-341 and 1100-341-BLK.

part number	description	price
1100-258	SE tip, screw-on black	



saliva ejector screen

Fits push-on or screw on style SE assemblies.

part number	description	price
0010-044	SE screen	



levers

Allows for on/off and flow control with use of finger tip. Both levers are autoclavable (325°F).

part number	description	price
1100-247	CV valve, black	
1100-248	CV valve, gray	
1100-311	SE, gray	
1100-311-BLK	SE valve, black	



plug

Stops suction and allows SE and HVE to be placed in universal holders after the SE/HVE have been removed for autoclaving or cleaning.

part number	description	price
1400-136	plug	



vacuum canister, 1/2" rear drain with locking lid

Features removable solids collector, two high velocity ports for 1/2" inside diameter tubing, one SE port for 3/8" inside diameter tubing, extra port plugs, bayonet type locking lid and rear drain. For use with telescoping arm systems.

part number	description	price
0020-101	solids collector, blue (package 10)	
1400-133	canister assembly, white	
1400-133-BLK	canister assembly, black	

vacuums

central vacuum



Features removable solids collector, two HVE, one SE port, optional $\frac{1}{2}$ " water drain port for $\frac{3}{8}$ " inside diameter tubing, optional water drain port for U.C., extra port plugs, mounting brackets and hardware.

part number	description	price
1400-111	canister assembly, gray	
1400-011	canister assembly, black	
1400-155	cap & o-ring only, black	
1400-055	cap & o-ring only, gray	
0015-033	o-ring only	
0020-100	solids collector, $1\frac{7}{8}$ outside diameter, green (package 10)	
1400-128	CV horizontal bracket, gray	
1400-128-BLK	CV horizontal bracket, black	
1400-129	CV vertical bracket, gray	
1400-129-BLK	CV vertical bracket, black	

vacuum canister, $\frac{1}{2}$ " drain with locking lid



Features are the same as above, except with bayonet type locking lid.

part number	description	price
1400-190	canister assembly, gray	
1400-090	canister assembly, black	
1400-193	cap & o-ring only, gray	
1400-093	cap & o-ring only, black	
0015-033	o-ring only	
0020-100	solids collector, $1\frac{7}{8}$ outside diameter, green (package 10)	
1400-128	CV horizontal bracket, gray	
1400-128-BLK	CV horizontal bracket, black	
1400-129	CV vertical bracket, gray	
1400-129-BLK	CV vertical bracket, black	

central vacuum canister plug



part number	description	price
1400-008	SE port, white	
1400-017	CV port, white	
1400-017-BLK	CV port, black	

vacuums



solids collector

For use with CV canisters. Green fits $1\frac{7}{8}$ " p/n 1400-111 standard and p/n 1400-190 (page 24) locking cannister assemblies. Blue fits $2\frac{1}{8}$ " p/n 1400-133. $\frac{5}{8}$ " rear drain cannister assembly.

part number	description	price
0020-100	green (package 10)	
0020-101	blue (package 10)	



solids collector with lid

For use with CV canisters.

Small- fits p/n 1400-111 (page 24) standard and p/n 1400-190 (page 24) locking cannister assembly.

Large- fits p/n 1400-133 (page 23), $\frac{5}{8}$ " rear drain canister assembly.

part number	description	price
1400-191	clear, small (each)	
1400-192	clear, large (each)	



post mount adapter

Allows either CV canister, p/n's 1400-111 (page 24) or 1400-190 (page 24), to be mounted to a 2" outside diameter post.

part number	description	price
0023-436	2" post mount adapter	



saliva ejector adapter

Adapts $\frac{3}{8}$ " outside diameter SE tubing to vacuum valve using $\frac{7}{16}$ " outside diameter tips.

part number	description	price
0016-114	$\frac{1}{4}$ " inside diameter x $\frac{7}{16}$ " inside diameter	



saliva ejector adapter

Adapts $\frac{3}{8}$ " outside diameter SE tubing to SE tip.

part number	description	price
0016-214	$\frac{1}{4}$ " inside diameter x $\frac{1}{4}$ " inside diameter	

vacuums



oral dry cup

Autoclavable polysulfone (beige). Fit all vacuum valves that accept $\frac{7}{16}$ " diameter tips.

part number	description	price
1400-010	oral dry cup	



standard HVE disposable tip

Made of lightweight polyethylene, $\frac{7}{16}$ " outside diameter x 6" long.

part number	description	price
1400-020	standard disposable hve tip (package 100)	



cuspidor strainer

Molded white plastic. Fits Forest cuspidor bowl.

part number	description	price
0010-625	cuspidor strainer (package 50)	
0010-627	cuspidor strainer (package 10)	



vacuum drain cuspidor kit

May be used to adapt all models of gravity drain cuspidor to CV system.

Includes:

- Vacuum drain valve
- Flush nipple
- Filter assembly
- Tubing
- Installation hardware

part number	description	price
1400-116	for $\frac{5}{8}$ " inside diameter drain	
1400-117	for $\frac{7}{8}$ " inside diameter drain	



air/water saliva ejector

Air/ water powered SE.

Includes:

- Vacuum generator (venturi)
- SE tip with tubing, 6"
- Toggle valve $\frac{1}{2}$ " inside diameter drain tubing, 9'
- Installation hardware

part number	description	price
1400-285	sterling tubing	
1400-085	gray tubing	
1400-885	black tubing	
1400-952	shadow tubing	

vacuums

saliva ejector hose assembly



SE assembly adapts to any vacuum valve using $\frac{7}{16}$ " outside diameter tips.

part number	description	price
1400-252	sterling tubing	
1400-052	gray tubing	
1400-852	black tubing	
1400-952	shadow tubing	

automatic dual air vacuum system



Air-operated oral evacuation system activated by automatic holders.

Includes:

- White standard vacuum lever valve with tubing, 5'
- White SE tip valve with tubing, 5'
- Two automatic holders
- Dual port venturi
- One quart (1 liter) waste container
- Air pilot valve

part number	description	price
1400-173	sterling tubing	
1400-174	gray tubin tubing	
1400-870	black tubing	
1400-970	shadow tubing	

manual dual air vacuum system



Air-operated oral evacuation system activated by manual toggle valve.

White standard vacuum lever valve with tubing, 5'

- White standard vacuum lever valve with tubing, 5'
- White SE tip valve with tubing, 5'
- Two automatic holders
- Dual port venturi
- One quart (1 liter) waste container
- Three-way toggle valve

part number	description	price
1400-183	sterling tubing	
1400-184	gray tubing	
1400-861	black tubing	
1400-961	shadow tubing	

holders



0010-235



0010-234-C



0010-835



0010-835

automatic holder for handpieces

Holder assembly with micro valve and actuator. Provides air signal to automatic control block for holdback when HP is placed into the holder.

Features:

- Micro valve and actuator
- On/off lockout toggle with red indicator
- Non-marring, plastic holder
- Aluminum mounting adapter for 1/2" diameter holder bar

part number	description	price
0010-235	with lockout toggle, gray	
0010-017	with lockout toggle, black	
0010-234	without lockout toggle, white	
0010-835-G	fixed mount holder bar, gray	
0010-835-SH	fixed mount holder bar, shadow	
0010-835-BLK	fixed mount holder bar, black	
0010-836-G	1/2" diameter holder bar, gray	
0010-836-SH	1/2" diameter holder bar, shadow	
0010-836-BLK	1/2" diameter holder bar, black	

7



0010-245



0010-094

automatic holder for accessories

Holder assembly with micro valve and actuator, provides air signal to an accessory, such as an air electric switch, while HP is out of the holder.

Features:

- Micro signal valve and actuator
- On/off lockout toggle with red indicator
- Non-marring, plastic holder
- Aluminum mounting adapter for 1/2" diameter holder bar

part number	description	price
0010-094	with lockout toggle, black	
0010-245	with lockout toggle, gray, with red indicator	
0010-846-G	1/2" diameter holder bar, gray	
0010-846-SH	1/2" diameter holder bar, shadow	
0010-846-BLK	1/2" diameter holder bar, black	
0010-847-G	fixed mount holder bar, gray	
0010-847-SH	fixed mount holder bar, shadow	
0010-847-BLK	fixed mount holder bar, black	

Images
Coming Soon :

0010-846 &
0010-847

holders

electric automatic holder



0010-255

Holder incorporates a low voltage, normally open or normally closed switch for activation of an electrical accessory when the HP is out of the holder. Normal state of the switch is when the SE or accessory is in holder.

Features:

- Underwriters Laboratories listed low voltage switch with 24" pigtail
- Lockout toggle with indicator
- Non-marring, plastic holder
- Aluminum mounting adapter for 1/2" diameter holder bar



0010-254

part number	description	price
-------------	-------------	-------

normally open:

0010-119	with lockout toggle, black
0010-254	without lockout toggle, gray
0010-255	with lockout toggle, gray
0010-856-G	1/2" diameter holder bar, gray
0010-856-SH	1/2" diameter holder bar,
0010-856-BLK	1/2" diameter holder bar, black
0010-857-G	fixed mount holder bar, gray
0010-857-SH	fixed mount holder bar, shadow
0010-857-BLK	fixed mount holder bar, black



0010-856



0010-857

normally closed:

0010-118	with lockout toggle, black
0010-264	without lockout toggle, gray
0010-265	with lockout toggle, gray

holders

vacuum holder, electric



0010-260

Holder incorporates a low voltage, normally open or normally closed switch for activation of an electrical accessory when the HP is out of the holder. Normal state of the switch is when the HVE is in holder.

Features:

- Underwriters Laboratories listed low voltage switch with 24" pigtail
- On/off lockout toggle with indicator
- Non-marring, plastic holder
- Aluminum mounting adapter for 1/2" diameter holder bar



0010-269



0010-858



0010-859

part number	description	price
-------------	-------------	-------

normally open:

0010-858-G	1/2" diameter holder bar, gray
0010-858-SH	1/2" diameter holder bar, shadow
0010-858-BLK	1/2" diameter holder bar, black
0010-859-G	fixed mount holder bar, gray
0010-859-SH	fixed mount holder bar, shadow
0010-859-BLK	fixed mount holder bar, black
0010-260	with lockout toggle, gray

normally closed:

0010-868-G	1/2" diameter holder bar, gray
0010-868-SH	1/2" diameter holder bar, shadow
0010-868-BLK	1/2" diameter holder bar, black
0010-269	without lockout toggle, gray
0010-270	with lockout toggle, gray



holders



0010-275

universal HVE holder

Holder assembly for use with HP, syringe, SE or HVE.

Features:

- Non-marring, plastic holder
- Aluminum mounting adapter for $\frac{1}{2}$ " diameter holder bar



0010-886



0010-885

part number	description	price
0010-275	universal, gray	
0010-025	universal, black	
0010-885-G	fixed mount holder, gray	
0010-885-SH	fixed mount holder, shadow	
0010-885-BLK	fixed mount holder, black	
0010-886-G	$\frac{1}{2}$ " diameter holder bar, gray	
0010-886-SH	$\frac{1}{2}$ " diameter holder bar, shadow	
0010-886-BLK	$\frac{1}{2}$ " diameter holder bar, black	



0010-285

HVE holder

Accepts HVE.

Features:

- Non-marring, plastic holder
- Aluminum mounting adapter for $\frac{1}{2}$ " diameter holder bar



0010-887



0010-888

part number	description	price
0010-285	standard, gray	
0010-227	standard, black	
0010-887-G	$\frac{1}{2}$ " diameter holder bar, gray	
0010-887-SH	$\frac{1}{2}$ " diameter holder bar, shadow	
0010-887-BLK	$\frac{1}{2}$ " diameter holder bar, black	
0010-888-G	fixed mount holder, gray	
0010-888-SH	fixed mount holder, shadow	
0010-888-BLK	fixed mount holder, black	

holders

replacement holder only

Non-marring, plastic holder.



part number	description	price
0010-212	auto, gray	
0010-018	auto, black	
0010-214	universal, gray	
0010-019	universal, black	
0010-216	vacuum HP, gray	
0010-219	vacuum HP, black	
0010-231	avs, gray	



0010-837

automatic holder, HVE

Non-marring, plastic holder.



0010-838

part number	description	price
0010-837-G	auto molded, gray	
0010-837-SH	auto molded, shadow	
0010-837-BLK	auto molded, black	
0010-838-G	molded notched, gray	
0010-838-SH	molded notched, shadow	
0010-838-BLK	molded notched, black	

avs handpiece holder

Features:

- Non-marring, plastic holder
- Aluminum mounting adapter for 1/2" diameter holder bar



part number	description	price
0010-276	avs, gray	

automatic holder micro valve

For use with automatic HP or automatic accessory holders.



part number	description	price
1100-350	HP holder micro valve normally closed	
1100-351	accessory holder valve normally open	

holders

holder bar mounted valve



Holder assembly with miniature valve and on/off toggle, provides manual air signal to accessories.

Features:

- Micro signal valve
- Lockout toggle with red indicator
- Aluminum mounting adapter for 1/2" diameter holder bar

part number	description	price
1100-285	holder bar mounted valve, gray	

tri holder



part number	description	price
0024-113	Delrin®, white	
0024-113-BLK	Delrin®, black	

quad holder



part number	description	price
0024-114	Delrin®, white	
0024-114-BLK	Delrin®, black	

swing mount holder




horizontal mount

Holder assembly for use with HP, syringe, SE or universal CV valve. Swivel mount allows holder to pivot for easy access.

Features:

- Delrin® holder bar
- Vertical or horizontal surface swing mount



vertical mount

part number	description	price
vertical surface mount:		
0024-115	tri: 3 universal	
0024-116	quad: 4 universal	
horizontal surface mount:		
0024-117	tri: 3 universal	
0024-118	quad: 4 universal	

holders

tri holder double HVE



part number	description	price
0024-125	Delrin®, tri holder double hve, white	
0024-125-BLK	Delrin®, tri holder double hve, black	

tri holder with HVE



part number	description	price
0024-126	Delrin®, tri holder with hve, white	
0024-126-BLK	Delrin®, tri holder with hve, black	

quad holder double HVE



part number	description	price
0024-127	Delrin®, quad holder double hve, white	
0024-127-BLK	Delrin®, quad holder double hve, black	

quad holder with HVE



part number	description	price
0024-128	Delrin®, quad holder with hve, white	
0024-128-BLK	Delrin®, quad holder with hve, black	

swing mount holder bar



vertical mount



horizontal mount

For use with automatic, universal and vacuum HP holders. Swivel mount allows holder bar to pivot for easy access.

Features:

- 1/2" diameter holder bar
- Vertical or horizontal surface swing mount
- Holds up to 4 holders

part number	description	price
0024-137	Delrin®, vertical surface mount	
0024-138	Delrin®, horizontal surface mount	

holders




horizontal mount

swing mount holder

Holder assembly for use with syringe or SE and two standard CV valves. Swivel mount allows holder to pivot for easy access.

Features:

- Delrin® holder bar
- Vertical or horizontal surface swing mount



vertical mount

part number	description	price
-------------	-------------	-------

vertical surface mount:

0024-141 tri: 1 universal, 2 standard

0024-142 quad: 2 universal, 2 standard

horizontal surface mount:

0024-143 tri: 1 universal, 2 standard

0024-144 quad: 2 universal, 2 standard




horizontal mount

swing mount holder

Holder assembly for use with syringe, SE and standard CV valve. Swivel mount allows holder to pivot for easy access.

Features:

- Delrin® holder bar
- Vertical or horizontal surface swing mount



vertical mount

part number	description	price
-------------	-------------	-------

vertical surface mount:

0024-145 tri: 2 universal, 1 standard

0024-146 quad: 3 universal, 1 standard

horizontal surface mount:

0024-147 tri: 3 universal, 1 standard

0024-148 quad: 4 universal, 1 standard



durr holder

part number	description	price
-------------	-------------	-------

0024-205 1091 int

0024-215 1591 int



holders

fixed mount holder bars

For use with automatic, universal and vacuum HP holders.



Features:

- Vertical surface mount
- $\frac{1}{2}$ " diameter

part number	description	price
0024-081	single $1\frac{1}{2}$ "	
0024-082	duo 4"	
0024-083	tri $6\frac{1}{4}$ "	
0024-084	quad $7\frac{3}{4}$ "	
0024-087	universal $14\frac{3}{8}$ "	
0024-088	small universal $10\frac{1}{2}$ "	

7

wrap around holder bars

$\frac{1}{2}$ " diameter holder bar for custom applications.



part number	description	price
0024-090	12" holder bar	
0024-091	16" holder bar	



0023-165

holder bar brackets

For use with wrap around holder bars. Includes set screws.



0023-166

part number	description	price
0023-165	$\frac{1}{2}$ " hole (1 set)	
0023-166	without hole (1 set)	

valves

toggle valve, on/off



Feautres:

- Positive on/ off control
- Brass construction
- Nickel plated
- Port size 10-32
- Panel mount with nuts, 15/32-32

part number	description	price
-------------	-------------	-------

2-way without exhaust:

1100-402	black	
1100-102	gray	

3-way with exhaust:

1100-103	black	
1100-403	gray	

toggle only:

0010-023	black	
0010-209	gray	

toggle valve, on/off, momentary



Feautres:

- Positive on/ off control
- Brass construction
- Nickel plated
- Port size 10-32
- Panel mount with nuts, 15/32-32

Note: No toggle only for black.

part number	description	price
-------------	-------------	-------

2-way without exhaust:

1100-125	black	
1100-425	gray	

3-way with exhaust:

1100-124	black	
1100-424	gray	

toggle only:

0010-208	gray	
-----------------	------	--

valves

push button, on/off, momentary



Feautres:

- Positive on/off control
- Brass construction
- Nickel plated
- Port size 10-32
- 15/32-32 panel mount with nuts

part number	description	price
2-way without exhaust, normally closed:		
1100-441	gray	
3-way with exhaust:		
1100-142	toggle, normally closed, gray	
1100-144	toggle, normally open, black	
1100-371	toggle, normally closed, black	

needle control valve



Feautres:

- Variable control for air and water
- Brass construction
- Nickel plated
- Port size 10-32
- 15/32-32 panel mount with nuts

Note: See below for knob.

part number	description	price
1100-100	flow without barbs	
1100-117	with long shaft	
1100-181	with style II knob & barbs, gray	
1100-372	with style II knob & barbs, black	

flow restrictor



Feautres:

- Variable restriction for air and water (pinch type)
- Accommodates tubing up to 1/4" outside diameter
- Brass construction
- Nickel plated
- Mounts directly on tubing

part number	description	price
1100-130	flow restrictor	

valves



style I



style II

plastic knob

Replacement knob for control valves with $\frac{1}{4}$ " outside diameter shaft.

part number	description	price
-------------	-------------	-------

style I:

0010-201

gray

style II:

0010-012

black

0010-202

gray



decorative trim nut

Decorative panel mount nut, 15/32-32.

part number	description	price
-------------	-------------	-------

0013-016

smooth nut, silver

valve, on/off: pilot air operated

Feautres:

- Positive on/off control
- Brass construction
- Nickel plated
- $\frac{1}{16}$ " barb fitting
- Port size 10-32
- Mounting hardware 15/32-32



part number	description	price
-------------	-------------	-------

1100-136

2-way on/off, without exhaust
normally closed

1100-137

3-way on/off with exhaust
normally closed

valve, on/off: pilot air operated

Feautres:

- Positive on/off control
- Brass construction
- Nickel plated
- $\frac{1}{16}$ " barb fitting
- Port size 10-32



part number	description	price
-------------	-------------	-------

1100-138

2-way on/off without exhaust,
normally open



valves

check valve

Allows air or water flow in one direction only.



Features:

- Brass construction
- Nickel plated $\frac{1}{16}$ " barb fitting
- Port size 10-32

part number	description	price
1100-105	check valve	

shuttle valve

Designed with two inlets and one outlet. Valve permits free flow from the high pressure side, while blocking the low pressure side.



Features:

- For use with air or water
- Inport (shuffled), 2
- Outport, 1
- Brass construction
- Nickel plated
- $\frac{1}{16}$ " x 10-32 barbs with sleeve

part number	description	price
1100-140	shuttle valve	

water relay valve

Designed for on/off control of HP water coolant.



Features:

- Non-retraction
- Brass construction
- Nickel plated
- $\frac{1}{16}$ " barb fitting
- Port size 10-32

part number	description	price
1100-039	water relay valve	

valves

water relay valve

Designed for on/off control of HP water coolant.

Features:

- Non-retraction
- Positive on/off control
- Brass construction
- Port size 10-32



part number	description	price
1100-145	2-way, without exhaust, normally closed	

routing valve

Features:

- Routing control for air or water
- Brass construction
- Nickel plated
- 15/32-32" panel mount with nuts
- Port size 10-32



part number	description	price
1100-104	black	
1100-404	gray	

valve mounting bracket

Used for mounting valves or switches. 1/2" inside diameter holes.



part number	description	price
0023-170	single	
0023-171	dual	
0023-282	dual	

automatic handpiece control block

Miniature assembly for selection and operation of drive air and air/water coolants for two to four air driven HP. Designed for use with automatic HP holders (see page 28).



part number	description	price
1100-175	duo-block (2 HP)	
1100-170	tri-block (3 HP)	
1100-270	quad-block (4 HP)	
0022-021	duo-block diameterphragm	
0022-020	tri-block diameterphragm	
0022-024	quad-block diameterphragm	

valves



1100-059



1100-365

handpiece logic valve

Designed for use with automatic HP holder.

part number	description	price
1100-059	pinch style, water cooland on/off (1 per HP)	
1100-365	flex valve (1 per HP)	

rectangular gauge

Features:

- 0-60psi
- 1/2" high x 1 3/4" wide
- Panel mount
- 10-32 port
- For use with air only



part number	description	price
0020-023	retangular gauge	

panel mount gauge

Features:

- 0-100psi
- 5/8" diameter
- 1/8" mpt with 10-32 port
- For use with air or water



part number	description	price
0020-024	panel mount gauge	

round gauge

Features:

- 1/8" MPT fitting with 10-32 port
- For use with air or water



part number	description	price
0020-021	1 1/2" diameter gauge (0-160 psi)	
0020-030	1" diameter gauge (0-60 psi)	
0020-033	1" diameter gauge (0-160 psi)	

valves

pneumatic indicator



Features:

- 25-120 psi operating range
- 10-32 port
- 15/32-32 panel mount with nut

part number	description	price
-------------	-------------	-------

green indicator:

0020-068	black bezel
0020-078	gray bezel

red indicator:

0020-069	black bezel
0020-079	gray bezel



umbilical cuff

Designed to adapt umbilical with spiral flex housing to utility center. Includes mounting hardware.

Note: See page 76 for custom umbilicals.

part number	description	price
-------------	-------------	-------

1115-090	for $\frac{3}{4}$ " inside diameter tubing
1115-091	for 1" inside diameter tubing



umbilical cuff

Designed to adapt 1 $\frac{1}{4}$ " inside diameter umbilical with spiral flex housing to utility center. Contains set of two cuffs.

Note: See page 76 for custom umbilicals.

part number	description	price
-------------	-------------	-------

1115-088	umbilical cuff (set of 2)
-----------------	---------------------------



umbilical cuff

Designed to adapt 2" inside diameter umbilical with spiral flex housing to utility center. Contains set of two cuffs.

Note: See page 76 for custom umbilicals.


part number	description	price
-------------	-------------	-------

1115-089	umbilical cuff (set of 2)
-----------------	---------------------------



valves

standard utility center



Features:

- Adjustable air and water pressure regulators with filters
- Master air and water shut-off valves
- Aluminum housing, powder coat finish
- Plastic cover, 4½" x 9" x 14"

part number	description	price
1115-001	complete utility center (accepts all umbilical cuffs except 2")	
1115-006	complete utility center (accepts 2" umbilical cuffs)	
1115-004	without housing & cover	
1115-020	housing & cover only, std	
1115-028	housing & cover only, 2"	
1115-022	cover only	

vacuum utility center




Features:

- Master air and water shut-off valves
- Aluminum housing, powder coat finish
- Plastic cover 4½" x 9" x 14"
- Accepts all umbilical cuffs except 2"

part number	description	price
1115-003	standard vacuum utility center	
1115-300	deluxe vacuum utility center	

complete deluxe utility center



deluxe utility center

Features:

- Right umbilical access
- Air pressure regulator with filter
- Master air shut-off valve



cover

part number	description	price
1115-005	factory option	
1115-005-G	forest grey	
1115-005-SH	shadow	
1115-005-BLK	black	

valves

deluxe utility center housing & cover

7 x 13 x 12⁵/₁₆ x 12¹¹/₁₆



part number	description		price
Rear umbilical access:			
1115-007	factory option		
1115-007-G	forest grey		
1115-007-SH	shadow		
1115-007-BLK	black		
Right umbilical access:			
1115-009	factory option		
1115-009-G	forest grey		
1115-009-SH	shadow		
1115-009-BLK	black		



regulator, air or water

Features:

- 0-125 psi operating range
- Brass body
- Adjustment knob: pull to unlock and adjust, push to lock
- Port size 1/8" ftp
- Fittings not included

part number	description		price
0020-020	regulator		

8

regulator with filter

Non-relieving regulator for air or water.

Features:

- 0-125 psi operating range
- Regulator with gauge
- Mounting bracket
- Filter with removable element, 35 micron filtration
- Inlet, 1/4" poly fitting
- Outlet, 10-32 position
- Barbs (1/8") with sleeves included



part number	description		price
0020-001	regulator with filter		



valves

regulator with filter and master shut-off

Non-relieving regulator for air or water service.

Features:

- 0-125 psi operating range
- Regulator with gauge
- Mounting bracket
- Filter with removable element, 35 micron filtration
- Adjustment knob: pull to unlock and adjust, push to lock
- Inlet, 1/4" poly fitting
- Outlet, 10-32 position



Note: Requires 1/2" MPT to 3/8" copper compression fitting (p/n 0016-036, page 63) to adapt to manual shut-off valve.

part number	description	price
0020-050	regulator with filter and master shut-off	

in-line filter assembly

Features:

- Nickel plated brass construction
- 35 micron filter disc
- Port size, 1/8" fpt
- Fittings 1/4" poly x 1/8" mpt
- 1" diameter x 1" long



part number	description	price
0020-039	filter assembly	
0020-040	filter only	
0020-041	replacement filter disc	

miniature adjustable regulator

Features:

For air or water service.

- 0-150 psi operating range
- Nickel plated brass construction
- Port size, 10-32
- Barbs with sleeves included
- 15/32-32 panel mount with nuts



part number	description	price
0020-025	mini regulator	
0020-043	mini regulator, in-line fixed, 40psi	

regulator mounting bracket

Used to mount regulator to wall or cabinet surface.



part number	description	price
0020-018	regulator mounting bracket	

valves

master shut-off valve



For air or water service.

Features:

- Nickel plated brass construction
- Inlet, 1/2" fpt
- Outlet, 1/4 poly fitting
- 1/16" x 10-32 barbs with seeves

Note: Require 1/2" mpt to 3/8" copper compression fitting (p/n 0016-036, page 63) to adapt to manual shut-off valve.

part number	description	price
1100-069	master shut-off valve	

master shut-off valve kit



For remote utility center installation

Features:

- Master shut-off valves (2), 1/2" FPT inlet
- Toggle valve
- Tubing
- Fittings

Note: Require 1/2" mpt to 3/8" copper compression fitting (p/n 0016-036, page 63) to adapt to manual shut-off valve.

part number	description	price
1100-079	master shut-off valve kit	

manual shut-off valve



Used at utility center to turn off air or water during installation or service.

Features:

- Chrome finish
- Inlet 1/2" ftp
- Outlet, 3/8" inside diameter compressin fitting

part number	description	price
1100-080	manual shut-off valve	

image
coming
soon

master shut-off valve

part number	description	price
1100-077	master shut-off vlave, air	
1100-078	master shut-off vlave, water	

accessories



self-contained water system with gauge, single bottle

Controlled dispensing system for antiseptic coolant or sterile water. May be used with HP controls and 3-way syringes.

Features:

- High density bottle, 1.5 liters
- Mini-air pressure regulator
- Interconnect tubing
- Mounting hardware

part number	description	price
1700-200	2" diameter post mount system	
1700-201	wall mount system	
1100-283	replacement pressure head with gauge	
0010-591	replacement water bottle/pressure head w/ gauge	
0011-011	replacement washer	



self-contained water system with gauge, dual bottle

Controlled dispensing system for antiseptic coolant or sterile water. Manual selection of one of the two bottles. May be used with HP controls and 3-way syringes.

Features:

- High density bottle, 1.5 liters (2)
- Mini-air pressure regulator
- Three-way on/off toggle valve
- Interconnect tubing
- Mounting hardware

part number	description	price
1700-226	2" diameter post mount system	
1700-221	wall mount system	
1700-220	beneath shelf mount system	
1100-281	replacement pressure head	
0010-593	replacement water bottle/pressure head w/ gauge	
0011-011	replacement washer	



high density bottle

For use with clean water systems.

Features:

- Opaque bottle
- High density polyethylene (HDPE)

part number	description	price
0010-364	high density 1 liter bottle	
0010-464	high density 1½ liter bottle	

accessories

oralite handpiece illumination system



A compact HP illumination system that installs easily to any HP control and is compatible with any standard power optic HP. An easy to replace krypton lamp provides brilliant, cool light for intra-oral use.

Note: 220 VAC is only offered in a euro style plug.

Features:

- ISO standard power optic HP terminal with stainless steel adapter and nut (5-hole Midwest)
- Asepsis style HP tubing, 7'
- Krypton lamp, 3.5 VAC, 700 mA
- Power transformer, output 3.6 VAC, 720 mA
- Air electric switch, normally open
- Installation hardware

part number	description	price
110 VAC (60 Hz):		
1700-018	straight sterling tubing	
1700-022	coiled sterling tubing	
1700-010	straight gray tubing	
1700-008	coiled gray tubing	
1700-002	straight black tubing	
1700-003	coiled black tubing	
1700-518	straight shadow tubing	
1700-522	coiled shadow tubing	
220 VAC (60 Hz):		
1700-048	straight sterling tubing	
1700-049	coiled sterling tubing	
1700-038	straight gray tubing	
1700-036	coiled gray tubing	
1700-039	straight black tubing	
1700-037	coiled black tubing	
1700-548	straight shadow tubing	
1700-549	coiled shadow tubing	

accessories

handpiece illumination tubing kit

Used to add extra power optic HP to oralite illumination system.



Features:

- ISO standard power optic HP terminal with stainless steel adapter and nut (5-hole Midwest)
- Asepsis style HP tubing, 7'
- Krypton lamp, 3.5 VAC, 720mA
- Air electric switch, normally open
- Installation hardware

part number	description	price
0017-694	straight sterling tubing	
0017-700	coiled sterling tubing	
0017-684	straight gray tubing	
0017-683	coiled gray tubing	
0017-682	straight black tubing	
0017-681	coiled black tubing	
0017-584	straight shadow tubing	
0017-500	coiled shadow tubing	

handpiece illumination replacement tubing

Used to replace power optic HP tubing.



Features:

- ISO standard power optic HP terminal with stainless steel adapter and nut (5-hole Midwest)
- Asepsis style HP tubing, 7'
- Krypton lamp, 3.5 VAC, 720 mA

part number	description	price
0017-692	straight sterling tubing	
0017-699	coiled sterling tubing	
0017-686	straight gray tubing	
0017-689	coiled gray tubing	
0017-675	straight black tubing	
0017-688	coiled black tubing	
0017-592	straight shadow tubing	
0017-599	coiled shadow tubing	

accessories

handpiece illumination system, 6-pin



A compact HP illumination system that installs easily to any HP control and is compatible with any ISO-C HP (6-pin).

Note: 220 VAC is only offered in a euro style plug. (ce 7120) (en 50075)

Features:

- ISO-C standard power optic HP terminal (6-Pin)
- Asepsis style HP tubing, 7'
- Power transformer, output 3.6 VAC, 720 mA
- Air electric switch, normally open
- Installation hardware

part number	description	price
110 VAC (60 Hz):		
1700-603	straight sterling tubing	
1700-600	straight gray tubing	
1700-800	straight black tubing	
1700-904	straight shadow tubing	
220 VAC (50 Hz):		
1700-604	straight sterling tubing	
1700-601	straight gray tubing	
1700-801	straight black tubing	
1700-902	straight shadow tubing	

handpiece illumination kit, 6-pin



Used to add extra power optic HP to above illumination system.

Features:

- ISO-C standard power optic HP terminal (6-pin)
- Asepsis style HP tubing, 7'
- Air electric switch, normally open
- Installation hardware

part number	description	price
0017-994	straight sterling tubing	
0017-984	straight gray tubing	
0017-982	straight black tubing	
0017-985	straight shadow tubing	

accessories

handpiece illumination 6-pin



Features ISO-C HP tubing, 7' with connector (6-Pin)

part number	description	price
0017-603	straight sterling tubing	
0017-600	straight gray tubing	
0017-903	straight black tubing	
0017-900	straight shadow tubing	

transformer



Note: 220 OVAC is only offered in a euro style plug.

part number	description	price
0014-115	3.6 VAC, 700 mA, 200 OVAC	

air actuated electric switch



Features:

- 120 vac, 5 amp switch
- 15-150 psi operating range
- 10-32 x 1/16" barb fitting
- 18 awg wire, 6"
- 5/32"-32 panel mount with nut

part number	description	price
0014-085	normally open	
0014-088	normally closed	

wire nut



Used for oralite HP illumination system wire connections in lieu of terminal strip, 16-22 gauge.

part number	description	price
0014-182	wire nut (package of 10)	

accessories

time delay kit



Time delay module for use with oralite HP illumination system. Easily installed kit retro-fits all existing oralite systems. One required per HP.

Features:

- 10 second delay after foot control is released
- Complete installation instructions

part number	description	price
0014-122	time delay kit	

terminal strip



Used for oralite HP illumination system p/n 1110-481 (page 55) wire connections.

part number	description	price
0014-064	6 circuit	
0014-066	4 circuit	

krypton lamp



For HP illumination system, 3.5VAC, 720mA.

part number	description	price
0014-067	krypton lamp	

accessories

storage organizer assembly

For use with mobile cabinets or assistant's carts. Provides two position rest for storage organizer (10° and 30°).

image
coming
soon

Features:

- Storage organizer with removable hinged cover
- Organizer holder assembly
- Universal mounting brackets, horizontal surface

part number	description	price
1110-350	4 position, 22 ³ / ₄ " width	
1110-351	3 position, 19 ³ / ₄ " width	

storage organizer with cover

For use with storage organizer holders.



Features:

- Storage organizer with removable hinged cover
- Storage dividers

part number	description	price
0010-164	hinge only	
0010-165	storage organizer with cover	
0010-343	full divider	
0010-344	half divider	

mobile stand

Features:

- Powder coat finish
- Adjustable height, 23" to 32"
- 2" diameter
- 15° angled upright
- Mounting plate, 6" x 6"
- Twin wheel casters



part number	description	price
4003-097	small stand, 18"w x 18"d	

accessories

h-type mobile stand



Features:

- Powder coat paint finish
- Adjustable height, 26" to 39"
- 2" diameter
- Vertical upright
- Mounting plate, 4½" x 6"
- Twin wheel casters

part number	description	price
4003-135	large, 24½ w x 18" d	

control post



Features:

- Powder coat paint finish
- 2" diameter
- 1½" inside diameter bronze bearing

part number	description	price
1110-480	control post, 24"	
1110-481	control post with slot, 24"	
1110-588	control post, 30"	
1110-592	control post with slot, 30"	

post cap



Provides finished end for above control posts when using 2" diameter post mount arms. Plastic.

part number	description	price
0010-140	post cap, control post	

accessories

post cap



Provides finished top for 2" diameter cuspidor posts, powder coat painted.

part number	description	price
0010-141	post cap, cuspidor post	


cart sleeve



Delrin® spacer for enlarging diameter of post to accept 2" diameter clamp.

part number	description	price
0010-195	cart sleeve, spacer, post enlarger	

post clamp



part number	description	price
1110-516	light post clamp only	

light post extension

Designed for used with any 2" diameter post.



Features:

- Powder coat paint finish
- 24" light post
- Post clamp
- Light adapter
- Power cord, 8'

part number	description	price
1110-517	for pelton & crane, 24"	
1110-518	for forest/castle, 24"	
1110-820	for forest/ castle, 35"	
1110-822	for pelton & crane, 35"	

accessories

light post



Features:

- Powder coat paint finish
- 2" outside diameter, 1.871 +/- .002 inner diameter

part number	description	price
1110-503	light post, 24"	
1110-504	light post, 40"	
1110-578	light post, 30"	
1110-643	light post, 44"	
1110-665	light post, 35"	

bushing



Designed to adapt light to above posts.

part number	description	price
1110-510	pelton & crane 1.705 +/- .001 inner diameter	
1110-582	forest/castle 1.378 +/- .002 inner diameter	

floor mounted



beneath chair

Features:

- Powder coat painted post, 2" outside diameter x 54" L
- Powder coat painted base
- 120V duplex outlet
- Powder cord
- Light adapters



beside chair

part number	description	price
beside chair:		
1115-094-F	for forest/castle	
1115-094-PC	for pelton & crane	
beneath chair:		
1115-098-F	for forest/castle	
1115-098-PC	for pelton & crane	

fittings



adjustable tie

part number	description	price
0016-167	4" length (package 10)	
0016-168	5½" length (package 10)	
0016-169	10" length (package 10)	



adjustable tie mount

Adhesive mount.

part number	description	price
0019-040	adjustable tie mount (package 10)	



clip

For use with 1/8" outside diameter tubing.

part number	description	price
0019-056	tube clip (package 10)	



sleeves

part number	description	price
0016-103	1/8" inside diameter (package 10)	
0016-104	1/4" inside diameter (package 10)	



uniclamp

part number	description	price
0016-278	1/8" inside diameter (package 10)	

fittings



plastic washer

part number	description	price
0016-264	for 6-32 fittings (package 100)	
0016-265	for 8-32 fittings (package 100)	
0016-266	for 10-32 fittings (package 100)	



plug

Includes washers. Plated.

part number	description	price
0016-165	hex plug, 10-32 (package 10)	



plug

part number	description	price
0016-073	socket head plug, 1/8" mpt (each)	



plug

part number	description	price
0016-074	socket head plug, 1/4" mpt (each)	



plug

Plated.

part number	description	price
0016-238	1/8" barb plug (package 10)	
0016-239	1/4" barb plug (package 10)	



connector

Includes washers.

part number	description	price
0016-262	coupler, male 10-32 connector (each)	

10



fittings



connector

Includes washers. Plated.

part number	description	price
0016-260	1/16" x 6-32 barb (package 10)	
0016-269	1/16" x 8-32 barb (package 10)	



connector

Includes washers. Plated.

part number	description	price
0016-102	1/16" x 10-32 (package 10)	



connector

Includes washers and sleeves. Plated.

part number	description	price
0016-101	barb, 1/8" x 10-32 (package 10)	



connector

part number	description	price
0016-296	1/16" x 1/8" mpt barb	



connector

part number	description	price
0016-297	1/8" mpt x 1/8" barb	



connector

part number	description	price
0016-210	10-32 x 1/8" & 1/16" stepped barb	



connector

part number	description	price
0016-252	1/8" in-line, twist to disconnect	

fittings

connector

Includes sleeves. Plastic.



part number	description	price
0016-151	1/16" x 10-32 barb elbow (package 10)	
0016-152	1/8" x 10-32 barb elbow (package 10)	

connector

Includes sleeves. Plastic.



part number	description	price
0016-155	1/16" barb tee (package 10)	
0016-156	1/8" barb tee (package 10)	
0016-157	1/8" x 1/16" x 1/8" barb tee (package 10)	

connector

Includes sleeves.



part number	description	price
0016-158	1/16" barb connector (package 10)	
0016-159	1/8" barb connector (package 10)	

connector

Includes sleeves. Plastic.



part number	description	price
0016-109	1/16" x 1/8" barb connector (package 10)	

connector

Includes sleeves. Plastic.



part number	description	price
0016-275	1/16" barb, "wye" (package 10)	

connector

Plastic.



part number	description	price
0016-433	1/16" barb cross	
0016-434	1/8" barb cross	

10



fittings



connector

part number	description	price
0016-140	6-port quick disconnect manifold assembly	
0016-154	8-port quick disconnect manifold assembly	
0022-023	6-port gasket	
0022-025	8-port gasket	



connector

Plated.

part number	description	price
0016-060	10-32, 4-way cross, female	



connector

Plated.

part number	description	price
0016-079	1/8" fpt, 4-way cross	



connector

part number	description	price
0016-080	1/8" fpt elbow	
0016-081	1/4" fpm elbow	



connector

part number	description	price
0016-129	1/8" mpt x 1/8" fpt elbow	



connector

part number	description	price
0016-100	1/8" mpt x 10-32 tee	



connector

part number	description	price
0016-153	1/8" mpt x 10-32 elbow	

fittings



connector

part number	description	price
0016-127	10-32 adjustable elbow	



connector

part number	description	price
0016-436	1/8" fpt x 10-32	



connector

Plated.

part number	description	price
0016-087	1/8" fpt x 10-32 adapter	



connector

Plated.

part number	description	price
0016-075	1/8" mpt x 1/4" fpt adapter	



connector

Plated.

part number	description	price
0016-071	1/4" mpt x 1/8" fpt reducer	
0016-300	1/2" mpt x 1/4" fpt reducer	



connector

part number	description	price
0016-049	10-32 to 1/4" poly	
0016-144	10-32 to 3/8" copper compression	



connector

part number	description	price
0016-029	1/8" mpt to 3/8" copper compression	
0016-036	1/2" mpt to 3/8" copper compression	



connector

part number	description	price
0016-077	1/8" mpt close nipple	

10



fittings



connector


part number	description	price
0016-083	1/8" mpt x 2" nipple	



connector

Plastic.

part number	description	price
0016-183	1/2" mpt x 1/2" barb connector	



connector

Plastic.

part number	description	price
0016-184	1/2" mpt x 1/2" barb elbow adapter	



connector

part number	description	price
0016-180	1/2" x 5/8" adapter	
0016-185	1/2" x 1/2" adapter	
0016-198	5/8" x 7/8" adapter	
0016-199	1/2" x 7/8" adapter	



connector

part number	description	price
0016-009	1/4" poly x 1/8" mpt	
0016-010	1/4" poly x 1/4" mpt	
0016-011	3/8" poly x 1/8" mpt	
0016-012	3/8" poly x 1/4" mpt	



connector

part number	description	price
0016-016	1/2" mpt x 1/2" barb elbow	
0016-017	3/8" x 3/8" poly union	
0016-018	3/8" x 1/4" poly reducing union	



connector

part number	description	price
0016-020	1/4" x 1/4" poly bulkhead union	

fittings



connector

part number	description	price
0016-027	1/4" poly x 1/4" mpt elbow	



connector

part number	description	price
0016-208	1/4" poly x 10-32 connector	



connector

part number	description	price
0016-014	1/4" poly tee	
0016-015	3/8" poly tee	
0016-019	3/8" x 3/8" x 1/4" poly tee	



connector

part number	description	price
0016-005	3/8" poly x 1/4" fpt	
0016-006	1/4" poly x 1/8" fpt	
0016-007	1/4" poly x 1/4" mpt	



connector

For use with manual shut-off valve (p/n 1100-080, see page 47).

part number	description	price
0016-299	1/4" poly x 3/8" copper compression	



connector

Includes poly sleeves.

part number	description	price
0016-122	1/4" poly nut (package 10)	
0016-139	3/8" poly nut (package 10)	



connector

part number	description	price
0016-123	1/4" poly sleeve (package 10)	
0016-124	3/8" poly sleeve (package 10)	

10



fittings



connector

part number	description	price
0016-273	1/4" poly plug (package 5)	
0016-274	3/8" poly plug (package 5)	



quick disconnect

Plated without shutoff.

part number	description	price
0016-040	1/4" outside diameter male x 1/4" poly	
0016-045	3/8" outside diameter male x 3/8" poly	



quick disconnect

Plated with shutoff.

part number	description	price
0016-041	1/4" outside diameter male x 1/4" poly, water	
0016-046	3/8" outside diameter male x 3/8" poly, air	



quick disconnect

Plated with shutoff.

part number	description	price
0016-042	1/4" inside diameter female x 1/4" poly, water	
0016-048	3/8" inside diameter female x 3/8" poly, air	



quick disconnect

Plated with shutoff.

part number	description	price
0016-043	1/4" inside diameter female x 1/8" mpt, water	
0016-044-BLK	1/4" inside diameter female x 1/4" mpt, water	
0016-047	3/8" inside diameter female x 1/4" mpt, air	



quick disconnect

Plated with shutoff.

part number	description	price
0016-290	1/4" inside diameter female x 1/8" mpt, water	

tubing

handpiece tubing, asepsis style (vinyl)



4-hole, straight with midwest terminal.

part number	description	price
0017-810	sterling	
0017-796	gray	
0017-398	black	
0017-510	shadow	

handpiece tubing, asepsis style (vinyl)



4-hole, bulk, per foot.

part number	description	price
0017-960	sterling	
0017-746	gray	
0017-340	black	
0017-560	shadow	

handpiece tubing, asepsis style (polyurethane)



3-hole, coiled without terminal.

part number	description	price
0017-959	sterling	
0017-729	gray	
0017-369	black	
0017-559	shadow	

handpiece tubing, asepsis style (polyurethane)



3-hole, coiled with borden terminal:

part number	description	price
0017-724	sterling	
0017-730	gray	
0017-370	black	
0017-524	shadow	

3-hole, coiled with midwest terminal:

part number	description	price
0017-725	sterling	
0017-731	gray	
0017-371	black	
0017-525	shadow	



tubing



handpiece tubing, asepsis style (polyurethane)

part number	description	price
-------------	-------------	-------

3-hole, straight with borden terminal:

0017-739	sterling	
-----------------	----------	--

0017-732	gray	
-----------------	------	--

0017-379	black	
-----------------	-------	--

0017-579	shadow	
-----------------	--------	--

3-hole, straight with midwest terminal:

0017-740	sterling	
-----------------	----------	--

0017-733	gray	
-----------------	------	--

0017-380	black	
-----------------	-------	--

0017-540	shadow	
-----------------	--------	--



handpiece tubing, asepsis style (polyurethane)

3-hole, bulk, per foot.

part number	description	price
-------------	-------------	-------

0017-958	sterling	
-----------------	----------	--

0017-734	gray	
-----------------	------	--

0017-381	black	
-----------------	-------	--

0017-558	shadow	
-----------------	--------	--



handpiece tubing, asepsis style (polyurethane)

4-hole, coiled without terminal.

part number	description	price
-------------	-------------	-------

0017-964	sterling	
-----------------	----------	--

0017-748	gray	
-----------------	------	--

0017-339	black	
-----------------	-------	--

0017-564	shadow	
-----------------	--------	--



handpiece tubing, asepsis style (polyurethane)

4-hole, coiled with midwest terminal.

part number	description	price
-------------	-------------	-------

0017-727	sterling	
-----------------	----------	--

0017 797	gray	
-----------------	------	--

0017-399	black	
-----------------	-------	--

0017-527	shadow	
-----------------	--------	--



tubing



handpiece tubing (polyurethane)

1-hole, bulk, per foot.

part number	description	price
0017-187	gray	



handpiece tubing (polyurethane)

Large twin hole, straight with midwest terminal.

part number	description	price
0017-289	sterling	
0017-287	gray	



handpiece tubing (polyurethane)

Large twin hole bulk, per foot.

part number	description	price
0017-957	sterling	
0017-297	gray	



handpiece tubing (polyurethane)

2-hole, coiled without terminal.

part number	description	price
0017-749	gray	



handpiece tubing (polyurethane)

2-hole, coiled with borden terminal.

part number	description	price
0017-750	gray	

tubing



handpiece tubing (polyurethane)

2-hole, straight with borden terminal.

part number	description	price
0017-751	gray	



handpiece tubing (polyurethane)

2-hole, bulk, per foot.

part number	description	price
0017-752	gray	



syringe tubing, asepsis style (polyurethane)

2-hole, straight, per foot.

part number	description	price
0017-935	sterling	
0017-735	gray	
0017-365	black	
0017-535	shadow	



syringe tubing, asepsis style (polyurethane)

3-hole, straight, per foot.

part number	description	price
0017-937	sterling	
0017-737	gray	
0017-367	black	
0017-567	shadow	

tubing



syringe tubing, asepsis style (polyurethane)

2-hole, coiled.

part number	description	price
0017-939	sterling	
0017-736	gray	
0017-366	black	
0017-539	shadow	



syringe tubing assembly, asepsis style, sheathed

2-hole, polyurethane, 7'.

part number	description	price
0017-865	sterling	
0017-765	gray	
0017-368	black	
0017-565	shadow	



heated syringe assembly, asepsis style

2-hole, straight, 7'.

part number	description	price
0017-122	sterling	
0017-118	gray	
0017-113	black	
0017-522	shadow	



heated syringe assembly, asepsis style

2-hole, coiled, 7'.

part number	description	price
0017-108	sterling	
0017-107	gray	
0017-103	black	
0017-508	shadow	

tubing



foot control tubing

part number	description	price
-------------	-------------	-------

2-hole, vinyl:

0017-945	sterling	
0017-702	gray	
0017-345	black	
0017-545	shadow	

2-hole, polyurethane:

0017-944	sterling	
0017-712	gray	
0017-344	black	
0017-544	shadow	



foot control tubing

3-hole polyurethane, per foot.

part number	description	price
-------------	-------------	-------

0017-713	gray	
-----------------	------	--



foot control tubing

part number	description	price
-------------	-------------	-------

4-hole, vinyl, per foot:

0017-946	sterling	
0017-704	gray	
0017-346	black	
0017-546	shadow	

4-hole, polyurethane, per foot:

0017-714	gray	
-----------------	------	--



foot control tubing, asepsis style (polyurethane)

4-hole, per foot.

part number	description	price
-------------	-------------	-------

0017-715	gray	
-----------------	------	--

tubing

foot control tubing, asepsis style (polyurethane)



5-hole, polyurethane, per foot.

part number	description	price
0017-705	gray	

interconnect tubing

part number	description	price
-------------	-------------	-------

$\frac{1}{16}$ " inside diameter, per foot:



0017-101	orange
0017-201	yellow
0017-301	black
0017-401	blue
0017-601	clear
0017-701	gray
0017-801	green
0017-901	red

$\frac{1}{8}$ " inside diameter, per foot:



0017-190	orange
0017-290	yellow
0017-390	black
0017-490	blue
0017-690	clear
0017-790	gray
0017-890	green
0017-990	red

tubing

smooth vacuum tubing



part number	description	price
$\frac{1}{2}$ " inside diameter:		
0017-933	sterling	
0017-908	gray	
0017-333	black	
0017-533	shadow	



$\frac{5}{8}$ " inside diameter:		
0017-934	sterling	
0017-910	gray	
0017-334	black	
0017-534	shadow	

$\frac{1}{2}$ " hex inside diameter:		
0017-932	sterling	
0017-912	gray	
0017-332	black	
0017-532	shadow	

$\frac{5}{8}$ " hex inside diameter:		
0017-929	sterling	
0017-913	gray	
0017-329	black	
0017-529	shadow	

11



saliva ejector tubing (vinyl)

$\frac{3}{8}$ " outside diameter:		
part number	description	price
0017-930	sterling	
0017-794	gray	
0017-330	black	
0017-530	shadow	

tubing



spiral flex tubing

Sold per foot.

part number	description	price
$\frac{1}{2}$ " inside diameter, per foot:		
0017-920	sterling	
0017-780	gray	
0017-320	black	
0017-520	shadow	
$\frac{5}{8}$ " inside diameter:		
0017-921	sterling	
0017-781	gray	
0017-321	black	
0017-521	shadow	
$\frac{3}{4}$ " inside diameter:		
0017-922	sterling	
0017-782	gray	
0017-327	black	
0017-536	shadow	
1" inside diameter:		
0017-923	sterling	
0017-783	gray	
0017-323	black	
0017-523	shadow	
$1\frac{1}{4}$ " inside diameter:		
0017-924	sterling	
0017-784	gray	
0017-324	black	
0017-528	shadow	
$1\frac{1}{2}$ " inside diameter:		
0017-925	sterling	
0017-785	gray	
0017-325	black	
0017-528	shadow	
$1\frac{3}{4}$ " inside diameter:		
0017-926	sterling	
0017-726	gray	
0017-326	black	
0017-526	shadow	
2" inside diameter:		
0017-787	gray	

tubing



braided tubing

part number	description	price
0017-100	1/4" outside diameter, gray	
0017-792	5/16" outside diameter, gray	
0017-793	3/8" outside diameter, gray	

braided tubing bundle



Features:

- Braided fabric housing
- 1/4" outside diameter inter connect (1)
- 1/8" outside diameter interconnect (5)
- 24 awg low voltage wire (2)

part number	description	price
0017-817	umbilical, gray	

braided tubing bundle



Features:

- Braided fabric housing
- 1/4" outside diameter interconnect (1)
- 1/8" outside diameter interconnect (5)

part number	description	price
0017-807	gray	

braided tubing bundle



Features:

- Braided fabric housing
- 1/8" outside diameter interconnect (8)

part number	description	price
0017-827	gray	

custom umbilical tubing



Custom umbilicals built to your specifications. Contact factory for quotation.

Note: See page 43 for umbilical cuffs.

service items

terminal and nuts

For HP tubing.



3-hole borden

part number	description	price
0016-201	s.s. nut with metal terminal	
0016-090	s.s. nut with plastic	
0016-145	metal terminal only	
0016-130	plastic terminal only	
0016-116	s.s. nut only	

terminal and nuts

For HP tubing.



4-hole midwest

part number	description	price
0016-202	s.s. nut with metal terminal	
0016-091	s.s. nut with plastic	
0016-146	metal terminal only	
0016-133	plastic terminal only	
0016-118	s.s. nut only	

autoclavable handpiece gasket (325°F)



borden



midwest, 4 hole



midwest 5 hole

part number	description	price
0022-007	midwest, 4 hole	
0022-008	borden, $\frac{2}{3}$ hole	
0022-009	midwest, 5 hole	

service items



0016-311

handpiece tubing terminal and nuts

part number	description	price
-------------	-------------	-------

0016-311	3-hole HP to 4-hole HP terminal and nut	
-----------------	--	--



0016-312

0016-312	4-hole HP to 2-hole HP terminal and nut	
-----------------	--	--



0016-313

0016-313	4-hole HP to 3-hole HP terminal and nut	
-----------------	--	--

sleeve installation tools

Used for installation of tubing sleeve clamps.



part number	description	price
-------------	-------------	-------

0010-048	standard 1/8" sleeve tool	
-----------------	---------------------------	--

0010-049	standard 1/4" sleeve tool	
-----------------	---------------------------	--

service tool

Adapts to 4-hole Midwest HP tubing to provide instant air source for cleaning of plugged tubings or testing of components.

Features:

- Nickel plated brass construction
- Universal barb fits 1/8" or 1/4" outside diameter tubing



part number	description	price
-------------	-------------	-------

0008-021	service tool	
-----------------	--------------	--

0016-210	replacement barb	
-----------------	------------------	--

service items

handpiece pressure gauge

Used to adjust drive air pressure at HP.

Features:

- Connects to 6-pin, 5-hole and 4-hole HP tubing assembly
- 0-160 PSI gauge
- Nickel plated brass construction



part number	description	price
0008-030	HP pressure gauge	

test gauge

Used for checking air or water pressures when troubleshooting.

Features:

- 0-160 PSI gauge
- Stepped barb, 10-32 x $\frac{1}{16}$ & $\frac{1}{8}$



part number	description	price
0008-045	test gauge	

handpiece water coolant non-retraction gauge

Used to determine if HP water coolant system meets standard for non-retraction.



part number	description	price
0008-040	water coolant non-retraction gauge	

o-ring lubricant

Used as lubricant for o-rings.



part number	description	price
0030-014	o-ring lubricant, water resistant, 2 oz.	

service items



LOCTITE® 242 threadlocker

Used to lock, seat, and protect threads. Hardens in 30-60 seconds, reaches full strength in 24 hours. Parts can be disassembled.

part number	description	price
0030-017	LOCTITE® 242 threadlocker	

forest allen wrench kit

Includes 1 each and holder:



- .05
- 1/16"
- 5/64"
- 3/32"
- 7/64"
- 1/8"
- 9/64"
- 5/32"
- 3/16"
- 7/32"
- 1/4"
- 5/16"

part number	description	price
0008-019	allen wrench kit	



service kit

Contains 120 o-rings in 12 popular sizes from 7/64" outside diameter to 1/2" outside diameter.

part number	description	price
0015-100	o-ring kit	

autoclavable o-rings

Viton.



0015-408




0015-413

part number	description	price
0015-408	for SE valve (#8 o-ring) (package 10)	
0015-413	for CV valve (#13 o-ring) (package 10)	

service items

o-ring kit refills

Package 10.



part number	description		price
0015-101	#1	1/32"	
0015-102	#2	3/64"	
0015-103	#3	1/16"	
0015-104	#4	5/64"	
0015-105	#5	3/32"	
0015-106	#6	1/8"	
0015-107	#7	5/32"	
0015-108	#8	3/16"	
0015-109	#9	7/32"	
0015-110	#10	1/4"	
0015-111	#11	5/16"	
0015-112	#12	3/8"	
0015-213	#13	7/16"	

barb fitting kit

Contains eight different fittings.

Features:

- 10-32 4-way cross female, (package 5)
- 10-32 female connector (package 5)
- 10-32 male coupler (package 5)
- 10-32 x 1/8" barbs and washers (package 10)
- 10-32 x 1/16" barbs and washers (package 10)
- 10-32 hex plugs and washers (package 10)
- 1/8" sleeve clamps (package 10)
- 1/16" sleeve clamps (package 10)



part number	description		price
0016-200	barb fitting kit		

water relay valve II kit

Includes o-rings, spring, cover and cap.




part number	description		price
1100-298	repair kit, gray cap		

service items

field repair kit

Contains 12 different repair kits.



Features:

1105-100	• FC drive air kit
1105-103	• FC chip blower valve kit (aluminum cap)
1105-104	• FC wet/dry toggle valve kit
1105-107	• FC signal relay kit
0022-021	• HP control, duo-black diameterphragm (1 each)
0022-020	• HP control, tri-block (1 each)
1100-298	• Water relay valve kit
1100-199	• Master shut-off kit
0022-024	• HP control, quad-block diameterphragm (1 each)
1400-049	• Flush cup kit
0016-268	• 6-32 barb and washer kit
0016-269	• 8-32 barb and washer kit

part number	description	price
0008-050	field repair kit	

master shut-off valve repair kit

Includes o-rings and return spring.



part number	description	price
1100-199	master shut-off valve kit	

automatic handpiece control block diameterphragm only



part number	description	price
0022-021	duo-block	
0022-020	tri-block	
0022-024	quad-block	

water flush cup kit


Includes valve stem, o-rings, spring and actuator ball.



part number	description	price
1400-049	water flush cup kit	

service items

foot control service kits (drive air valve kit)



Includes plunger, poppet, o-rings and springs.

part number	description	price
1105-073	poppet with insert only	
1105-100	complete kit	

chip blower valve kit



Includes o-rings, keeper pin, spring and cap.

part number	description	price
1105-103	for valve with aluminum cap	

wet/dry valve kit, toggle



Includes plunger, o-rings and spring.

part number	description	price
1105-104	wet/dry kit, toggle	

signal replay repair kit



Includes plunger, o-rings, spring, cover and screws.

part number	description	price
1105-107	signal relay repair kit	

index

Part #	Page #	Price	Part #	Page #	Price	Part #	Page #	Price
0008-019	80		0010-269	30		0010-868-SH	30	
0008-021	78		0010-270	30		0010-885-BLK	31	
0008-030	79		0010-275	31		0010-885-G	31	
0008-040	79		0010-276	32		0010-885-SH	31	
0008-045	79		0010-285	31		0010-886-BLK	31	
0008-050	82		0010-343	54		0010-886-G	31	
0010-012	39		0010-344	54		0010-886-SH	31	
0010-017	28		0010-364	48		0010-887-BLK	31	
0010-018	32		0010-447	84		0010-887-G	31	
0010-019	32		0010-449	84		0010-887-SH	31	
0010-023	37		0010-464	48		0010-888-BLK	31	
0010-025	31		0010-591	48		0010-888-G	31	
0010-025	46		0010-593	48		0010-888-SH	31	
0010-044	23		0010-625	26		0011-011	48	
0010-048	78		0010-627	26		0011-011	48	
0010-049	78		0010-835-BLK	28		0013-016	39	
0010-091	9		0010-835-G	28		0014-004	14	
0010-094	28		0010-835-SH	28		0014-064	53	
0010-118	29		0010-836-BLK	28		0014-065	84	
0010-119	29		0010-836-G	28		0014-066	53	
0010-140	55		0010-836-SH	28		0014-067	53	
0010-141	56		0010-837-BLK	32		0014-085	52	
0010-164	54		0010-837-G	32		0014-088	52	
0010-165	54		0010-837-SH	32		0014-099	3	
0010-191	9		0010-838-BLK	32		0014-107	3	
0010-195	56		0010-838-G	32		0014-115	52	
0010-200	9		0010-838-SH	32		0014-117	3	
0010-200-B	9		0010-846-BLK	28		0014-122	53	
0010-201	39		0010-846-G	28		0014-143	2	
0010-202	39		0010-846-SH	28		0014-148	2	
0010-208	37		0010-847-BLK	28		0014-182	52	
0010-209	37		0010-847-G	28		0014-254	3	
0010-212	32		0010-847-SH	28		0014-300	2	
0010-214	32		0010-856-BLK	29		0014-301	3	
0010-216	32		0010-856-G	29		0014-314	3	
0010-219	32		0010-856-SH	29		0014-315	3	
0010-227	31		0010-857-BLK	29		0014-316	3	
0010-231	32		0010-857-G	29		0015-033	24	
0010-234	28		0010-857-SH	29		0015-033	24	
0010-235	28		0010-858-BLK	30		0015-089	16	
0010-245	28		0010-858-G	30		0015-096	14	
0010-254	29		0010-858-SH	30		0015-099	16	
0010-255	29		0010-859-BLK	30		0015-100	80	
0010-260	30		0010-859-G	30		0015-101	81	
0010-264	29		0010-859-SH	30		0015-102	81	
0010-265	29		0010-868-BLK	30		0015-103	81	
0010-266	59		0010-868-G	30		0015-104	81	



index

Part #	Page #	Price	Part #	Page #	Price	Part #	Page #	Price
0015-105	81		0016-083	64		0016-202	77	
0015-106	81		0016-087	63		0016-208	65	
0015-107	81		0016-090	77		0016-210	60	
0015-108	81		0016-091	77		0016-210	78	
0015-109	81		0016-092	16		0016-214	25	
0015-110	81		0016-100	62		0016-238	59	
0015-111	81		0016-101	60		0016-239	59	
0015-112	81		0016-102	60		0016-252	60	
0015-213	81		0016-103	58		0016-260	60	
0015-408	80		0016-104	58		0016-260	82	
0015-413	80		0016-109	61		0016-262	59	
0016-005	65		0016-114	25		0016-264	59	
0016-006	65		0016-116	77		0016-265	59	
0016-007	65		0016-118	77		0016-269	60	
0016-009	64		0016-122	65		0016-269	82	
0016-010	64		0016-123	65		0016-273	66	
0016-011	64		0016-124	65		0016-274	66	
0016-012	64		0016-127	63		0016-275	61	
0016-014	65		0016-129	62		0016-278	58	
0016-015	65		0016-130	77		0016-290	66	
0016-016	64		0016-133	77		0016-296	60	
0016-017	64		0016-139	65		0016-297	60	
0016-018	64		0016-140	62		0016-299	65	
0016-019	65		0016-144	63		0016-300	63	
0016-020	64		0016-145	77		0016-311	78	
0016-027	65		0016-146	77		0016-312	78	
0016-029	63		0016-151	61		0016-313	78	
0016-036	63		0016-152	61		0016-433	61	
0016-040	66		0016-153	62		0016-434	61	
0016-041	66		0016-154	62		0016-436	63	
0016-042	66		0016-155	61		0017-100	76	
0016-043	66		0016-156	61		0017-101	73	
0016-044-BLK	66		0016-157	61		0017-103	71	
0016-045	66		0016-158	61		0017-107	71	
0016-046	66		0016-159	61		0017-108	71	
0016-047	66		0016-165	59		0017-113	71	
0016-048	66		0016-167	58		0017-118	71	
0016-049	63		0016-168	58		0017-122	71	
0016-060	62		0016-169	58		0017-187	69	
0016-071	63		0016-180	64		0017-190	73	
0016-073	59		0016-183	64		0017-201	73	
0016-074	59		0016-184	64		0017-287	69	
0016-075	63		0016-185	64		0017-289	69	
0016-077	63		0016-198	64		0017-290	73	
0016-079	62		0016-199	64		0017-297	69	
0016-080	62		0016-200	81		0017-301	73	
0016-081	62		0016-201	77		0017-320	75	



index

Part #	Page #	Price	Part #	Page #	Price	Part #	Page #	Price
0017-321	75		0017-533	74		0017-729	67	
0017-323	75		0017-534	74		0017-730	67	
0017-324	75		0017-535	70		0017-731	67	
0017-325	75		0017-536	75		0017-732	68	
0017-326	75		0017-539	71		0017-733	68	
0017-327	75		0017-540	68		0017-734	68	
0017-329	74		0017-544	72		0017-735	70	
0017-330	74		0017-545	72		0017-736	71	
0017-332	74		0017-546	72		0017-737	70	
0017-333	74		0017-558	68		0017-739	68	
0017-334	74		0017-559	67		0017-740	68	
0017-339	68		0017-560	67		0017-746	67	
0017-340	67		0017-564	68		0017-748	68	
0017-344	72		0017-565	71		0017-749	69	
0017-345	72		0017-567	70		0017-750	69	
0017-346	72		0017-579	68		0017-751	70	
0017-365	70		0017-584	50		0017-752	70	
0017-366	71		0017-592	50		0017-765	71	
0017-367	70		0017-599	50		0017-780	75	
0017-368	71		0017-600	52		0017-781	75	
0017-369	67		0017-601	73		0017-782	75	
0017-370	67		0017-603	52		0017-783	75	
0017-371	67		0017-675	50		0017-784	75	
0017-379	68		0017-681	50		0017-785	75	
0017-380	68		0017-682	50		0017-787	75	
0017-381	68		0017-683	50		0017-790	73	
0017-390	73		0017-684	50		0017-792	76	
0017-398	67		0017-686	50		0017-793	76	
0017-399	68		0017-688	50		0017-794	74	
0017-401	73		0017-689	50		0017-796	67	
0017-490	73		0017-690	73		0017-797	68	
0017-500	50		0017-692	50		0017-801	73	
0017-508	71		0017-694	50		0017-807	76	
0017-510	67		0017-699	50		0017-810	67	
0017-520	75		0017-700	50		0017-817	76	
0017-521	75		0017-701	73		0017-827	76	
0017-522	71		0017-702	72		0017-865	71	
0017-523	75		0017-704	72		0017-890	73	
0017-524	67		0017-705	73		0017-900	52	
0017-525	67		0017-712	72		0017-901	73	
0017-526	75		0017-713	72		0017-903	52	
0017-527	68		0017-714	72		0017-908	74	
0017-528	75		0017-715	72		0017-910	74	
0017-528	75		0017-724	67		0017-912	74	
0017-529	74		0017-725	67		0017-913	74	
0017-530	74		0017-726	75		0017-920	75	
0017-532	74		0017-727	68		0017-921	75	



index

Part #	Page #	Price	Part #	Page #	Price	Part #	Page #	Price
0017-922	75		0020-100	25		0024-126-BLK	34	
0017-923	75		0020-101	25		0024-127	34	
0017-924	75		0020-101	23		0024-127-BLK	34	
0017-925	75		0022-007	77		0024-128	34	
0017-926	75		0022-008	77		0024-128-BLK	34	
0017-929	74		0022-009	77		0024-133	8	
0017-930	74		0022-020	41		0024-134	8	
0017-932	74		0022-020	82		0024-135	8	
0017-933	74		0022-020	82		0024-136	8	
0017-934	74		0022-021	41		0024-137	34	
0017-935	70		0022-021	82		0024-138	34	
0017-937	70		0022-021	82		0024-141	35	
0017-939	71		0022-023	62		0024-142	35	
0017-944	72		0022-024	41		0024-143	35	
0017-945	72		0022-024	82		0024-144	35	
0017-946	72		0022-024	82		0024-145	35	
0017-957	69		0022-025	62		0024-146	35	
0017-958	68		0023-165	36		0024-147	35	
0017-959	67		0023-166	36		0024-148	35	
0017-960	67		0023-170	41		0024-205	35	
0017-964	68		0023-171	41		0024-215	35	
0017-982	51		0023-282	41		0030-014	79	
0017-984	51		0023-406	10		0030-017	80	
0017-985	51		0023-436	25		1100-039	40	
0017-990	73		0023-665	19		1100-059	42	
0017-994	51		0024-066	7		1100-069	47	
0019-040	58		0024-081	36		1100-077	47	
0019-056	58		0024-082	36		1100-078	47	
0020-001	45		0024-083	36		1100-079	47	
0020-018	46		0024-084	36		1100-080	47	
0020-020	45		0024-087	36		1100-100	38	
0020-021	42		0024-088	36		1100-102	37	
0020-023	42		0024-090	36		1100-103	37	
0020-024	42		0024-091	36		1100-104	41	
0020-030	42		0024-113	33		1100-105	40	
0020-033	42		0024-113-BLK	33		1100-106	85	
0020-039	46		0024-114	33		1100-107	85	
0020-040	46		0024-114-BLK	33		1100-108	22	
0020-041	46		0024-115	33		1100-111	85	
0020-043	46		0024-116	33		1100-115	84	
0020-050	46		0024-117	33		1100-117	38	
0020-068	43		0024-118	33		1100-124	37	
0020-069	43		0024-123	7		1100-125	37	
0020-078	43		0024-124	7		1100-130	38	
0020-079	43		0024-125	34		1100-136	39	
0020-100	24		0024-125-BLK	34		1100-137	39	
0020-100	24		0024-126	34		1100-138	39	



index

Part #	Page #	Price	Part #	Page #	Price	Part #	Page #	Price
1100-140	40		1105-006-BLK	19		1107-057	16	
1100-142	38		1105-006-G	19		1107-062	14	
1100-144	38		1105-073	83		1107-063	15	
1100-145	41		1105-079	19		1107-070	15	
1100-170	41		1105-080	19		1107-086	11	
1100-175	41		1105-082	19		1107-096	12	
1100-181	38		1105-100	82		1107-185	11	
1100-199	82		1105-100	83		1107-186	11	
1100-199	82		1105-103	82		1107-195	11	
1100-230-BLK	85		1105-103	83		1107-196	11	
1100-230-WHT	85		1105-104	82		1107-196-BLK	11	
1100-231-BLK	85		1105-104	83		1107-196-SH	11	
1100-231-WHT	85		1105-107	82		1107-203	11	
1100-247	23		1105-107	83		1107-207	11	
1100-248	23		1105-113	17		1107-208	11	
1100-258	22		1105-114	17		1107-213	11	
1100-261	84		1105-123	17		1107-217	11	
1100-262	85		1105-124	17		1107-218	11	
1100-264	85		1105-133	17		1107-223	11	
1100-267	85		1105-134	17		1107-227	11	
1100-267-WHT	85		1105-143	17		1107-228	11	
1100-270	41		1105-144	17		1107-233	11	
1100-281	48		1105-173	18		1107-237	11	
1100-283	48		1105-174	18		1107-238	11	
1100-285	33		1105-183	18		1107-243	12	
1100-298	81		1105-184	18		1107-244	12	
1100-298	82		1105-223	18		1107-245	12	
1100-308	22		1105-224	18		1107-246	12	
1100-311	23		1105-233	18		1107-286	11	
1100-311-BLK	23		1105-234	18		1107-295	12	
1100-341	22		1105-356	17		1107-296	12	
1100-341-BLK	22		1105-357	17		1107-330	15	
1100-342	22		1105-358	17		1107-333	15	
1100-343	22		1105-359	17		1107-403	14	
1100-350	32		1105-360	18		1107-407	14	
1100-351	32		1105-361	18		1107-408	14	
1100-365	42		1105-362	18		1107-413	14	
1100-371	38		1105-363	18		1107-417	14	
1100-372	38		1105-576	17		1107-418	13	
1100-402	37		1105-577	17		1107-433	14	
1100-403	37		1105-578	17		1107-437	14	
1100-404	41		1105-579	17		1107-443	14	
1100-424	37		1105-580	18		1107-444	14	
1100-425	37		1105-581	18		1107-447	14	
1100-441	38		1105-582	18		1107-448	14	
1105-005-BLK	19		1105-583	18		1107-463	14	
1105-005-G	19		1107-020	15		1107-467	14	



index

Part #	Page #	Price	Part #	Page #	Price	Part #	Page #	Price
1107-468	14		1110-578	57		1115-007-G	45	
1107-691	15		1110-582	57		1115-007-SH	45	
1107-693	15		1110-588	55		1115-009	45	
1107-694	16		1110-592	55		1115-009-BLK	45	
1107-695	14		1110-603	9		1115-009-G	45	
1107-721	15		1110-604	9		1115-009-SH	45	
1107-770	15		1110-606	9		1115-020	44	
1107-906	14		1110-608	9		1115-022	44	
1107-908	11		1110-643	57		1115-028	44	
1107-916	14		1110-649	2		1115-088	43	
1107-918	11		1110-652	8		1115-089	43	
1107-928	11		1110-653	8		1115-090	43	
1107-938	11		1110-654	8		1115-091	43	
1107-942	14		1110-655	8		1115-094-F	57	
1107-946	12		1110-656	2		1115-094-PC	57	
1107-947	14		1110-657	2		1115-098-F	57	
1107-969	14		1110-665	57		1115-098-PC	57	
1107-986	11		1110-744	8		1115-300	44	
1107-995	12		1110-760	5		1300-030	4	
1110-130	84		1110-761	5		1400-008	24	
1110-130-G	84		1110-762	5		1400-010	26	
1110-131	84		1110-763	5		1400-011	24	
1110-350	54		1110-764	5		1400-017	24	
1110-351	54		1110-765	5		1400-017-BLK	24	
1110-376	6		1110-768	5		1400-020	26	
1110-377	10		1110-769	5		1400-049	82	
1110-378	10		1110-770	5		1400-049	82	
1110-380	10		1110-771	5		1400-052	27	
1110-386	6		1110-808	2		1400-055	24	
1110-446	7		1110-815	3		1400-074	21	
1110-447	7		1110-816	2		1400-077	21	
1110-450	7		1110-820	56		1400-078	21	
1110-459	7		1110-822	56		1400-082	21	
1110-464	6		1110-898	5		1400-085	26	
1110-480	55		1110-899	5		1400-090	24	
1110-481	55		1110-970	2		1400-093	24	
1110-503	57		1110-973	84		1400-111	24	
1110-504	57		1115-001	44		1400-116	26	
1110-510	57		1115-003	44		1400-117	26	
1110-516	56		1115-004	44		1400-128	24	
1110-517	56		1115-005	44		1400-128	24	
1110-518	56		1115-005-BLK	44		1400-128-BLK	24	
1110-525	6		1115-005-G	44		1400-128-BLK	24	
1110-526	6		1115-005-SH	44		1400-129	24	
1110-527	6		1115-006	44		1400-129	24	
1110-536	6		1115-007	45		1400-129-BLK	24	
1110-537	6		1115-007-BLK	45		1400-129-BLK	24	



index

Part #	Page #	Price	Part #	Page #	Price	Part #	Page #	Price
1400-133	23		1400-805	20		1700-200	48	
1400-133-BLK	23		1400-806	20		1700-201	48	
1400-136	23		1400-807	20		1700-220	48	
1400-155	24		1400-852	27		1700-221	48	
1400-163	21		1400-861	27		1700-226	48	
1400-164	21		1400-866	21		1700-518	49	
1400-167	21		1400-870	27		1700-522	49	
1400-173	27		1400-871	21		1700-548	49	
1400-174	27		1400-872	21		1700-549	49	
1400-177	21		1400-873	21		1700-600	51	
1400-183	27		1400-874	21		1700-601	51	
1400-184	27		1400-876	21		1700-603	51	
1400-190	24		1400-877	21		1700-604	51	
1400-191	25		1400-878	21		1700-800	51	
1400-192	25		1400-885	26		1700-801	51	
1400-193	24		1400-900	20		1700-902	51	
1400-252	27		1400-901	20		1700-904	51	
1400-263	21		1400-902	20		2000-010	4	
1400-264	21		1400-903	20		2000-011	4	
1400-267	21		1400-904	20		4003-097	54	
1400-274	21		1400-905	20		4003-135	55	
1400-276	21		1400-906	20				
1400-277	21		1400-907	20				
1400-278	21		1400-952	26				
1400-282	21		1400-952	27				
1400-285	26		1400-961	27				
1400-500	20		1400-970	27				
1400-501	20		1400-971	21				
1400-502	20		1400-972	21				
1400-503	20		1400-973	21				
1400-504	20		1400-974	21				
1400-505	20		1400-975	21				
1400-506	20		1400-976	21				
1400-507	20		1400-977	21				
1400-600	20		1400-978	21				
1400-601	20		1700-002	49				
1400-602	20		1700-003	49				
1400-603	20		1700-008	49				
1400-604	20		1700-010	49				
1400-605	20		1700-018	49				
1400-606	20		1700-022	49				
1400-607	20		1700-036	49				
1400-800	20		1700-037	49				
1400-801	20		1700-038	49				
1400-802	20		1700-039	49				
1400-803	20		1700-048	49				
1400-804	20		1700-049	49				

