




V&KL dental rotor

Catalogue 2025



Dental Spare Parts



Quality Products, Fast Delivery, Technical Support ...

... the Advantages That Make Your Job Easier!



V&KL dental rotor is a trading company supplying dental replacement rotors since 1995 to all European countries and worldwide. In almost 30 years of our successful existence, the product assortment has been extended with various replacement parts which relate to and are necessary for the repairing of dental handpieces.

Our aim is to be the best specialised store in supporting dental handpiece repair-shops. From the very beginning we try to supply customers with highest possible quality aftermarket products at reasonable prices. A large stock of products allows us to ship out orders the very same day we receive an order so our customers can continuously work with minimal stock or even without it.

Working with us you won't only get a great support with products we offer but also with a willingness to hear and fulfil your special needs. We like to share various technical specifications and our knowledge about dental rotary instruments.

We are flexible at negotiating terms of business relations, extended payment terms and prices because we like to be a part of your success. Take a look at the catalogue and see our offer on the following pages. We are sure you will find something resourceful for your business!

Technical Features of V&KL Highspeed Cartridges

Rotors are supplied with a guide bushing made of **HARDENED STEEL**. The guide bushing is **SEPARATE** from a chuck just like the original rotors. Thanks to this advanced manufacturing technique, V&KL rotors perform like original rotors, costing less! Compared to other aftermarket replacement rotors, V&KL rotors allow more precise work. It is a reliable product that will surely gain the trust of your customers.

- ✓ **PRECISE WORK**
- ✓ **NO VIBRATION**
- ✓ **LONGER LIFETIME**



Every rotor is marked with a serial No.


Balancing: 150µg
Accuracy of bur in rotation: MAX 0,01mm

Shaft assembly: Laser Welded
Shaft hardness (Ø 3.175±0.001): 52 HRC

4 - lines clamping mechanism;
Bur extracting force: 40-50N;
Bur guide hardness: 60HRC;
Opening dimension: Ø 0/1,6+0,005;
The clamping accuracy after 12 months of exploitation: 0,005 - 0,010

Push button opening force: 22N - 25N

High speed Ceramic Bearings with torlon cage,
appropriate for autoclave, with Radial play or
Angular contact configuration



How to read this catalogue

Legend of Symbols

Products presented in this catalogue are divided into several groups and categories. Each category has its own symbol.

- Turbines, Rotors, Cartridges
- Back Caps
- Back Cap Wrenches
- O-Rings
- Wave Washers
- Bulbs
- Bearings
- Electrical Rotors
- Carbon Brushes
- Fiber Optic
- Intermediate Shafts
- Gears

Legend of Colours

For easier recognition of products among various brands within groups and categories, each brand has its own colour.

- BA international
- Bien Air
- Castellini
- Fona
- Kavo
- Midwest
- Morita
- NSK
- Shick
- Siemens
- Sirona
- Star
- W&H


Trademarks

All brands and models that are mentioned in this catalogue and any V&KL marketing material are registered trade names. No implication is made that any of our products are original products or that any of the referenced companies endorse our products. This clause applies to the following brands: KaVo, Star, Midwest, NSK, W&H, Sirona, Siemens, Morita, Castellini, Bien Air, Mk-Dent, B.A. International, Fona, Shick.


Head of Lists and Tables Explained




Rotors / Turbines / Cartridges Explained



Non-turbine Bearings Explained



Turbine Bearings Explained



Glossary

Preferred over OEM	the term denotes the reference number from Original Equipment Manufacturer and suggests products are interchangeable but never the same.
SKU	the term denotes V&KL reference number of the product.
Related products	parts that fit the product and can be ordered separately, as an individual product.
RPM	is a unit to measure the speed of bearings (Revolutions Per Minute). The unit describes how many turns are completed in one minute around fixed axis.
Spindle & Impeller	or a combo. The term refers to a cartridge without bearings. If a particular turbine model does not have a spindle & impeller option, it means that this option is not available to order due to technical reasons.
Rotor / Turbine / Cartridge	the terms are used interchangeably in this catalogue and mean a rotor package equipped with bearings and all necessary O-rings or wave washers. It is air driven.
Non-turbine Bearing	the term refers to bearings that rotate up to 250K RPM and are used for contra angle, attachment handpieces and various low speed instruments.
Turbine Bearing	the term refers bearings that rotate from 250K up to 450K RPM and are used in air driven turbine handpieces.
40K - 100K	40.000 to 100.000 RPM. A bearing that rotates at speeds up to 100,000 RPM is considered a lowspeed bearing.
150K - 250K	150.000 to 250.000 RPM. A bearing that rotates at speeds up to 250.000 RPM is regarded a highspeed bearing.
250K - 450K	250.000 to 450.00 RPM. A bearing that rotates at speeds in this range is regarded a turbine bearing and it is propelled by air.
Bearing dimension	measure in mm is presented as: inner ring x outer ring x width / unique dimension.
Smooth Bearing	a bearing's outer ring is smooth and plain without any grooves nor flanges.
Flanged Bearing	a bearing's outer ring has a flange – a flat extended rim which reinforces a bearing's strength and makes the fitting of the bearing in its position easier.
Grooved / Stepped Bearing	a bearing's outer ring is not smooth. Both terms are interchangeably used in this catalogue.
Ceramic Bearing	the term refers to the material of the bearing balls. Balls coated with a layer of ceramic (silicon nitride) have a smoother surface, which causes less noise and friction, thus extending the life of the bearing.
Integral Shield Bearing	the bearing shield is also its outer ring. These bearings have the longest lifetime, as the bearing is integrally closed and thus completely protected from debris and dirt.
Angular Bearing	have asymmetrical inner and outer ring raceways which allows an angle contact that sustains significant loads in radial and axial direction at the same time. Because they withstand these "combined loads" well they are more reliable, last longer, and reach higher speeds than radial bearings.
Radial Bearing	or deep groove ball bearing has symmetrical inner and outer ring raceways. This type of configuration makes the bearing's cage less mechanically stable and thus its lifetime is shorter in comparison to angular bearings.
Closed/closed	bearing is closed on both sides, either with a metallic shield or rubber seal.
Closed/Open	bearing is closed on one side, either with a metallic shield or rubber seal, and open on the other side.
Open/Open	bearing is open on both sides.
Torlon	the term refers to a cage material. Torlon has improved temperature and chemical resistance in comparison to phenolic but lubricant absorption is lower, making torlon more sensitive to insufficient or lack of lubrication.
Head Gear	an assembly in the head housing of an attachment handpiece that holds a bur and is by means of gears driven by a motor.
Intermediate Shaft	a part consisting of a long shaft and two gears, which is located in the neck of an attachment handpiece and meshes with a head gear rotor.

Content

Highspeed Handpiece Parts



Turbines / Rotors / Cartridges pages 12 to 29

V&KL offers a wide range of rotors of all world-famous dental brands - from the oldest to the latest models within each brand. Rotors are available with ceramic bearings or without bearings. When buying rotors with bearings all the correspondent O-rings, wavy washers and slingers come along the package for free.

B.A. International	page 13
Bien Air	page 15
Castellini	page 16
Fona	page 16
Kavo	page 17
Midwest	page 19
Morita	page 20
NSK	page 20
Siemens	page 22
Sirona	page 23
Star	page 27
W&H	page 27

Back Caps pages 30 to 36

Back caps do not break as often as rotors but still, sometimes the real problem lies in a worn out back cap and not necessarily in a rotor's chuck. In the aim to offer a back cap for each replacement rotor in a portfolio, new back caps are constantly added to our offer.

Back Caps Wrenches pages 37 to 40

A repair work is much easier if you have the right tools – especially the wrenches. V&KL supports you with a wide range of back cap wrenches to help you disassemble various instruments.

O-rings and Wave Washers pages 41 to 43

O-Rings and wave washers provide a firm sitting of a rotor inside a handpiece's head while at the same time they help to minimise vibrations. Because O-rings are made of rubber, they become flattened sooner or later. Sometimes only the replacement of O-rings can prolong a life span of a rotor. V&KL rotor packages always come equipped with O-Rings and Wavy Washers, where applicable. Should you need them extra they are also available as a single product.

Bulbs pages 44 to 47

You can choose between LED or XENON bulbs.

Contra Angle Parts



Parts listed by an instrument brand and model pages 50 to 75

In this section the replacement parts (head gears, intermediate shafts and gears, lower shaft gears, bearings, fiber optic rods) for a particular model of a brand are gathered together.

Parts for 1:5 Bien Air models	page 51
Parts for 1:1 Bien Air models	page 52
Parts for 1:5 Kavo models	page 53
Parts for 1:1 Kavo models	page 59
Parts for 1:5 NSK models	page 62
Parts for 1:1 NSK models	page 64
Parts for 1:5 Sirona models	page 65
Parts for 1:1 Sirona models	page 67
Parts for 1:5 W&H models	page 68
Parts for 1:1 W&H models	page 72

Content

Parts listed by a category pages 76 to 80

In this section all the replacement parts available within a category are listed together.

1:5 Head Gears for Contra Angles	page 77
1:1 Head Gears for Contra Angles	page 78
1:5 Intermediate Shafts for Contra Angles	page 79
1:1 Intermediate Shafts for Contra Angles	page 80
Light Guides / Fiber Optic Rods	page 80

Micro Motors Parts



Electrical Rotors and Sets pages 84 to 87

For Bien Air and Kavo Micro Motorss we offer rotors, bearing sets and carbon brushes. Rotors for Bien Air Micro Motorss are NOT equipped with bearings whereas rotors for Kavo Micro Motors come with assembled bearings.

Rotors, Bearing Sets and Carbon Brushes for Bien Air	page 85
Rotors, Bearing Sets and Carbon Brushes for Kavo	page 86

Bearings



Non-turbine Lowspeed Bearings 40K – 100K RPM pages 90 to 97

Bearings for various lowspeed instruments. The majority of bearings within this group are manufactured by NMB Minebea (Japan).

Bien Air	page 91
Kavo	page 91
NSK	page 94
Shick	page 95
Sirona	page 96
W&H	page 97

Non-turbine Highspeed Bearings 150K – 250K RPM pages 98 to 100


A wide range of bearings that come handy when repairing highspeed contra angles are included in this section. The majority of bearings within this group are manufactured by NMB Minebea (Japan) or myonic (Germany).

Bien Air	page 99
Kavo	page 99
NSK	page 100
Sirona	page 100
W&H	page 100

Turbine Bearings 250K – 450K RPM pages 101 to 102

Ball bearings for air driven highspeed rotors / turbines / cartridges are offered in Angular and Radial (Deep Groove) configuration. Angular ball bearings provide lower noise level operation than radial ball bearings and at the same time give more power to a turbine and generally prolong a life span of a rotor. Angular and radial bearings are available with Ceramic balls only. All bearings within this group are manufactured by myonic (Germany).

Highspeed Handpiece Parts





Highspeed Handpiece Parts



Turbines Rotors Cartridges

	Ultimate 650/660 (new)	SKU: 30.076 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10330PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1)
	Ultimate 670 (old)	SKU: 30.000 Spindle & Impeller SKU: 30.006 W/ Ceramic Radial Bearings + O-Rings SKU: 30.008 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 42004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33109 - O-Ring (6,36 × 0,82)
	Ultimate 670 (new)	SKU: 30.076 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10330PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1)
	Ultimate 675/685/635FM	SKU: 30.010 Spindle & Impeller SKU: 30.016 W/ Ceramic Radial Bearings + O-Rings SKU: 30.018 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
	Ultimate 678/758	SKU: 30.020 Spindle & Impeller SKU: 30.026 W/ Ceramic Radial Bearings + O-Rings SKU: 30.028 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
	Ultimate 688 after 2015	SKU: 30.078 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22113 - Back Cap SKU: 33101 - O-Ring (6 × 1)
	Ultimate 695 (old) before 2016	SKU: 30.046N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
	Ultimate 695 (new) after 2016	SKU: 30.068 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 33101 - O-Ring (6 × 1)
	Ultimate 695 FB/FM fixed base	SKU: 30.066 W/ Ceramic Radial Bearings + O-Rings + Wave Washer	Related products: SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1)



	Ultimate 755 (old)	SKU: 30.030 Spindle & Impeller SKU: 30.036 W/ Ceramic Radial Bearings + O-Rings SKU: 30.038 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33101 - O-Ring (6 × 1)
	Ultimate 755 (new)	SKU: 30.068 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 33101 - O-Ring (6 × 1)
	Optima 535	SKU: 30.010 Spindle & Impeller SKU: 30.016 W/ Ceramic Radial Bearings + O-Rings SKU: 30.018 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
	Optima 525	SKU: 30.056 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer



	Black Pearl TLC/ Croma/Boralina Original P/N: 209.40.12, 1500164	SKU: 24.000 Spindle & Impeller SKU: 24.006 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 24003 - Back Cap SKU: 44004 - Back Cap Wrench SKU: 10110PT - Ceramic Radial Bearing SKU: 33100 - O-Ring (8 × 1)
	Bora/ Tornado X Original P/N: 209.34.30, 1500181	SKU: 24.010 Spindle & Impeller SKU: 24.016 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 24013 - Back Cap SKU: 44014 - Back Cap Wrench SKU: 10110PT - Ceramic Radial Bearing SKU: 33100 - O-Ring (8 × 1)
	Prestige/Prestilina Original P/N: 209.31.51, 1500485	SKU: 24.026 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 33118 - O-Ring (6,3 × 1) SKU: 33110 - O-Ring (3,5 × 1)
	TDS783/TDS890 Original P/N: 208.55.01, 1500922	SKU: 24.036 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10110PT - Ceramic Radial Bearing SKU: 33100 - O-Ring (8 × 1) SKU: 33110 - O-Ring (3,5 × 1)
	Tornado Original P/N: 1502252	SKU: 24.046 Spindle & Impeller SKU: 24.048 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 24023 - Back Cap SKU: 44044 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33118 - O-Ring (6,3 × 1) SKU: 33108 - Wave Washer
	Tornado S Original P/N: 1501969	SKU: 24.056 W/ Ceramic Radial Bearings + O-Rings SKU: 24.059 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 10810PT - Ceramic Radial Bearing SKU: 1087N - Ceramic Angular Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer
NEW 	Bora 2 Original P/N: 1503003	SKU: 24.066 W/ Ceramic Radial Bearings + O-Rings SKU: 24.069 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 24033 - Back Cap SKU: 44034 - Back Cap Wrench SKU: 10830PT - Ceramic Radia Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer



	High Power M / Clean Air 2000 / Clean Light Original P/N: L0000874	SKU: 25.000 Spindle & Impeller SKU: 25.006 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10210PT - Ceramic Radial Bearing SKU: 33118 - O-Ring (6,3 × 1)
	Titanium Gold 2 / Hi Power Light 2 Original P/N: L0001370	SKU: 25.010 Spindle & Impeller SKU: 25.016 W/ Ceramic Radial Bearings + O-Rings SKU: 25.018 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing
	Silent Power / 4L / 4LK / 2L / 2 Original P/N: L0001857	SKU: 25.020 Spindle & Impeller SKU: 25.026 W/ Ceramic Radial Bearings + O-Rings SKU: 25.028 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 25003 - Back Cap SKU: 45004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33118 - O-Ring (6,3 × 1)
	Silent Power EVO / 4L / 4LK / 2L / 2 Original P/N: 95529230	SKU: 25.030 Spindle & Impeller SKU: 25.036 W/ Ceramic Radial Bearings + O-Rings SKU: 25.038 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 25013 - Back Cap SKU: 45014 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearingw SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33118 - O-Ring (6,3 × 1) SKU: 33108 - Wave Washer



	6080MR/8060SR Original P/N: 6080MR / 8060SR	SKU: 31.000 Spindle & Impeller SKU: 31.006 W/ Ceramic Radial Bearings + O-Rings SKU: 31.008 W/ Ceramic Angular Integral Shield Bearing + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
	8080SR Original P/N: 8080SR	SKU: 31.010 Spindle & Impeller SKU: 31.016 W/ Ceramic Radial Bearings + O-Rings SKU: 31.018 W/ Ceramic Angular Integral Shield Bearing + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33101 - O-Ring (6 × 1)
	Tern Advance / Colibri Original P/N: NA	SKU: 31.028 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 33101 - O-Ring (6 × 1)



	605/615 Original P/N: 1.008.4164	SKU: 21.178 Spindle & Impeller SKU: 21.176 W/ Ceramic Radial Bearings + O-Rings SKU: 21.179 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 21173 - Back Cap SKU: 10810PT - Ceramic Radial Bearing SKU: 1032M - Ceramic Angular Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33108 - Wave Washer
	609/619 Original P/N: 1.008.4163	SKU: 21.000 Spindle & Impeller SKU: 21.006 W/ Ceramic Radial Bearings + O-Rings SKU: 21.009 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 21093 - Back Cap SKU: 41094 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1084M - Ceramic Angular Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33108 - Wave Washer
	625/630/640 Original P/N: 0.553.5851	SKU: 21.010 Spindle & Impeller SKU: 21.016 W/ Ceramic Radial Bearings + O-Rings SKU: 21.019 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 210134 - Back Cap SKU: 21083 - Back Cap SKU: 41014 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 1035M - Ceramic Angular Stepped Bearing SKU: 33101 - O-Ring (6 × 1)
	635 Original P/N: 0.553.8581	SKU: 21.036 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 21033 - Back Cap SKU: 41034 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 10610PT - Ceramic Radial Stepped Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33115 - O-Ring (6 × 1,06)
	625/636/646 Original P/N: 0.553.8911	SKU: 21.040 Spindle & Impeller SKU: 21.046 W/ Ceramic Radial Bearings + O-Rings SKU: 21.049 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 21083 - Back Cap SKU: 210134 - Back Cap SKU: 41094 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 1035M - Ceramic Angular Stepped Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33115 - O-Ring (6 × 1,06)
	637/635 BR Original P/N: 0.553.8661	SKU: 21.056 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 21033 - Back Cap SKU: 41034 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 10610PT - Ceramic Radial Stepped Bearing SKU: 33115 - O-Ring (6 × 1,06)
	632/639/642/645 Original P/N: 0.553.6911, 0.553.7471, 0.553.6131	SKU: 21.060 Spindle & Impeller SKU: 21.066 W/ Ceramic Radial Bearings + O-Rings SKU: 21.069 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 21063 - Back Cap SKU: 41014 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 1035M - Ceramic Angular Stepped Bearing SKU: 33101 - O-Ring (6 × 1)
	647/649 Original P/N: 0.553.8541	SKU: 21.076 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 41094 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 33101 - O-Ring (6 × 1)



	650 Original P/N: 0.553.8231	SKU: 21.086 W/ Ceramic Radial Bearings + O-Rings SKU: 21.089 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 21083 - Back Cap SKU: 41094 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 1034M - Ceramic Angular Stepped Bearing SKU: 1035M - Ceramic Angular Stepped Bearing SKU: 33115 - O-Ring (6 × 1,06)
	655/660 Original P/N: 1.000.2078	SKU: 21.090 Spindle & Impeller SKU: 21.096 W/ Ceramic Radial Bearings + O-Rings SKU: 21.099 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 21093 - Back Cap SKU: 41094 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer
	675/677 Original P/N: 1.007.9457	SKU: 21.190 Spindle & Impeller SKU: 21.199 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 1036M - Ceramic Angular Stepped Bearing SKU: 1032M - Ceramic Angular Stepped Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer
	679/680 Original P/N: 1.007.9313	SKU: 21.100 Spindle & Impeller SKU: 21.106 W/ Ceramic Radial Bearings + O-Rings SKU: 21.109 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 21103 - Back Cap SKU: 10830PT - Ceramic Radial Bearing SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 1036M - Ceramic Angular Stepped Bearing SKU: 1084M - Ceramic Angular Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer SKU: 33135 - Wave Washer
	4500/5000 Original P/N: 1.006.4564	SKU: 21.110 Spindle & Impeller SKU: 21.119 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 1036M - Ceramic Angular Stepped Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33135 - Wave Washer
	6000 Original P/N: 1.001.3888	SKU: 21.120 Spindle & Impeller SKU: 21.126 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 41094 - Back Cap Wrench SKU: 10610PT - Ceramic Radial Stepped Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer SKU: 33135 - Wave Washer
	8000/6500 BR Original P/N: 1.003.6884	SKU: 21.166 W/ Ceramic Radial Bearings + O-Rings SKU: 21.169 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 21163 - Back Cap SKU: 41094 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 10830PT - Ceramic Radial Bearing SKU: 1084M - Ceramic Angular Bearing SKU: 1032M - Ceramic Angular Stepped Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer
	6500 B Original P/N: 1.003.9664	SKU: 21.136 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 41094 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer SKU: 33135 - Wave Washer



	7000 Original P/N: 1.0004772	SKU: 21.140 Spindle & Impeller SKU: 21.156 W/ Ceramic Radial Bearings + O-Rings SKU: 21.149 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 21093 - Back Cap SKU: 41094 - Back Cap Wrench SKU: 1032M - Ceramic Angular Smooth Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer
	8700 / Mini Original P/N: 3.001.9915	SKU: 21.209 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 21163 - Back Cap SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 1036M - Ceramic Angular Stepped Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33108 - Wave Washer SKU: 33135 - Wave Washer
	8900/9000 Original P/N: 2.000.2266	SKU: 21.180 Spindle & Impeller SKU: 21.186 W/ Ceramic Radial Bearings + O-Rings SKU: 21.189 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 41094 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1084M - Ceramic Angular Bearing SKU: 33115 - O-Ring (6 × 1,06) SKU: 33108 - Wave Washer



	Tradition Lever Chuck Original P/N: 780118	SKU: 29.000 Spindle & Impeller SKU: 29.006 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10210PT - Ceramic Radial Flanged Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33149 -O-Ring (6 × 0,8)
	Tradition Push Button/XGT Original P/N: 790118	SKU: 29.010 Spindle & Impeller SKU: 29.016 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10210PT - Ceramic Radial Flanged Ball Bear
	Stylus Original P/N: 790510	SKU: 29.020 Spindle & Impeller SKU: 29.029 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 1090M - Ceramic Angular Shield Bearing
	Stylus 540/541 Original P/N: 792310	SKU: 29.030 Spindle & Impeller SKU: 29.039 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 1090M - Ceramic Angular Shield Bearing
	Quiet Air Lever Chuck Original P/N: 484300	SKU: 29.040 Spindle & Impeller SKU: 29.046 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10210PT - Ceramic Radial Flanged Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33108 - Wave Washer



	Quiet Air Push Button Original P/N: NA	SKU: 29.050 Spindle & Impeller SKU: 29.056 Ceramic Radial Bearings + O-Rings	Related products: SKU: 10210PT - Ceramic Radial Flanged Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33108 - Wave Washer
---	--	---	---



	TwinPower Standard Original P/N: PAR-4HEX-O N5010862	SKU: 32.006 W/ Ceramic Radial Bearings + O-Rings SKU: 32.009 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 1087N - Ceramic Angular Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33149 - O-Ring (6 × 0,8)
	TwinPower Torque Original P/N: PAR-4HX-O N5010863	SKU: 32.016 W/ Ceramic Radial Bearings + O-Rings SKU: 32.019 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 1087N - Ceramic Angular Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
	TwinPower Mini Original P/N: PAR-4HUEX-O N5011147	SKU: 32.026 W/ Ceramic Radial Bearings + O-Rings SKU: 32.029 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 1087N - Ceramic Angular Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33109 - O-Ring (6,36 × 0,82)



	Ti Max X700/ NL9000T Original P/N: TIX-TU03 Torque	SKU: 26.006 W/ Ceramic Radial Bearings + O-Rings SKU: 26.009 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 26003 - Back Cap SKU: 46004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1084M - Ceramic Angular Bearing SKU: 33136 - O-Ring (6,3 × 0,9) SKU: 33108 - Wave Washer
	Ti Max X600/ X650 Original P/N: TIX-SU03 Standard	SKU: 26.016 Ceramic Radial Bearings + O-Rings SKU: 26.019 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 26013 - Back Cap SKU: 46014 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33136 - O-Ring (6,3 × 0,9) SKU: 33108 - Wave Washer



	Ti Max X500/ NL9000M Original P/N: TIX-MU03 Mini	SKU: 26.026 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 26023 - Back Cap SKU: 46024 - Back Cap Wrench SKU: 10610PT - Ceramic Radial Stepped Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33136 - O-Ring (6,3 × 0,9) SKU: 33108 - Wave Washer
	S Max M600/ PanaMax Plus/ M900 Original P/N: SX-SU03 Standard	SKU: 26.036 W/ Ceramic Radial Bearings + O-Rings SKU: 26.039 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 26033 - Back Cap SKU: 46034 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1084M - Ceramic Angular Bearing
	Ti Max A600/ NL95 S Original P/N: Ti-SU03 Standard	SKU: 26.056 W/ Ceramic Radial Bearings + O-Rings SKU: 26.059 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 26043 - Back Cap SKU: 46054 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 1087N - Ceramic Angular Bearing
	Ti Max A700 L/ NL95 T Original P/N: Ti-TU03 Torque	SKU: 26.066 W/ Ceramic Radial Bearings + O-Rings SKU: 26.068 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 46054 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing
	Mach Lite A300 S/ Mach2 2S / Phatelus Machlite 2S Original P/N: NMC-SU03 Standard	SKU: 26.076 W/ Ceramic Radial Bearings + O-Rings SKU: 26.079 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 26073 - Back Cap SKU: 46074 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1084M - Ceramic Angular Bearing
	Mach Lite A400 T/ Mach2 2T Original P/N: NMC-TU03 Torque	SKU: 26.086 W/ Ceramic Radial Bearings + O-Rings SKU: 26.089 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 46074 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1084M - Ceramic Angular Bearing
	S Max M500/ M800/ PanaMax Plus M Original P/N: SX-MU03 Mini	SKU: 26.106 W/ Ceramic Radial Bearings + O-Rings SKU: 26.109 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 26103 - Back Cap SKU: 46034 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 1087M - Ceramic Angular Bearing
	Ti Max Z900 Original P/N: TIZ-SU03 Standard	SKU: 26.116 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 26113 - Back Cap SKU: 46114 - Back Cap Wrench SKU: 10210PT - Ceramic Radial Flanged Bearing SKU: 10810PT - Ceramic Radial Bearing SKU: 33136 - O-Ring (6,3 × 0,9)
	Ti Max Z800 & Back Cap Original P/N: TIZ-MU03 Mini	SKU: 26.129 W/ Ceramic Angular Bearings + O-Rings + Back Cap	Related products: SKU: 26123 - Back Cap SKU: 10810PT - Ceramic Radial Bearing SKU: 33136 - O-Ring (6,3 × 0,9) SKU: 33108 - Wave Washer



	Pana Air Original P/N: NPA-SU03 Standard	SKU: 26.136 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10810PT - Ceramic Radial Bearing
	Pana Max 2S Original P/N: PAX2-SU03 Standard	SKU: 26.146 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10810PT - Ceramic Radial Bearing



	3000/4000S Original P/N: NA	SKU: 20.000 Spindle & Impeller SKU: 20.006 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 20003 - Back Cap SKU: 10610PT - Ceramic Radial Stepped Bearing SKU: 33101 - O-Ring (6 × 1)
	TM1 Mini Original P/N: 5951595	SKU: 20.010 Spindle & Impeller SKU: 20.016 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10610PT - Ceramic Radial Stepped Bearing SKU: 33101 - O-Ring (6 × 1)
	TS1/TS2 Original P/N: 5953054	SKU: 20.020 Spindle & Impeller SKU: 20.026 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 20003 - Back Cap SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 33101 - O-Ring (6 × 1)
	TC3 Original P/N: 4175217	SKU: 22.000 Spindle & Impeller SKU: 22.006 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22003 - Back Cap SKU: 42014 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 33101 - O-Ring (6 × 1)



T1 Control	 T1 Control (SN<65.000) Original P/N: 4175217	SKU: 22.000 Spindle & Impeller SKU: 22.006 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22003 - Back Cap SKU: 42014 - Back Cap Wrench SKU: 10310PT - Ceramic Radial Stepped Bearing SKU: 33101 - O-Ring (6 × 1)
T1 Control	 T1 Control (SN>300.000) Original P/N: 5868299	SKU: 22.010 Spindle & Impeller SKU: 22.016 W/ Ceramic Radial Bearings + O-Rings SKU: 22.018 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22013 - Back Cap SKU: 42014 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
T1 Control	 T1 Control (SN>600.000) (SN>700.000) Original P/N: 6432608	SKU: 22.096N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22073 - Back Cap SKU: 42024 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T1 Control	 T1 Control (SN>800.000) Original P/N: 6710367	SKU: 22.186N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22133 - Back Cap SKU: 42014 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T1 Boost	 T1 Boost (SN>600.000) Original P/N: 6432624	SKU: 22.126N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22063 - Back Cap SKU: 42034 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T1 Boost	 T1 Boost CBM (SN>700.000) Original P/N: 6547694	SKU: 22.136N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22093 - Back Cap SKU: 42034 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T1 Boost	 T1 Boost (SN>800.000) Original P/N: 6710409	SKU: 22.206N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22083 - Back Cap SKU: 42054 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T1 mini	 T1 mini (SN>502.000) Original P/N: 5868315	SKU: 22.040 Spindle & Impeller SKU: 22.046 W/ Ceramic Radial Bearings + O-Rings SKU: 22.048 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22043 - Back Cap SKU: 42004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)



T1 mini		T1 mini <i>(SN>600.000)</i> Original P/N: 6432640	SKU: 22.108 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 22113 - Back Cap SKU: 42074 - Back Cap Wrench SKU: 33101 - O-Ring (6 × 1)
T1 mini		T1 mini CBM <i>(SN>700.000)</i> Original P/N: 6548106	SKU: 22.156N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22103 - Back Cap SKU: 42044 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T2 Control		T2 Control <i>(SN<600.000)</i> Original P/N: 5868299	SKU: 22.010 Spindle & Impeller SKU: 22.016 W/ Ceramic Radial Bearings + O-Rings SKU: 22.018 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22013 - Back Cap SKU: 42014 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
T2 Control		T2 Control <i>(SN>600.000)</i> <i>(SN>700.000)</i> Original P/N: 6432590	SKU: 22.096N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22073 - Back Cap SKU: 42024 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T2 Control		T2 Control <i>(SN>800.000)</i> Original P/N: 6710425	SKU: 22.196N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22133 - Back Cap SKU: 42014 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T2 Boost		T2 Boost <i>(SN>600.000)</i> Original P/N: 6432616	SKU: 22.066N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22063 - Back Cap SKU: 42034 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T2 Boost		T2 Boost CBM <i>(SN>700.000)</i> Original P/N: 6548122	SKU: 22.076N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22093 - Back Cap SKU: 42034 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T2 Boost		T2 Boost <i>(SN>800.000)</i> Original P/N: 6710433	SKU: 22.146N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22083 - Back Cap SKU: 42054 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer



T2 Racer		T2 Racer <i>(SN<600.000 smaller than)</i> Original P/N: 5868307	SKU: 22.020 Spindle & Impeller SKU: 22.026 W/ Ceramic Radial Bearings + O-Rings SKU: 22.028 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22023 - Back Cap SKU: 42004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
T2 mini		T2 mini <i>(SN<502.000)</i> Original P/N: 5868315	SKU: 22.040 Spindle & Impeller SKU: 22.046 W/ Ceramic Radial Bearings + O-Rings SKU: 22.048 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22043 - Back Cap SKU: 42004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
T2 mini		T2 mini <i>(SN>600.000)</i> Original P/N: 6432632	SKU: 22.166N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22103 - Back Cap SKU: 42044 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T2 mini		T2 mini <i>(SN>800.000)</i> Original P/N: 6710441	SKU: 22.176N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22123 - Back Cap SKU: 42124 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
T3 Boost		T3 Boost <i>(SN>600.000)</i> Original P/N: 6432616	SKU: 22.066N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22063 - Back Cap SKU: 42034 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T3 Boost		T3 Boost CBM <i>(SN>700.000)</i> Original P/N: 6548122	SKU: 22.076N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22093 - Back Cap SKU: 42034 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T3 Boost		T3 Boost <i>(SN>800.000)</i> Original P/N: 6710433	SKU: 22.146N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22083 - Back Cap SKU: 42054 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T3 Racer		T3 Racer <i>(SN<600.000 smaller than) / P Racer</i> Original P/N: 5868323	SKU: 22.020 Spindle & Impeller SKU: 22.026 W/ Ceramic Radial Bearings + O-Rings SKU: 22.028 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22023 - Back Cap SKU: 42004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)



T3 mini		T3 mini (SN<502.000 smaller than) Original P/N: 5868331	SKU: 22.040 Spindle & Impeller SKU: 22.046 W/ Ceramic Radial Bearings + O-Rings SKU: 22.048 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22043 - Back Cap SKU: 42004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
T3 mini		T3 mini (SN>600.000) Original P/N: 6432632	SKU: 22.166N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22113 - Back Cap SKU: 42074 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 33101 - O-Ring (6 × 1) SKU: 33139 - Wave Washer
T3 mini		T3 mini (SN>800.000) Original P/N: 6710441	SKU: 22.176N W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22123 - Back Cap SKU: 42124 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
SiroBoost		SiroBoost Original P/N: 6334598	SKU: 22.050 Spindle & Impeller SKU: 22.056 W/ Ceramic Radial Bearings + O-Rings SKU: 22.058 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22053 - Back Cap SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33101 - O-Ring (6 x1)
SiroPure Control		SiroPure Control Original P/N: 6084979	SKU: 22.010 Spindle & Impeller SKU: 22.016 W/ Ceramic Radial Bearings + O-Rings SKU: 22.018 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22013 - Back Cap SKU: 42014 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
SiroPure Racer		SiroPure Racer Original P/N: 5868299	SKU: 22.010 Spindle & Impeller SKU: 22.016 W/ Ceramic Radial Bearings + O-Rings SKU: 22.018 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 22013 - Back Cap SKU: 42014 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82)
TMC mini		TMC mini (SN<65.000) Original P/N: 4178591	SKU: 22.030 Spindle & Impeller SKU: 22.036 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 22033 - Back Cap SKU: 10330PT - Ceramic Radial Stepped Bearing SKU: 33101 - O-Ring (6 × 1)




	430/Solara/Vista Original P/N: 263739, 263968	SKU: 27.000 Spindle & Impeller SKU: 27.006 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10710PT - Ceramic Radial Flanged/Unique Bearing
	395/396/695/696 Original P/N: 02137000, 02765500	SKU: 23.000 Spindle & Impeller SKU: 23.006 W/ Ceramic Radial Bearings + O-Rings SKU: 23.009 W/ Ceramic Angular Bearings + O-Rings	Related products: SKU: 23003 - Back Cap SKU: 43004 - Back Cap Wrench SKU: 10810PT - Ceramic Radial Bearing SKU: 1087N - Ceramic Angular Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	398/795/895/898 Original P/N: 02760000	SKU: 23.010 Spindle & Impeller SKU: 23.016 W/ Ceramic Radial Bearings + O-Rings SKU: 23.018 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23013 - Back Cap SKU: 43004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	796/896 Original P/N: 02765500	SKU: 23.020 Spindle & Impeller SKU: 23.026 W/ Ceramic Radial Bearings + O-Rings SKU: 23.028 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Trend TC-95 Original P/N: 04396680	SKU: 23.030 Spindle & Impeller SKU: 23.036 W/ Ceramic Radial Bearings + O-Rings SKU: 23.038 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 0270100 - Back Cap SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Synea TA-98 Original P/N: 04396600, 04386000	SKU: 23.040 Spindle & Impeller SKU: 23.046 W/ Ceramic Radial Bearings + O-Rings SKU: 23.048 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23043 - Back Cap SKU: 43044 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer



	Synea TA-98 LED/C LED Original P/N: 05306900, 05307000	SKU: 23.050 Spindle & Impeller SKU: 23.056 W/ Ceramic Radial Bearings + O-Rings SKU: 23.058 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23053 - Back Cap SKU: 43064 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Synea TA-97 LED/C LED (2007) / TK-97 Original P/N: 05307200, 07272600	SKU: 23.060 Spindle & Impeller SKU: 23.066 W/ Ceramic Radial Bearings + O-Rings SKU: 23.068 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23063 - Back Cap SKU: 43074 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Synea TA-96 Original P/N: 04396700, 04386100	SKU: 23.070 Spindle & Impeller SKU: 23.076 W/ Ceramic Radial Bearings + O-Rings SKU: 23.078 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23073 - Back Cap SKU: 43074 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Alegra TE-95 / Synea TA-98 Original P/N: 04457300	SKU: 23.080 Spindle & Impeller SKU: 23.086 W/ Ceramic Radial Bearings + O-Rings SKU: 23.088 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23083 - Back Cap till 2009 SKU: 23103 - Back Cap from 2009 SKU: 43004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Alegra TE-98 Original P/N: 05958900	SKU: 23.090 Spindle & Impeller SKU: 23.096 W/ Ceramic Radial Bearings + O-Rings SKU: 23.098 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23083 - Back Cap 2003 SKU: 23103 - Back Cap 2009 SKU: 43004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Alegra TE-95 (2009)/TE-98 LED Original P/N: 05958900	SKU: 23.100 Spindle & Impeller SKU: 23.106 W/ Ceramic Radial Bearings + O-Rings SKU: 23.108 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23083 - Back Cap SKU: 23103 - Back Cap SKU: 43004 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Synea Vision TK-98 Original P/N: 07272800	SKU: 23.110 Spindle & Impeller SKU: 23.116 W/ Ceramic Radial Bearings + O-Rings SKU: 23.118 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23053 - Back Cap SKU: 43064 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer



	Synea Fusion TG-98 Original P/N: 07005900	SKU: 23.120 Spindle & Impeller SKU: 23.126 W/ Ceramic Radial Bearings + O-Rings SKU: 23.128 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 23053 - Back Cap SKU: 43064 - Back Cap Wrench SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Synea Vision TK-100 Original P/N: 07272900	SKU: 23.130 Spindle & Impeller SKU: 23.136 W/ Ceramic Radial Bearings + O-Rings SKU: 23.138 W/ Ceramic Angular Integral Shield Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 1086M - Ceramic Angular Integral Shield Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer
	Alegra TE-95 RM / TE-98 (after 2014) Original P/N: 07495401	SKU: 23.146 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 23083 - Back Cap Alegra TE-95 RM SKU: 23103 - Back Cap Alegra TE-98 SKU: 43004 - Back Cap Wrench SKU: 1131M - Ceramic Radial Stepped Double Flanged Bearing
	Alegra TE-97 (after 2014) Original P/N: 07234101	SKU: 23.156 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 1131M - Ceramic Radial Stepped Double Flanged Bearing SKU: 33149 - O-Ring (6 × 0,8)
NEW 	Alegra TE-97 / LED (before 2014) Original P/N: 05959100	SKU: 23.160 Spindle & Impeller SKU: 23.166 W/ Ceramic Radial Bearings + O-Rings	Related products: SKU: 10830PT - Ceramic Radial Bearing SKU: 33109 - O-Ring (6,36 × 0,82) SKU: 33108 - Wave Washer



Highspeed Handpiece Parts



Back Caps

BIEN AIR / Back Caps

Highspeed Handpiece Parts



Black Pearl (with Washers)

SKU: **24003**
Original P/N: 130.07.00, 130.07.01, 130.07.02, 130.07.03

Related Wrench: SKU: **44004**



Tornado

SKU: **24023**
Original P/N: 1502650

Related Wrench: SKU: **44024**



Bora (with Washers)

SKU: **24013**
Original P/N: 130.10.86, 130.10.85, 150.02.34, 130.45.39

Related Wrench: SKU: **44014**

NEW



Bora 2

SKU: **24033**
Original P/N: 1502999

Related Wrench: SKU: **44034**

CASTELLINI / Back Caps

Highspeed Handpiece Parts



Silent Power

SKU: **25003**
Original P/N: F699F40R0

Related Wrench: SKU: **45004**



Silent Power EVO

SKU: **25013**
Original P/N: 95529229

Related Wrench: SKU: **45014**

KAVO / Back Caps

Highspeed Handpiece Parts



635/637

SKU: **21033**
Original P/N: 0.553.8491

Related Wrench: SKU: **41034**



680

SKU: **21103**
Original P/N: 1.007.8903

Related Wrench: SKU: **41094**



632/639/642

SKU: **21063**
Original P/N: 0.553.5642

Related Wrench: SKU: **41014**



614/625/630

SKU: **210134**
Original P/N: 0.553.1782

Related Wrench: SKU: **41014**

KAVO / Back Caps
Highspeed Handpiece Parts



640/636/646/650

SKU: **21083**
Original P/N: 0.553.6511, 0.553.8742,
0.553.806

Related Wrench: SKU: **41094**



8000/9000

SKU: **21163**
Original P/N: 1.003.5525, 3.000.1874

Related Wrench: SKU: **41094**



609/619/655/660/7000

SKU: **21093**
Original P/N: 1.007.0854, 1.000.2076,
1.000.4770

Related Wrench: SKU: **41094**



605/615

SKU: **21173**
Original P/N: 1.007.7277

NSK / Back Caps
Highspeed Handpiece Parts



Ti Max X700

SKU: **26003**
Original P/N: TiX-TU04

Related Wrench: SKU: **46004**



Ti Max A600

SKU: **26043**
Original P/N: Ti-SU04

Related Wrench: SKU: **46054**



Ti Max X600

SKU: **26013**
Original P/N: TiX-SU04

Related Wrench: SKU: **46014**



Mach Lite A300

SKU: **26073**
Original P/N: NCH-MU04

Related Wrench: SKU: **46074**



Ti Max X500

SKU: **26023**
Original P/N: TiX-MU04

Related Wrench: SKU: **46024**



Ti Max Z900

SKU: **26113**
Original P/N: TiZ-SU04

Related Wrench: SKU: **46114**



S Max M500/M800

SKU: **26103**
Original P/N: SX-MU04

Related Wrench: SKU: **46034**



Ti Max Z800

SKU: **26123**
Original P/N: TiZ-MU04



S Max M600

SKU: **26033**
Original P/N: SX-SU04

Related Wrench: SKU: **46034**

SIEMENS / Back Caps
Highspeed Handpiece Parts



4000S/TS1/TS2

SKU: **20003**
Original P/N: NA



TC3

SKU: **22033**
Original P/N: 3327595

SIRONA / Back Caps
Highspeed Handpiece Parts



T1 Control



T1 Control < 65.000

SKU: **22003**
Original P/N: 3327595

Related Wrench: SKU: **42014**

T1 Boost



T1 Boost > 700.000

SKU: **22093**
Original P/N: 6548155

Related Wrench: SKU: **42034**

T1 Control



**T1 Control < 600.000
(SN lower than 600.000)**

SKU: **22013**
Original P/N: 6000884

Related Wrench: SKU: **42004**

T1 Boost



T1 Boost > 800.000

SKU: **22083**
Original P/N: 6725712

Related Wrench: SKU: **42054**

T1 Control



T1 Control > 600.000

SKU: **22073**
Original P/N: 6465798

Related Wrench: SKU: **42024**

T1 mini



**T1 mini < 600.000
(SN lower than 600.000)**

SKU: **22043**
Original P/N: 5828871

Related Wrench: SKU: **42024**

T1 Control



T1 Control > 700.000

SKU: **22073**
Original P/N: 6548148

Related Wrench: SKU: **42024**

T1 mini



T1 mini > 600.000

SKU: **22113**
Original P/N: 6465814

Related Wrench: SKU: **42074**

NEW

T1 Control



T1 Control > 800.000

SKU: **22133**
Original P/N: 6725704

Related Wrench: SKU: **42134**

NEW

T1 mini



T1 mini > 700.000

SKU: **22103**
Original P/N: 6548163

Related Wrench: SKU: **42044**

T1 Boost



T1 Boost > 600.000

SKU: **22063**
Original P/N: 6465806

Related Wrench: SKU: **42034**

NEW

T1 mini



T1 mini > 800.000

SKU: **22123**
Original P/N: 6725720

Related Wrench: SKU: **42124**



T2 Control		T2 Control < 600.000 (SN lower than 600.000) SKU: 22013 Original P/N: 6000884 Related Wrench: SKU: 42004	NEW		T2 mini > 700.000 SKU: 22103 Original P/N: 6548163 Related Wrench: SKU: 42044
T2 Control		T2 Control > 600.000 SKU: 22073 Original P/N: 6465798 Related Wrench: SKU: 42024	NEW		T2 mini > 800.000 SKU: 22123 Original P/N: 6725720 Related Wrench: SKU: 42124
T2 Control		T2 Control > 700.000 SKU: 22073 Original P/N: 6465798 Related Wrench: SKU: 42024			T2 Racer SKU: 22023 Original P/N: 6000892 Related Wrench: SKU: 42004
NEW		T2 Control > 800.000 SKU: 22133 Original P/N: 6725704 Related Wrench: SKU: 42134			T3 Boost > 600.000 SKU: 22063 Original P/N: 6465806 Related Wrench: SKU: 42034
		T2 Boost > 600.000 SKU: 22063 Original P/N: 6465806 Related Wrench: SKU: 42034			T3 Boost > 700.000 SKU: 22093 Original P/N: 6548155 Related Wrench: SKU: 42034
		T2 Boost > 700.000 SKU: 22093 Original P/N: 6548155 Related Wrench: SKU: 42034			T3 Boost > 800.000 SKU: 22083 Original P/N: 6725712 Related Wrench: SKU: 42054
		T2 Boost > 800.000 SKU: 22083 Original P/N: 6725712 Related Wrench: SKU: 42054			T3 mini < 600.000 (SN lower than 600.000) SKU: 22043 Original P/N: 6000876 Related Wrench: SKU: 42024
		T2 mini < 600.000 (SN lower than 600.000) SKU: 22043 Original P/N: 6000876 Related Wrench: SKU: 42024	NEW		T3 mini > 600.000 SKU: 22113 Original P/N: 6465814 Related Wrench: SKU: 42074
NEW		T2 mini > 600.000 SKU: 22113 Original P/N: 6465814 Related Wrench: SKU: 42074	NEW		T3 mini > 700.000 SKU: 22103 Original P/N: 6548163 Related Wrench: SKU: 42044



NEW		T3 mini > 800.000 SKU: 22123 Original P/N: 6725720 Related Wrench: SKU: 42124	SiroPure Control		SKU: 22013 Original P/N: 6000884 Related Wrench: SKU: 42004
		T3 Racer < 600.000 (SN lower than 600.000) SKU: 22023 Original P/N: 6000892 Related Wrench: SKU: 42004	SiroPure mini		SKU: 22043 Original P/N: 6000876 Related Wrench: SKU: 42024
		T4 Racer < 600.000 (SN lower than 600.000) SKU: 22023 Original P/N: 6000892 Related Wrench: SKU: 42004	TMC mini		SKU: 22033 Original P/N: 6000918
SiroBoost		SiroBoost SKU: 22053 Original P/N: 6401439 Related Wrench: SKU: 42004			



	395 SKU: 23003 Original P/N: NA Related Wrench: SKU: 43004		Synea TA98LED/TA98 after 2007/ Vision TK98/Fusion TG98 SKU: 23053 Original P/N: 06446300, 06571200 Related Wrench: SKU: 43064
	398/795/895/898 SKU: 23013 Original P/N: 2701005 Related Wrench: SKU: 43004		Synea TA97/TA97LED after 2007/ TG97 SKU: 23063 Original P/N: 06446400 Related Wrench: SKU: 43064
	Synea TA98 before 2007 SKU: 23043 Original P/N: 3503000 Related Wrench: SKU: 43004		Synea TA96 before 2007 SKU: 23073 Original P/N: 03503900 Related Wrench: SKU: 43004



Alegria TE95/TE98 before 2009

SKU: **23083**
Original P/N: **04232400**
Related Wrench: SKU: **43004**



Alegria TE95/TE98/TE98LED after 2009


SKU: **23103**
Original P/N: **05489800**
Related Wrench: SKU: **43004**

NEW



Trend TC95RM/TC80

SKU: **02701000**
Original P/N: **02701000**



Highspeed Handpiece Parts



Back Cap Wrenches

BIEN AIR / Back Cap Wrenches
Highspeed Handpiece Parts



	Wrench for Bora Outer Ring SKU: 44014		Wrench for Back Cap for Tornado SKU: 44044
NEW 	Wrench for Bora Inner Ring SKU: 44014-1		Wrench for Back Cap for Bora 2 SKU: 44034
	Wrench for Croma Outer Ring SKU: 44004	NEW 	Wrench for Prestige Outer Ring SKU: 44024
NEW 	Wrench for Croma Inner Ring SKU: 44014-1	NEW 	Wrench for Prestige Inner Ring SKU: 44024-1

CASTELLINI / Back Cap Wrenches
Highspeed Handpiece Parts



	Wrench for Back Cap for Silent Power SKU: 45004		Wrench for Back Cap for Silent Power EVO SKU: 45014
---	--	---	--

KAVO / Back Cap Wrenches
Highspeed Handpiece Parts



	Wrench for Back Cap for 630/639/640/642/645 SKU: 41014		Kavo Universal Back Cap Wrench SKU: 41094 4-point Wrench
	Wrench for Back Cap for 635/637 SKU: 41034		Wrench for Back Cap for 650 SKU: 41035

NSK / Back Cap Wrenches
Highspeed Handpiece Parts



	Wrench for Back Cap for Ti Max X700 SKU: 46004		Wrench for Back Cap for Ti Max A600/A700 SKU: 46054
	Wrench for Back Cap for Ti Max X600 SKU: 46014		Wrench for Back Cap for Mach/A300L/A400L/Pana Air SKU: 46074
	Wrench for Back Cap for Ti Max X500 SKU: 46024		Wrench for Back Cap for Ti Max A500 SKU: 46094
	Wrench for Back Cap for S-Max M600/M500/Pana Max Plus SKU: 46034		Wrench for Back Cap for Ti Max Z900 SKU: 46114


SIRONA / Back Cap Wrenches
Highspeed Handpiece Parts



	Wrench for Back Cap for T1/T2/TC3 Control SKU: 42014 3-point Wrench		Wrench for Back Cap for T1/T2/T3 Boost > 800.000 SKU: 42054
	Sirona Universal Back Cap Wrench SKU: 42004		Wrench for Back Cap for T1/T2/T3 Mini > 700.000 SKU: 42044
	Wrench for Back Cap for T1/T2 Control > 600.000 SKU: 42024		Wrench for Back Cap for T1/T2/T3 mini > 600.000 SKU: 42074
	Wrench for Back Cap for T1/T2 Control > 800.000 SKU: 42134		Wrench for Back Cap for T1/T2/T3 mini > 800.000 SKU: 42124
	Wrench for Back Cap for T1/T2/T3 Boost > 600.000 T1/T2/T3 Boost > 700.000 SKU: 42034 3-point Wrench		

W&H / Back Cap Wrenches
Highspeed Handpiece Parts

	W&H Universal Back Cap Wrench SKU: 43004 2-point Wrench		Wrench for Back Cap for Synea TA97/TA97LED/TA98LED SKU: 43064
	Wrench for Back Cap for Synea TA98 SKU: 43044 2-point Wrench		Wrench for Back Cap for Synea TA96 SKU: 43074 2-point Wrench



Highspeed Handpiece Parts



**O-Rings &
 Wave Washers**

BIEN AIR / O-Rings and Wave Washers
Highspeed Handpiece Parts



O-ring for Bien Air Bora/Croma/Boralina
8 × 1

SKU: **33100**



O-ring for Bien Air Prestige/Prestilina/TD783/TDS890 – smaller
3,5 × 1

SKU: **33110**



O-ring for Bien Air Prestige/Tornado,
6,3 × 1

SKU: **33118**

KAVO / O-Rings and Wave Washers
Highspeed Handpiece Parts



O-ring for Kavo Rotors – older models
6 × 1

SKU: **33101**
Original P/N: **1.002.1149**



O-ring for Kavo Rotors – newer
models 6 × 1,06

SKU: **33115**
Original P/N: **0.200.6172**



Wave Washer for Kavo rotors
4,9 × 6,3 × 0,5 × 0,1

SKU: **33108**



Wave Washer for Kavo rotors
4,7 × 6,0 × 0,5 × 0,09

SKU: **33135**

NSK / O-Rings and Wave Washers
Highspeed Handpiece Parts



O-ring for NSK Ti Max rotors
6,3 × 0,9

SKU: **33136**

SIRONA / O-Rings and Wave Washers
Highspeed Handpiece Parts



O-ring for Sirona rotors
6 × 1

SKU: **33101**

STAR / O-Rings and Wave Washers
Highspeed Handpiece Parts



O-ring for Star 430
6 × 0,8

SKU: **33103**

W&H / O-Rings and Wave Washers
Highspeed Handpiece Parts



O-ring for W&H rotors
6,36 × 0,82

SKU: **33109**



Wave Washer for W&H rotors
4,9 × 6,3 × 0,5 × 0,1

SKU: **33108**



O-ring for W&H rotors
6 × 0,8

SKU: **33149**



Highspeed Handpiece Parts



Bulbs

BIEN AIR / Bulbs
Highspeed Handpiece Parts



	LED Bulb HS Handpieces and Quick Connectors SKU: 67000LED Original P/N: 1500000/204.33.07		Xenon Bulb HS Handpieces and Quick Connectors SKU: 77000 Original P/N: 1500000/204.33.07
	LED Bulb Motors SKU: 67002LED Original P/N: 1500007/211.60.16		Xenon Bulb Motors SKU: 77002 Original P/N: 1500007/211.60.16

KAVO / Bulbs
Highspeed Handpiece Parts



	LED Bulb Quick Connectors SKU: 67007LEDT		Xenon Bulb Quick Connectors and Motors SKU: 77001 Original P/N: 1.002.2928
	LED Bulb Motors SKU: 67001LEDM Original P/N: 1.007.5372		

NSK / Bulbs
Highspeed Handpiece Parts



	LED Bulb Quick Connectors SKU: 67006LED Original P/N: Y900529		Xenon Bulb Quick Connectors SKU: 77006 Original P/N: Y900529
	LED PTL-CL-LED Quick Connectors SKU: 67008 Original P/N: 10000198		

SIRONA / Bulbs
Highspeed Handpiece Parts




	LED Bulb Quick Connectors (brown) SKU: 67003LED Original P/N: 5940291		Xenon Bulb Quick Connectors SKU: 77003 Original P/N: 5940291
	LED Bulb Motors (green) SKU: 67004LED Original P/N: 6314558		Xenon Bulb Motors SKU: 77004 Original P/N: 5945407


W&H / Bulbs
Highspeed Handpiece Parts



	LED Bulb HS Handpieces SKU: 67005LED Original P/N: 02207600/02207610		Xenon Bulb HS Handpieces SKU: 77005 Original P/N: 02207600/02207610
--	---	--	--



Contra Angle Parts





Contra Angle Parts

Parts Listed by Brand and Model

for red **Bien Air CA 1:5 L**

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83864 Original P/N: 1501933	
	Set of Bearings for Head Gear SKU: 11025 Original P/N: 02649400	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 83877 Original P/N: 1500057	
	Front Gear for Intermediate Shaft SKU: 83876 Original P/N: NA	
	Set of Gears and Bearings for Intermediate Shaft SKU: 83882 Set consists: Front Gear + Bearings (701.01.75)	


for red **Bien Air CA EVO 15**

Head Gear / Cartridge ▼		
NEW 	Head Gear / Cartridge SKU: 84046 Original P/N: 1502645	
	Set of Bearings for Head Gear SKU: 11006 Original P/N: 5899863	

for blue **Bien Air CA 1:1**

Head Gear / Cartridge ▼		
NEW		Head Gear / Cartridge SKU: 84060 Original P/N: 1501085 
Back Cap ▼		
NEW		Back Cap SKU: 84160 Original P/N: 1500550 
Intermediate Shaft ▼		
NEW		Intermediate Shaft SKU: 84062 Original P/N: 1501083 

for red **Kavo 25LP/A / 200XD/R**


Head Gear / Cartridge ▼		
		Head Gear / Cartridge SKU: 83840 Original P/N: 1.001.3520, 1.002.6135 
		Set of Bearings for Head Gear SKU: 11001 Original P/N: 1.001.2836, 1.001.3167 
Back Cap ▼		
NEW		Back Cap SKU: 84065 Original P/N: 1.001.3283 
Intermediate Shaft ▼		
		Intermediate Shaft SKU: 83856 Original P/N: 1.001.3515 
		Front Gear for Intermediate Shaft SKU: 83871 Original P/N: NA 
		Rear Gear for Intermediate Shaft SKU: 83872 Original P/N: NA 
		Set of Bearings for Intermediate Shaft SKU: 1001 2834 Original P/N: 1.001.2834 
		Set of Gears and Bearings for Intermediate Shaft SKU: 83880 Original P/N: NA 
Insert/Intermediate Drive ▼		
		Insert/Intermediate Drive SKU: 83875 Original P/N: 1.001.3519 

... to be continued on the next page

for red **Kavo 25LP/A / 200XD/R**

Lower Shaft ▼		
	Gear for Lower Shaft SKU: 1001 3071 Original P/N: 1.001.3071	
Fiber Optic ▼		
	Fiber Optic Rod SKU: 1001 3320 Original P/N: 1.001.3320	

for red **Kavo 24LN / 25LH/A/C**

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83821 Original P/N: 0.553.7341, 0.553.5761, 1.002.4748	
	Set of Bearings for Head Gear SKU: 11002 Original P/N: 1.001.2836	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 83852 Original P/N: 0.553.7201	
	Front Gear for Intermediate Shaft SKU: 83873 Original P/N: NA	
	Rear Gear for Intermediate Shaft SKU: 83874 Original P/N: NA	
	Set of Gears and Bearings for Intermediate Shaft SKU: 83881 Original P/N: NA	
	Set of Bearings for Intermediate Shaft SKU: 220 0730 Original P/N: 0.220.0730	
Fiber Optic ▼		
	Fiber Optic Rod for 24LN SKU: 553 6421 Original P/N: 0.553.6421	
	Fiber Optic Rod for 25LH/LHA SKU: 553 0283 Original P/N: 0.553.0283	

for red **Kavo LUX M25L**

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 84040 Original P/N: 1.001.0129	
	Set of Bearings for Head Gear SKU: 11026 Original P/N: 1.010.9766, 1.010.9767	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 83856 Original P/N: 1.001.3515	
	Front Gear for Intermediate Shaft SKU: 83871 Original P/N: NA	
	Rear Gear for Intermediate Shaft SKU: 83872 Original P/N: NA	
	Set of Bearings for Intermediate Shaft SKU: 1001 2834 Original P/N: 1.001.2834	
Insert/Intermediate Drive ▼		
	Insert/Intermediate Drive SKU: 83875 Original P/N: 1.001.3519	
Lower Shaft ▼		
	Gear for Lower Shaft SKU: 1001 3071 Original P/N: 1.001.3071	
Fiber Optic ▼		
	Fiber Optic Rod SKU: 1001 3320 Original P/N: 1.001.3320	

for red **Kavo LUX M05L Mini**

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 84041 Original P/N: 1.011.1374	
	Set of Bearings for Head Gear SKU: 11026 Original P/N: 1.010.9767, 1.010.9766	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 84055 Original P/N: 1.011.3206	
Insert/Intermediate Drive ▼		
	Insert/Intermediate Drive SKU: 83875 Original P/N: 1.001.3519	
Lower Shaft ▼		
	Gear for Lower Shaft SKU: 1001 3071 Original P/N: 1.001.3071	

for red **Kavo E25L/C**

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83865 Original P/N: 1.008.4040, 1.014.5555	
	Set of Bearings for Head Gear SKU: 11024 Original P/N: 1.008.4038, 1.008.4039	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 83897 Original P/N: 1.007.9916, 1.014.5557	
Tools ▼		
	Back Cap Wrench SKU: 47004	

for blue **Kavo 68LH/LU / E20L/C / 2068LHC/CHC**

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83822 Original P/N: 0.540.3661	
	Set of Bearings for Head Gear SKU: 11004 Original P/N: 0.220.0051, 0.220.2005	
	Head Gear Fitting Set SKU: 83813 Original P/N: 0.540.4072, 0.540.4102, 0.540.4522	
Back Cap ▼		
	Back Cap SKU: 83892 Original P/N: 0.540.3741	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 83837 Original P/N: 0.540.1151	
Tools ▼		
	Back Cap Wrench SKU: 47001 Original P/N: 0.411.2371	

for blue **Kavo 68LDN / 2068D**

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83823 Original P/N: 0.540.2491	
	Set of Bearings for Head Gear SKU: 11003 Original P/N: 0.220.0105, 0.220.0713	
	Head Gear Fitting Set SKU: 83813 Original P/N: 0.540.4072, 0.540.4102, 0.540.4522	
Back Cap ▼		
NEW 	Back Cap SKU: 84100 Original P/N: 0.540.2501	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 83837 Original P/N: 0.540.1151	
Tools ▼		
	Back Cap Wrench SKU: 47001 Original P/N: 0.411.2371	

for green **Kavo 67LH / E15L/C / 2967CHC**

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83827 Original P/N: 0.540.3671	
	Set of Bearings for Head Gear SKU: 11004 Original P/N: 0.220.0051, 0.220.2005	
	Head Gear Fitting Set SKU: 83813 Original P/N: 0.540.4072, 0.540.4102, 0.540.4522	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 83825 Original P/N: 0.540.1161	

for **Kavo CL3/C3** Surgical Instrument

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83868 Original P/N: 1.011.5595	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 84050 Original P/N: 1.003.4207	

for red **NSK S-Max M95L / Ti-Max X95L / Ti95**

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83883 Original P/N: C600C015	
Back Cap ▼		
NEW 	Back Cap SKU: 84056 Original P/N: C600B050/10000664	
Intermediate Shaft ▼		
	Intermediate Shaft for S-Max M95L SKU: 83895 Original P/N: 10000559	
	Intermediate Shaft for Ti-Max X95L SKU: 83886 Original P/N: C600C22A	
Tools ▼		
	Back Cap Wrench SKU: 47003 Original P/N: DT-350	

for red **NSK Z95** after SN 0BCX0302 (new style)

Head Gear / Cartridge ▼		
NEW 	Head Gear / Cartridge SKU: 84049 Original P/N: C103A4001	
	Set of Bearings for Head Gear SKU: 11029 Original P/N: NA NOTE: Bearing size front: 3,8 x 6,6 x1,74 Bearing size rear: 4 x 7 x 2	
Back Cap ▼		
NEW 	Back Cap SKU: 84044 Original P/N: C103A050	
Intermediate Shaft ▼		
NEW 	Intermediate Shaft SKU: 84070 Original P/N: C103A022002	

for red **NSK Z95** before SN 0BCX0302 (old style)

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 84045 Original P/N: C103A017	
	Set of Bearings for Head Gear SKU: 103A015B Original P/N: NA NOTE: Bearing size (both) 3,8 x 6,6 x1,74	
Intermediate Shaft ▼		
NEW 	Intermediate Shaft SKU: 84070 Original P/N: C103A022002	

for blue **NSK Ti-Max X15/X25 / S-Max M15/M25**

Head Gear / Cartridge



Head Gear / Cartridge

SKU: 83898
Original P/N: C601015



Intermediate Shaft

NEW



Intermediate Shaft

SKU: 84080
Original P/N: 10000574, C601022



for blue **NSK Ti-Max 15L/25/Ti25**

Head Gear / Cartridge



Head Gear / Cartridge

SKU: 84053
Original P/N: C875060



for blue **NSK Ti-Max Z25L**

Head Gear / Cartridge

NEW



Head Gear / Cartridge

SKU: 84087
Original P/N: C1038015



Back Cap

NEW



Back Cap

SKU: 84088
Original P/N: C1038050001



SIRONA/ Contra Angle Parts

red 1:5

for red **Sirona A200/L / C200/L / S200/L**
(T13, $\Phi 9,7$)

Head Gear / Cartridge



Head Gear / Cartridge

SKU: 83860
Original P/N: 6175587, 6405687

NOTE: Head size $\Phi 9,7$; 13 teeth



Set of Bearings for Head Gear



Set of Bearings for Head Gear

SKU: 11005
Original P/N: 5890855



Intermediate Shaft



Intermediate Shaft

SKU: 83861
Original P/N: 6175603



for red **Sirona C200 / C200L**
(old - T11, $\Phi 9,7$)

Head Gear / Cartridge



Head Gear / Cartridge

SKU: 83866
Original P/N: 6000991

NOTE: Head size $\Phi 9,7$; 11 teeth



Set of Bearings for Head Gear



Set of Bearings for Head Gear

SKU: 11005
Original P/N: 58908556



Intermediate Shaft



Intermediate Shaft

SKU: 83861
Original P/N: 6175603



64  www.vkl.biz / info@vkl.biz

info@vkl.biz / www.vkl.biz  65

for red **Sirona T1 Line C200L / S200L Mini**
(T13, Φ8,7)

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83869T Original P/N: 6175611, 6562586, 6405679 NOTE: Head size Φ8,7; 13 teeth	
	Set of Bearings for Head Gear SKU: 11006 Original P/N: 5899863	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 83861 Original P/N: 6175603	

for blue **Sirona T1 C40L / T2 A40 / T3 E40** (Φ9,7)

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83867 Original P/N: 6000983 NOTE: Front bearing dimension 3x7x2	
	Set of Bearings for Head Gear SKU: 11028 Original P/N: 4181157, 220 0105	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 84043 Original P/N: 6000959	


for blue **Sirona T1 C40L / T2 A40 Mini** (Φ8,7)

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 84042 Original P/N: 6405695 NOTE: Front bearing dimension 3x6x2	
	Set of Bearings for Head Gear SKU: 11027 Original P/N: 5887174, 6360593	
Back Cap ▼		
	Back Cap SKU: 6049329 Original P/N: 6049329	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 84043 Original P/N: 6000959	

for red **W&H WA99 A/LT**
(before 2007 - 2 red rings)

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83841T Original P/N: 04517200	
	Set of Bearings for Head Gear SKU: 11011 Original P/N: 04486000	
Intermediate Shaft ▼		
	Set of Bearings for Intermediate Shaft SKU: 11012 Original P/N: 01579600	
Lower Shaft ▼		
	Gear for Lower Shaft SKU: 3815800 Original P/N: 03815800	
	Set of Bearings for Lower Shaft SKU: 11013 Original P/N: 01565800	
Fiber Optic ▼		
	Fiber Optic Rod SKU: 4058500 Original P/N: 04058500	

for red **W&H WA99 A/LT**
(after 2007 - 1 red ring)

Head Gear / Cartridge ▼		
	Head Gear / Cartridge SKU: 83842T Original P/N: 05220300	
	Set of Bearings for Head Gear SKU: 11014 Original P/N: 05798700, 05257200	
Intermediate Shaft ▼		
	Intermediate Shaft SKU: 83847 Original P/N: 06420900	
	Set of Bearings for Intermediate Shaft SKU: 11012 Original P/N: 01579600	
Lower Shaft ▼		
	Gear for Lower Shaft SKU: 3815800 Original P/N: 03815800	
	Set of Bearings for Lower Shaft SKU: 11013 Original P/N: 01565800	
Fiber Optic ▼		
	Fiber Optic Rod SKU: 4058500 Original P/N: 04058500	

for red **W&H WG99 A/LT** before SN 007801 (old)

Head Gear / Cartridge ▼



Head Gear / Cartridge

SKU: 83859
Original P/N: 06993900



Intermediate Shaft ▼



Intermediate Shaft

SKU: 83847
Original P/N: 06420900



Fiber Optic ▼



Fiber Optic Rod

SKU: 4058500
Original P/N: 04058500



for red **W&H WG99 A/LT / WK99 A/LT** after SN 007802(new)

Head Gear / Cartridge ▼



Head Gear / Cartridge

SKU: 83893T
Original P/N: 07473400



Set of Bearings for Head Gear ▼



Set of Bearings for Head Gear

SKU: 05257200
Original P/N: 05257200



Intermediate Shaft ▼



Intermediate Shaft

SKU: 83887
Original P/N: 06579200



Fiber Optic ▼



Fiber Optic Rod

SKU: 6584600
Original P/N: 06584600



for red **W&H WE99 LED G**

Head Gear / Cartridge ▼



Head Gear / Cartridge

SKU: 83853
Original P/N: 05846100



for blue **W&H WA56** / green **WA66**
(before 2007 - 2 blue/green rings)

Head Gear / Cartridge ▼		
	Set of Bearings for Head Gear SKU: 11007 Original P/N: 01589600	
Intermediate Shaft ▼		
	Set of Bearings for Intermediate Shaft SKU: 11008 Original P/N: 01579300	
Lower Shaft ▼		
	Set of Bearings for Lower Shaft SKU: 11009 Original P/N: 01589400, 01577000	
Fiber Optic ▼		
	Fiber Optic Rod SKU: 4058500 Original P/N: 04058500	

for blue **W&H WA56** / green **WA66**
(after 2007 - 1 blue/ green ring)

Head Gear / Cartridge ▼		
	Set of Bearings for Head Gear SKU: 11010 Original P/N: 05039400	
Intermediate Shaft ▼		
	Intermediate Shaft for WA56 SKU: 83849 Original P/N: 05233200	
	Intermediate Shaft for WA66 SKU: 05273000 Original P/N: 05273000	
	Set of Bearings for Intermediate Shaft SKU: 11008 Original P/N: 01579300	
Lower Shaft ▼		
	Set of Bearings for Lower Shaft SKU: 11009 Original P/N: 01589400, 01577000	
Fiber Optic ▼		
	Fiber Optic Rod SKU: 4058500 Original P/N: 04058500	

for blue **W&H WE56A/LED G**

Head Gear / Cartridge ▼



Head Gear / Cartridge

SKU: 83843
Original P/N: 05844400



for blue **W&H WE56E** / green **WE66/E/LED G**
(after 2007 - 1 blue/green ring)

Head Gear / Cartridge ▼



Head Gear / Cartridge

SKU: 83858
Original P/N: 05843300



for blue **W&H WG56 A/LT** / green **WG66 A/LT** (old)

Head Gear / Cartridge ▼



Head Gear / Cartridge

SKU: 83857
Original P/N: 06994300

NOTE: Front (bushing) side gears.
Front and rear bearings are different



Intermediate Shaft ▼



Intermediate Shaft for WG56 LT (SN before 14634)
Intermediate Shaft for WG56 A (SN before 02462)

SKU: 83849
Original P/N: 05233200





Intermediate Shaft for WG66 LT (SN before 2950)
Intermediate Shaft for WG66 A (SN before 1185)

SKU: 05273000
Original P/N: 05273000



Fiber Optic ▼



Fiber Optic Rod

SKU: 4058500
Original P/N: 04058500



for blue **W&H WG56 A/LT** (new)

Head Gear / Cartridge ▼



Head Gear / Cartridge

SKU: 83896
Original P/N: 07605900

NOTE: Rear side gears.
Front and rear bearings are same.



Intermediate Shaft ▼



Intermediate Shaft for WG56 LT (SN after 14634)
Intermediate Shaft for WG56 A (SN after 02462)

SKU: 0679600
Original P/N: 0679600



Fiber Optic ▼



Fiber Optic Rod

SKU: 6584600
Original P/N: 06584600





Contra Angle Parts

Parts Listed by Category

 1:5 Head Gears for Contra Angles

red 1:5

Bien Air		<i>Bien Air CA 1:5 L</i> SKU: 83864 Original P/N: 1501933	NEW		<i>NSK Z95 after SN OBCX0302 (new)</i> SKU: 84049 Original P/N: C103A4001
Bien Air	NEW 	<i>Bien Air EVO 15</i> SKU: 84046 Original P/N: 1502645	Sirona		<i>Sirona A200/C200L/E200/S200L; 13Teeth; ⌀9,7</i> SKU: 83860 Original P/N: 6175587, 6405687
B.A. Int.		<i>B.A. Int. BA200LTS/BA250LT; 13Teeth; ⌀9,7</i> SKU: 83860 Original P/N: 6000991	Sirona		<i>Sirona C200/L; old; 11Teeth; ⌀9,7</i> SKU: 83866 Original P/N: 6000991
Kavo		<i>Kavo 25LPA/200XD/XDR</i> SKU: 83840 Original P/N: 1.001.3520, 1.002.6135	Sirona		<i>Sirona S200L/C200L Mini; 13Teeth; ⌀8,7</i> SKU: 83869T Original P/N: 6175611, 6405679, 6562586
Kavo		<i>Kavo 24LN/25LH/LN/LHA/LHC</i> SKU: 83821 Original P/N: 0.553.7341, 0553.5761	W&H		<i>W&H WA99A/LT before 2007</i> SKU: 83841T NOTE: 2 red rings Original P/N: 04517200
Kavo		<i>Kavo E25L/ E25C</i> SKU: 83865 Original P/N: 1.008.4040	W&H		<i>W&H WA99A/LT after 2007</i> SKU: 83842T NOTE: 1 red ring Original P/N: 05220300
Kavo		<i>Kavo LUX M25L</i> SKU: 84040 Original P/N: 1.001.0129	W&H		<i>W&H WG99 till SN 7801 / WG99A till SN 1818</i> SKU: 83859 Original P/N: 06993900
Kavo		<i>Kavo LUX M05L Mini</i> SKU: 84041 Original P/N: 1.011.1374	W&H		<i>W&H WG99 after SN7801 / WG99A after SN 1818</i> SKU: 83893T Original P/N: 07473400
NSK		<i>NSK M95L/X95L</i> SKU: 83883 Original P/N: C600C015, C600C015001	W&H		<i>W&H WE99 LED G</i> SKU: 83853 Original P/N: 05846100
NSK		<i>NSK Z95 before SN OBCX0302 (old)</i> SKU: 84045 Original P/N: C1034017			



1:1 Head Gears for Contra Angles


blue/green 1:1

NEW Bien Air		Bien Air CA 1:1 SKU: 84060 Original P/N: 1501085	Sirona		Sirona T1 C40L / T2 A40L / T3 E40L; Ø9,7 SKU: 83867 Original P/N: 6000983
Kavo		Kavo 68LH/LU/ E20L/ 2068LHC/ CHC SKU: 83822 Original P/N: 0.540.3661	Sirona		Sirona T1 C40 / T1 S40 / T2 A40 Mini; Ø8,7 SKU: 84042 Original P/N: 6405695
Kavo		Kavo 68LDN/ 2068D SKU: 83823 Original P/N: 0.540.2491	W&H		W&H WA56/66 after 2007 / WK56/66 SKU: 83845 NOTE: 1 blue ring Original P/N: 05431100
Kavo		Kavo 67LH/ E15L/C / 0767LHC / 2967CHC SKU: 83827 Original P/N: 0.540.3671	W&H		W&H WE56A/LED G SKU: 83843 NOTE: rear gears Original P/N: 05844400
NSK		NSK Ti-Max X15/X25/ S-Max M15/25 SKU: 83898 NOTE: 16 teeth Original P/N: C601015	W&H		W&H WE56A/T / WE66E / LED G SKU: 83858 NOTE: 1 blue ring Original P/N: 05843300
NEW NSK		NSK Ti-Max Z25L SKU: 84087 Original P/N: C1038015	W&H		W&H WG56 A/LT / 66 A/TL SKU: 83857 NOTE: front gears Original P/N: 06994300
NSK		NSK Ti-Max Ti15L/25/Ti25 SKU: 84053 Original P/N: C875060	W&H		W&H WG56 A/LT after SN 014634 SKU: 83896 Original P/N: 07605900



1:5 Intermediate Shafts for Contra Angles

red 1:5

Bien Air		Bien Air CA 1:5 L SKU: 83877 Original P/N: 01500057	NSK		NSK Ti-Max X95L SKU: 83886 Original P/N: C600C022A
Kavo		Kavo 25LPA/M25L SKU: 83856 Original P/N: 1.001.3515	NSK		NSK S-Max M95L SKU: 83895 Original P/N: 10000559
Kavo		Kavo 25LH/25LN SKU: 83852 Original P/N: 0.553.7201	NEW NSK		NSK Ti-Max Z95L (old and new) SKU: 84070 Original P/N: C1034022002
Kavo		Kavo E25L/ E25C SKU: 83897 Original P/N: 1.007.9916	Sirona		Sirona A200/C200/S200 SKU: 83861 Original P/N: 6175603
Kavo		Kavo M05L Mini SKU: 84055 Original P/N: 1.011.3206	W&H		W&H WA99 A/LT after 2007 1 red ring / WG99 A/LT before SN 7801 SKU: 83847 Original P/N: 06420900
Kavo		Insert Drive Kavo M05L/ M25L/ 25LP SKU: 83875 Original P/N: 1.001.3519	W&H		W&H WK99 all models / WG99 A/LT after SN 7802 SKU: 83887 Original P/N: 06579200



1:1 Intermediate Shafts for Contra Angles


blue/green 1:1

Bien Air	NEW 	Bien Air 1:1 SKU: 84062 Original P/N: 1501083	NEW 	NSK S-Max M25L SKU: 84080 Original P/N: 10000574
Kavo		Kavo 68LH/68LD/E20L SKU: 83837 Original P/N: 0.540.1151	Sirona	 Sirona A40L/C40L; Ø9,7; Ø8,7 SKU: 84043 Original P/N: 6000959
Kavo		Kavo 67LH/E15L/0767/2967 SKU: 83825 Original P/N: 0.540.1161	W&H	 W&H WA56 after 2007 / WG56/ WK56 SKU: 83849 NOTE: 1 blue ring Original P/N: 05233200
Kavo		Kavo Intra Lux CL3 / S201 SKU: 84050 NOTE: 20:1 Original P/N: 1.003.4207	W&H	 W&H WG66 A/LT old SKU: 05273000 Original P/N: 05273000






Light Guides / Fiber Optics Rod

Kavo		Kavo 24LN/LS SKU: 553 6421 Original P/N: 0.553.6421	NSK		NSK Z SERIES SKU: 1034026 Original P/N: C1034026
Kavo		Kavo 25LH/LN/LCR/LCA/LHC SKU: 553 0283 Original P/N: 0.553.0283/ 0.553.7281	W&H		W&H WA56/62/66/67/86/99 / WG56/66/67/69/99 / WK56/66/86 SKU: 4058500 Original P/N: 04058500
Kavo		Kavo 25LP/LPA/LPR/ 200XD/XDR/ XDA SKU: 1001 3320 Original P/N: 1.001.3320	W&H		W&H WK99/93 / WG99 SN > 7802 SKU: 6584600 Original P/N: 06584600



Micro Motors Parts





Micro Motors Parts



Electrical Rotors and Sets

BIEN AIR / Electrical Rotors

Micro Motors Parts

	Electrical Rotor Motor MC2 IR SKU: 83807 Original P/N: 1500932		Electrical Rotor Motor MC3 LK/IR SKU: 83809 Original P/N: 211.75.29
	Electrical Rotor Motor MC2 40 LK SKU: 83808 Original P/N: 211.61.16, 1500931		

BIEN AIR / Bearing Sets

Micro Motors Parts

	Bearing Set for MC2 Motor SKU: 11018 Original P/N: 701.01.14, 701.01.23		Bearing Set for MC3 Motor SKU: 11019 Original P/N: 701.01.70, 701.01.71
--	--	--	--

BIEN AIR / Carbon Brushes

Micro Motors Parts

	MC2 Motor SKU: 93814 Original P/N: 211.45.06		MC3 Motor SKU: 93815T Original P/N: 211.75.08
---	---	---	--



Electrical Rotor Motor K186

SKU: **83802**
Original P/N: **0.676.0211**



Electrical Rotor Motor K193

SKU: **83805**
Original P/N: **0.676.0991**



Bearing Set for K9 Motor

SKU: **11015**
Original P/N: **0.220.0157, 0.220.0158**



Bearing Set for K10/K11/K12 Motor

SKU: **11016**
Original P/N: **0.220.0778, 0.220.0779**



*Bearing Set for K186/K190/K192/
K193/K196 Motor*

SKU: **11017**
Original P/N: **0.220.0089, 0.220.0136**



Bearing Set for K4/K5 Motor

SKU: **11020**
Original P/N: **0.220.0058, 0.220.0056**



Bearing Set for K12 attachment

SKU: **11021**
Original P/N: **1.000.3586, 0.220.0044**



Bearing Set for K200 Motor

SKU: **11023**
Original P/N: **0.220.0209, 0.220.7185**



K10 Motor

SKU: **93810**
Original P/N: **0.674.2481**




K9/K186/K187/K188/K189 Motor

SKU: **93811**
Original P/N: **0.676.0241**




*K190/K192/K193/K196/K197/K200/
K201 Motor*

SKU: **93812**
Original P/N: **0.676.0881**



Bearings





Bearings



Non-turbine Lowspeed Bearings

RPM from 40K to 100K

BIEN AIR / Bearings (RPM from 40K to 100K)



	3,175 × 7,938 × 3,571 mm <i>Smooth Closed/Closed</i> SKU: 701.01.14 Original P/N: 701.01.14		4 × 11 × 4 mm <i>Smooth Closed/Closed</i> SKU: 701.01.71 Original P/N: 701.01.71
	2,38 × 7,938 × 3,571 mm <i>Flanged Closed/Closed</i> SKU: 701.01.23 Original P/N: 701.01.23		2,35 × 4,763 × 2,35 mm <i>Smooth Closed/Closed</i> SKU: 701.01.75 Original P/N: 701.01.75
	3 × 10 × 4 mm <i>Smooth Closed/Closed</i> SKU: 701.01.65 Original P/N: 701.01.65		4,763 × 9,525 × 3,175 mm <i>Smooth Closed/Closed</i> SKU: 701.01.93 Original P/N: 701.01.93
	2,38 × 7,938 × 3,571 mm <i>Smooth Closed/Closed</i> SKU: 701.01.70 Original P/N: 701.01.70		3,967 × 7,938 × 3,175 mm <i>Smooth Closed/Closed</i> SKU: 701.01.94 Original P/N: 701.01.94

KAVO / Bearings (RPM from 40K to 100K)



	3 × 10 × 4 mm <i>Smooth Closed/Closed</i> SKU: 220 0005 Original P/N: 0.220.0005		5 × 8 × 2 mm <i>Smooth Closed/Closed</i> SKU: 220 0023 Original P/N: 0.220.0023
	2,38 × 4,763 × 1,588 mm <i>Smooth Closed/Open</i> SKU: 220 0012 Original P/N: 0.220.0012		7,859 × 15,898 × 5,98 mm <i>Smooth Closed/Closed</i> SKU: 220 0044 Original P/N: 0.220.0044
	3 × 7 × 2 mm <i>Smooth Open/Open</i> SKU: 220 0020 Original P/N: 0.220.0020		7,859 × 15,898 × 7,5 mm <i>Smooth Closed/Closed</i> SKU: 220 0045 Original P/N: 0.220.0045
	3 × 7 × 2 mm <i>Smooth Open/Open</i> SKU: 220 0021 Original P/N: 0.220.0021		3 × 7 × 2,5 / 3,35 mm <i>Ext. Inner Ring Closed/Closed</i> SKU: 220 0051 Original P/N: 0.220.0051



	7 × 14,1 × 7,5 mm Groove Closed/Closed SKU: 220 0055 Original P/N: 0.220.0055
	10 × 19 × 7 mm Smooth Closed/Closed SKU: 220 0056 Original P/N: 0.220.0056
	6,95 × 14,10 × 7,50 mm K4-K5 Set Groove Closed/Closed SKU: 220 0058 Original P/N: 0.220.0058
	4 × 13 × 5 mm Smooth Open/Open SKU: 220 0068 Original P/N: 0.220.0068
	4 × 13 × 5 mm Smooth Closed/Open SKU: 220 0069 Original P/N: 0.220.0069
	4 × 11 × 4 mm Smooth Closed/Closed SKU: 220 0089 Original P/N: 0.220.0089
	3 × 8 × 3 mm Smooth Closed/Closed SKU: 220 0097 Original P/N: 0.220.0097
	3 × 8 × 3 mm Smooth Closed/Closed SKU: 220 0098 Original P/N: 0.220.0098
	3 × 7 × 3 mm Smooth Closed/Closed SKU: 220 0102 Original P/N: 0.220.0102
	4 × 9 × 4 mm Smooth Closed/Closed SKU: 220 0103 Original P/N: 0.220.0103
	3 × 6 × 2 mm Smooth Closed/Closed SKU: 220 0104 Original P/N: 0.220.0104
	4 × 7 × 2 mm Smooth Open/Open SKU: 220 0105 Original P/N: 0.220.0105
	5 × 8 × 2 mm Smooth Open/Open SKU: 220 0106 Original P/N: 0.220.0106
	7 × 14 × 5 mm Smooth Closed/Closed SKU: 220 0112 Original P/N: 0.220.0112
	8 × 16 × 6 mm Smooth Closed/Closed SKU: 220 0113 Original P/N: 0.220.0113
	5 × 11 × 5 mm Smooth Closed/Closed SKU: 220 0126 Original P/N: 0.220.0126
	2 × 6 × 2,3 mm Smooth Closed/Open SKU: 220 0136 Original P/N: 0.220.0136
	4 × 11 × 4 mm Smooth Closed/Open SKU: 220 0138 Original P/N: 0.220.0138
	4 × 11 × 4 mm Smooth Open/Open SKU: 220 0139 Original P/N: 0.220.0139
	2 × 6 × 2,3 mm Smooth Closed/Closed SKU: 220 0144 Original P/N: 0.220.0144



	8 × 14 × 4 mm Smooth Closed/Closed SKU: 220 0148 Original P/N: 0.220.0148
	4 × 11 × 4 mm Smooth Closed/Closed SKU: 220 0157 Original P/N: 0.220.0157
	2 × 6 × 2,3 mm Smooth Closed/Closed SKU: 220 0158 Original P/N: 0.220.0158
	3 × 7 × 2 mm Smooth Open/Open SKU: 220 0161 Original P/N: 0.220.0161
	4 × 9 × 2,5 mm Smooth Open/Open SKU: 220 0171 Original P/N: 0.220.0171
	4 × 11 × 4 mm Smooth Closed/Open SKU: 220 0174 Original P/N: 0.220.0174
	3 × 7 × 2 mm Smooth Open/Open SKU: 220 0175 Original P/N: 0.220.0175
	3 × 8 × 3 mm Smooth Closed/Closed SKU: 220 0218 Original P/N: 0.220.0218
	5 × 11 × 5 mm Smooth Closed/Closed SKU: 220 0219 Original P/N: 0.220.0219
	5 × 11 × 4 mm Smooth Closed/Closed SKU: 220 0223 Original P/N: 0.220.0223
	3 × 7 × 3 mm Smooth Closed/Closed SKU: 220 0713 Original P/N: 0.220.0713
	2,35 × 5 × 1,5 mm Smooth Closed/Open SKU: 220 0730 Original P/N: 0.220.0730
	3 × 8 × 4 mm Smooth Closed/Closed SKU: 220 0778 Original P/N: 0.220.0778
	4 × 11 × 4 mm Smooth Closed/Closed SKU: 220 0779 Original P/N: 0.220.0779
	4 × 9 × 4 mm Smooth Closed/Closed SKU: 220 1301 Original P/N: 0.220.1301
	5 × 8 × 3 mm Smooth Closed/Closed SKU: 220 1302 Original P/N: 0.220.1302
	5 × 9,525 × 3,175 mm Smooth Closed/Closed SKU: 220 1306 Original P/N: 0.220.1306
	3 × 7 × 2,5 mm Smooth Closed/Closed SKU: 220 2004 Original P/N: 0.220.2004
	4 × 7 × 2 mm Smooth Closed/Closed SKU: 220 2005 Original P/N: 0.220.2005
	8 × 22 × 7 mm Smooth Closed/Closed SKU: 247 6056 Original P/N: 0.247.6056

KAVO / Bearings (RPM from 40K to 100K)



	8 × 22 × 7 mm Smooth Closed/Closed SKU: 247 6057 Original P/N: 0.247.6057		3 × 8 × 3 mm Smooth Closed/Closed SKU: 1003 5818 Original P/N: 1.003.5818
	K11 Set Smooth Closed/Closed SKU: 674 3841 Original P/N: 0.674.3841		2 × 5 × 1,5 mm Smooth Open/Open SKU: 1003 5383 Original P/N: 1.003.5383
	6,95 × 14,10 × 7,5 mm K9 Set Groove Closed/Closed SKU: 1000 3586 Original P/N: 1.000.3586		5 × 8 × 3 mm Smooth Closed/Closed SKU: 1003 6023 Original P/N: 1.003.6023
	5 × 8 × 3 mm Smooth Closed/Closed SKU: 1000 7185 Original P/N: 1.000.7185		6,95 × 14,10 × 7,5 mm 4944 Set Groove Closed/Closed SKU: 1004 0203 Original P/N: 1.004.0203
	2 × 5 × 1,5 mm Smooth Closed/Open SKU: 1001 2834 Original P/N: 1.001.2834		7 × 14 × 5 mm Smooth Closed/Closed SKU: 1004 0207 Original P/N: 1.004.0207

NSK / Bearings (RPM from 40K to 100K)



	4 × 7 × 2 mm Smooth Open/Open SKU: B10111-0704 Original P/N: B10111-0704		6 × 13 × 5 mm Smooth Closed/Closed SKU: B11223-1306 Original P/N: B11223-1306
	3 × 7 × 2 mm Smooth Closed/Open SKU: B10313-0703 Original P/N: B10313-0703		4 × 9 × 4 mm Smooth Closed/Closed SKU: B11253-0904 Original P/N: B11253-0904
	5 × 10 × 3 mm Smooth Closed/Open SKU: B10321-1005 Original P/N: B10321-1005		4 × 7 × 2,5 mm Smooth Closed/Closed SKU: B12213-0704 Original P/N: B12213-0704

NSK / Bearings (RPM from 40K to 100K)



	4 × 7 × 2,5 mm Smooth Closed/Closed SKU: B13313-0703 Original P/N: B13313-0703		2,38 × 7,938 × 3,571 mm Smooth Closed/Closed SKU: B22227-7923 Original P/N: B22227-7923
	2,38 × 4,763 × 1,588 mm Smooth Closed/Open SKU: B19157-4723 Original P/N: B19157-4723		5 × 10 × 4 mm Smooth Closed/Closed SKU: B11223-1005 Original P/N: B11223-1005
	6,35 × 12,7 × 4,763 mm Smooth Closed/Closed SKU: B20215-2763 Original P/N: B20215-2763		6 × 13 × 5 mm Smooth Closed/Closed SKU: B11243-1306 Original P/N: B11243-1306
	6,35 × 12,7 × 4,763 mm Smooth Closed/Closed SKU: B20245-2763 Original P/N: B20245-2763		

SHICK / Bearings (RPM from 40K to 100K)



	8 × 16 × 6 mm Smooth Closed/Closed SKU: 4120 Original P/N: 4120		10,15 × 19 × 7 mm Smooth Closed/Closed SKU: 7120 Original P/N: 7120
	4 × 9 × 4 mm Smooth Closed/Closed SKU: 5212 Original P/N: 5212		5 × 9,525 × 3,75 / 4 mm Ext. Inner Ring Closed/Closed SKU: 7212 Original P/N: 7212
	7 × 14 × 5 mm Smooth Closed/Closed SKU: 611/1+2 Original P/N: 611/1+2		6 × 13 × 5 mm Smooth Closed/Closed SKU: 9160 Original P/N: 9160
	8 × 15,875 × 6 mm Smooth Closed/Closed SKU: 7111 Original P/N: 7111		5 × 11 × 5 mm Smooth Closed/Closed SKU: 9162 Original P/N: 9162



	5 × 8 × 2 mm Smooth Closed/Closed SKU: 4181108 Original P/N: 4181108		3 × 6 × 2 mm Smooth Closed/Closed SKU: 6360593 Original P/N: 6360593
	3 × 8 × 3 mm Smooth Open/Open SKU: 4181132 Original P/N: 4181132		2 × 6 × 2,3 mm Flanged Closed/Closed SKU: 7036130 Original P/N: 7036130
	3 × 7 × 2 mm Smooth Open/Open SKU: 4181157 Original P/N: 4181157		4 × 9 × 4 mm Flanged Closed/Closed SKU: 7040801 Original P/N: 7040801
	3 × 10 × 4 mm Smooth Open/Open SKU: 4181199 Original P/N: 4181199		6 × 10,35 × 2,98 mm Smooth Closed/Open SKU: 7044688, 6360601 Original P/N: 7044688, 6360601
	4 × 7 × 1,6 mm Smooth Open/Open SKU: 5887174 Original P/N: 5887174		



	3 × 8 × 4 mm Smooth Closed/Closed SKU: 01206200 Original P/N: 01206200		3,175 × 7,983 × 3,571 mm Smooth Closed/Closed SKU: 01850600, 04096800 Original P/N: 01850600, 04096800
	5 × 8 × 2 mm Smooth Open/Open SKU: 01381900 Original P/N: 01381900		5 × 11 × 5 mm Smooth Closed/Closed SKU: 02685500 Original P/N: 02685500
	5 × 8 × 2 mm Smooth Closed/Open SKU: 01565800 Original P/N: 01565800		7 × 14 × 5 mm Smooth Closed/Closed SKU: 03202200 Original P/N: 03202200
	3 × 10 × 4 mm Smooth Closed/Closed SKU: 01577000 Original P/N: 01577000		8 × 16 × 6 mm Smooth Closed/Closed SKU: 03202300 Original P/N: 03202300
	2,35 × 5,5 × 2 mm Smooth Open/Open SKU: 01579300, 04733900 Original P/N: 01579300, 04733900		4 × 7 × 2 mm Smooth Closed/Open SKU: 04027300 Original P/N: 04027300
	2,35 × 5,5 × 2 mm Smooth Closed/Open SKU: 01579600 Original P/N: 01579600		3 × 6 × 2 mm Smooth Closed/Open SKU: 04033500, 05765800 Original P/N: 04033500, 05765800
	3 × 7 × 2 mm Smooth Open/Open SKU: 01589400 Original P/N: 01589400		2,35 × 5,5 × 2 mm Smooth Open/Open SKU: 04733900 Original P/N: 04733900
	4 × 7 × 2 mm Smooth Open/Open SKU: 01589600 Original P/N: 01589600		4 × 7 × 1,6 mm Smooth Open/Open SKU: 05039400 Original P/N: 05039400
	5 × 8 × 3 mm Smooth Closed/Closed SKU: 01672300, 05294400 Original P/N: 01672300, 05294400		



Bearings



Non-turbine Highspeed Bearings

RPM from 150K to 250K

BIEN AIR / Bearings (RPM from 150K to 250K)



4 × 7 × 2 mm

Smooth **Closed/Closed**

SKU: **701.01.68**
Original P/N: **701.01.68**

KAVO / Bearings (RPM from 150K to 250K)



4 × 7 × 2 mm

Smooth **Closed/Closed**

SKU: **220 0734**
Original P/N: **0.220.0734**



4 × 7 / 1,6 × 2 mm

Smooth **Closed/Open**

SKU: **1008 4039**
Original P/N: **1.008.4039**



4 × 7 × 2 mm

Smooth **Closed/Closed**

SKU: **1001 2836**
Original P/N: **1.001.2836**



3,5 × 6,5 × 2 mm

Smooth **Closed/Closed**

SKU: **1011 0127**
Original P/N: **1.010.9767**



3,175 × 6,350 × 2,779 / 1,910 mm

Extended Inner Ring **Closed/Open**

SKU: **1001 3167**
Original P/N: **1.001.3167**



3,175 × 5,8 × 2,28/1,42 mm

Ext. Inner Ring **Closed/Closed**

SKU: **1011 0128**
Original P/N: **1.010.9766**



3,175 × 6,350 × 2,779 / 1,910 mm

Extended Inner Ring **Closed/Open**

SKU: **1008 4038**
Original P/N: **1.008.4038**

NSK / Bearings (RPM from 150K to 250K)



3,8 × 6,6 × 1,74 mm
Smooth **Closed/Closed**
SKU: **1034015B**
Original P/N: **NA**

SIRONA / Bearings (RPM from 150K to 250K)



3,175 × 6 × 2 mm
Smooth **Closed/Open**
SKU: **5899863**
Original P/N: **5899863**



3,175 × 7 × 2 mm
Smooth **Closed/Open**
SKU: **5890855, 4181173**
Original P/N: **5890855, 4181173**

W&H / Bearings (RPM from 150K to 250K)



4 × 7 × 2 mm
Smooth **Closed/Open**
SKU: **02649400R**
Original P/N: **02649400**



3,4 × 6,5 × 1,6 mm
Smooth **Closed/Open**
SKU: **05257200**
Original P/N: **05257200**



4 × 7 × 2 mm
Smooth **Closed/Closed**
SKU: **04025900**
Original P/N: **04025900**



3,4 × 6,5 × 1,6 mm
Bearing with a Cup
SKU: **05798700**
Original P/N: **05798700**



4 × 7 × 2 mm
Smooth **Closed/Open**
SKU: **04486000R**
Original P/N: **04486000**

Bearings



Turbine Bearings
RPM from 250K to 450K



	3,175 × 7,938 × 3,571 mm <i>Ceramic Radial Smooth Torlon</i>	SKU: 10110PT myonic		3,175 × 6,350 × 2,381 / 2,779 mm <i>Ceramic Angular Smooth Torlon</i>	SKU: 1032M myonic
	3,175 × 6,350 × 2,779 mm <i>Ceramic Radial Flanged Torlon</i>	SKU: 10210PT myonic		3,175 × 6,350 × 4,159 / 2,779 mm <i>Ceramic Angular Stepped Torlon</i>	SKU: 1034M myonic
	3,175 × 6,350 × 2,779 mm <i>Ceramic Radial Stepped Torlon</i>	SKU: 10310PT myonic		3,175 × 6,350 × 2,779 mm <i>Ceramic Angular Stepped Torlon</i>	SKU: 1035M myonic
	3,175 × 6,350 × 2,381 / 2,779 mm <i>Ceramic Radial Smooth Torlon</i>	SKU: 10320PT myonic		3,175 × 6,450 × 2,381 / 2,779 mm <i>Ceramic Angular Stepped Torlon</i>	SKU: 1036M myonic
	3,175 × 6,350 × 2,381 / 2,779 mm <i>Ceramic Radial Stepped Torlon</i>	SKU: 10330PT myonic		3,175 × 6,350 × 2,779 mm <i>Ceramic Angular Smooth Torlon</i>	SKU: 1084M myonic
	3,175 × 6,350 × 2,380 mm <i>Ceramic Radial Stepped Torlon</i>	SKU: 10610PT myonic		3,175 × 6,350 × 2,779 mm <i>Ceramic Angular Smooth Integral Shield Torlon</i>	SKU: 1086M myonic
	3,175 × 6,350 × 2,779 mm <i>Ceramic Radial Flanged/Unique Torlon</i>	SKU: 10710PT myonic		3,175 × 6,350 × 2,779 mm <i>Ceramic Radial Stepped Double Flanged Torlon</i>	SKU: 1131M
	3,175 × 6,350 × 2,380 mm <i>Ceramic Radial Smooth Torlon</i>	SKU: 10810PT myonic		3,175 × 6,350 × 2,380 mm <i>Ceramic Angular Smooth Torlon</i>	SKU: 1087N
	3,175 × 6,350 × 2,779 mm <i>Ceramic Radial Smooth Torlon</i>	SKU: 10830PT myonic			



www.vkl.biz

Visit our website – your one stop shop for dental replacement parts

- On the website you can find
- new products available
 - on-line shop
 - info about the first order discount
 - special order discounts
 - regular discounts
 - info about who we are
 - info about warranty
 - info about shipping policy
 - info about returns & exchanges
 - useful tips and advice
 - digital catalogue



www.vkl.biz

V&KL dental rotor d.o.o.
Ulica gledališča BTC 6
SI-1000 Ljubljana, Slovenia
tel: +386 1 361 61 57
info@vkl.biz
www.vkl.biz

